Research Article

Genome-Wide Linkage and Association Scans for Quantitative Trait Loci of Serum Lactate Dehydrogenase—The Framingham Heart Study

Jing-Ping Lin,¹ Gang Zheng,¹ Jungnam Joo,¹ and L. Adrienne Cupples²

¹ Office of Biostatistics Research, Division of Cardiovascular Sciences, National Heart, Lung and Blood Institute, NIH 6701 Rockledge Dr. Suite 9196, Bethesda, MD 20892-7913, USA

² Department of Biostatistics, Boston University School of Public Health, MA 02118, USA

Correspondence should be addressed to Jing-Ping Lin, linj@nhlbi.nih.gov

Received 29 January 2010; Revised 11 May 2010; Accepted 18 July 2010

Academic Editor: Federico Innocenti

Copyright © 2010 Jing-Ping Lin et al. This is an open access article distributed under the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Serum lactate dehydrogenase (LDH) is used in diagnosing many diseases and is significantly determined by genetic factors. Three genes coding for LDH isoenzymes were mapped to chromosome 11q15 and 12p12. We used 330 Framingham Heart Study largest families for microsatellite linkage scan and 100K SNPs association scan to determine quantitative trait loci of LDH level. We estimated the heritability at 41%. Our genome-wide linkage analysis yielded several chromosomal regions, other than 11q and 12p, with LOD scores between 1 and 2.5. None of the 100K SNPs with a *P*-value < 10^{-4} in our genome-wide association study was close to the chromosomal regions where the LDH genes reside. Our study demonstrated a strong genetic effect on the variation of LDH levels. There may not be a single gene with a large effect, instead may be several genes with small effects in controlling the variation of serum LDH. Those genes may be located on chromosomal regions that differ from where the genes encoding LDH isoenzymes reside.

1. Introduction

LDH is a cytoplasmatic enzyme found in almost all body tissues. A small amount of the enzyme exists in the blood. When cells are damaged or destroyed, they release LDH into the bloodstream and cause blood levels to rise. LDH is used in diagnosing heart, liver, muscle, lung diseases, and testicular cancer. Its sensitivity and specificity in helping to diagnose acute myocardial infarction are very high [1].

The enzyme is composed of 4 peptide chains of 2 types, A and B. A is encoded by LDH A gene and predominant in muscle, and B by LDH B gene predominant in the heart. Five different combinations of A and B compose five LDH isoenzymes: LDH1 (A4), LDH2 (A3B), LDH3 (A2B2), LDH4 (AB3), and LDH5 (B4). While all of the five isoenzymes are represented in total serum LDH; LDH2 makes up the greatest percentage. In addition, there is homotetrameric LDH C isoenzyme, encoded by the LDH C gene, found only in mature testes and spermatozoa in humans. The LDH A and C genes are syntenic on human chromosome 11q15 [2, 3] while the LDH B gene maps to chromosome 12p12 [4].

Family studies have shown that serum LDH variation is significantly determined by genetic factors with heritability estimate as 50% [5]. So far, no genome-wide linkage or association studies on serum levels of LDH have been reported. We carried out a 10 cm genome-wide linkage analysis and an Affymetrix Human GeneChip 100 K genomewide association study for quantitative trait loci of serum LDH level in a community-based Caucasian cohort, the Framingham Heart Study.

2. Subjects and Methods

The Framingham Heart Study, a population-based study, began in 1948 with the recruitment of 5209 residents aged 28–62 years (mean age 44.1) from Framingham, Massachusetts [6]. The participants have undergone biennial

examinations since the study began. In 1971, the Framingham Offspring Study [7] was started, in part, to evaluate the genetic components of cardiovascular disease etiology. In total, 5124 subjects aged 5-70 years (mean age 36.3) including the offspring of the original cohort and the spouses of the offspring were recruited. The offspring cohort has been examined every four years (except the first two examinations with eight years intervening). Within the study, the 330 largest extended families were selected and genotyped with a 10 cm density microsatellite marker and the Affymetrix Human GeneChip 100 K genome-wide scans. The number of subjects with microsatellite genotyped was 1702 and with the 100 K SNPs genotyped was 1343. A recent study reported no evidence of major population substructure in the Framingham Heart Study [8]. LDH was measured during the first (1971-1975) and second (1980-1983) examinations of the offspring. Our analysis was limited to the first examination, since far fewer individuals have serum LDH measured in the second examination.

Total LDH was measured by kinetic method using Dow UV LDH-Reagent Kit. Weight was measured with the subject in light clothing and shoes off. Laboratory measurements were made on 12 hours fasting venous blood samples that were collected in tubes containing 0.1% EDTA. Lipid determinations were performed at the Framingham Heart Study laboratory, which participates in the Standardization Program of the Centers for Disease Control. All subjects provided informed consent prior to each clinic visit, and the examination protocol was approved by the Institutional Review Board at Boston Medical Center (Boston, MA). The clinical and laboratory methods have been detailed elsewhere [5].

2.1. Genotyping Methods. For the microsatellite genotyping, genomic DNA was isolated from nucleated blood cells. DNA samples were sent to the Marshfield Mammalian Genotyping Service (http://research.marshfieldclinic.org/genetics/). At an average of 10 cm density, 399 microsatellite markers (Screening Set 9) [9] covered the genome with an average marker heterozygosity of 0.77. The genotyping data were cleaned in two steps. First, the sib-kin program in Aspex (ftp://lahmed.standord.edu/pub/aspex/index.html) was used to verify family relationships based on all markers available. Second, the GENTEST program, as a precursor of INFER, created by Southwest Foundation for Biomedical Research (http://www.sfbr.org/sfbr/public/software/software.html), was used to identify and eliminate additional genotype

inconsistencies. When inconsistencies were found, the genotyping values in all members of the nuclear family were set to missing.

For the Affymetrix 100 K SNPs genotyping, the details are provided in [10]. In summary, SNPs on the Affymetrix Human GeneChip 100 K (n = 112,990 autosomal SNPs) were genotyped in a sample of 1343 individuals from 330 families. SNPs with minor allele frequency <10% or call rate <80% or Hardy-Weinberg equilibrium *P*-value <.001 were excluded, leaving 70591 SNPs for analysis, criteria set by the Framingham Heart Study for the Affymetrix 100 K SNPs association genome-wide scans [11]. The 100 K SNPs data are publicly available through dbGaP (http://www.ncbi.nlm.nih.gov/gap) under Framingham SHARe project.

2.2. Statistical Methods

2.2.1. Genome-Wide Linkage Analysis. Variation in LDH from known factors was identified and removed by regression modeling incorporated in SOLAR [12–14], to enhance the ability of linkage analysis to detect genetically determined variation using a maximum-likelihood-based variance decomposition method. The covariates selected (P < .05) and incorporated into both the heritability estimation and the linkage analyses were age, sex, height, weight, high density lipoprotein cholesterol (HDL-C), and alkaline phosphatase (ALP).

An estimate of heritability was obtained using the variance-component method. Heritability is the proportion of total phenotypic variation due to additive genetic effects, after removing the variation attributable to covariates. The variance component method analysis was also used for the linkage analysis of LDH, adjusted for known covariates, using the random microsatellite DNA markers covering the entire genome. Marker allele frequencies were estimated from the study participants and then used to estimate the proportion of a single-point alleles shared identical by descent (IBD) among all relative pairs. Multipoint IBDs were estimated based on the singlepoint estimates by a regression approach [13]. A likelihood ratio test was used to evaluate linkage by comparing a purely polygenic model (without consideration of genetic marker information) to a model that incorporates IBD information at the marker. The LOD score was the log (base 10) of the ratio of the likelihoods of two models, purely polygenic versus one that also included IBD information at the marker.

Since the variance-component method is based on the assumption of a multivariate normal distribution, violations of this assumption may result in inaccurate results [14–16]. We found that LDH had high kurtosis and thus used an LOD score adjustment method implemented in SOLAR to ensure more reliable results. A fully informative marker linked to the trait studied was simulated. The IBD information for this marker was calculated, and then linkage analysis to the trait was performed. The LOD score adjustment process [15, 16] regresses the observed LOD scores from simulation (10,000 replicates) on the expected LOD scores for a multivariate normal trait to obtain a correction factor for the LOD scores from the analyses of the observed trait. This method provides robust LOD scores for data with nonnormal distributions.

2.2.2. Genome-Wide Association Study. One individual from each family was randomly selected for the Hardy-Weinberg equilibrium (HWE) test using a chi-squared statistic with one degree of freedom. For each SNP, we modeled the logtransformed trait value adjusted for the same covariates as used in the linkage analysis. Trait values within each pedigree were assumed to be correlated, while they were independent between pedigrees. A special case of linear mixed effects

Variable	Mean \pm SD	Range		
LDH (mg/dl)	155.9 ± 34.7	55.0-600.0		
Age (years)	32.7 ± 10.5	10–64		
ALP* (mg/dl)	58.3 ± 17.8	16.0–98.0		
Height (inchs)	66.2 ± 4.0	52.3-77.5		
Weight (pound)	156.6 ± 34.6	54.0-281.0		
HDL-C [#] (mg/dl)	50.4 ± 14.0	16.0-118.0		

TABLE 1: Characteristics of the 1603 individuals with LDH and covariates used in the linkage analysis.

*alkaline phosphatase.

#high density lipoprotein cholesterol.

model (LME) in SAS was used to fit the model with the covariates and calculate the residual for each individual. The following model was used for family i, $Y = X\beta + \epsilon$, where Yis a vector of log-transformed LDHs of subjects in the family, and X contains covariates to adjust (age, sex, height, weight, HDL-C, and ALP), which were treated as fixed effects. The covariance matrix for Y in family i was unstructured. Once the estimate of β was obtained, denoted as $\hat{\beta}$, we calculated the residual for person j in family i with covariates X_{ij} as $\epsilon_{ii}^* = Y_{ii} - X_{ii} \hat{\beta}$. Then, in the second stage, we modeled the association between the residual and genetic effect by LME given by $\epsilon^* = \mu + \beta c(G) + \epsilon$, where the covariance matrix of ϵ^* was unstructured within each family and diagonal between families, and c(G) was a coding for the genetic models underlying the disease locus, which were defined as recessive, additive, or dominant. For each genetic model, a P-value was obtained. The smallest of the three P-values was used to rank all SNPs. In genome-wide association studies, some SNPs associated with common diseases have little additive effect while there is strong deviation from Hardy-Weinberg equilibrium, which indicates strong recessive or dominant effects [17, 18]. These SNPs would not be detected using the P-value only based on the additive model. Using the minimum P-value to rank SNPs will be able to detect SNPs with all three genetic models. Since the minimum P-value is no longer a true p-value, we used the following threshold values: if the minimum *P*-value is less than 5×10^{-7} , the association is regarded as strong, while if the minimum Pvalue is between 5×10^{-7} and 1×10^{-5} , the association is moderate. No correction for multiple testing was given for the minimum *P*-value.

3. Results

3.1. Genome-Wide Linkage Analysis. The total number of individuals with measured LDH and with all covariates used for the heritability estimates and linkage analysis in offspring was 1603. The basic characteristics of the clinical covariates of these individuals (about 50% male) are displayed in Table 1. Of the 1603 individuals, 1260 had microsatellite marker genotypes.

The skewness and kurtosis of LDH was 1.87 ± 0.04 and 13.92 ± 0.07 , respectively. The heritability estimate for LDH, after adjusting for the covariates, was $40.9 \pm 6.4\%$. The

TABLE 2: Chromosomal regions in the genome scan with multipoint LODs ≥ 1.0 .

Chromosome	cm (Mbp)	LOD score		
3	55 (29)	1.03		
4	195 (186)	2.01		
7	51 (33)	1.08		
8	37 (18)	2.25		
9	164 (137)	2.41		
19	91 (59)	1.24		
20	96 (57)	1.28		

proportion of variation due to all covariates included in the model was approximately 5.8%.

From the multipoint linkage analysis of LDH, several maximum LOD scores between 1 and 2.5 were observed on different chromosomes (Table 2). The highest LOD score was 2.41 on chromosome 9.

3.2. Genome-Wide Association Study. After log-transformed, only one individual's LDH value, 1.74, was out of the range of 5SD from the mean, 1.77–2.64. Since the value of 1.74 was so close to the lower bound of 1.77, we did not exclude this individual from the association study and reanalyzed the data. In the genome-wide association study, no SNP showed strong association, while there were two SNPs, rs9318892 and rs10518949, which showed moderate association. The top ranked SNP on chromosome 13 had a minimum P-value of 7.1×10^{-7} , while the second one on chromosome 15 had a mini p = 1.6×10^{-7} . There were 32 SNPs having a minimum $P < 10^{-4}$, Table 3. All SNPs within 60 kb of LDHA gene (rs4237721, rs2279900, rs8018, rs3781640) and LDHB gene (rs10505873, rs1677106, rs1677104, rs1030452, rs1012354) have a *P*-value > .20.

4. Discussion

In our study, the heritability of LDH was estimated at 41% indicating that a substantial portion of the variation in serum LDH was attributable to additive genetic factors. This is consistent with a previous finding [5]. In the linkage study, we identified several chromosomal regions with LOD scores between 1 and 2.5. None of them are located on the chromosomes where the genes encoding the LDH isoenzymes reside. No obvious candidate genes were found in those chromosomal regions. In the genome-wide association study, no SNP with genome-wide significance was observed. Furthermore, testing three models per SNP in the genome-wise association study adds burden to the multiple test and this is not reflected in the P-values in Table 3. Among the top 32 SNPs with a *P*-values $<10^{-4}$, none of them were near the LDH structure genes, and none were within the regions of linkage peaks. Similarly as in the linkage analysis, no obvious candidate genes were found on those chromosomal regions.

TABLE 3: Top 32 SNPs associated with LDH levels with a minimum P-value < 10^{-4} by LME test.

SNP Chromosome Physical position cm P recessive P additive P dominant Minume P rs9318892 13 82117408 64.375 7.08385E-07 4.61436E-06 0.0159389 7.0844E-07 rs10518949 15 69415390 71.581 0.134216811 9.45936E-06 1.557E-06 1.557F-06 rs727614 15 51032409 46.056 0.00013203 1.28881E-05 0.00141981 0.085699 2.0672E-05 rs1010154 13 29188260 18.909 2.06735E-05 0.00141981 0.085699 2.047E-05 rs716039 5 158083224 161.32 2.48992E-05 0.00172614 0.0504592 2.949E-05 rs190x451 10 5347766 7.028 0.08510484 7.70997E-05 2.999E-05 3.158E-05 rs190x525 18 24919668 54.364 0.014370952 0.001911018 3.286E-05 3.286E-05 rs1803037 13 13422380 48.33 0.05721322 0.00198124 4.393E-05 <								
rs9318892 13 82117408 64.375 7.08385E-07 4.61436E-06 0.0159389 7.084E-07 rs10518949 15 69415390 7.1581 0.134216811 9.45936E-06 1.557E-06 1.557E-06 rs727614 15 51032409 46.056 0.000132039 1.2881E-05 0.000433 1.288E-05 rs1010154 13 29188260 18.909 2.06735E-05 0.0114981 0.0885699 2.067E-05 rs7716399 5 158083224 161.32 2.48992E-05 0.001732614 0.0504592 2.499E-05 rs1700451 10 55347706 75.028 0.085104884 7.7097E-05 2.999E-05 2.999E-05 rs1722248 12 2225840 41.55 3.05871E-05 0.00051173 0.1632718 3.058E-05 rs156305 10 81162789 10.389 0.036085304 3.15764E-05 6.928E-05 3.286E-05 rs15780402 2 117190012 12.61.13 3.3709E-05 7.68122E-05 0.01918111 3.377E-05 rs2041980 17 63989613 89.283 4.36211E-05 <	SNP	Chromosome	Physical position	cm	P recessive	P additive	P dominant	Minume P
rs10518949156941539071.5810.1342168119.45936E-061.557E-061.557E-06rs727614155103240946.0560.0001320391.28881E-050.00040331.2889-05rs1953433143251544934.6381.65253E-050.0011498910.00851092.0677E-05rs101015413219825018.9092.06735E-050.0011498910.08855092.0677E-05rs7163995158083224161.322.48992E-050.0017326140.05045922.499E-05rs1900451105534770675.0280.0851048847.70997E-052.999E-052.999E-05rs15730501083162789103.3890.0360853043.1576E-066.928E-053.158E-05rs157305174134827467131.8510.298266750.0005710183.286E-053.458E-05rs16395174134829467131.8510.2022766190.0231164.362E-05rs845813112162238043.6211E-050.0006726863.458E-053.458E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs9306664134844616131.864.8235E-050.0001818160.36573474.882E-05rs93076664134844616131.860.829192640.000892125.233E-055.306E-05rs403803543201269.410.919473444.990754E-050.00118160.36673606.412E-05rs403838<	rs9318892	13	82117408	64.375	7.08385 <i>E</i> –07	4.61436 <i>E</i> 06	0.0159389	7.084 <i>E</i> -07
rs727614155103240946.0560.0001320391.28881E-050.00040331.289E-05rs1953433143251544934.6381.65233E-050.011049810.10251131.653E-05rs101014132918826018.9092.06735E-050.0017326140.008850922.049E-05rs100451105534770675.0280.0051048447.7097E-052.999E-052.999E-05rs120248122225834041.553.05871E-050.0056157330.16327183.059E-05rs15763051083162789103.3890.0360853043.15764E-056.928E-053.158E-05rs1505255182491966854.3640.0143709520.0006726863.458E-053.458E-05rs758402821719012126.1133.3770-057.68122E-050.01981113.377E-05rs14935174134829467131.8510.2982666750.0006726863.458E-053.458E-05rs803037915823819147.894.6323E-050.0011070.1095424.362E-05rs803037915823819147.894.53395E-050.0001101670.1095424.539E-05rs1050466387903123494.610.012473464.99754E-050.00180114.998E-05rs1050466387903123494.610.012473464.99754E-050.00018014.998E-05rs136303543201269.410.910493010.009816065.306E-055.306E-05rs136303<	rs10518949	15	69415390	71.581	0.134216811	9.45936 <i>E</i> 06	1.557 <i>E</i> -06	1.557 <i>E</i> -06
rs1953433143251544934.6381.65253E-050.0110498910.10251131.653E-05rs1010154132918826018.9092.06735E-050.0011419810.08856992.067E-05rs77163995158083224161.222.4892E-050.0017326140.00545922.49E-05rs702045110534770675.0280.085104847.70997E-052.99E-052.99E-05rs7522248122225834041.553.05871E-050.0056157330.16327183.059E-05rs150502525182491966854.3640.0143709520.0005110183.286E-053.286E-05rs75840282117190012126.1133.3770E-057.68122E-050.01981113.377E-05rs4935174134829467131.8510.0292666750.000672683.458E-053.458E-05rs403037915823819178.924.5895E-050.00116070.01995424.539E-05rs803037915823819178.924.5895E-050.000116070.01995424.539E-05rs105046638793123474.824.5825E-050.0001818160.36573474.882E-05rs10504663137829536138.6850.8829192640.000891225.233E-055.306E-05rs1050463137829536138.6850.8829192640.000816965.306E-055.306E-05rs1050463137829536138.6850.8829192640.000780737.489E-050.0001816rs2362620141348	rs727614	15	51032409	46.056	0.000132039	1.28881 <i>E</i> -05	0.0004033	1.289 <i>E</i> -05
rs1010154 13 29188260 18.909 2.06735E-05 0.001041981 0.0885699 2.067E-05 rs7716399 5 158083224 161.32 2.48992E-05 0.001732614 0.0504592 2.999E-05 rs1700451 10 55347706 75.028 0.085104884 7.7097E-05 2.999E-05 2.999E-05 rs722248 12 22258340 41.55 3.05871E-05 0.005615733 0.1632718 3.058716 rs1576305 10 83162789 103.389 0.036085304 3.15764E-05 6.928E-05 3.158E-05 rs15784028 2 117190012 126.113 3.37709E-05 7.68122E-05 0.0198111 3.377E-05 rs840379 17 63989613 89.283 4.3621E-05 0.000195121 3.493E-05 3.493E-05 rs803079 15 82381914 78.92 4.53895E-05 0.00011607 0.019542 4.539E-05 rs9307696 4 134844616 131.86 4.88235E-05 0.0001801 4.998E-05 rs93076	rs1953433	14	32515449	34.638	1.65253 <i>E</i> -05	0.011049891	0.1025113	1.653 <i>E</i> -05
rs77163995158083224161.322.48992E-050.0017326140.05045922.49E-05rs1900451105534770675.0280.0851048847.70997E-052.999E-052.999E-05rs722248122225834041.553.05871E-050.0056157330.16327183.059E-05rs15763051083162789103.3890.0360853043.15764E-056.928E-053.158E-05rs75840282117190012126.1133.37709E-057.68122E-050.01981113.377E-05rs7840282117190012126.1133.37709E-057.68122E-050.01981113.377E-05rs8581312162238048.530.0572213220.0001951213.493E-053.493E-05rs8030379158238191478.924.36211E-050.0027976190.20231164.362E-05rs93076964134844616131.864.88235E-050.00116070.01995244.539E-05rs93076964134829546131.860.882192640.0008921225.233E-055.306E-05rs93076964134829546131.850.882192640.0008921225.233E-055.306E-05rs15054663879012349.410.919493010.0009816965.306E-055.306E-05rs23626916.1722349.41.80.91243460.01246360.06736096.472E-05rs1050406311.642159813.510.000367937.08984E-050.00045627.09E-05rs105090	rs1010154	13	29188260	18.909	2.06735 <i>E</i> -05	0.001041981	0.0885699	2.067 <i>E</i> -05
rs1900451105534770675.0280.0851048847.70997E-052.999E-052.999E-05rs722248122225834041.553.05871E-050.0056157330.16327183.059F-05rs15763051083162789103.3890.0360853043.15764E-056.928E-053.158E-05rs75840282117190012126.1133.37709E-057.68122E-050.01981113.377E-05rs14935174134829467131.8510.2928666750.0006726863.458E-053.493E-05rs8030379158238191478.924.53895E-050.000116070.019951214.362E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs9306664134849616131.860.8829192640.000891120.36573474.882E-05rs9307696413482956138.650.8829192640.000891225.233E-055.306E-05rs930666387903123494.610.012473464.99754E-050.00018014.998E-05rs93216116137829536138.650.8829192640.000892125.233E-055.306E-05rs4238934134829546131.8516.18458E-050.00018014.499E-05rs4238316.177223494.186.47241E-050.000126160.60746096.472E-05rs4238316.167223494.186.47241E-050.000126160.0045627.09E-05rs458331	rs7716399	5	158083224	161.32	2.48992 <i>E</i> -05	0.001732614	0.0504592	2.49 <i>E</i> -05
rs722248122225834041.553.05871E-050.0056157330.16327183.059E-05rs15763051083162789103.3890.0360853043.15764E-056.928E-053.158E-05rs10502525182491966854.3640.0143709520.0005110183.286E-053.286E-05rs7584028211719012126.1133.37709E-057.68122E-050.01981113.377E-05rs14935174134829467131.8510.2982666750.0006726863.458E-053.458E-05rs2041980176398961389.2834.36211E-050.0279766190.20231164.362E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs903066634134844616131.864.88235E-050.001181160.36573474.882E-05rs403037543201269.410.0124734464.99754E-050.00018014.998E-05rs403803543201269.410.9190493010.0008911665.306E-055.336E-05rs1389384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs10521232171348052932.756.91427E-050.00124630.06756096.472E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs10509003105543762875.1450.001369127.09645E-050.00038547.096E-05 <td< td=""><td>rs1900451</td><td>10</td><td>55347706</td><td>75.028</td><td>0.085104884</td><td>7.70997<i>E</i>–05</td><td>2.999<i>E</i>-05</td><td>2.999<i>E</i>-05</td></td<>	rs1900451	10	55347706	75.028	0.085104884	7.70997 <i>E</i> –05	2.999 <i>E</i> -05	2.999 <i>E</i> -05
rs15763051083162789103.3890.0360853043.15764E-056.928E-053.158E-05rs10502525182491966854.3640.0143709520.0005110183.286E-053.286E-05rs75840282117190012126.1133.37709E-057.68122E-050.01981113.377E-05rs14935174134829467131.8510.298266750.0001510183.438E-053.458E-05rs85851312162238048.530.0572213220.0001951213.493E-053.439E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs8030379158238191478.924.53895E-050.000118160.36573474.882E-05rs93076964134844616131.864.88235E-050.00118160.36573474.882E-05rs93076964137829536138.6850.8829192640.0008921225.233E-055.233E-05rs93076964134829467131.8516.18438E-050.0023860730.49995846.185E-05rs93076964134829546131.8516.18438E-050.002860730.49995846.185E-05rs93076964134829546131.8516.18438E-050.002860730.49995846.185E-05rs1389384134829546131.8516.18438E-050.002860730.49995846.185E-05rs139232171348052932.756.7421E-050.00101012510.03494116.914E-05<	rs722248	12	22258340	41.55	3.05871 <i>E</i> -05	0.005615733	0.1632718	3.059 <i>E</i> -05
rs10502525182491966854.3640.0143709520.0005110183.286E-053.286E-05rs75840282117190012126.1133.37709E-057.68122E-050.01981113.377E-05rs14935174134829467131.8510.2982666750.0006726863.458E-053.458E-05rs8581312162238048.530.0572213220.0001951213.493E-053.493E-05rs2041980176398961389.2834.36211E-050.0279766190.20231164.362E-05rs8030379158238191478.924.53895E-050.00116070.01995424.539E-05rs93076964134844616131.864.88235E-050.00118160.36573474.882E-05rs93076964137829536138.6850.8821912640.0008921225.233E-055.233E-05rs93216116137829536138.6850.8821912640.0008921225.233E-055.306E-05rs43803543201269.410.9190493010.0009816965.306E-055.306E-05rs7230422334535351.3180.3350523640.0175672776.41E-056.41E-05rs25265916177223494.1866.47241E-050.0001012510.03494116.914E-05rs205203105543762875.1450.0003679037.8984E-050.00045627.09E-05rs205903105543762875.1450.001861217.09645E-050.00035847.09E-05rs2059003	rs1576305	10	83162789	103.389	0.036085304	3.15764 <i>E</i> 05	6.928 <i>E</i> –05	3.158 <i>E</i> -05
rs75840282117190012126.1133.37709E-057.68122E-050.01981113.377E-05rs14935174134829467131.8510.2982666750.0006726863.458E-053.458E-05rs8581312162238048.530.0572213220.0001951213.493E-053.493E-05rs2041980176398961389.2834.36211E-050.0279766190.20231164.362E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs93076964134844616131.864.88235E-050.0001818160.36573474.882E-05rs93076964137829536138.6850.8829192464.99754E-050.00018104.998E-05rs93216116137829536138.6850.8829192460.0008921225.233E-055.306E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs453834134829546131.8516.18458E-050.001246360.06736096.472E-05rs236526916177223494.1866.47241E-050.0012046360.06736096.472E-05rs10509003105543762875.1450.001367037.08984E-050.00085847.09E-05rs125244461204640851223.6970.02912397.371E-050.00013977.372E-05rs134082218195427107.560.2336071020.0008019417.815E-05rs13406821	rs10502525	18	24919668	54.364	0.014370952	0.000511018	3.286 <i>E</i> -05	3.286 <i>E</i> -05
rs14935174134829467131.8510.2982666750.0006726863.458E-053.458E-05rs88581312162238048.530.0572213220.0001951213.493E-053.493E-05rs2041980176398961389.2834.36211E-050.0279766190.20231164.362E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs93076964134844616131.864.88235E-050.0001818160.36573474.882E-05rs1050466387903123494.610.0124734464.99754E-050.00018014.998E-05rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs62562916177223494.186.47241E-050.0001016360.06736096.72E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs10521232171348052932.627.09645E-050.000818447.096E-05rs1050903105543762875.1450.001836127.09645E-050.000818447.096E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs1406821	rs7584028	2	117190012	126.113	3.37709 <i>E</i> -05	7.68122 <i>E</i> –05	0.0198111	3.377 <i>E</i> -05
rs88581312162238048.530.0572213220.0001951213.493E-053.493E-05rs2041980176398961389.2834.36211E-050.0279766190.20231164.362E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs93076964134844616131.864.88235E-050.0001181160.36573474.882E-05rs93076964134844616131.864.88235E-050.0010818160.36573474.882E-05rs1050466387903123494.610.0124734464.99754E-050.00018014.998E-05rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs2304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.000112510.03494116.914E-05rs10521232171348052932.756.91427E-050.000112510.03494116.914E-05rs25458331165421598185.510.003679037.8984E-050.00038547.90E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs134	rs1493517	4	134829467	131.851	0.298266675	0.000672686	3.458 <i>E</i> -05	3.458 <i>E</i> -05
rs2041980176398961389.2834.36211E-050.0279766190.20231164.362E-05rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs93076964134844616131.864.88235E-050.0001181160.36573474.882E-05rs1050466387903123494.610.0124734464.99754E-050.00018014.998E-05rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs2304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.000112510.03494116.914E-05rs10521232171348052932.756.91427E-050.000112510.03494116.914E-05rs25458331165421598185.510.003679037.08984E-050.00085847.096E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs1340682180195427107.560.233607100.0008019417.815E-057.815E-05rs1340682180195427107.560.233607120.0006963520.61936367.889E-05rs13	rs885813	1	21622380	48.53	0.057221322	0.000195121	3.493 <i>E</i> -05	3.493 <i>E</i> -05
rs8030379158238191478.924.53895E-050.000116070.01995424.539E-05rs93076964134844616131.864.88235E-050.0010818160.36573474.8822E-05rs1050466387903123494.610.0124734464.99754E-050.00018014.998E-05rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs1389384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs236526916177223494.1866.47241E-050.0012046360.06736096.472E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs5458331165421598185.510.003679037.08984E-050.00035847.096E-05rs139003105543762875.1450.0123937.3719E-050.00013977.372E-05rs1340682180195427107.560.233607120.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs4939561688430532.627.88943E-050.000204568.577E-058.577E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs213588	rs2041980	17	63989613	89.283	4.36211 <i>E</i> -05	0.027976619	0.2023116	4.362 <i>E</i> -05
rs93076964134844616131.864.88235E-050.0010818160.36573474.882E-05rs1050466387903123494.610.0124734464.99754E-050.00018014.998E-05rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.0001012510.03494116.914E-05rs10521232171348052932.756.91427E-050.0001012510.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00038847.096E-05rs12340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs2061999158237374478.920.0279382910.00024568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs213794815514646480.5650.1127764420.0003541019.779E-059.779E-05<	rs8030379	15	82381914	78.92	4.53895 <i>E</i> -05	0.00011607	0.0199542	4.539 <i>E</i> -05
rs1050466387903123494.610.0124734464.99754E-050.00018014.998E-05rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.000112510.03494116.914E-05rs10521232171348052932.756.91427E-050.000112510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0123937.3719E-050.00013977.372E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.000294568.577E-058.577E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs9307696	4	134844616	131.86	4.88235 <i>E</i> -05	0.001081816	0.3657347	4.882 <i>E</i> -05
rs93216116137829536138.6850.8829192640.0008921225.233E-055.233E-05rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.00012046360.06736096.472E-05rs10521232171348052932.756.91427E-050.00011012510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.00026558.577E-058.577E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.00021160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05 <td>rs10504663</td> <td>8</td> <td>79031234</td> <td>94.61</td> <td>0.012473446</td> <td>4.99754<i>E</i>05</td> <td>0.0001801</td> <td>4.998<i>E</i>-05</td>	rs10504663	8	79031234	94.61	0.012473446	4.99754 <i>E</i> 05	0.0001801	4.998 <i>E</i> -05
rs463803543201269.410.9190493010.0009816965.306E-055.306E-05rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.0001012510.03494116.914E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs13406821204640851223.6970.0299123937.3719E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008619417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs9321611	6	137829536	138.685	0.882919264	0.000892122	5.233E-05	5.233 <i>E</i> -05
rs13893384134829546131.8516.18458E-050.0023860730.49995846.185E-05rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.0012046360.06736096.472E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs13406821204640851223.6970.0299123937.3719E-050.00013977.372E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.00021160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs463803	5	4320126	9.41	0.919049301	0.000981696	5.306E-05	5.306 <i>E</i> -05
rs72304223345335351.3180.3350523640.0175672776.41E-056.41E-05rs236526916177223494.1866.47241E-050.0012046360.06736096.472E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs125744461204640851223.6970.0299123937.3719E-050.00013977.372E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.0003541019.779E-059.779E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs1389338	4	134829546	131.851	6.18458 <i>E</i> -05	0.002386073	0.4999584	6.185 <i>E</i> –05
rs236526916177223494.1866.47241E-050.0012046360.06736096.472E-05rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00013977.372E-05rs22744461204640851223.6970.0299123937.3719E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs723042	2	33453353	51.318	0.335052364	0.017567277	6.41 <i>E</i> –05	6.41 <i>E</i> -05
rs10521232171348052932.756.91427E-050.0001012510.03494116.914E-05rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00085847.096E-05rs22744461204640851223.6970.0299123937.3719E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.00006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.0003541019.779E-059.779E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs2365269	1	61772234	94.186	6.47241 <i>E</i> -05	0.001204636	0.0673609	6.472 <i>E</i> –05
rs5458331165421598185.510.0003679037.08984E-050.0045627.09E-05rs10509003105543762875.1450.0018369127.09645E-050.00085847.096E-05rs22744461204640851223.6970.0299123937.3719E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.000204568.577E-058.577E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000321160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs10521232	17	13480529	32.75	6.91427 <i>E</i> -05	0.000101251	0.0349411	6.914 <i>E</i> –05
rs10509003105543762875.1450.0018369127.09645E-050.00085847.096E-05rs22744461204640851223.6970.0299123937.3719E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.0003541019.779E-059.779E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs545833	1	165421598	185.51	0.000367903	7.08984 <i>E</i> 05	0.004562	7.09 <i>E</i> –05
rs22744461204640851223.6970.0299123937.3719E-050.00013977.372E-05rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000321160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs10509003	10	55437628	75.145	0.001836912	7.09645 <i>E</i> 05	0.0008584	7.096 <i>E</i> -05
rs1340682180195427107.560.2336071020.0008019417.815E-057.815E-05rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs2274446	1	204640851	223.697	0.029912393	7.3719 <i>E</i> -05	0.0001397	7.372 <i>E</i> –05
rs105177024160646045157.997.84578E-050.1671934460.27642487.846E-05rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs1340682	1	80195427	107.56	0.233607102	0.000801941	7.815 <i>E</i> -05	7.815 <i>E</i> -05
rs64939561688430532.627.88943E-050.0006963520.61936367.889E-05rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs10517702	4	160646045	157.99	7.84578 <i>E</i> –05	0.167193446	0.2764248	7.846 <i>E</i> -05
rs2061999158237374478.920.0279382910.000204568.577E-058.577E-05rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs649395	6	16884305	32.62	7.88943 <i>E</i> -05	0.000696352	0.6193636	7.889 <i>E</i> –05
rs2135880158238061778.929.11851E-050.000221160.02681949.119E-05rs231794815514646480.5650.1127764420.0003541019.779E-059.779E-05	rs2061999	15	82373744	78.92	0.027938291	0.00020456	8.577 <i>E</i> -05	8.577 <i>E</i> -05
rs2317948 1 55146464 80.565 0.112776442 0.000354101 9.779E-05 9.779E-05	rs2135880	15	82380617	78.92	9.11851 <i>E</i> 05	0.00022116	0.0268194	9.119 <i>E</i> –05
	rs2317948	1	55146464	80.565	0.112776442	0.000354101	9.779 <i>E</i> –05	9.779 <i>E</i> 05

A power study for linkage on the 330 Framingham families with traits measured only in the Offspring Cohort, similar to our study, was performed using SOLAR. The power was estimated to be 97%, 84%, and 62% to detect a QTL heritability of 30% using a LOD score cut off of 1, 2, or 3, respectively, as significant (http://www.framinghamheartstudy.org/). The results of the simulation studies imply that this study sample only has sufficient power to detect large QTL effects. A power study for association on the same study population under additive genetic model demonstrated 97%, 100%, 94%, and 77% power for SNP minor allele frequencies of 5%, 10%, 20%, and 30%, assuming 10% QTL heritability, the allele frequency for the QTL to be 0.10, and the QTL and the marker are in linkage disequilibrium (D' = 1). The power was calculated using PBAT, a software for the family-based

association test (the same Framingham website mentioned above). Usually LME is expected to have higher power than PBAT, since only heterozygous parents are informative for a family-based study. The results of the simulation studies imply that this study population may have good power for linkage study. However, since the Affymetrix 100 K GeneChip has a mean marker distance of 26 kb, in most situations, the genomic coverage may not be high enough to reach the power mentioned above for association study unless there is a major gene with a larger effect size in linkage disequilibrium with the adjacent SNPs.

Although there is a strong genetic effect on the variation of serum LDH levels, the results of both linkage and association genome scans did not show evidence of a major gene effect. Instead, there may be many genes with small effects in controlling the variation of serum LDH levels. A limitation of this study is that our cohort is mostly composed of Caucasians. Therefore, caution is advised in extrapolating our results to other ethnicities.

Acknowledgments

The authors are very grateful to Nancy Geller for her helpful comments on the data analysis and critical review of the manuscript. This work was supported by the National Heart, Lung, and Blood Institute's Framingham Heart Study (Contract no. N01-HC-25195).

References

- G. Vasudevan, D. W. Mercer, and M. A. Varat, "Lactic dehydrogenase isoenzyme determination in the diagnosis of acute myocardial infarction," *Circulation*, vol. 57, no. 6, pp. 1055–1057, 1978.
- [2] S. S.-L. Li, M. Luedemann, F. S. Sharief, T. Takano, and L. L. Deaven, "Mapping of human lactate dehydrogenase-A, -B, and -C genes and their related sequences: the gene for LDHC is located with that for LDHA on chromosome 11," *Cytogenetics and Cell Genetics*, vol. 48, no. 1, pp. 16–18, 1988.
- [3] C. L. Markert, J. B. Shaklee, and G. S. Whitt, "Evolution of a gene. Multiple genes for LDH isozymes provide a model of the evolution of gene structure, function, and regulation," *Science*, vol. 189, no. 4197, pp. 102–114, 1975.
- [4] M. O. Rethore, C. Junien, and G. Malpuech, "Localization of the gene of the glyceraldehyde 3 phosphate dehydrogenase on the distal segment of the short arm of the chromosome 12," *Annales de Genetique*, vol. 19, no. 2, pp. 140–142, 1976.
- [5] L. Bathum, H. C. Petersen, J.-U. Rosholm, P. H. Petersen, J. Vaupel, and K. Christensen, "Evidence for a substantial genetic influence on biochemical liver function tests: results from a population-based Danish twin study," *Clinical Chemistry*, vol. 47, no. 1, pp. 81–87, 2001.
- [6] T. R. Dawber, The Framingham Heart Study: The Epidemiology of atherosclerotic Diseases, Harvard University Press, Cambridge, Mass, USA, 1980.
- [7] W. B. Kannel, M. Feinleib, and P. M. McNamara, "An investigation of coronary heart disease in families. The Framingham offspring study," *American Journal of Epidemiology*, vol. 110, no. 3, pp. 281–290, 1979.
- [8] J. B. Wilk, A. K. Manning, J. Dupuis, et al., "No evidence of major population substructure in the Framingham Heart Study," *Genetic Epidemiology*, vol. 29, no. 3, 286 pages, 2005.
- [9] B. Yuan, D. Vaske, J. L. Weber, J. Beck, and V. C. Sheffield, "Improved set of short-tandem-repeat polymorphisms for screening the human genome," *American Journal of Human Genetics*, vol. 60, no. 2, pp. 459–460, 1997.
- [10] A. Herbert, N. P. Gerry, M. B. McQueen et al., "A common genetic variant is associated with adult and childhood obesity," *Science*, vol. 312, no. 5771, pp. 279–283, 2006.
- [11] L. A. Cupples, H. T. Arruda, E. J. Benjamin, et al., "The Framingham Heart Study 100K SNP genome-wide association study resource: overview of 17 phenotype working group reports," *BMC Medical Genetics*, vol. 8, supplement 1, article S1, 2007.
- [12] L. Almasy and J. Blangero, "Multipoint quantitative-trait linkage analysis in general pedigrees," *American Journal of Human Genetics*, vol. 62, no. 5, pp. 1198–1211, 1998.

- [13] J. Blangero and L. Almasy, "Multipoint oligogenic linkage analysis of quantitative traits," *Genetic Epidemiology*, vol. 14, no. 6, pp. 959–964, 1997.
- [14] D. B. Allison, M. C. Neale, R. Zannolli, N. J. Schork, C. I. Amos, and J. Blangero, "Testing the robustness of the likelihood-ratio test in a variance- component quantitativetrait loci-mapping procedure," *American Journal of Human Genetics*, vol. 65, no. 2, pp. 531–544, 1999.
- [15] J. Blangero, J. T Williams, and L. Almasy, "Robust LOD scores for variance component-based linkage analysis," *Genetic Epidemiology*, vol. 19, supplement 1, pp. S8–S14, 2000.
- [16] J. Blangero, J. T. Williams, and L. Almasy, "12 Variance component methods for detecting complex trait loci," *Advances in Genetics*, vol. 42, pp. 151–181, 2001.
- [17] P. R. Burton, D. G. Clayton, L. R. Cardon, et al., "Genomewide association study of 14,000 cases of seven common diseases and 3,000 shared controls," *Nature*, vol. 447, no. 7145, pp. 661–678, 2007.
- [18] G. Zheng, J. Joo, D. Zaykin, C. O. Wu, and N. L. Geller, "Robust tests in genome-wide scans under incomplete linkage disequilibrium," *Statistical Science*, vol. 24, pp. 503–516, 2009.