

MEETING ABSTRACT

Open Access

# Transient transfection of insect Sf-9 cells in TubeSpin<sup>®</sup> bioreactor 50 tubes

Xiao Shen<sup>1†</sup>, Patrik O Michel<sup>1†</sup>, Qiuling Xie<sup>1,2</sup>, David L Hacker<sup>1</sup>, Florian M Wurm<sup>1\*</sup>

From 22nd European Society for Animal Cell Technology (ESACT) Meeting on Cell Based Technologies Vienna, Austria. 15-18 May 2011

## Background

Sf-9 cells, derived from *Spodoptera frugiperda*, are widely used for recombinant protein production using the baculovirus expression vector system (BEVS). However, this results in a productive viral infection and cell lysis. Therefore, a non-lytic, plasmid-based expression system for suspension Sf-9 cells would be a valuable alternative to the BEVS for rapid, scalable, and high-yielding recombinant protein production and for the generation of stable Sf-9 cell lines [4]. In this work, we present a simple, efficient and cost-effective plasmid-based method for transient expression of recombinant proteins in Sf-9 cells cultivated in serum-free suspension mode in a high-throughput culture system, TubeSpin<sup>®</sup> bioreactor 50 tubes (TubeSpins).

## Materials and methods

Sf-9 cells were maintained in suspension in TubeSpins (TPP, Trasadingen, Switzerland) at 28 °C [5]. The human tumor necrosis factor receptor-Fc fusion protein (TNFR-Fc) gene was cloned into pIEx10 (Novagen, Merck, Darmstadt, Germany) to generate pIEx-TNFR-Fc. The cells in exponential growth phase were inoculated in fresh Sf900 II medium (Invitrogen, Carlsbad, CA) one day prior to transfection. The next day, the cells were transfected with 1.5 µg pTNFR-Fc and 2.25 µg linear 25 kDa polyethylenimine (PEI, Polysciences, Warrington, PA) per 10<sup>6</sup> cells. The DNA and PEI were first mixed in de-ionized water at room temperature for 10 min prior to addition to the culture. At the time of transfection the cell density was 20 x 10<sup>6</sup> cells per mL. Cells were subsequently diluted to 4 x 10<sup>6</sup> cells per mL

with fresh Sf900 II media to allow for growth. The culture was maintained at 28 °C in an incubator shaker for 7 d [5]. The TNFR-Fc concentration in the medium was determined by ELISA as described [6].

## Results

Sf-9 cells were transfected in TubeSpins with pIEx-TNFR-Fc. By 7 d post-transfection, the TNFR-Fc concentration reached 42 mg/L (Figure 1). In a separate transfection with a plasmid expression the enhanced green fluorescent protein (GFP) gene, 58 % of the cells were GFP-positive at 5 d post-transfection (data not shown).

## Conclusions

This study validates the use of PEI for transient expression in suspension Sf-9 cell cultures. This system was used to produce both intracellular (GFP) and secreted (TNFR-Fc) proteins. The results show that PEI-mediated transient transfection is a fast and efficient alternative to BEVS for high-yielding protein expression in Sf-9 cells.

## Acknowledgements

This work has been supported by the CTI Innovation Promotion Agency of the Swiss Federal Department of Economic Affairs, by grants from the Guangdong Provincial Department of Science and Technology (2008B030301349), the MOE of China (211 Grant), and the Academy of Finland (decision no. 135820). TPP (Trasadingen, Switzerland) is acknowledged for providing TubeSpin<sup>®</sup> bioreactor 50 tubes.

## Author details

<sup>1</sup>Laboratory of Cellular Biotechnology, Faculty of Life Sciences, Ecole Polytechnique Fédérale de Lausanne, CH-1015 Lausanne, Switzerland.


<sup>2</sup>Institute of Bioengineering, Jinan University, Guangzhou 510632, Guangdong, P. R. China.

Published: 22 November 2011

\* Correspondence: florian.wurm@epfl.ch

† Contributed equally

<sup>1</sup>Laboratory of Cellular Biotechnology, Faculty of Life Sciences, Ecole Polytechnique Fédérale de Lausanne, CH-1015 Lausanne, Switzerland  
Full list of author information is available at the end of the article



**Figure 1** Sf-9 cells were transfected with pTNFR-Fc (■) plasmid DNA and PEI. The TNFR-Fc concentration in the medium was determined by ELISA at the times indicated.

#### References

1. De Jesus M, Wurm FM: Manufacturing recombinant proteins in kg-ton quantities using animal cells in bioreactors. *Eur J Pharm Biopharm* 2011, **78**:184-188.
2. Harrison RL, Jarvis DL: Transforming lepidopteran insect cells for continuous recombinant protein expression. *Methods Mol Biol* 2007, **388**:299-316.
3. Jarvis DL: Baculovirus-insect cell expression systems. *Methods Enzymol* 2009, **463**:191-222.
4. Li E, Brown SL, Dolman CS, Brown GB, Nemerow GR: Production of functional antibodies generated in a nonlytic insect cell expression system. *Prot Expr Purif* 2001, **21**:121-128.
5. Xie Q, Michel PO, Baldi L, Hacker DL, Zhang X, Wurm FM: TubeSpin bioreactor 50 for the high-density cultivation of Sf-9 insect cells in suspension. *Biotechnol Lett* 2011, **33**:897-902.
6. Oberbek A, Matasci M, Hacker DL, Wurm FM: Generation of stable, high-producing CHO cell lines by lentiviral vector-mediated gene transfer in serum-free suspension culture. *Biotechnol Bioeng* 2011, **108**:600-610.

doi:10.1186/1753-6561-5-S8-P37

Cite this article as: Shen *et al.*: Transient transfection of insect Sf-9 cells in TubeSpin® bioreactor 50 tubes. *BMC Proceedings* 2011 **5**(Suppl 8):P37.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at  
[www.biomedcentral.com/submit](http://www.biomedcentral.com/submit)

 BioMed Central