
iScience



Article

Human iPS cells engender corneal epithelial stem cells with holoclone-forming capabilities



Shinya Watanabe, Ryuhei Hayashi, Yuzuru Sasamoto, ..., Andrew J. Quantock, Natasha Y. Frank, Kohji Nishida

ryuhei.hayashi@ophthal.med. osaka-u.ac.jp (R.H.) knishida@ophthal.med. osaka-u.ac.jp (K.N.)

Highlights

Human iPS cell-derived corneal epithelia contain ABCB5-positive stem cells

The ABCB5-positive cells possess holocloneforming capabilities

An antibody-mediated ABCB5 blockade suppresses holoclone formation

Holoclone-forming stem cells are present in a human iPS cell-derived tissue construct

Watanabe et al., iScience 24, 102688 June 25, 2021 © 2021 The Authors. https://doi.org/10.1016/ j.isci.2021.102688

iScience

Article

Human iPS cells engender corneal epithelial stem cells with holoclone-forming capabilities

CellPress OPEN ACCESS

Shinya Watanabe,¹ Ryuhei Hayashi,^{1,2,3,12,*} Yuzuru Sasamoto,⁴ Motokazu Tsujikawa,⁵ Bruce R. Ksander,⁶ Markus H. Frank,^{7,8,9} Andrew J. Quantock,¹⁰ Natasha Y. Frank,^{4,8,11} and Kohii Nishida^{1,3,*}

SUMMARY

Human induced pluripotent stem cells (hiPSCs) can generate a multiplicity of organoids, yet no compelling evidence currently exists as to whether or not these contain tissue-specific, holoclone-forming stem cells. Here, we show that a subpopulation of cells in a hiPSC-derived corneal epithelial cell sheet is positive for ABCB5 (ATP-binding cassette, sub-family B, member 5), a functional marker of adult corneal epithelial stem cells. These cells possess remarkable holocloneforming capabilities, which can be suppressed by an antibody-mediated ABCB5 blockade. The cell sheets are generated from ABCB5⁺ hiPSCs that first emerge in 2D eye-like organoids around six weeks of differentiation and display corneal epithelial immunostaining characteristics and gene expression patterns, including sustained expression of ABCB5. The findings highlight the translational potential of ABCB5-enriched, hiPSC-derived corneal epithelial cell sheets to recover vision in stem cell-deficient human eyes and represent the first report of holoclone-forming stem cells being directly identified in an hiPSCderived organoid.

INTRODUCTION

Previously, we reported how human induced pluripotent stem cells (hiPSCs) in defined culture conditions can generate a two-dimensional eye-like organoid, referred to as an SEAM (self-formed, ectodermal, autonomous, multi-zone) (Hayashi et al., 2016, 2017, 2018). In a number of respects, SEAM formation mimics whole-eye development in that the innermost central area (zone 1) forms first, followed by the emergence of three more radially distant, concentric cell populations: zones 2–4 (Figures 1A and 1B). Specifically, cells in zone 1 are akin to presumptive neuroectoderm; zone 2 cells are developmental analogs of neural eye tissues comprising the neuroretina, neural crest, and retinal pigment epithelium (RPE); cells in zone 3 have characteristics of corneal and conjunctival epithelial cells, while zone 4 cells resemble general surface ectodermal cells that will likely differentiate into epidermal keratinocytes. Further investigation of cells in SEAM zone 3 indicated that they could be sorted by flow cytometry and that those that were positive for the cell surface makers, SSEA4 and CD104, and negative for CD200 (i.e. SSEA4⁺/CD104⁺/CD200⁻) possessed corneal epithelial characteristics and could be fashioned into functional corneal epithelial cell sheets (Hayashi et al., 2016, 2017, 2018). Our investigations have also identified ABCB5 as a functional marker of adult corneal epithelial stem cells (Frank et al., 2003; Ksander et al., 2014; Sasamoto et al., 2018, 2020).

Vision loss arising from a corneal epithelial stem cell deficiency invariably requires surgical intervention if a patient's sight is to be restored. Contemporary treatments often involve the *ex vivo* expansion of corneal epithelial stem cells obtained from the stem cell niche at the edge of the cornea (a region of the eye known as the limbus) to form a functional corneal epithelial cell sheet that can be grafted onto the eye (Pellegrini et al., 1997; Lindberg et al., 1993; Ang et al., 2007; Schwab et al., 2000). It has been suggested, based on retrospective immunostaining of expanded corneal epithelial cell sheets with the proposed stem cell marker, p63, that a successful surgical outcome correlates with the number of stem cells likely to be present in the expanded sheet (Rama et al., 2010). However, establishing whether or not holoclone-forming stem cells are present in an expanded corneal epithelial cell sheet—or in any other hiPSC-derived organoid for

¹Department of Ophthalmology, Osaka University Graduate School of Medicine, Suita, Osaka 565-0871, Japan

²Department of Stem Cells and Applied Medicine, Osaka University Graduate School of Medicine, Suita, Osaka 565-0871, Japan

³Integrated Frontier Research for Medical Science Division, Institute for Open and Transdisciplinary Research Initiatives, Osaka University, Suita, Osaka 565-0871, Japan

⁴Division of Genetics, Brigham and Women's Hospital, Harvard Medical School, Boston, MA 02115, USA

⁵Department of Biomedical Informatics, Osaka University Graduate School of Medicine, Division of Health Sciences, Suita, Osaka 565-0871, Japan

⁶Massachusetts Eye and Ear, Schepens Eye Research Institute, Harvard Medical School, Boston, MA 02114, USA

⁷Transplant Research Program, Division of Nephrology, Boston Children's Hospital, Harvard Medical School, Boston, MA 02115, USA

⁸Harvard Stem Cell Institute, Harvard University, Boston, MA 02138, USA

⁹School of Medical and Health Sciences, Edith Cowan University, Perth, WA 6027, Australia

¹⁰School of Optometry and Vision Sciences, Cardiff University, Cardiff CF24 4HQ, Wales, UK

¹¹Department of Medicine, VA Boston Healthcare System, Boston, MA 02130, USA ¹²Lead Contact

*Correspondence: ryuhei.hayashi@ophthal. med.osaka-u.ac.jp (R.H.), knishida@ophthal.med. osaka-u.ac.jp (K.N.) https://doi.org/10.1016/j.isci. 2021.102688












Figure 1. ABCB5⁺ corneal epithelial cells can be induced and isolated from hiPSCs

(A) Strategy to generate and sort hiPSC-derived cells. CDM, corneal differentiation medium; CEM, corneal epithelium maintenance medium; DM, differentiation medium; d, day(s); w, week(s).

(B) SEAM of hiPSCs mimics ocular development. CNS, central nervous system; NR, neuroretina; RPE, retinal pigment epithelium; NC, neural crest; LE; lens; CE, corneal epithelium; EK, epidermal keratinocyte.

(C) ABCB5 expression in SEAM zones 1-4 during differentiation (each n = 5 independent experiments). Error bars are standard deviation.

(D) Flow cytometric analysis of the differentiated SEAM. $0.71 \pm 0.62\%$ of hiPSC-derived corneal epithelial cells (SSEA4⁺/CD104⁺/CD200⁻) were ABCB5⁺. Data shown are representative of 58 independent cell sorting experiments.

that matter—has not been achieved. Here, we combine the expansion of hiPSCs into SEAMs with cell sorting based on ABCB5 specificity to investigate this issue.

RESULTS

ABCB5⁺ corneal epithelial cells can be induced and isolated from hiPSCs

ABCB5 mRNA was detected at 6 weeks in SEAM zone 2 and at 7 weeks in zone 3 and was maintained until 12 weeks of differentiation in both zones (Figure 1C). In contrast, no expression was seen in zone 1 or zone 4 throughout the whole cultivation period. Flow cytometry revealed that $0.71 \pm 0.62\%$ of SSEA4⁺/CD104⁺/ CD200⁻ corneal epithelial cells derived from zone 3 were ABCB5⁺ (Figure 1D). Of note, a similar percentage of ABCB5⁺ cells were found in corneal limbal epithelial cells isolated from a human cadaveric donor eye (Sasamoto et al., 2020). These findings revealed that ABCB5⁺ corneal epithelial cells can be induced and isolated from hiPSCs via an SEAM and that the localization of *ABCB5* expression in an SEAM mirrors that in human eyes, where it is found in the RPE and limbal epithelium (Ksander et al., 2014; Chen et al., 2005). We also observed that *ABCB5* expression first appears midway through SEAM formation, which would be consistent with the notion that ABCB5 is associated with adult stem cell activity during human ocular differentiation.

hiPSC-derived ABCB5⁺ corneal epithelial cells have self-renewal potential

To ascertain the self-renewal potential of ABCB5⁺ and ABCB5⁻ cells isolated from SEAM-derived corneal epithelial cells, we performed colony-forming assays and holoclone assays (Figures S1A and S1B). These indicated that the colony-forming efficiency was greater and the colony size was larger for ABCB5⁺ cells when compared to ABCB5⁻ cells (13.5 \pm 3.8% vs. 7.5 \pm 2.4%, P < 0.05; 1.56 \pm 0.65 mm vs. 0.79 \pm 0.43 mm, P < 0.05, respectively) (Figure 2A). Moreover, holoclone formation was a characteristic of ABCB5⁺ cells but not ABCB5⁻ cells (13.8 \pm 4.0% vs. 0.0 \pm 0.0%, P < 0.05) (Figure 2B), and even though the number of ABCB5⁺ cells in the SEAM-derived corneal epithelial cells was few, the holoclone-forming capability resided exclusively in this subpopulation. A gene expression analysis by quantitative real-time reverse-transcriptase PCR (qRT-PCR) further showed that ⊿Np63, a widely used marker of corneal epithelial stem cells (Rama et al., 2010; Pellegrini et al., 2001), was enriched in ABCB5⁺ cell-derived colonies compared to ABCB5⁻ cell-derived colonies (5.17 \pm 1.77% vs. 4.30 \pm 1.69%, P < 0.05) (Figure 2C). The expression of KRT12 and PAX6, known markers of differentiated corneal epithelial cells, was also lower in ABCB5⁺ cell-derived colonies (4.05 \pm 12.41% vs. 11.32 \pm 24.31%, P < 0.05; 0.96 \pm 2.30% vs. 2.51 \pm 3.63%, P < 0.05, for ABCB5⁺ and ABCB5⁻, respectively) (Figure 2C). These findings reflect the undifferentiated status of SEAM-derived ABCB5⁺ corneal epithelial cells and their holoclone-forming capacity.

ABCB5 blockade inhibits self-renewal of hiPSC-derived ABCB5⁺ corneal epithelial cells

To further investigate the connection between ABCB5 and SEAM-derived corneal epithelial cell self-renewal, colony-forming assays and holoclone assays were performed with and without the inclusion of a specific monoclonal antibody (mAb)-mediated ABCB5 blockade (Ksander et al., 2014; Wilson et al., 2014; Lutz et al., 2016; Schatton et al., 2015) (Figure 3A). This revealed that a consecutive 10-day treatment of SEAM-derived ABCB5⁺ cells with 50 µg/mL of anti-ABCB5 blocking mAb (clone 3C2-1D12) suppressed colony formation to levels that were significantly below those of the control, non-mAb-treated group (1.6 ± 0.7% vs. 12.2 ± 5.6%, P < 0.05) (Figure 3B). Holoclone formation, moreover, was found to be completely inhibited in the anti-ABCB5 mAb-treated group compared with the control group (0.0 ± 0.0% vs. 9.0 ± 5.4%, P < 0.05) (Figure 3C), which is supportive of a defining role for ABCB5 in the self-renewal potential of hiPSC-derived corneal epithelial stem cells.








Figure 2. hiPSC-derived ABCB5⁺ corneal epithelial cells have self-renewal potential

(A) Representative colony-forming assays of seven independent experiments for ABCB5⁺/ABCB5⁻ hiPSC-derived corneal epithelial cells plus colony-forming efficiency and colony diameter (left panel, n = 7 independent experiments; right panel, ABCB5⁺; n = 172 colonies: ABCB5⁻; n = 93 colonies from 4 independent experiments).
(B) Representative images of a holoclone and a non-holoclone and percentage of holoclones (n = 4 independent experiments). See also Figure S1.

(C), Gene expression analysis for stem cell marker $\Delta Np63$ and differentiated corneal epithelium-related markers PAX6 and KRT12 in ABCB5⁺/ABCB5⁻ colonies (ABCB5⁺; n = 131 colonies: ABCB5⁻; n = 55 colonies from 6 independent experiments). *p < 0.05, Wilcoxon signed-rank test. Error bars are standard deviation.

Corneal epithelial cell sheets can be fabricated from hiPSC-derived ABCB5⁺/ABCB5⁻ cells

It was discovered that ex vivo expanded corneal epithelial cell sheets could be fabricated from both ABCB5⁺ and ABCB5⁻ cells obtained from SEAM zone 3, although the ABCB5⁺ cell-derived sheets were more highly stratified and significantly thicker (10.41 \pm 0.88 µm vs. 5.65 \pm 0.81 µm, P < 0.05) (Figure 4A). The corneal epithelial marker, KRT12, the ocular cell marker, PAX6, and the epithelial/surface ectodermal marker, p63, were all identified in the ABCB5⁺ and ABCB5⁻ cell-derived sheets (Figure 4B). This represents a typical corneal epithelial immunostaining pattern, the validity of which was reinforced by the results of a gene expression analysis of the cell sheets (Figure 4C). The fact that the pattern of gene expression changes seen in the hiPSC-derived colonies, on the one hand, and in the expanded corneal epithelial cell sheets (21–28 days cultivation) compared to cells in the colonies








Figure 3. ABCB5 blockade inhibits self-renewal in hiPSC-derived ABCB5⁺ corneal epithelial cells (A) Schematic representation of mAb-mediated ABCB5 blocking followed by colony-forming assay (CFA) and holoclone assay (HA) for the ABCB5⁺ hiPSC-derived corneal epithelial cells (iCECs). Sorted ABCB5⁺ iCECs were treated with anti-ABCB5 blocking mAb or isotype control IgG for 10 consecutive days from the initial day of sorting. (B) Representative images of CFA after treatment with anti-ABCB5 blocking mAb or isotype control IgG (n = 4 independent experiments) and colony-forming efficiency after the treatments (n = 4 independent experiments). (C) Representative images of a holoclone and a non-holoclone after antibody treatment and percentage of holoclones after treatments with anti-ABCB5 blocking mAb or isotype control IgG (n = 4 independent experiments). *p < 0.05, Wilcoxon signed-rank test. Error bars are standard deviation.

(10–12 days cultivation). Indeed, the culture medium used for cell sheet fabrication in the current study is the same as that employed previously to generate mature corneal epithelial cell sheets (Miyashita et al., 2013; Hayashi et al., 2016, 2017). Alternative culture conditions more suited to stem cell cultivation would presumably lessen the maturation propensity of cells in the expanded cell sheets. Immunostaining also








Figure 4. Corneal epithelial cell sheets can be fabricated from hiPSC-derived ABCB5⁺/ABCB5⁻ cells

(A) Hematoxylin and eosin (H&E) staining of cell sheets fabricated from ABCB5⁺/ABCB5⁻ hiPSC-derived corneal epithelial cells (iCECs). Scale bar, 50 μ m. Comparative thickness of cell sheets fabricated from hiPSC-derived ABCB5⁺/ABCB5⁻ cells. (n = 10 measurements from 1 experiment). *p < 0.05, Wilcoxon signed-rank test. Error bars are standard deviation. (B) Immunostaining for ABCB5 and ocular epithelial markers (green) in ABCB5⁺/ABCB5⁻ iCEC-derived sheets. ABCB5, stem cell marker; KRT12, corneal epithelial marker, PAX6, ocular cell marker; p63, epithelial/surface ectodermal markers; KRT13, conjunctival epithelial marker; MITF, retinal pigment epithelium (RPE) marker. Nuclei, red. Scale bar, 50 μ m. Dotted lines indicate border between cell sheets and culture inserts.

(C) Gene expression analysis for corneal epithelial stem cell markers and corneal epithelium-related markers in ABCB5⁺/ ABCB5⁻ iCEC-derived sheets (n = 3 independent experiments). Error bars are standard deviation.

disclosed that neither the conjunctiva-specific keratin, KRT13, nor the RPE marker, MITF, was detected in the ABCB5⁺ or ABCB5⁻ cell-derived sheets (Figure 4B), suggesting that no contamination with non-corneal epithelial cells occurred during the fabrication of the hiPSC-derived corneal epithelial cell sheets.





DISCUSSION

Here, we report the successful isolation of ABCB5⁺ corneal epithelial cells from hiPSCs via an SEAM. Collectively, our data indicate that ABCB5⁺ corneal epithelial cells possess functional stem cell characteristics, namely a high holoclone-forming propensity accompanied by the expression of a corneal epithelial stem cell marker, *ANp63*, and the capability to differentiate and form the corneal epithelium (Rama et al., 2010; Pellegrini et al., 1999, 2001). The identification of ABCB5⁺ stem cells in an SEAM-derived corneal epithelium represents the first definitive report of holoclone-forming human stem cells in an organoid originating from hiPSCs. This has theoretical and practical relevance in the fields of stem cell biology and regenerative medicine. In particular, the results suggest that ABCB5⁺ SEAM-derived corneal epithelial stem cells are similar in nature to the population of adult corneal epithelial stem cells in human eyes, which can differentiate into transient amplifying cells and drive corneal epithelial turnover and repair (Lehrer et al., 1998). Translational relevance lies in the fact that although surgeries to regain vision in eyes blinded by a corneal epithelial stem cell deficiency already exist (Pellegrini et al., 1997; Lindberg et al., 1993; Ang et al., 2007; Schwab et al., 2000), the corneal epithelial cell sheets that are used in contemporary graft surgery are fabricated from corneal epithelial cells obtained from a small biopsy of corneal epithelial cells taken from the limbal stem cell niche of a healthy eye or from the limbus of postmortem donor tissue sourced from an eye bank. The desire to avoid invasive biopsies of healthy eyes, aligned to the worldwide shortage of donated tissue (Sasamoto et al., 2020; Gain et al., 2016), makes the prospect of generating hiPSC-derived corneal epithelial cell sheets for transplant surgery an attractive one. This is especially true when we remember that hiPSC-derived graft tissue containing corneal epithelial stem cells can be applied for patients with bilateral disease whose limbal stem cell niche has been destroyed in both eyes and for whom an autograft using a biopsy of healthy would otherwise not be an option. Previously, we discovered that transplantation of hiPSC/SEAM-derived corneal epithelial cell sheets onto experimental ocular surface wounds in rabbit eyes could recover ocular surface function (Hayashi et al., 2016), and we have recently started to perform a small series of first-in-human corneal epithelial graft surgeries using hiPSC/SEAM-derived corneal epithelial cell sheets as part of a clinical trial (UMIN000036539) that has real potential to revolutionize the treatment of corneal blindness caused by a corneal epithelial stem cell deficiency. The current discovery that the fabricated sheets do, indeed, contain holoclone-forming corneal epithelial stem cells represents a major advance in our knowledge, which provides us with confidence in the concept of the translational approach and demonstrates that functional stem cells can be formed in organoids derived from hiPSCs.

Limitations of the study

We demonstrated that hiPSC-derived ABCB5⁺ corneal epithelial stem cells possess self-renewal potency within corneal epithelial cell sheets formed *in vitro*. We reason that the hiPSC-derived ABCB5⁺ cells will likely contribute to the successful repair of an injured corneal surface after epithelial cell sheet graft surgery, but this has not been proven. *In vivo* experiments using experimental ocular surface wounds in animal eyes would be required to corroborate this supposition.

STAR***METHODS**

Detailed methods are provided in the online version of this paper and include the following:

- KEY RESOURCES TABLE
- **RESOURCE AVAILABILITY**
 - O Lead contact
 - Materials availability
 - Data and code availability
- EXPERIMENTAL MODEL AND SUBJECT DETAILS
 - Cell lines
- METHOD DETAILS
 - O Differentiation of hiPSCs using the SEAM method
 - O Quantitative real-time reverse-transcriptase PCR
 - Flow cytometry and cell sorting
 - O Colony-forming assay (CFA) and holoclone assay (HA)
 - O mAb-mediated ABCB5 blockade
 - O Fabrication of ex vivo expanded cell sheets
 - O Hematoxylin and eosin (H&E) staining





- O Immunofluorescence staining
- QUANTIFICATION AND STATISTICAL ANALYSIS
 - O Statistical analysis

SUPPLEMENTAL INFORMATION

Supplemental information can be found online at https://doi.org/10.1016/j.isci.2021.102688.

ACKNOWLEDGMENTS

We thank T Katayama, S Inoue, Y. Yasukawa, M Morita, Y Yamate, S Shibata, Y Kobayashi, and Y Ishikawa of Osaka University for technical assistance. This work was supported, in part, by a project for the realization of regenerative medicine of the Japan Agency for Medical Research and Development (AMED), a Grants-in-Aid for Scientific Research from the Ministry of Education, Culture, Sports, Science and Technology of Japan, an AMED-CREST grant, an Osaka Eye Bank Research Grant 2017, and an International Joint Research Promotion Program of Osaka University.

AUTHOR CONTRIBUTIONS

S.W., R.H., and Y.S. designed the research; S.W. performed the *in vitro* experiments and acquired the data; M.H.F. and N.Y.F. provided reagents (3C2-1D12); S.W. and R.H. analyzed the data and wrote the respective methods and results; Y.S., M.T., B.R.K, M.H.F, A.J.Q., N.Y.F., and K.N. supervised the project; and S.W., R.H., A.J.Q., and K.N. wrote the paper.

DECLARATION OF INTERESTS

M.H.F., B.R.K., and N.Y.F. are inventors or co-inventors of US and international patents assigned to Brigham and Women's Hospital, Boston Children's Hospital, the Massachusetts Eye and Ear Infirmary, and/or the VA Boston Healthcare System, Boston, MA, licensed to TICEBA GmbH (Heidelberg, Germany) and RHEACELL GmbH & Co. KG (Heidelberg, Germany). M.H.F. serves as a scientific advisor for TICEBA GmbH and RHEACELL GmbH & Co. KG.

Received: January 11, 2021 Revised: April 28, 2021 Accepted: June 2, 2021 Published: June 25, 2021

REFERENCES

Ang, L.P., Sotozono, C., Koizumi, N., Suzuki, T., Inatomi, T., and Kinoshita, S. (2007). A comparison between cultivated and conventional limbal stem cell transplantation for Stevens-Johnson syndrome. Am. J. Ophthalmol. 143, 178–180.

Barrandon, Y., and Green, H. (1987). Three clonal types of keratinocyte with different capacities for multiplication. Proc. Natl. Acad. Sci. U S A *84*, 2302–2306.

Chen, K.G., Szakács, G., Annereau, J., Rouzaud, F., Liang, X., Valencia, J.C., Nagineni, C.N., Hooks, J.J., Hearing, V.J., and Gottesman, M.M. (2005). Principal expression of two mRNA isoforms (ABCB 5alpha and ABCB 5beta) of the ATP-binding cassette transporter gene ABCB 5 in melanoma cells and melanocytes. Pigment Cell Res. 18, 102–112.

Frank, N.Y., Pendse, S.S., Lapchak, P.H., Margaryan, A., Shlain, D., Doeing, C., Sayegh, M.H., and Frank, M.H. (2003). Regulation of progenitor cell fusion by ABCB5 P-glycoprotein, a novel human ATP-binding cassettle transporter. J. Biol. Chem. 278, 47156–47165. Gain, P., Jullienne, R., He, Z., Aldossary, M., Acquart, S., Cognasse, F., and Thuret, G. (2016). Global survey of corneal transplantation and eye banking. JAMA Ophthalmol. *134*, 167–173.

Hayashi, R., Ishikawa, Y., Sasamoto, Y., Katori, R., Nomura, N., Ichikawa, T., Araki, S., Soma, T., Kawasaki, S., Sekiguchi, K., et al. (2016). Coordinated ocular development from human iPS cells and recovery of corneal function. Nature 531, 376–380.

Hayashi, R., Ishikawa, Y., Katori, R., Sasamoto, Y., Taniwaki, Y., Takayanagi, H., Tsujikawa, M., Sekiguchi, K., Quantock, A.J., and Nishida, K. (2017). Coordinated generation of multiple ocular-like cell lineages and fabrication of functional corneal epithelial cell sheets from human iPS cells. Nat. Protoc. *12*, 683–696.

Hayashi, R., Ishikawa, Y., Katayama, T., Quantock, A.J., and Nishida, K. (2018). CD200 facilitates the isolation of corneal epithelial cells derived from human pluripotent stem cells. Sci. Rep. *8*, 16550.

Ksander, B.R., Kolovou, P.E., Wilson, B.J., Saab, K.R., Guo, Q., Ma, J., McGuire, S.P., Gregory, M.S., Vincent, W.J.B., Perez, V.L., et al. (2014). ABCB5 is a limbal stem cell gene required for corneal development and repair. Nature 511, 353–357.

Lehrer, M.S., Sun, T.T., and Lavker, R.M. (1998). Strategies of epithelial repair: modulation of stem cell and transit amplifying cell proliferation. J. Cell Sci. 111, 2867–2875.

Lindberg, K., Brown, M.E., Chaves, H.V., Kenyon, K.R., and Rheinwald, J.G. (1993). In vitro propagation of human ocular surface epithelial cells for transplantation. Invest. Ophthalmol. Vis. Sci. *34*, 2672–2679.

Lutz, N.W., Banerjee, P., Wilson, B.J., Ma, J., Cozzone, P.J., and Frank, M.H. (2016). Expression of cell-surface marker ABCB5 causes characteristic modifications of glucose, amino acid and phospholipid metabolism in the G3361 melanoma-initiating cell line. PLoS One 11, e0161803.

Miyashita, H., Yokoo, S., Yoshida, S., Kawakita, T., Yamagami, S., Tsubota, K., and Shimmura, S. (2013). Long-term maintenance of limbal epithelial progenitor cells using rho kinase inhibitor and keratinocyte growth factor. Stem Cells Transl. Med. 2, 758–765.

iScience Article



Miyazaki, T., Futaki, S., Suemori, H., Taniguchi, Y., Yamada, M., Kawasaki, M., Hayashi, M., Kumagai, H., Nakatsuji, N., Sekiguchi, K., et al. (2012). Laminin E8 fragments support efficient adhesion and expansion of dissociated human pluripotent stem cells. Nat. Commun *3*, 1236.

Nakagawa, M., Taniguchi, Y., Senda, S., Takizawa, N., Ichisaka, T., Asano, K., Morizane, A., Doi, D., Takahashi, J., Nishizawa, M., et al. (2014). A novel efficient feeder-free culture system for the derivation of human induced pluripotent stem cells. Sci. Rep. *4*, 3594.

Pellegrini, G., Traverso, C.E., Franzi, A.T., Zingirian, M., Cancedda, R., and Luca, M.D. (1997). Long-term restoration of damaged corneal surfaces with autologous cultivated corneal epithelium. Lancet *349*, 990–993.

Pellegrini, G., Golisano, O., Paterna, P., Lambiase, A., Bonini, S., Rama, P., and Luca, M.D. (1999). Location and clonal analysis of stem cells and their differentiated progeny in the human ocular surface. J. Cell Biol. 145, 769–782.

Pellegrini, G., Dellambra, E., Golisano, O., Martinelli, E., Fantozzi, I., Bondanza, S., Ponzin, D., McKeon, F., and Luca, M.D. (2001). p63 identifies keratinocyte stem cells. Proc. Natl. Acad. Sci. U S A *98*, 3156–3161.

Rama, P., Matuska, S., Paganoni, G., Spinelli, A., Luca, M.D., and Pellegrini, G. (2010). Limbal stemcell therapy and long-term corneal regeneration. N. Engl. J. Med. 363, 147–155.

Sasamoto, Y., Ksander, B.R., Frank, M.H., and Frank, N.Y. (2018). Repairing the corneal epithelium using limbal stem cells or alternative cell-based therapies. Expert Opin. Biol. Ther. 18, 505–513.

Sasamoto, Y., Sasamoto, N., Tran, J., Mishra, A., Ksander, B.R., Frank, M.H., and Frank, N.Y. (2020). Investigation of factors associated with ABCB5positive limbal stem cell isolation yields from human donors. Ocul. Surf 18, 114–120. Schatton, T., Yang, J., Kleffel, S., Uehara, M., Barthel, S.R., Schlapbach, C., Zhan, Q., Dudeney, S., Mueller, H., Lee, N., et al. (2015). ABCB5 identifies immunoregulatory dermal cells. Cell Rep. *12*, 1564–1574.

Schwab, I.R., Reyes, M., and Isseroff, R.R. (2000). Successful transplantation of bioengineered tissue replacements in patients with ocular surface disease. Cornea *19*, 421–426.

Takahashi, K., Tanabe, K., Ohnuki, M., Narita, M., Ichisaka, T., Tomoda, K., and Yamanaka, S. (2007). Induction of pluripotent stem cells from adult human fibroblasts by defined factors. Cell 131, 861–872.

Wilson, B.J., Saab, K.R., Ma, J., Schatton, T., Putz, P., Zhan, Q., Murphy, G.F., Gasser, M., Waaga-Gasser, A.M., Frank, N.Y., et al. (2014). ABCB5 maintains melanoma-initiating cells through a proinflammatory cytokine signaling circuit. Cancer Res. 74, 4196–4207.





STAR*METHODS

KEY RESOURCES TABLE

Mouse anti-ABCBSKrander et al., 2014N/AGost anti-KRT12Santa Cruz BiotechnologyCattso: 17098, RID:AB_639266Mouse anti-KRT13AbcamCattabiol112; RID:AB_30267Mouse anti-MIFExalpha BiologicalsCattaAc21202; RID:AB_30267ConjugateCattaAc21202; RID:AB_141607CattaAc21202; RID:AB_141607ConjugateDonkey anti-mouse IgG Alexa Fluor 488Life TechnologiesCattaAc21206; RID:AB_253102ConjugateDonkey anti-tabbit IgG Alexa Fluor 488Life TechnologiesCattaAc1206; RID:AB_253102ConjugateBiologondCatta327808; RID:AB_2534102ConjugateBiologondCatta327808; RID:AB_205274Mouse anti-SCD0 PE-Cy7 conjugateBiologondCatta22780; RRID:AB_2129146Mouse anti-CD0 PE-Cy7 conjugateBiologondCatta2208; RID:AB_2129146Mouse anti-CD0 PE-Cy7 conjugateBiologondCatta2200;Chemicals, peptides, and recombinant proteinsCattaAc1200;CattaB22002Chemicals, peptides, and recombinant proteinsCattaB22002CattaB22002Chemicals, peptides, and recombinant proteinsCattaB22002CattaB22002DMEM Filty Blucose, provate, no glutamine)Life TechnologiesCattB1120-003Chemicals, peptides, and recombinantLife TechnologiesCattB1120-013Chemicals, peptides, and recombinantLife TechnologiesCattB1120-021Maris STI (MST1E8)NippiCattaB22028Chemicals, peptides, and recombinantLife TechnologiesCattB120012Chemicals, peptides, and RECOCattB120-001Y2	REAGENT or RESOURCE	SOURCE	IDENTIFIER
Norse anti-p63Santa Cruz BiotechnologyCatise-2431 (RD-AB_62007)Mause anti-ABCB5Kander et al., 2014N/AGaat anti-KT12Santa Cruz BiotechnologyCatise-17098, RD-AB_6320267Mouse anti-KT13AbcamCatise-17098, RD-AB_6320267Mouse anti-MTFEsalpha BiologicaisCatise-27098, RD-AB_14107Donkey anti-mouse IgG Alexa Fluor 488Life TechnologiesCatise-2102, RD-AB_141708Donkey anti-rabbit IgG Alexa Fluor 488Life TechnologiesCatise-2102, RD-AB_141708Donkey anti-goat IgG Alexa Fluor 488Life TechnologiesCatise-2100, RD-AB_2534102ConjugateBiologendCatise-2100, RD-AB_109204Mouse anti-SEA4 FITC conjugateBiologendCatis230410, RND-AB_109204Mouse anti-SEA4 FITC conjugateBiologendCatis230410, RND-AB_1219146Mouse anti-SEA4 FITC conjugateMiponotoCatis230410, RND-AB_1219146Mouse anti-SEA4 FITC conjugateBiologiesCatis1700, Catis240202Mouse anti-SEA4 FITCMiponotoCatis230410, RND-AB_1219146Mouse anti-SEA4 FITCMiponotoCatis230410, RND-AB_230410Matri-S111 (MS1TEB)NipoiCatis23021Mouse anti-SEA4 FITCMiponotoCatis1700, Catis240202Mouse An	Antibodies		
Mouse anti-ABCBSKander et al. 2014N/AGat anti-KRT12Santa Cruz Biote/mologyCat/sci 7098; RID:AB_639265Mouse anti-KRT13AbcamCat/sci 71098; RID:AB_639265Mouse anti-MITFExalpha BiologicalsCat/sci 7098; RID:AB_10109Donkey anti-mouse IgG Alexa Fluor 488Life TechnologiesCat/sci 7098; RID:AB_141007ConjugateCat/sci 7008; RID:AB_141007Cat/sci 7008; RID:AB_141008Donkey anti-rabelit IgG Alexa Fluor 488Life TechnologiesCat/sci 7008; RID:AB_141708ConjugateDiokey anti-goat IgG Alexa Fluor 488Life TechnologiesCat/sci 7008; RID:AB_1231002Mouse anti-SCI 4FIC conjugateBiologendCat/sci 7008; RID:AB_1232024Mouse anti-CD100 FE conjugateBiologendCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateBiologendCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateBiologendCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateAjinomotoCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateAjinomotoCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateJinomotoCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateJinomotoCat/sci 7008; RID:AB_1212914Cherlicki, poptifes, and recombinant proteinsCat/sci 7008; RID:AB_1212914Mouse anti-CD200 FE-Cy7 conjugateLife TechnologiesCat/sci 7004; RID:AB_120031Cherlicki, poptifes, and recombinant proteinsCat/sci 7004; RID:AB_120031Cherlicki, poptifes, and recombinant proteinsCa	Rabbit anti-PAX6	Covance Research Products Inc	Cat#PRB-278P; RRID:AB_291612
Gast anti-KR112Santa Cruz BiotechnologyCatilac-17098, RBD.AB_639264Mouse anti-KR113AbcamCatilac 2108, RBD.AB_030267Mouse anti-MITFExalpha BiologialsCattlA 21002, RBD.AB_141607Donkey anti-mouse IgG Alexa Fluor 488Life TechnologiesCattlA 21206, RBD.AB_14167Donkey anti-rabbit IgG Alexa Fluor 488Life TechnologiesCattlA 21206, RBD.AB_111708ConjugateCattlA 21206, RBD.AB_141708CattlA 21206, RBD.AB_231102Donkey anti-gast IgG Alexa Fluor 488Life TechnologiesCattlA 21206, RBD.AB_231102ConjugateBiologendCattlA 21206, RBD.AB_2129146Mouse anti-SCHA FITC conjugateBiologendCattlA 21206, RBD.AB_2129146Mouse anti-CD104 FE conjugateBiologendCattlA 2008, RBD.AB_2129146Mouse anti-CD104 FE conjugateBiologendCattlA 2008Mouse anti-CD200 FE-Cy1 conjugateBiologendCattlA 2008Mouse anti-CD200 FE-Cy1 conjugateAjnomotoCattlA 2003Chemicals, peptides, and recombinant proteinsCattlA 2003CattlA 2003Mouse CD20 FE-Cy1 conjugateAjnomotoCattlA 2003Chemicals (CG Ef and FGF)CHLILEC Advanced Call SystemsCattlA 2003Chemicals (CG Ef and FGF)CHLILEC Advanced Call SystemsCattlA 2003DMEM/F12 (1:1)Life TechnologiesCattlA 2004Mark 5-12 Nutrient MixLife TechnologiesCattlA 2004Y-27632WakoCattlA 2004Rocombinant Human KGFWakoCattlA 2004Rocombinant Human KGFWakoCattlA 2004<	Mouse anti-p63	Santa Cruz Biotechnology	Cat#sc-8431; RRID:AB_628091
Abcam Catter Mouse anti-KRT13 Abcam Catter Catt	Mouse anti-ABCB5	Ksander et al., 2014	N/A
Mouse anti-MITF Exalpha Biologicals Cat#X2378M Denkey anti-mouse IgG Alexa Fluor 488 Life Technologies Cat#A-2120; RRID:AB_141708 Donkey anti-rabbit IgG Alexa Fluor 488 Life Technologies Cat#A-11055; RRID:AB_2534102 Conjugate Donkey anti-rabbit IgG Alexa Fluor 488 Life Technologies Cat#A-11055; RRID:AB_2534102 Conjugate BioLegend Cat#A-21206; RRID:AB_2129146 Mouse anti-CD104 PE conjugate BioLegend Cat#A-21206; RRID:AB_2129146 Mouse anti-CD104 PE conjugate BioLegend Cat#A-21206; RRID:AB_2129146 Mouse anti-CD104 PE conjugate Cat#A-21206; RRID:AB_2129146 Mouse anti-CD104 PE conjugate BioLegend Cat#A-21206; RRID:AB_2129146 Cat#A-21206; RRID:AB_2129146 Mouse anti-CD104 PE conjugate BioLegend Cat#A-21206; RRID:AB_2129146 Cat#A-21206; RRID:AB_2129146 Mouse anti-CD200 PE Cy7 conjugate BioLegend Cat#A-21206; RRID:AB_212916 Cat#A-21206; RRID:AB_212916 Chemicky, peptides, and recombinant proteins File Technologies Cat#1700-001 Cat#A-21206; RRID:AB_212916 Mouse anti-CD104 FE conjugate Life Technologies Cat#1110-035 Cat#1130-032 Monthing flug luccse, pyruwate, no	Goat anti-KRT12	Santa Cruz Biotechnology	Cat#sc-17098; RRID:AB_639266
Denkay anti-mouse IgG Alava Fluor 488 Lif Technologies Cat#A-21202; RRID:AB_141607 Donkay anti-rabbit IgG Alexa Fluor 488 Life Technologies Cat#A-21206; RRID:AB_141708 Donkay anti-rabbit IgG Alexa Fluor 488 Life Technologies Cat#A-21206; RRID:AB_2331102 Donkay anti-rabbit IgG Alexa Fluor 488 Life Technologies Cat#A-21206; RRID:AB_2331102 Donkay anti-rabbit IgG Alexa Fluor 488 Life Technologies Cat#A-21206; RRID:AB_2331102 Mouse anti-CD104 PE conjugate Biolegend Cat#327808; RRID:AB_2129146 Mouse anti-CD200 PE Cy7 conjugate BD PharMingen Cat#Ac330 Mouse anti-CD104 PE conjugate Ajnomoto Cat#Ac330 Chemicab, peptides, and recombinant proteins Sammiti Cat#Ac330 Sammiti Ajnomoto Cat#Ac330 MDMEM/T2 (1:1) Life Technologies Cat#1320-033 DMEM/M12 (1:1) Life Technologies Cat#1120-033 DMEM/T3 (2:1) Life Technologies Cat#1320-033 DMEM/T3 (2:1) Life Technologies Cat#1120-0313 DMEM/T3 (2:1) Life Technologies Cat#1120-0313 DMEM/T3 (2:1) Life Technologies	Mouse anti-KRT13	Abcam	Cat#ab16112; RRID:AB_302267
conjugate Cat#A-21206, RRID:AB, 141708 Donkey anti-abbit IgG Alexa Fluor 488 Life Technologies Cat#A-11055, RRID:AB, 141708 Donkey anti-goat IgG Alexa Fluor 488 Life Technologies Cat#330410; RRID:AB, 2034102 Mouse anti-CD200 PE-Cy7 conjugate BioLegend Cat#32080; RRID:AB, 2129146 Mouse anti-CD200 PE-Cy7 conjugate BD PharMingen Cat#624052 Chemicals, peptdes, and recombinant proteins StemFit Ajnomoto Cat#624052 Mixtris-S11 (USS1168) Nipp I Cat#624052 Cat#624052 GMEM Life Technologies Cat#A11033 Cat#624052 Mixtris-S11 (USS1168) Nipp I Cat#624052 Cat#624052 GMEM Life Technologies Cat#70310 Cat#62031 GMEM (LISS1168) Nipp I Cat#62031 Cat#7170-035 Chemicel Work (EGF and FGF) CELLITEC Advanced Cell Systems Cat#7130-031 DMEM (high glucose, pyruvate, no glutamine) Life Technologies Cat#1120-0313-021 Ham's F12 Nutrient Mix Life Technologies Cat#11200813-021 Non-ossential anino acids Life Technologies Cat#1780-022	Mouse anti-MITF	Exalpha Biologicals	Cat#X2398M
conjugateDankey anti-goat IgG Alexa Fluor 488Life TechnologiesCat#A11055; RRIDxAB_2534102Mouse anti-SEA4 FITC conjugateBioLegendCat#330410; RRIDxAB_1089204Mouse anti-CD104 PE conjugateBioLegendCat#327808; RRIDxAB_2129146Mouse anti-CD200 PE-Cy7 conjugateBD PharMingenCat#62052Chemicals, peptides, and recombinant proteinsTTStem FitiMinomotoCat#62012GMEMLife TechnologiesCat#1110-035GNEMLife TechnologiesCat#1110-035GNEMLife TechnologiesCat#1110-035DMEM/F12(1:1)Life TechnologiesCat#1031-021Marin SF112 Nutriert MixLife TechnologiesCat#1105062Y-27532WakoCat#1105062Sociam proventsLife TechnologiesCat#110405062Y-27532WakoCat#11005062Sociam proventsLife TechnologiesCat#1104004Sociam p	, ,	Life Technologies	Cat#A-21202; RRID:AB_141607
conjugateNumber of the second sec	Donkey anti-rabbit IgG Alexa Fluor 488 conjugate	Life Technologies	Cat#A-21206; RRID:AB_141708
Mouse anti-CD104 PE conjugateBioLegendCat#327808; RRID-AB_2129146Mouse anti-CD200 PE-Cy7 conjugateBD PharMingenCat#624052Chemicals, peptides, and recombinant proteinsCat#624052StemFitAjinomotoCat#4K03NMidatix-S11 (LN511E8)NippiCat#692012GMEMCat#052012Cat#110-035Chemicals, peptides, and recombinant proteinsCat#110-035StemFitMide Te chonologiesCat#1110-035Cht-Prime (w/o; EGF and FGF)CELL/TEC Advanced Cell SystemsCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033Nore System (FSR)Ufe TechnologiesCat#1140-050Recombinant Human KGFWakoCat#1140-050Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050PhizerCat#195-15791Cat#252301492-MercaptoethanolKakoCat#130-0145Phosphate buffer solutionWako	Donkey anti-goat IgG Alexa Fluor 488 conjugate	Life Technologies	Cat#A-11055; RRID:AB_2534102
Mouse anti-CD200 PE-Cy7 conjugateBD PharMingenCat#624052Chemicals, peptides, and recombinant proteinsStem FitAjinomotoCat#AK03NStem FitAjinomotoCat#892012GMEMLife TechnologiesCat#892012GMEMLife TechnologiesCat#11710-035Cnt-Prime (w/o; EGF and FGF)CELLnTEC Advanced Cell SystemsCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM fnigh glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033DMEM frigh glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033DMEM frigh glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033DMEM frigh glucose, pyruvate, no glutamine)Life TechnologiesCat#11200813B27 supplementLife TechnologiesCat#11320-033Rocokut serum replacement (KSR)Life TechnologiesCat#11320-030Non-essential amino acidsLife TechnologiesCat#11320-033Penicillin-streptomycinLife TechnologiesCat#12030149L-GlutamineLife TechnologiesCat#25030149L-GlutamineLife TechnologiesCat#25030149L-GlutamineMakoCat#25030149L-GlutamineWakoCat#3240-07951MonchrioghycerolWako<		BioLegend	Cat#330410; RRID:AB_1089204
Chemicals, peptides, and recombinant proteins StemFit Ajinomoto Cat#AK03N iMatrix-511 (LN511E8) Nippi Cat#892012 GMEM Life Technologies Cat#11710-035 Cnt-Prime (w/o; EGF and FGF) CELLITEC Advanced Cell Systems Cat#11320-033 DMEM/F12 (1:1) Life Technologies Cat#11320-033 DMEM (high glucose, prywate, no glutamine) Life Technologies Cat#11320-033 DMEM (high glucose, prywate, no glutamine) Life Technologies Cat#1031-021 Ham's F-12 Nutrient Mix Life Technologies Cat#11765062 Y-27632 Wako Cat#1172-00813 B27 supplement Life Technologies Cat#1120-00813 B27 supplement Life Technologies Cat#11300-070 Non-assential amino acids Life Technologies Cat#11300-070 Non-assential amino acids Life Technologies Cat#11300-070 Penicillin-streptomycin Life Technologies Cat#11300-070 Non-assential amino acids Life Technologies Cat#1130-070 Penicillin-streptomycin Life Technologies Cat#11300-070 Valeocatrisone sodium succinate Phizer Cat#2	Mouse anti-CD104 PE conjugate	BioLegend	Cat#327808; RRID:AB_2129146
StemFitAjinomotoCat#AK03NiMatrix-511 (LN511E8)NippiCat#892012GMEMLife TechnologiesCat#11710-035Cnt-Prime (w/o; EGF and FGF)CELLnTEC Advanced Cell SystemsCat#CnT-PR-EFDMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#1165062Y-27632WakoCat#11765062Y-27632WakoCat#112-00813B27 supplementLife TechnologiesCat#112-00813Knockout serum replacement (KSR)Life TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Lydrocortisone sodium succinatePhizerCat#25030149L-GlutamineLife TechnologiesCat#195-157912-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#1750-023MonothioglycerolWakoCat#163-201452-MercaptoethanolLife TechnologiesCat#163-2014510% formaldehyde neutral buffer solution<	Mouse anti-CD200 PE-Cy7 conjugate	BD PharMingen	Cat#624052
Matrix-S11 (LNS11E8)NippiCat#892012GMEMLife TechnologiesCat#11710-035GNEMCat#CnT-PR-EFDMEM(F12 (1:1)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033PMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11320-033PMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11360-021Paris F-12 Nutrient MixLife TechnologiesCat#112600813Recombinant Human KGFWakoCat#1120-031B27 supplementLife TechnologiesCat#11360-070Knockout serum replacement (KSR)Life TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1540-122LGlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolLife TechnologiesCat#163-201452-MercaptoethanolWakoCat#163-201452-MercaptoethanolWakoCat#164-079512-Mercaptoethanol	Chemicals, peptides, and recombinant proteins		
GMEMLife TechnologiesCat#11710-035Cnt-Prime (w/o; EGF and FGF)CELLnTEC Advanced Cell SystemsCat#CnT-PR-EFDMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#11031-021Ham's F-12 Nutrient MixLife TechnologiesCat#11765062Y-27632WakoCat#1120813B27 supplementLife TechnologiesCat#11200813B27 supplementLife TechnologiesCat#11360-070Knockout serum replacement (KSR)Life TechnologiesCat#1140050Sodium pyruvateLife TechnologiesCat#1140-050Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#162-0144% (wt/ot) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#13-02110% formaldehyde neutral buffer solutionWakoCat#13-0214MethanolWakoCat#162-01641MethanolWakoCat#13-021MethanolWakoCat#13-01826MethanolWakoCat#13-01826MethanolWakoCat#13-01826MethanolWakoCat#13-01826MethanolWakoCat#13-01826MethanolWakoCat#13-01826MethanolWakoCat#142065-010Methanol<	StemFit	Ajinomoto	Cat#AK03N
Act-Prime (w/c); EGF and FGF)CELLnTEC Advanced Cell SystemsCat#CnT-PR-EFDMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#1031-021Ham's F-12 Nutrient MixLife TechnologiesCat#11765062Y-27632WakoCat#112-00813B27 supplementLife TechnologiesCat#112-00813Knockout serum replacement (KSR)Life TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#2282L-GlutamineLife TechnologiesCat#21985-023MonothioglycerolWakoCat#21985-023MonothioglycerolWakoCat#21985-023MonothioglycerolWakoCat#21985-023MonothioglycerolWakoCat#132-015Yé (w/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#131-01826MothioglycerolWakoCat#131-01826MethanolWakoCat#31-01826Hoechst 33342WakoCat#313-01826Hoechst 33342Life TechnologiesCat#346-07951TrypLE ExpressLife TechnologiesCat#12-0161Mitomycin C (Medicine)Kyowa Hakko KirinN/A	iMatrix-511 (LN511E8)	Nippi	Cat#892012
DMEM/F12 (1:1)Life TechnologiesCat#11320-033DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#1031-021Ham's F-12 Nutrient MixLife TechnologiesCat#11765062Y-27632WakoCat#034-24024Recombinant Human KGFWakoCat#112-00813B27 supplementLife TechnologiesCat#110828-028Knockout serum replacement (KSR)Life TechnologiesCat#1140-050Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#252L-GlutamineLife TechnologiesCat#220301492-MercaptoethanolLife TechnologiesCat#203149MonothioglycerolWakoCat#195-15791MonothioglycerolWakoCat#131-01826MonothioglycerolWakoCat#131-01826Hoechst 33342WakoCat#31-01826Hoechst 33342Life TechnologiesCat#31-01826Hordycin C (Medicine)Kyowa Hakko KirinN/A	GMEM	Life Technologies	Cat#11710-035
DMEM (high glucose, pyruvate, no glutamine)Life TechnologiesCat#10313-021Ham's F-12 Nutrient MixLife TechnologiesCat#11765062Y-27632WakoCat#034-24024Recombinant Human KGFWakoCat#12-00813B27 supplementLife TechnologiesCat#10828-028Knockout serum replacement (KSR)Life TechnologiesCat#1140-050Non-essential amino acidsLife TechnologiesCat#1140-050Non-essential amino acidsLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#20301492-MercaptoethanolLife TechnologiesCat#21985-023NonothioglycerolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#163-2014510% formaldehyde neutral buffer solutionWakoCat#31-01826HethanolWakoCat#31-01826Hoechst 33342WakoCat#31-01826TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A	Cnt-Prime (w/o; EGF and FGF)	CELLnTEC Advanced Cell Systems	Cat#CnT-PR-EF
Ham's F-12 Nutrient MixLife TechnologiesCat#11765062Y-27632WakoCat#034-24024Recombinant Human KGFWakoCat#112-00813B27 supplementLife TechnologiesCat#17504-044Knockout serum replacement (KSR)Life TechnologiesCat#11860-070Sodium pyruvateLife TechnologiesCat#11140-050Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat#1062-01641MorothioglycerolWakoCat#1062-01661MethanolWakoCat#131-01826MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A	DMEM/F12 (1:1)	Life Technologies	Cat#11320-033
Y27632WakoCat#034-24024Recombinant Human KGFWakoCat#112-00813B27 supplementLife TechnologiesCat#17504-044Knockout serum replacement (KSR)Life TechnologiesCat#10828-028Sodium pyruvateLife TechnologiesCat#1140-050Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#250301492-MercaptoethanolLife TechnologiesCat#250301492-MercaptoethanolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#163-2014510% formaldehyde neutral buffer solutionWakoCat#131-01826MethanolWakoCat#11-01826MethanolWakoCat#11-01826TypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A	DMEM (high glucose, pyruvate, no glutamine)	Life Technologies	Cat#10313-021
Recombinant Human KGFWakoCat#112-00813B27 supplementLife TechnologiesCat#17504-044Knockout serum replacement (KSR)Life TechnologiesCat#10828-028Sodium pyruvateLife TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat#195-157914% (w/ol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#163-2014510% formaldehyde neutral buffer solutionWakoCat#131-01826MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TryLE ExpressKiowa Makko KirinN/A	Ham's F-12 Nutrient Mix	Life Technologies	Cat#11765062
B27 supplementLife TechnologiesCat#17504-044Knockout serum replacement (KSR)Life TechnologiesCat#10828-028Sodium pyruvateLife TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#163-2014510% formaldehyde neutral buffer solutionWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A	Y-27632	Wako	Cat#034-24024
Knockout serum replacement (KSR)Life TechnologiesCat#10828-028Sodium pyruvateLife TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#1140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#872452L-GlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#195-023MonothioglycerolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#31-0182610% formaldehyde neutral buffer solutionWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010MorothioglycerolWako KirrinN/A	Recombinant Human KGF	Wako	Cat#112-00813
Sodium pyruvateLife TechnologiesCat#11360-070Non-essential amino acidsLife TechnologiesCat#11140-050Penicillin-streptomycinLife TechnologiesCat#1540-122Hydrocortisone sodium succinatePhizerCat#872452L-GlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#104-20166110% formaldehyde neutral buffer solutionWakoCat#31-01826MothanolWakoCat#346-07951Hoechst 33342WakoCat#12605-010TrypLE ExpressLife TechnologiesCat#12605-010Monothicicine)Waya Hakko KirinN/A	B27 supplement	Life Technologies	Cat#17504-044
Non-essential amino acidsLife TechnologiesCat#11140-050Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#872452L-GlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat# 195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat# 163-2014510% formaldehyde neutral buffer solutionWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Minor, C (Medicine)Kyowa Hakko KirinN/A	Knockout serum replacement (KSR)	Life Technologies	Cat#10828-028
Penicillin-streptomycinLife TechnologiesCat#15140-122Hydrocortisone sodium succinatePhizerCat#872452L-GlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA)WakoCat# 163-20145phosphate buffer solutionWakoCat#163-2014510% formaldehyde neutral buffer solutionWakoCat#131-01826MethanolWakoCat#131-01826Hoechst 33342WakoCat#12605-010TrypLE ExpressLife TechnologiesCat#12605-010MonotinoglycerolKyowa Hakko KirinN/A	Sodium pyruvate	Life Technologies	Cat#11360-070
Hydrocortisone sodium succinatePhizerCat#872452L-GlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#195-023MonothioglycerolWakoCat#195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat#163-2014510% formaldehyde neutral buffer solutionWakoCat#062-01661MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Minorycin C (Medicine)Kyowa Hakko KirinWako	Non-essential amino acids	Life Technologies	Cat#11140-050
L-GlutamineLife TechnologiesCat#250301492-MercaptoethanolLife TechnologiesCat#21985-023MonothioglycerolWakoCat# 195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat# 163-2014510% formaldehyde neutral buffer solutionWakoCat#062-01661MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010MethanolKyowa Hakko KirinN/A	Penicillin-streptomycin	Life Technologies	Cat#15140-122
2-Mercaptoethanol Life Technologies Cat#21985-023 Monothioglycerol Wako Cat 195-15791 4% (wt/vol) Paraformaldehyde (PFA) Wako Cat 163-20145 phosphate buffer solution Wako Cat#062-01661 Methanol Wako Cat#062-01661 Methanol Wako Cat#131-01826 Hoechst 33342 Wako Kirin Cat#346-07951 TrypLE Express Life Technologies Cat#12605-010	Hydrocortisone sodium succinate	Phizer	Cat#872452
MonothioglycerolWakoCat# 195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat# 163-2014510% formaldehyde neutral buffer solutionWakoCat#062-01661MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinM/A	L-Glutamine	Life Technologies	Cat#25030149
MonothioglycerolWakoCat# 195-157914% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solutionWakoCat# 163-2014510% formaldehyde neutral buffer solutionWakoCat#062-01661MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinWako	2-Mercaptoethanol	Life Technologies	Cat#21985-023
phosphate buffer solution Wako Cat#062-01661 Methanol Wako Cat#131-01826 Hoechst 33342 Wako Cat#346-07951 TrypLE Express Life Technologies Cat#12605-010 Mitomycin C (Medicine) Kyowa Hakko Kirin N/A	Monothioglycerol		Cat# 195-15791
MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A	4% (wt/vol) Paraformaldehyde (PFA) phosphate buffer solution	Wako	Cat# 163-20145
MethanolWakoCat#131-01826Hoechst 33342WakoCat#346-07951TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A	10% formaldehyde neutral buffer solution	Wako	Cat#062-01661
TrypLE ExpressLife TechnologiesCat#12605-010Mitomycin C (Medicine)Kyowa Hakko KirinN/A		Wako	Cat#131-01826
Mitomycin C (Medicine) Kyowa Hakko Kirin N/A	Hoechst 33342	Wako	Cat#346-07951
Mitomycin C (Medicine) Kyowa Hakko Kirin N/A	TrypLE Express	Life Technologies	Cat#12605-010
· · ·	Mitomycin C (Medicine)	•	
	FBS	Life Technologies	Cat#12483-020

(Continued on next page)

iScience

Article



Continued		
REAGENT or RESOURCE	SOURCE	IDENTIFIER
3,3′,5-Triiodo-l-thyronine sodium salt	Sigma-Aldrich	Cat#T5516
Cholera toxin	List Biological Laboratory	Cat#100B
Insulin transferrin selenium	Life Technologies	Cat#41400045
Rhodamine B	Wako	Cat#180-00132
Tissue-Tek OCT compound	Sakura finetek Japan	Cat#4583
D-PBS	Wako	Cat#045-29795
TBS(Tris-Buffered Saline)powder	TaKaRa Bio	Cat#T903
Triton X-100	Sigma-Aldrich	Cat#T8787-100ML
Normal donkey serum	Jackson ImmunoResearch	Cat#017-000-121
Experimental models: cell lines		
Human: iPS cell line 201B7	RIKEN Bio Resource Cente	Cat#HPS0063; RRID:CVCL_A324
Mouse: 3T3-J2 cell line	H. Green, Harvard Medical School, Boston, MA	N/A
Oligonucleotides		
Taqman probe ABCB5 (Hs02889060_m1)	Life Technologies	Cat#4331182
Taqman probe KRT12 (Hs00165015_m1)	Life Technologies	Cat#4331182
Taqman probe PAX6 (Hs00240871_m1)	Life Technologies	Cat#4331182
Taqman probe ΔNp63 (Hs00978339_m1)	Life Technologies	Cat#4331182
Taqman probe GAPDH (Hs99999905_m1)	Life Technologies	Cat#4331182
Software and algorithms		
BD FACSDiva Software v8.0.3	BD Biosciences	RRID:SCR_001456
JMP software program v14.0.0	SAS institute	RRID:SCR_014242
ImageJ v1.51k	National Institutes of Health	RRID:SCR_003070

RESOURCE AVAILABILITY

Lead contact

Further information and requests for resources and reagents should be directed to, and will be fulfilled upon reasonable request, by the lead contact, Ryuhei Hayashi (ryuhei.hayashi@ophthal.med.osaka-u.ac.jp).

Materials availability

This study did not generate new unique reagents.

Data and code availability

The authors declare that all relevant data are available within the article and its Supplemental information files or from the corresponding author upon reasonable request. This study did not generate new code.

EXPERIMENTAL MODEL AND SUBJECT DETAILS

Cell lines

hiPSCs (cell line 201B7) were purchased from the RIKEN Bio Resource Center (Tsukuba, Japan) and cultivated in a feeder-free system in the presence of laminin-511 E8 fragments (LN511E8; Nippi, Tokyo, Japan) and serum-free medium (StemFit; Ajinomoto, Tokyo, Japan) (Hayashi et al., 2016, 2017; Takahashi et al., 2007; Nakagawa et al., 2014; Miyazaki et al., 2012).

METHOD DETAILS

Differentiation of hiPSCs using the SEAM method

The fabrication of SEAMs from hiPSCs was performed according to our previously reported method (Hayashi et al., 2016, 2017). Briefly, hiPSCs were cultivated for 12-14 weeks on LN511E8-coated ($0.5 \ \mu g/cm^2$) plates in serum-free differentiation medium (DM), corneal differentiation medium (CDM), and corneal epithelium maintenance medium (CEM) (Hayashi et al., 2016, 2017) (Figure 1A). DM consists of GMEM





(Life Technologies, Carlsbad, CA) supplemented with 10% knockout serum replacement (KSR; Life Technologies), 1 mM sodium pyruvate (Life Technologies), 0.1 mM non-essential amino acids (Life Technologies), 2 mM L-glutamine (Life Technologies), 1% penicillin-streptomycin solution (Life Technologies) and 55 μ M 2-mercaptoethanol (Life Technologies). CDM is diluted 1:1 in DM and Cnt-PR (without EGF or FGF2) (CELLnTEC Advanced Cell Systems, Bern, Switzerland) supplemented with 10 ng/mL keratinocyte growth factor (KGF) (Wako, Osaka, Japan), 10 μ M Y-27632 (Wako) and 1% penicillin-streptomycin (Life Technologies). CEM consists of DMEM:F-12 medium (1:1) (Life Technologies) supplemented with 2% B-27 (Thermo Fisher Scientific, Waltham, MA), 10 ng/mL KGF (Wako), 10 μ M Y-27632 (Wako), and 1% penicillin-streptomycin (Life Technologies).

Quantitative real-time reverse-transcriptase PCR

qRT-PCR was performed using the ABI Prism 7500 Fast Sequence Detection System (Life Technologies) in accordance with the manufacturer's instructions. Targets were amplified with TaqMan probe (Life Technologies) assays against *ABCB5* (Hs02889060_m1), *KRT12* (Hs00165015_m1), *PAX6* (Hs00240871_m1), $\Delta Np63$ (Hs00978339_m1), and *GAPDH* (Hs99999905_m1). The thermocycling program comprised an initial cycle at 95°C for 20 s, followed by 45 cycles at 95°C for 3 s and 60°C for 30 s.

Flow cytometry and cell sorting

Flow cytometry was performed according to a previously published method (Hayashi et al., 2018). Briefly, hiPSC-derived SEAMs were dissociated using TrypLE Express (Thermo Fisher Scientific) and re-suspended in ice-cold keratinocyte culture medium (KCM) (Hayashi et al., 2016, 2017, 2018) supplemented with 5% FBS (Japan Bio Serum, Hiroshima, Japan). Harvested cells were then stained with 5 μ g/1×10⁶ cells of mouse anti-ABCB5 mAb (clone 3C2-1D12) (Ksander et al., 2014; Wilson et al., 2014; Lutz et al., 2016; Schatton et al., 2015) for 30 min on ice. Following three washes with PBS they were then stained with 1 μ g/1×10⁶ cells of Alexa Fluor 647-conjugated secondary antibody (Thermo Fisher Scientific) for 30 min on ice. Cells were subsequently stained with 0.5 μ g/1×10⁶ cells of FITC-conjugated anti-CD104 (ITGB4; 58XB4, Biolegend, San Diego, CA), 0.375 μ g/1×10⁶ cells of FITC-conjugated anti-SSEA-4 (MC813-70, Biolegend), and 0.25 μ g/1×10⁶ cells of PE-Cy7-conjugated anti-CD200 (OX-104, BD Biosciences, San Diego, CA) antibodies for 1 hr on ice, followed by three washes with PBS. Cell sorting was performed with a FACSAria II instrument (BD Biosciences), and the data analyzed using BD FACSDiva Software v8.0.3 (BD Biosciences).

Colony-forming assay (CFA) and holoclone assay (HA)

For the CFA investigation ABCB5⁺/ABCB5⁻ SEAM-derived corneal epithelial cells that had been isolated by cell sorting were seeded in 6-well plates onto mitomycin C (MMC)-treated 3T3-J2 feeder cells (kindly provided by Prof H. Green, Harvard Medical School, Boston, MA) in corneal epithelium maturation medium (CMM; i.e.10% KCM medium containing 10% FBS, 10 ng/mL KGF and 10 μM Y-27632) at a density of 500 cells/well. Following cultivation for 8-10 days, cells were fixed with 10% formaldehyde neutral buffer solution (Wako) and stained with rhodamine B (Wako). Colony-forming efficiency was ascertained by counting colonies under a dissecting microscope, and colony diameters were measured using ImageJ v1.51k (National Institutes of Health, Bethesda, MD). For the HAs, ABCB5⁺/ABCB5⁻ SEAM-derived corneal epithelial cells were cultivated in 6-well plates on 3T3-J2 feeder cells in CMM at a density of 500 cells/well. After 10-13 days, single colonies were picked up under a dissecting microscope. After dissociating the colonies using TrypLE Express (Thermo Fisher Scientific), the dissociated cells were again seeded on 3T3-J2 feeder cells and cultivated in CMM for 8-9 days. Colonies were then fixed with 10% formaldehyde neutral buffer solution (Wako) and stained with rhodamine B (Wako), and viewed under a microscope to classified as holoclones or non-holoclones based on a modification of a previously described method (Hayashi et al., 2016; Barrandon and Green, 1987). See also Figure S1.

mAb-mediated ABCB5 blockade

Sorted ABCB5⁺ hiPSC-derived corneal epithelial cells were treated with 50 μ g/mL of anti-ABCB5 blocking mAb (clone 3C2-1D12) (Ksander et al., 2014; Wilson et al., 2014; Lutz et al., 2016; Schatton et al., 2015) for 10 consecutive days from the initial day of sorting. The mAb-mediated ABCB5 blockade was followed by CFAs and HAs for the treated ABCB5⁺ hiPSC-derived corneal epithelial cells (Figure 3A).

iScience Article



Fabrication of ex vivo expanded cell sheets

Fabrication of corneal epithelial cell sheets was accomplished according to previous reports (Hayashi et al., 2016, 2017). Briefly, ABCB5⁺/ABCB5⁻ SEAM-derived corneal epithelial cells that had been isolated by cell sorting were seeded on LN511E8-coated dishes at a density of 3,000-9,600 cells/well (96-well plates or cell culture inserts) and cultivated in CEM for 21–28 days.

Hematoxylin and eosin (H&E) staining

Corneal epithelial cell sheets generated from hiPSC-derived ABCB5⁺/ABCB5⁻ cells were embedded in OCT compound and frozen at -80° C. Thin (10 μ m) sections were cut, fixed with 10% formaldehyde neutral buffer solution (Wako), washed with distilled water and stained with hematoxylin and eosin (H&E). Following deparaffinization and hydration, they were imaged on a NanoZoomer-XR C12000 (Hamamatsu Photonics, Hamamatsu, Japan) microscope. The thickness of the H&E stained cell sheets was measured using ImageJ v1.51k (National Institutes of Health) at ten equally spaced positions in individual sheets and the results were averaged (Figure 4A).

Immunofluorescence staining

Immunofluorescence staining of ABCB5 and KRT12 was performed on unfixed material, whereas staining of PAX6 and p63 was performed on tissue that had been fixed in cold methanol (Wako) and washed, three times for 10 min each, with Tris-buffered saline (TBS, TaKaRa Bio, Shiga, Japan). Staining of KRT13 and MITF was conducted on material that had been fixed in 4% paraformaldehyde (PFA, Wako) and washed with TBS three times for 10 min each time. All sections were incubated in TBS containing 5% donkey serum and 0.3% Triton X-100 (Sigma-Aldrich, St. Louis, MO) for 1 hr to block non-specific reactions. Subsequently they were incubated overnight at 4°C with primary antibodies against ABCB5 (3C2-1D12) (2.5 μ g/mL; Ksander et al., 2014; Wilson et al., 2014; Lutz et al., 2016; Schatton et al., 2015), KRT12 (1:200; N-16; Santa Cruz Biotechnology, Santa Cruz, CA), PAX6 (1:1200; PRB-278P; Covance Research Products, Denver, PA), p63 (1:200; 4A4; Santa Cruz Biotechnology), KRT13 (1:200; AE8; Abcam, Cambridge, MA), and MITF (1:1000; C5; Exalpha Biologicals, Watertown, MA). After three more washes with TBS for 10 min each, samples were incubated with AF488-conjugated secondary antibodies (1:200; Life Technologies) for 1 hr at room temperature. Counterstaining was performed with Hoechst 33,342 (1:200; Life Technologies) prior to confocal fluorescence microscopy (FV3000, Olympus Corporation, Tokyo, Japan).

QUANTIFICATION AND STATISTICAL ANALYSIS

Statistical analysis

The data are expressed as means \pm standard deviation (s.d.). The statistical analyses were performed using the Wilcoxon signed-rank test. All of the statistical analyses were performed using the JMP software program v14.0.0 (SAS institute Inc., Cary, NC). No statistical methods were used to predetermine sample size.