Stem cell-based 3D brain organoids for mimicking, investigating, and challenging Alzheimer's diseases

Federica Cordella, Carlo Brighi, Alessandro Soloperto^{*}, Silvia Di Angelantonio^{*}

Alzheimer's disease (AD) is a progressive and irreversible neurodegenerative disorder that causes a decline of cognitive functions and a deterioration of behavioral and social performances. According to current estimations, AD is considered the prevalent cause of dementia, accounting for 60% and 80% of cases every year. The etiology of the disease is related to both genetic and environmental factors, making AD a cultifactorial disorder, with a major classification as late-onset AD and earlyonset AD.

AD's clinical manifestations, including memory impairment and altered cognitive functions, are caused by synaptic and neuronal loss together with a prominent inflammatory state in the entorhinal cortex, hippocampus, and neocortex (Livingston et al., 2020).

Despite the large number of theories that attempt to explain AD-related pathomorphological changes in the brain, evidences coming from *in vitro* and *in vivo* AD models have highlighted a cytotoxic effect of amyloid β (A β) fragments and hyperphosphorylated TAU on neurotransmission, axonal transport, signaling cascades and immune response, leading to synaptic loss and dysfunction (Rajmohan and Reddy, 2017).

Mutations in the amyloid precursor protein (APP) or in its processing enzyme, γ -secretase (subunit-PSEN1 (Presenilin 1) and subunit-PSEN2 (Presenilin 2)) correlates with the familial form of AD, and lead to increased β -site cleavage of APP and production of longer, aggregationprone variants of A β peptides. The major constituents of amyloid plaques are two abnormally folded products of APP metabolism: A β_{40} and A β_{42} , which show a higher rate of oligomerization and insolubility.

Another distinctive AD hallmark is represented by the presence of neurofibrillary tangles. These structures are primarily composed of paired helical filaments of hyperphosphorylated TAU proteins, sometimes mixed with straight filaments (Pizzarelli et al., 2020). TAU protein, encoded by the MAPT gene, represents the major microtubule associated protein (MAP) of a normal adult brain. This protein interacts with tubulin enhancing its polymerization into microtubules and helping their stability. The alternative splicing of TAU pre-mRNA originates 6 different isoforms: these proteins differ in containing three (3R) or four (4R) microtubule binding repeats in the carboxyl terminal and one (1N), two (2N) or zero (ON) amino terminal inserts. Thus, the alternative splicing of TAU is developmentally regulated, resulting in the production of two major groups: the 3R TAUs and 4R TAUs. The 2N4R TAU represents the most abundant TAU isoform inside the adult brain, while ON3R is only expressed in the fetal human brain.

Like the other MAP proteins, TAU is recognized as a phosphoprotein and its biological activity is regulated by its phosphorylation level. The presence of several missense mutations on the MAPT gene, such as G272V, P301L, V337M, and R406W, leads to TAU hyperphosphorylation and subsequent fibrillization. Moreover, phospho-TAU per se inhibits the assembly of tubulin into microtubules instead of favoring its stability and promotes the sequestration of both healthy TAU and the other MAP proteins (Chang et al., 2008). Neuropathological study showed that in the AD brains, Tau spreads among several brain areas, possibly transferred by the neuronal network (Braak and Braak, 1991). This phenomenon, known as Tau propagation, can be considered a relevant starting point to understand the molecular mechanisms of AD progression, helping to identify new targets and new therapeutic strategies. Indeed, Tau misfolding is associated not only with AD but also with other neurodegenerative pathologies known as Tauopathies, such as fronttemporal dementia or chronic traumatic encephalopathy.

Stem cell-based strategy for investigating neurodegeneration *in vitro*: Given the late phenotypic onset of AD and its poor prognosis, there is an urgent need to develop novel in vivo and in vitro models that better recapitulate human AD, to dissect molecular pathological mechanisms, identify novel therapeutic targets, and test drug candidates. Over the past few years, a combination of transgenic animal models has been employed to elucidate the causative link between structural and functional brain changes and genetic mutations characterizing AD. Although several mouse models are able to recapitulate some AD key hallmarks, such as protein aggregation and the timing by which the pathology takes place, most of them recapitulate only the A β plaques phenotype by overexpressing human genes commonly associated with familial forms of AD, without taking in account Tau pathology. However, Alzheimer's condition is defined as the concomitant manifestation of AB plaques deposition and neurofibrillary tau tangles which together exacerbates a progressive neuroinflammatory state and leads to a massive neuronal loss.

In this scenario, induced pluripotent stem cells (iPSCs) represent an important bridge to overcome these gaps, fostering the establishment of an in vitro human model that might display both AB and Tau phenotype with a relevant timeframe as that observed in humans. The discovery that terminally differentiated cells, like fibroblast, can be reverted into a pluripotent condition with an exogenous cocktail of transcriptional factors such as OCT4, SOX2, c-MYC and KLF4 (Takahashi and Yamanaka, 2006), has been a groundbreaking leap for the development of human-based in vitro disease models. Indeed, the improvement in human somatic cells reprogramming into iPSCs, and in protocols for iPSCs differentiation into brain cells is paving the road for a deeper analysis of development and functions of the human brain in physiological and pathological conditions. Moreover, the availability of patientderived iPSCs, displaying the genome and the molecular phenotype of the affected individuals, allows to go beyond the limitation of transgenic mouse models, and to develop a more faithful disease model to investigate the pathology at a cellular level.

Patient-specific iPSC-derived neurons recapitulate some pathological features of the Alzheimer pathology, such as the presence of soluble A β accumulation or insoluble A β aggregation, key hallmarks of the AD early stages. However, in

2D cultures, even harboring the most aggressive mutations of familiar AD, only low levels of AB species have been detected, without remarkable extracellular β-amyloid aggregation, probably because of the lack of interstitial compartment. Despite the major improvements made over the past years in 2D cell culture techniques, they still remain oversimplified platforms presenting many limitations including the absence of the in vivo-like cytoarchitectural organization and the brain's spatial synaptic connections, which instead are necessary to faithfully replicate Alzheimer pathophysiology (Liu et al., 2018).

Three-dimensional revolution in cell

biology: The number of innovative, clinically relevant constructs for disease modeling is rapidly growing with a variety of 3D *in vitro* platforms, including organoid cultures, 3D bioprinting, and 3D cellular scaffolding approaches. These advanced technologies have been successfully applied to replicate *in vitro* the pathogenesis of several human neurodegenerative diseases such as AD, Parkinson's disease, and amyotrophic lateral sclerosis, contributing to dissecting the molecular pathways involved in the disease.

Application of iPSC-derived neuronal and glial cells in the generation of brain organoids represent a step forward in order to fill the gap between 2D cultures and animal models (Brighi et al., 2020). Using these constructs, defined by selforganized and self-patterning features, it is possible to generate a complex 3D environment, suitable for the investigation of pathologies characterized by the presence of protein aggregates, such as AD. In current protocols, inhibitors of the dual SMAD signaling are used to enhance neural fate choice, mimicking the complex signaling conditions which characterize the embryonic state. Specifically, after initial neural commitment, the anterior/ dorsalization of the neuroepithelium is promoted by high levels of WNT signaling. This process recapitulates the in vivo corticogenesis through which glia and neurons of the 6 cortical layers are generated (Sloan et al., 2018). The heterogeneous but developmentally regulated differentiation of iPSCs into brain cells within 3D structures allow their direct interaction with the extracellular matrix, the formation of cell-to-cell connections and a more consistent gradients of neurogenic growth factors, leading eventually to complex and functional networks (Yakoub, 2019). Notably, despite in 2D cultures, the proteins secreted by the cells, like AB peptides, are removed every time the medium is changed, thus preventing protein aggregation, in organoids, the presence of a 3D extracellular matrix promotes the interstitial accumulation in microdomains of extracellularly released proteins. It has been speculated that iPSC-derived brain organoids with an AD genetic background may recapitulate the extracellular protein aggregation and the neuronal degeneration evidenced in human post-mortem brain tissues. In this line, patient-derived cortical organoids harboring presenilin-1 M146V mutation showed Aβ42 overexpression and altered matrix remodeling protein expression compared to healthy control derived organoids (Yan et al., 2018). Furthermore, cortical organoids generated from familial AD patient-derived iPSCs spontaneously displayed pathological AD characteristics, including accumulation of structures highly comparable with amyloid plaques and neurofibrillary tangles (Gonzalez et al., 2018).

Even though brain organoids generation represent an outstanding technology, attention is required when applied in neurodegenerative diseases research. As a matter of fact, despite cortical organoids revealed the presence of astrocytes and oligodendrocytes together with an patterned neuronal population, the absence of microglia and a vascular network, both derived from the mesodermal lineage, still represents a major drawback that needs to be rectified. Microglia, the brain's resident macrophages, express both neuroprotective and immunomodulatory roles and play important functions during neurodevelopment by modulating neuronal connection through synaptic pruning, as well as during adulthood, patrolling brain parenchyma, triggering prompt pro- or anti-inflammatory responses. Given that, the establishment of an organoid model that includes the presence of this cell subtype is paramount for interfacing correctly with the disease. However, most of the current in vitro 3D brain models lack microglia, thus being unable to provide an accurate framework for elucidating the neurodegenerative mechanisms implicated in the disease manifestation and in the pathologyspreading from cell to cell throughout the entire brain.

Thus, modelling AD in a dish still represents a big challenge, due to the complex etiology of the disease. A promising step to increase our knowledge, to identify druggable targets, and to foster drug development will be the integration in cortical organoids of iPSC derived immune and vascular cells harboring AD related mutations. In this scenario, the opportunity to exploit genome editing technologies such as CRISPR-CAS9, zinc finger nucleases, and transcription activator-like effector nucleases for genetically engineer healthy iPSCs with a panel of AD-relevant mutations might reveal new biochemical pathways involved in the pathogenesis of the disease and offer novel therapeutic strategies in tackling the pathology (Figure 1).

For example, human iPSC engineered ad hoc to express TAU-P301L mutation could be used as bricks to build up specific brain organoids that better characterize the development of the TAU-associated neurodegenerative processes, looking at the cause-effect temporal relationship between A_β plaques, neurofibrillary tangles and synaptic degeneration. Indeed, TAU-P301L mutation is frequently used in animal models to recapitulate pathological features of human AD, frontotemporal dementia, and parkinsonism. This mutation represents the most common and more aggressive MAPT mutation leading to the formation of twisted ribbon-like filaments of 4R-TAU isoforms in the temporal cortex, hippocampus, and substantia nigra, and to a functional TAU impairment in promoting microtubules assembly. In this manner, a more compliant framework could be generated by dissecting all cellular and molecular events occurring in different cell types in response to specific genetic alteration, posing, from individual evidence, valuable resources to integrate key molecular elements of AD in a stagedependent manner. Furthermore, the use of iPSC derived neurons expressing genetically encoded calcium, chloride, or voltage sensors (as GCaMP, Clop-Hensor or Archon1) will enable a functional high throughput analysis of synaptic and neuronal degeneration induced by TAU misfolding.

In a more holistic approach, to better model AD using an *in vitro* 3D system, efforts should be made to implement the organoid model considering the interplay between neuron, astrocytes, microglia, and the neurovascular system, the so called neurovascular unit, that *in vivo* ensures

Perspective




Figure 1 | Potential applications of cortical organoids for fighting Alzheimer's disease.

Three-dimensional cortical organoids represent a powerful tool to study tissue architecture and matrix remodeling in response to extracellular protein aggregation with better mimicking brain tissue respect to two-dimensional cultures. High throughput transcriptomics may foster the discovery of novel pathways involved in the disease, shedding light on human-specific phenotypic features otherwise not revealed by humanized animal models. Current advances in culturing and differentiating induced pluripotent stem cells (iPSCs) allow to integrate cortical organoids with cellular types not originating from the neuroectodermal line, such as microglia, paving the road to high-throughput study of the neuro-immune profiles. Lastly, new fast-scanning and high-performance optical approaches allow fine dissection of the functional communications among cells within cortical organoids, unveiling changes in neural network dynamics during the Alzheimer's disease associated neurodegenerative process.

the correct formation and function and the proper immune response of the brain. This can be achieved combining brain organoid with iPSC derived vascular structures and microglial cells though the hematopoietic differentiation, in microfluidic devices (multi-organ-on chip) or in bio-printed constructs. The combination of cells, with the same genetic background with or without TAU mutation will provide a powerful tool to identify the causative cascade of AD pathology.

In conclusion, human-derived brain organoids might represent an innovative in vitro model which considers for the very first time the whole complexity of the human brain otherwise not fully replicated within 2D cultures or animal models. iPSC-based organoids provide a unique technology for deepening our understanding of the neuron-gliamicroglia interaction and crosstalk in both physiological and pathological conditions. In this regard, 3D brain organoids might help to determine how different brain cells play key roles in the establishment of neurodegenerative mechanisms that underlie AD and other dementia-like phenotypes, expanding our knowledge in their pathogenesis and progression. Lastly, the use of genome editing approaches to introduce AD mutations and/or patientderived iPSC could reveal new molecular targets paving the way to establish more

compliant therapeutic treatments likely reducing the use of animal models in drug screening.

Alessandro Soloperto and Silvia Di Angelantonio contributed equally to this article.

Federica Cordella, Carlo Brighi, Alessandro Soloperto^{*}, Silvia Di Angelantonio

Center for Life Nanoscience, Istituto Italiano di Tecnologia, Rome, Italy (Cordella F, Brighi C, Soloperto A, Di Angelantonio S) Department of Physiology and Pharmacology, Sapienza University, Rome, Italy (Cordella F, Brighi C, Di Angelantonio S)

*Correspondence to: Alessandro Soloperto, PhD, alessandro.soloperto@iit.it; Silvia Di Angelantonio, PhD, silvia.diangelantonio@uniroma1.it. https://orcid.org/0000-0002-8142-6059 (Alessandro Soloperto) https://orcid.org/0000-0003-1434-3648 (Silvia Di Angelantonio) Date of submission: December 2, 2020 Date of decision: January 15, 2021 Date of acceptance: March 4, 2021

Date of acceptance: March 4, 2021 Date of web publication: July 8, 2021

https://doi.org/10.4103/1673-5374.317976

How to cite this article: Cordella F, Brighi C, Soloperto A, Di Angelantonio S (2022) Stem cell-based 3D brain organoids for mimicking, investigating, and challenging Alzheimer's diseases. Neural Regen Res 17(2):330-332. Copyright license agreement: The Copyright License Agreement has been signed by all authors before publication.

Plagiarism check: Checked twice by iThenticate.

Peer review: Externally peer reviewed. Open access statement: This is an open access journal, and articles are distributed under the terms of the Creative Commons Attribution-NonCommercial-ShareAlike 4.0 License, which allows others to remix, tweak, and build upon the work non-commercially, as long as appropriate credit is given and the new creations are licensed under the identical terms.

Open peer reviewer: Liang Yuan, Tufts University,

Additional file: Open peer review report 1.

References

- Braak H, Braak E (1991) Demonstration of amyloid deposits and neurofibrillary changes in whole brain sections. Brain Pathol 1:213-216.
- Brighi C, Cordella F, Chiriatti L, Soloperto A, Di Angelantonio S (2020) Retinal and brain organoids: bridging the gap between in vivo physiology and in vitro micro-physiology for the study of Alzheimer's diseases. Front Neurosci 14:655.
- Chang E, Kim S, Yin H, Nagaraja HN, Kuret J (2008) Pathogenic missense MAPT mutations differentially modulate tau aggregation propensity at nucleation and extension steps. J Neurochem 107:1113-1123.
- Gonzalez C, Armijo E, Bravo-Alegria J, Becerra-Calixto A, Mays CE, Soto C (2018) Modeling amyloid beta and tau pathology in human cerebral organoids. Mol Psychiatry 23:2363-2374.
- Liu C, Oikonomopoulos A, Sayed N, Wu JC (2018) Modeling human diseases with induced pluripotent stem cells: From 2D to 3D and beyond. Development 145:dev156166.
- Livingston G, Huntley J, Sommerlad A, Ames D, Ballard C, Banerjee S, Brayne C, Burns A, Cohen-Mansfield J, Cooper C, Costafreda SG, Dias A, Fox N, Gitlin LN, Howard R, Kales HC, Kivimäki M, Larson EB, Ogunniyi A, Orgeta V, et al (2020) Dementia prevention, intervention, and care: 2020 report of the Lancet Commission. Lancet 396:413-446.
- Pizzarelli R, Pediconi N, Di Angelantonio S (2020) Molecular imaging of Tau protein: new insights and future directions. Front Mol Neurosci 13:586169.
- Rajmohan R, Reddy PH (2017) Amyloid-beta and phosphorylated Tau accumulations cause abnormalities at synapses of Alzheimer's disease neurons. J Alzheimers Dis 57:975-999.
- Sloan SA, Andersen J, Paşca AM, Birey F, Paşca SP (2018) Generation and assembly of human brain region–specific three-dimensional cultures. Nat Protoc 13:2062-2085.
- Takahashi K, Yamanaka S (2006) Induction of pluripotent stem cells from mouse embryonic and adult fibroblast cultures by defined factors. Cell 126:663-76.
- Yakoub AM (2019) Cerebral organoids exhibit mature neurons and astrocytes and recapitulate electrophysiological activity of the human brain. Neural Regen Res 14:757-761.
- Yan Y, Song L, Bejoy J, Zhao J, Kanekiyo T, Bu G, Zhou Y, Li Y (2018) Modeling neurodegenerative microenvironment using cortical organoids derived from human stem cells. Tissue Eng Part A 24:1125-1137.

P-Reviewer: Yuan L; C-Editors: Zhao M, Qiu Y; T-Editor: Jia Y