

Complete Genome Sequence of a Thermophilic Hydrogenotrophic Methanogen, *Methanothermobacter* sp. Strain CaT2

Tomoyuki Kosaka,^{a,b} Hidehiro Toh,^c Atsushi Toyoda^d

Department of Biological Chemistry, Faculty of Agriculture, Yamaguchi University, Yamaguchi, Japan^a; PRESTO, Japan Science and Technology Agency (JST), Kawaguchi, Saitama, Japan^b; Medical Institute of Bioregulation, Kyushu University, Higashi-ku, Fukuoka, Japan^c; Center for Information Biology, National Institute of Genetics, Mishima, Shizuoka, Japan^d

We isolated a thermophilic hydrogenotrophic methanogen, *Methanothermobacter* sp. strain CaT2, which is able to aggregate and utilize formate. Here, we report the complete genome sequence of this organism.

Received 29 July 2013 Accepted 31 July 2013 Published 29 August 2013

Citation Kosaka T, Toh H, Toyoda A. 2013. Complete genome sequence of a thermophilic hydrogenotrophic methanogen, *Methanothermobacter* sp. strain CaT2. *Genome Announc.* 1(4):e00672-13. doi:10.1128/genomeA.00672-13.

Copyright © 2013 Kosaka et al. This is an open-access article distributed under the terms of the [Creative Commons Attribution 3.0 Unported license](https://creativecommons.org/licenses/by/3.0/).

Address correspondence to Tomoyuki Kosaka, tkosaka@yamaguchi-u.ac.jp.

Hydrogenotrophic methanogens, which are obligate autotrophs, grow on H₂ and CO₂, some of which utilize formate, and are ubiquitous in a number of anaerobic environments (1). *Methanothermobacter* species are thermophilic hydrogenotrophic methanogens, occur largely in thermophilic methanogenic environments, and grow in a high temperature range, from 40 to 70°C, and in a pH range from 6.0 to 8.0 (1). Here, we report the complete genome sequence of a thermophilic and hydrogen- and formate-utilizing methanogen, *Methanothermobacter* sp. strain CaT2, isolated in our laboratory. CaT2 shows self-aggregation. CaT2 has been deposited in Deutsche Sammlung von Mikroorganismen und Zellkulturen GmbH (DSMZ) and NITE Biological Resource Center (NBRC) under the accession no. DSM 24414 and NBRC 107770, respectively.

The genome sequencing of CaT2 was performed by a 454 Life Sciences GS-FLX sequencer (Roche). Sequence reads were assembled with Newbler. A fosmid library (768 clones) was constructed. End sequencing of this fosmid library was performed by ABI 3730 capillary sequencers (Applied Biosystems). The 454 contig data and the fosmid end sequence data were assembled using the Phred-Phrap-Consed systems (2–4). Gap closing and resequencing of low-quality regions in the assembly data were performed by the nested deletion method (5), PCR, primer walking, use of shattered insert libraries (6), and direct sequencing of fosmid clones.

The genome sequence of CaT2 consists of a circular chromosome of 1,720,003 bp and an 11,015-bp plasmid. The plasmid, pCaT2, is nearly identical to the plasmid pFZ1 (11,014 bp, accession no. X68367) of *Methanothermobacter thermotrophicus* strain Z-245 (7). The chromosome and the plasmid contain 1,749 and 11 predicted protein-coding genes, respectively. We compared the genome of CaT2 with those of the nonaggregating hydrogenotrophic methanogens *M. thermotrophicus* ΔH (8) and *Methanothermobacter marburgensis* Marburg (9). Genome alignments among the CaT2, ΔH, and Marburg strains. Sixty-seven genes in CaT2 were absent in both the ΔH and Marburg genomes, 35 (52%) of which encode hypothetical or conserved hypothetical

proteins. The 67 genes include the *fdhABC* genes encoding formate dehydrogenase (MTCT_1438 and MTCT_1439) and formate transporter (MTCT_1440), which were not conserved in the ΔH and Marburg genomes; this is supported by a previous report that ΔH and Marburg are not able to grow on formate (1).

Nucleotide sequence accession numbers. The sequence data for the CaT2 genome have been deposited in DDBJ/GenBank/EMBL under the accession no. [AP011952](https://www.ncbi.nlm.nih.gov/nuclink/AB011952) (chromosome) and [AP011953](https://www.ncbi.nlm.nih.gov/nuclink/AB011953) (plasmid).

ACKNOWLEDGMENTS

We thank the technical staff of the Comparative Genomics Laboratory in the National Institute of Genetics.

This research was supported by the PRESTO project, Japan Science and Technology Agency (JST).

REFERENCES

1. Wasserfallen A, Nöling J, Pfister P, Reeve J, Conway de Macario E. 2000. Phylogenetic analysis of 18 thermophilic *Methanobacterium* isolates supports the proposals to create a new genus, *Methanothermobacter* gen. nov., and to reclassify several isolates in three species, *Methanothermobacter thermotrophicus* comb. nov., *Methanothermobacter wolfeii* comb. nov., and *Methanothermobacter marburgensis* sp. nov. *Int. J. Syst. Evol. Microbiol.* 50(Pt 1):43–53. doi:10.1099/00207713-50-1-43.
2. Ewing B, Hillier L, Wendl MC, Green P. 1998. Base-calling of automated sequencer traces using Phred. I. Accuracy assessment. *Genome Res.* 8:175–185.
3. Ewing B, Green P. 1998. Base-calling of automated sequencer traces using Phred. II. Error probabilities. *Genome Res.* 8:186–194.
4. Gordon D, Abajian C, Green P. 1998. Consed: a graphical tool for sequence finishing. *Genome Res.* 8:195–202.
5. Hattori M, Tsukahara F, Furuhashi Y, Tanahashi H, Hirose M, Saito M, Tsukuni S, Sakaki Y. 1997. A novel method for making nested deletions and its application for sequencing of a 300 kb region of human APP locus. *Nucleic Acids Res.* 25:1802–1808.
6. McMurray AA, Sulston JE, Quail MA. 1998. Short-insert libraries as a method of problem solving in genome sequencing. *Genome Res.* 8:562–566.
7. Nöling J, van Eeden FJ, Eggen RI, de Vos WM. 1992. Modular organization of related archaeal plasmids encoding different restriction-modification systems in *Methanobacterium thermoformicum*. *Nucleic Acids Res.* 20:6501–6507.

8. Smith DR, Doucette-Stamm LA, Deloughery C, Lee H, Dubois J, Al-dredge T, Bashirzadeh R, Blakely D, Cook R, Gilbert K, Harrison D, Hoang L, Keagle P, Lumm W, Pothier B, Qiu D, Spadafora R, Vicaire R, Wang Y, Wierzbowski J, Gibson R, Jiwani N, Caruso A, Bush D, Safer H, Patwell D, Prabhakar S, McDougall S, Shimer G, Goyal A, Pietrokovski S, Church GM, Daniels CJ, Mao J-I, Rice P, Nölling J, Reeve JN. 1997. Complete genome sequence of *Methanobacterium thermoautotrophicum* deltaH: functional analysis and comparative genomics. *J. Bacteriol.* 179: 7135–7155.
9. Kaster AK, Goenrich M, Seedorf H, Liesegang H, Wollherr A, Gottschalk G, Thauer RK. 2011. More than 200 genes required for methane formation from H₂ and CO₂ and energy conservation are present in *Methanothermobacter marburgensis* and *Methanothermobacter thermoautotrophicus*. *Archaea* 2011:973848. doi:[10.1155/2011/973848](https://doi.org/10.1155/2011/973848).