Original Article

Radiat Oncol J 2015;33(1):42-49 http://dx.doi.org/10.3857/roj.2015.33.1.42 pISSN 2234-1900 · eISSN 2234-3156



New methods for optical distance indicator and gantry angle quality control tests in medical linear accelerators: image processing by using a 3D phantom

Mahdi Heravian Shandiz, MSc¹, Ghorban Safaeian Layen, PhD², Kazem Anvari, PhD³, Mohammadmahdi Khalilzadeh, PhD¹

¹Department of Biomedical Engineering, College of Engineering, Mashhad Branch, Islamic Azad University, Mashhad, Iran; ²Department of Radiologic Technology, School of Paramedical Science, Mashhad University of Medical Science, Mashhad, Iran; ³Cancer Research Center, School of Medicine, Mashhad University of Medical Sciences, Mashhad, Iran

Purpose: In order to keep the acceptable level of the radiation oncology linear accelerators, it is necessary to apply a reliable quality assurance (OA) program.

Materials and Methods: The QA protocols, published by authoritative organizations, such as the American Association of Physicists in Medicine (AAPM), determine the quality control (QC) tests which should be performed on the medical linear accelerators and the threshold levels for each test. The purpose of this study is to increase the accuracy and precision of the selected QC tests in order to increase the quality of treatment and also increase the speed of the tests to convince the crowded centers to start a reliable QA program. A new method has been developed for two of the QC tests; optical distance indicator (QDI) QC test as a daily test and gantry angle QC test as a monthly test. This method uses an image processing approach utilizing the snapshots taken by the CCD camera to measure the source to surface distance (SSD) and gantry angle.

Results: The new method of ODI QC test has an accuracy of 99.95% with a standard deviation of 0.061 cm and the new method for gantry angle QC has a precision of 0.43°. The automated proposed method which is used for both ODI and gantry angle QC tests, contains highly accurate and precise results which are objective and the human-caused errors have no effect on the results. **Conclusion:** The results show that they are in the acceptable range for both of the QC tests, according to AAPM task group 142.

Keywords: Gantry angle, Linear accelerator, Optical distance indicator (ODI), Quality control, Source to surface distance (SSD)

Introduction

Quality control (QC) tests, recommended by institutions, such as the American Association of Physicists in Medicine (AAPM), are one of the important parts of the radiation oncology department's overall quality assurance (QA) program [1]. It is fundamentally important that the prescribed dose delivers

accurately and precisely to the target volume. Geometric and mechanical accuracy of the medical linear accelerators (LINACs) should be tested and verified on a regular basis [2].

The source to surface distance (SSD) QC test of a LINAC is measured by using an optical device for measuring the SSD called optical distance indicator (ODI) which needs to have a daily QC test [3,4]. ODI is a device which shows the SSD on

Received 14 August 2014, Revised 26 October 2014, Accepted 22 January 2015.

Correspondence: Ghorban Safaeian Layen, PhD, Department of Radiologic Technology, School of Paramedical Science, Mashhad University of Medical Science, Mashhad, Iran. Tel: +98-51-38-846716, Fax: +98-51-36070058, E-mail: safaeiangh@mums.ac.ir

This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (http://creativecommons.org/licenses/by-nc/3.0/) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

www.e-roj.org



the surface using an optical graded ruler. In other words, ODI shows the current SSD by projecting the number on the couch.

Currently a mechanical distance indicator is used to verify the ODI [5]. Welsh et al. [2] used a tape measure with a digital screen to show the SSD for the ODI verification. Geyer [6] used a laser based tool to overcome the subjective human errors caused due to measuring the SSD by means of manual methods.

Usually spirit levels are used to perform the gantry angle QC test in the radiation oncology centers [7]. Welsh et al. [2] used a digital goniometer to find out the gantry angle. In order to improve the accuracy and precision of ODI QC test, as one of the important daily quality assurance tests, a new method has been developed which uses image processing.

The proposed method is also able to carry out the monthly gantry angle QC test as one of the other important mechanical QC tests without the need of any extra hardware.

The system is composed of a low-cost CCD camera, a simple phantom and software. The in-house software uses a distance measurement and angle detection method based on the number of pixel variation of images by referencing to four selected points in the image [8].

This method offers a quick setup and automatic measurement of the SSD and gantry angle. Since the values are automatically and directly saved into the database, any human-caused errors while reading the results could be eliminated. In the software, the real values are compared to the target values, and if the difference is acceptable by AAPM task group 142, the tests would be announced as it is passed.

The AAPM task group 40 and task group 142 recommend that up to 2 mm deviation between the real SSD and measured SSD and up to 1° between the real gantry angle and the measured angle are acceptable [3,4]

The purpose of this paper is to improve the accuracy and precision of ODI and gantry angle quality control tests by eliminating the error causes in other methods. In order to achieve this goal, the proposed method has been tested in two radiotherapy centers in Mashhad, Imam Reza (AS) Hospital and Reza (AS) Radiation Oncology Center.

Materials and Methods

1. System development

In order to measure the SSD and the gantry angle, our system includes a CCD camera, a 3D phantom and software. The SSD and gantry angle were calculated by using a pixel variation based algorithm. This system was established and tested on

two Linacs in two radiation oncology centers in Mashhad, a Siemens Primus in Reza (AS) radiation oncology center and an ELEKTA Precise in Imam Reza (AS) hospital of Mashhad.

1) CCD camera: A 2-megapixel CCD camera was used to take photos from different distances of the surface and different angles around the zero angle. A CCD camera was used because of the better fill factor, uniformity and system noise compared to CMOS cameras [9]. The detailed information and photo of the camera could be found in Table 1 and Fig. 1, respectively.

2) 3D phantom: The 3D phantom was made of six acrylic sheets ($180 \times 180 \times 0.5 \text{ mm}^3$) which create a cubic phantom together as illustrated in Fig. 2. Four circular fiducial marks, each of them with a diameter of 10 mm, has been engraved on the corners of an imaginary square with 150 mm sides on the surface of the 3D phantom. The image processing algorithm used these fiducial marks, as its references, to calculate the distance and angle.

3) Distance calculation algorithm: Hsu et al. [8] introduced

Table 1. Detailed information of the CCD camera

Sensor type	1/1.8" progressive scan CCD
Pixels (H×V)	$1,280 \times 1,024$
Pixel size (μm²)	4.65×4.65
Sensing area (mm²)	5.95×4.76
Pixel depth	8-bit
Frame rate (fps)	15



Fig. 1. CCD camera used for the proposed method.

a distance measurement algorithm based on pixel variations of images which was used for detecting source to surface distance in this study. Fig. 3 shows a schematic outline of the CCD camera capturing images at h_1 and h_2 heights between the lens and the surface in which h_s shows the gap between the CCD sensor and the lens, $D(h_1)$ and $D(h_2)$ are the maximum horizontal length in centimetres available in the photos at each height, formed by the field of view of the camera at photographing distances h_1 and h_2 respectively, I is the distance




Fig. 2. 3D phantom used to perform the quality control tests.

between two adjacent reference points and $N(h_1)$ and $N(h_2)$ are the pixel counts of / at the two photographing distances. According to the triangles similarity theorem, Eqs. (1)–(4) are extracted:

$$\frac{l}{N(h1)} = \frac{D(h1)}{Nmax} \tag{1}$$

$$\frac{l}{N(h2)} = \frac{D(h2)}{Nmax} \tag{2}$$

$$\frac{D(h1)}{D(h2)} = \frac{N(h2)}{N(h1)} \tag{3}$$

And finally:

$$\frac{N(h2)}{N(h1)} = \frac{D(h1)}{D(h2)} + \frac{h1 + hs}{h2 + hs} \tag{4}$$

Therefore, by having the height in one distance (h_1) which is 100 cm by default, any other heights (h_2) could be calculated by using Eq. (4).

BC in Fig. 4 shows the side view of the top side of the 3D phantom in h_1 which is the standard height. A and B indicate the lens and the CCD sensor of the camera, respectively, and B represents the X-ray source. Since the gantry angle is B0°, ABC1 is an isosceles triangle.

In an unknown height as demonstrated in Fig. 5, h_2 could be calculated in A'B'C' after figuring out A'B' and A'C' according to the Eqs. (5)–(9):




Fig. 3. Schematic outline of the CCD camera capturing images at h_1 and h_2 heights between the lens and the surface.





Fig. 4. Side view of the gantry, camera and the phantom at standard SSD = 100 cm. SSD, source to surface distance.



$$A'C' = \left(\frac{N4}{L'4}\right)(AC + hs) - hs \tag{6}$$

$$A'C'^{2} = A'B'^{2} + L^{2} - 2(A'B')(L)\cos B'$$
(7)

$$B' = \cos^{-1} \left(\frac{A'B'^2 + L^2 - A'C'^2}{2(A'B')(L)} \right)$$
 (8)

$$A'D' = h2 = \sqrt{A'B'^2 + \left(\frac{L}{2}\right)^2 - 2(A'B')\left(\frac{L}{2}\right)\cos B'}$$
 (9)

And finally the SSD which is composed of h_2 and the offset of the source to camera lens (XA or XA) distance could be calculated.

$$h_{2Real} = XA + h_2 \tag{10}$$

4) Angle detection algorithm: The gantry angle could be calculated at SSD = 100 cm in a wide range of about $\pm 150^{\circ}$. In Fig. 6, A is the camera lens at 0° , while A is the camera lens at an unknown angle. The unknown angle of the gantry is C_2 . A is calculated using the triangles theorem for ABC and since ABC is an isosceles triangle, C could also be figured out as shown in Eq. (11)–(13).

$$A = \cos^{-1} \left(\frac{AB^2 + AC^2 - L^2}{2(AB)(AC)} \right) \tag{11}$$

$$A + 2C = 180 (12)$$




Fig. 5. Side view of the gantry, camera and the phantom at an unknown source to surface distance.




Fig. 6. Side view of the camera lens and the phantom in 0° and an unknown degree for the gantry angle.

$$C = 90 - \frac{A}{2} \tag{13}$$

Since A C and AB distances are calculated by using the previous section algorithm, C_1 is also calculated by means of triangle theorem, so having C and C_1 , both available, the gantry angle (C_2) is figured out by using Eq. (14) and (15):

$$C1 = \cos^{-1} \left(\frac{A' C^2 + L^2 - A' B^2}{2(A' C)(L)} \right) \tag{14}$$

$$C_2 = C - C_1 \tag{15}$$

5) Computer software: The SSD and gantry angle calculating source codes are written in MATLAB 7.12. The steps which the software takes for SSD and gantry angle detection from the camera detection and snapshot to the reference points extraction and SSD and gantry angle calculation are




Fig. 7. Steps which the software takes to calculate source to surface distance and gantry angle.



Fig. 8. Quality assurance system setup with the camera on the gantry and the 3D phantom on the couch.

mentioned in Fig. 7.

2. System performance

The first step to start the QC tests is to place the camera board in the accessory tray of the gantry and connect the cable to the laptop. The in-house software, developed for this system, should be started after placing the 3D phantom on the couch. Fig. 8 indicates the CCD camera and the 3D phantom while performing the QC tests.

The position of the CCD camera should always be calibrated and fixed towards the 3D phantom. As it is not easy to have the CCD camera in the same position every time, we have added a software position alignment and calibration. As it is illustrated in Fig. 9, there is a cross on the screen, dividing it into four sections. These sections help to locate the phantom




Fig. 9. Camera calibration tab in the software. The millimetre numbers in front of the arrows demonstrate the shifts needed in the direction of the arrow.

approximately in the middle of the screen. In order to adjust the location of the phantom exactly in the middle, some arrows would appear on the screen, guiding the user to move the phantom towards X and Y directions. In other words, there are four arrows showing up, down, right, and left. If there is any misalignment, a number would appear in front of one or more arrows, indicating the millimetres that the phantom should be moved in that direction.

The system is not so sensitive to the environment light and it could work if there are some lights on in the room, but since the best result is acquired when the environment has a higher brightness, it is strongly recommended to keep all the lights turned on.

After adjusting the phantom in the correct place, 'SSD QC'



button on the main tab of the software would perform the SSD QC test automatically and would show the result while the 'gantry angle QC' button would execute the gantry angle QC test and display the conclusion.

In order to evaluate the accuracy and precision of the proposed method, the ODI QC test was performed at 95, 105, 110, and 115 cm for the SSD and the gantry angle QC test was carried out at 0° , 2° , 5° , 80° , 90° , and 150° .

Results

1. ODI QC test

The ODI QC test has been performed on an ELEKTA Precise and a SIEMENS Primus LINACs. The measurements for ODI QC test were made at 95, 105, 110 and 115 cm for the SSD and in all the heights the gantry angle was at the zero degree.

Table 2 shows the results of the ODI QC test for the ELEKTA Precise LINAC and compares them to the calibration rods measurement results as a reference method. The test had been repeated 5 times for each SSD to ponder the precision. Table 3 indicate the same tests for the SIEMENS Primus LINAC. Since the system is completely symmetric, the measurement and results of 0° to 180° are exactly the same as -180° to 0°; therefore, in order to keep the tables legible and as brief as possible, we have only evaluated 0° to 180°.

Tables 2 and 3 reveal that the accuracy of the proposed method is 99.96% with a standard deviation of 0.084.

2. Gantry angle QC test

The gantry angle QC test has been done on the same two LINACs at the constant 100 cm SSD and six different gantry angles, 0°, 2°, 5°, 80°, 90°, and 150°. As an ODI QC test, it has

Table 2. ODI quality control results for ELEKTA Precise comparing the proposed method to calibration rods method

SSD (ODI)	SSD (calibration rods)	Average SSD (proposed method)	Error SSD2-SSD3	Standard deviation (proposed method)	% Accuracy (proposed method)	p-value ^{a)}
95	95.01	95.04	0.03	0.071	99.97	0.344
105	105.01	105.04	0.02	0.033	99.98	0.152
110	110.01	109.99	0.02	0.043	99.98	0.450
115	115.01	114.96	0.05	0.049	99.96	0.077

ODI, optical distance indicator; SSD, source to surface distance.

Table 3. ODI quality control results for SIEMENS Primus comparing the proposed method to calibration rods method

SSD (ODI)	SSD (calibration rods)	Average SSD (proposed method)	Error SSD2-SSD3	Standard deviation (proposed method)	% Accuracy (proposed method)	p-value ^{a)}
95	95.02	95.04	0.02	0.084	99.98	0.560
105	105.02	105.05	0.03	0.045	99.97	0.171
110	110.02	109.99	0.03	0.061	99.97	0.363
115	115.02	114.98	0.04	0.049	99.97	0.126

ODI, optical distance indicator; SSD, source to surface distance.

Table 4. Gantry angle quality control results for ELEKTA Precise with 5 repeats for each angle

SSD (cm)	Gantry angle (°)	Average (°)	Accuracy (%)	Standard deviation	p-value ^{a)}
100	0	0.06	99.98	0.142	0.371
100	2	1.91	99.97	0.172	0.317
100	5	4.94	99.98	0.327	0.676
100	80	80.06	99.98	0.236	0.588
100	90	90.06	99.98	0.191	0.520
100	150	150.03	99.99	0.462	0.885

ODI, optical distance indicator; SSD, source to surface distance.

^{a)}Paired t-test (reference data – our data).

^{a)}Paired t-test (reference data – our data).

^{a)}Paired t-test (reference data – our data).

Table 5. Gantry angle quality control results for SIEMENS Primus with 5 repeats for each angle

SSD (cm)	Gantry angle (°)	Average (°)	Accuracy (%)	Standard deviation	p-value ^{a)}
100	0	0.01	99.99	0.039	0.670
100	2	2.17	99.95	0.308	0.275
100	5	5.10	99.97	0.430	0.644
100	80	79.97	99.99	0.201	0.755
100	90	89.91	99.97	0.240	0.469
100	150	149.89	99.97	0.291	0.460

ODI, optical distance indicator; SSD, source to surface distance.

been repeated 5 times for each angle to consider the precision. The results for the two LINACs are indicated in Tables 4 and 5.

According to Tables 4 and 5, the accuracy of the proposed method is 99.95% with a standard deviation of 0.462.

Discussion and Conclusion

The proposed method is used both for ODI QC and gantry angle QC tests.

The simple QA method has high accuracy and precision and consists of a portable camera inserted in the accessory tray and a simple phantom. The proposed method offers a simple way to replace the subjective traditional used methods for ODI and gantry angle QC tests. Currently most of the radiation oncology centers use subjective methods for performing the QC tests, so the accuracy and precision of these conventional and subjective methods cannot be measured, but the proposed methods are totally objective and have precise accuracy and precision which is 99.96% accuracy and 0.084 standard deviation for the SSD QC test and 99.95% accuracy and 0.462 standard deviation for the gantry angle QC test. One thing that makes radiation oncology centers reluctant to follow the complete QC tests is that the tests performed in a traditional way are time-consuming and the centers do not have enough time to spend on them. The objective result of the proposed method takes less than 2 minutes to fix and adjust the camera and phantom and only a click for performing each test, so the total consumed time would not exceed 3 minutes for both of the tests and there would be no problem such as lack of time left for these QC tests.

In conclusion, for the ODI QC test, the accuracy is 99.95% with a standard deviation of 0.061 cm, so it is completely acceptable by AAPM task group 142 which defines a tolerance level of 2 mm for this test [3]. The gantry angle QC test is performed with the same hardware of the ODI QC test and no mentionable part has been added except a few lines of MATLAB code. The result of this test has a precision of 0.43° which is in the acceptable range, defined by AAPM task group 142 that should be less than 1° [3].

Conflict of Interest

No potential conflict of interest relevant to this article was reported.


Acknowledgments

The authors would like to thank Shahram Bayani, Mohammad Hasan Laki, Mohammadreza Alagheband Hoseini, and Amirhosein Ziaie for their helpful tips and cooperation in the tests.

References

- 1. Welsh KT, Reinstein LE, Wlodarczyk RA. Automated image based quality assurance for medical linear accelerators. In: Schlegel W, Bortfeld T, editors. The use of computers in radiation therapy. Heidelberg: Springer; 2000. p. 368-70.
- 2. Welsh KT, Wlodarczyk RA, Reinstein LE. A new geometric and mechanical verification device for medical LINACs. J Appl Clin Med Phys 2002;3:154-61.
- 3. Klein EE, Hanley J, Bayouth J, et al. Task group 142 report: quality assurance of medical accelerators. Med Phys 2009;36: 4197-212.
- 4. Kutcher GJ, Coia L, Gillin M, et al. Comprehensive QA for radiation oncology: report of AAPM Radiation Therapy Committee Task Group 40. Med Phys 1994;21:581-618.
- 5. Mayles WP. Physics aspects of quality control in radiotherapy. New York: Institute of Physics and Engineering in Medicine; 1999.
- 6. Geyer P. Source-skin-distance measurements using an offaxis laser device. In: AAPM 45th Annual Meeting; 2003 Aug 10-14; San Diego, CA.
- 7. Chang L, Ho SY, Wu JM, Yu CY, Sung CC. Technical innovation

^{a)}Paired t-test (reference data – our data).





- to calibrate the gantry angle indicators of linear accelerators. J Appl Clin Med Phys 2001;2:54-8.
- 8. Hsu CC, Lu MC, Wang WY, Lu YY. Distance measurement based
- on pixel variation of CCD images. ISA Trans 2009;48:389-95.
- 9. Litwiller D. CCD vs. CMOS: facts and fiction. Photonics Spectra 2001;35:154-8.