



Original Article

Implementation of a single quad MS detector in routine QC analysis of peptide drugs[☆]Matthias D'Hondt, Bert Gevaert, Evelien Wynendaele, Bart De Spiegeleer^{*}

Drug Quality and Registration (DruQuaR) Group, Faculty of Pharmaceutical Sciences, Ghent University, Ottergemsesteenweg 460, 9000 Ghent, Belgium

ARTICLE INFO

Article history:

Received 24 November 2014

Received in revised form

21 September 2015

Accepted 22 September 2015

Available online 26 September 2015

Keywords:

Quality control

Peptide

Impurity profiling

Single quad MS detector

ABSTRACT

A newly developed single quad mass spectrometry (MS) detector was coupled to a ultra-high performance liquid chromatography (UPLC) system and implemented in the routine quality control (QC) and impurity analysis of four therapeutic peptides, namely bleomycin sulfate, tyrothricin, vancomycin HCl and bacitracin, which were selected given their multi-component drug nature and their closely structurally related impurity profiles. The QC and impurity profiling results obtained using the ultra-high performance liquid chromatography ultraviolet/mass spectrometry (UPLC-UV/MS) detection system were analyzed against the results obtained using traditional high performance liquid chromatography-ultraviolet detection (HPLC-UV) methods derived from pharmacopoeial methods. In general, the used stationary phases of sub-2 μm particle (UPLC) technology resulted in lower limits of detection and higher resolution separations, which resulted in more detected impurities and shorter overall run times contrasting the traditional HPLC columns. Moreover, online coupling with a single quad MS detector allowed direct peak identification of the main compounds as well as small impurities, hereby increasing the information content without the need of reference standards.

© 2015 Xi'an Jiaotong University. Production and hosting by Elsevier B.V. All rights reserved. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

1. Introduction

Over the past few decades, mass spectrometry (MS) has become an essential tool in fundamental disease research [1–4]. Mass spectrometry-based proteomics, metabolomics, peptidomics, glycomics, phosphoproteomics and lipidomics are used to quantitatively and qualitatively differentiate protein, peptide, sugar and lipid structures of patients and healthy volunteers, in order to identify new biomarkers or seek treatment opportunities [5]. Due to the complexity of these human samples, as well as the quest to obtain greater sensitivity, a wide variety of high-end mass spectrometry instruments have been developed in recent years. An excellent review written by Domon and Aebersold provides a performance overview of different types of MS used in proteomic field. For example, triple quadrupole instruments were used for quantification of proteins and ultra-high resolution Fourier transform ion cyclotron resonance (FT-ICR) for their identification [3]. MS has also been applied in drug development such as drug metabolism studies [6,7]. The factors that limit the widespread use of high-end MS instruments are high purchasing costs and expensive maintenance. Furthermore, highly skilled operators for method optimization and subsequent data interpretation are needed.

Factors like those mentioned above have limited the use of MS

in routine pharmaceutical quality control (QC) environments. The performance demands on MS instrumentation in routine pharmaceutical setting are less challenging than those in the “omics”-field, as matrix effects are less evident and drug concentrations are much higher. However, simplified operation and maintenance within a good manufacturing practice (GMP) environment are paramount, which leads to the need for downsized and lower-end MS systems. Recently, a number of small single quad MS detectors that are directly compatible with the existing high performance liquid chromatography (HPLC)/ultra-high performance liquid chromatography (UPLC) equipment and software platforms used in QC, have been developed and introduced into the pharmaceutical setting [8–10]. The availability of these new MS detectors has created the potential to bring routine use of MS to the QC laboratory and pharmaceutical QC to a higher level. The present work considered a new ultra-high performance liquid chromatography ultraviolet/mass spectrometry detection (UPLC-UV/MS) equipment set-up in relation to the traditional high performance liquid chromatography-ultraviolet detection (HPLC-UV) set-up for the QC and impurity profiling of complex therapeutic peptides.

2. Materials and methods

2.1. Materials

Tyrothricin, bacitracin, propionic acid, sodium pentane sulfate,

[☆]Peer review under responsibility of Xi'an Jiaotong University.

^{*}Corresponding author.

E-mail address: Bart.DeSpiegeleer@UGent.be (B. De Spiegeleer).

Table 1
Overview of selected peptides.

Name	Major compounds ^a	Formula	Molecular weight ^b	Sample solvent ^c	Ph. Eur. concentration ^c		
Bleomycin sulfate	Bleomycin A ₂	C ₅₅ H ₈₄ N ₁₇ O ₂₁ S ₃	1414.52	H ₂ O	Test solution: 0.5 mg/mL ^d		
	Bleomycin B ₂	C ₅₅ H ₈₄ N ₂₀ O ₂₁ S ₂	1424.56				
Tyrothricin	Gramicidin A ₁	C ₉₉ H ₁₄₀ N ₂₀ O ₁₇	1881.07	15 volumes H ₂ O	Test solution: 1.0 mg/mL		
	Gramicidin A ₂	C ₁₀₀ H ₁₄₂ N ₂₀ O ₁₇	1895.09	85 volumes MeOH	Ref. solution B: 0.02 mg/mL ^d		
	Gramicidin C ₁	C ₉₇ H ₁₃₉ N ₁₉ O ₁₈	1858.05				
	Gramicidin C ₂	C ₉₈ H ₁₄₁ N ₁₉ O ₁₈	1872.07				
	Tyrocidin A	C ₆₆ H ₈₈ N ₁₃ O ₁₃	1270.66				
	Tyrocidin B	C ₆₈ H ₈₉ N ₁₄ O ₁₃	1309.67				
	Tyrocidin C	C ₇₀ H ₉₀ N ₁₅ O ₁₃	1348.68				
Vancomycin HCl	Tyrocidin D	C ₇₂ H ₉₁ N ₁₆ O ₁₂	1371.70	10 volumes THF 70 volumes ACN 920 volumes H ₂ O 40 volumes ACN	Test solution A: 2.0 mg/mL Test solution B: 0.08 mg/mL Test solution C: 0.002 mg/mL ^d		
	Tyrocidin E	C ₆₆ H ₈₈ N ₁₃ O ₁₂	1254.67				
	Vancomycin B	C ₆₆ H ₇₅ Cl ₂ N ₉ O ₂₄	1447.43			400 volumes H ₂ O	Test solution: 2.0 mg/mL
	Bacitracin	Bacitracin A	C ₆₆ H ₁₀₃ N ₁₇ O ₁₆ S				
		Bacitracin B ₁	C ₆₅ H ₁₀₁ N ₁₇ O ₁₆ S			1407.73	520 volumes MeOH
	Bacitracin B ₂	C ₆₅ H ₁₀₁ N ₁₇ O ₁₆ S	1407.73				
	Bacitracin B ₃	C ₆₅ H ₁₀₁ N ₁₇ O ₁₆ S	1407.73				

^a As listed in Ph. Eur. monograph.^b Monoisotopic molecular weight.^c Based upon Ph. Eur. 8.0 monographs, buffer solution substituted with H₂O.^d For disregard limit purposes.**Table 2**
Chromatographic conditions of selected peptides (Ph. Eur.).

Compound	Mobile phase A	Mobile phase B	Run time (min)	Flow rate (mL/min)	Column temperature (°C)	Injection volume (μL)	Quant. wavelength (nm)
Bleomycin sulfate	MeOH	Pentanesulfonate buffer	100	1.2	Room	20	254
Tyrothricin	75 MeOH 25 sulfate buffer	–	60	1.2	60	25	280
Vancomycin HCl	70 ACN 10 THF	290 ACN 10 THF	35	1	Room	20	280
	920 triethylamine 40 ACN	700 triethylamine	75	1	Room	100	254
Bacitracin	300 H ₂ O 520 MeOH 100 phosphate buffer	–	75	1	Room	100	254

Table 3
Different gradient systems used for UPLC-PDA/MS analysis.

Sample	Time (min)	Mobile phase A (% v/v)	Mobile phase B (% v/v)
Bleomycin sulfate	0	100	0
	10	70	30
	12.50	70	30
	14	100	0
	20	100	0
Tyrothricin	0	60	40
	10.50	10	90
	11.86	10	90
	12.60	60	40
	16.80	60	40
Vancomycin HCl	0	100	0
	10	90	10
	20	70	30
	22	100	0
	30	100	0
Bacitracin	0	95	5
	10.50	50	50
	11.86	50	50
	12.60	95	5
	16.80	95	5

sodium sulfate and triethylamine were bought from Sigma-Aldrich (Diegem, Belgium), whereas vancomycin HCl was purchased from Bufa (Ijsselstein, the Netherlands). These peptides were ordered

without explicit pharmacopoeial quality requirements, in order to increase the probability of containing related impurities. Bleomycin sulfate was obtained from Sanofi Aventis (Brussel, Belgium). Acetonitrile, methanol and tetrahydrofuran of HPLC grade were acquired from Fisher Scientific (Aalst, Belgium), whereas acetonitrile, methanol and trifluoroacetic acid (TFA) of UPLC grade were bought from Biosolve (Valkenswaard, the Netherlands). Water was purified using an Arium 611 purification system (Sartorius, Göttingen, Germany) yielding $\geq 18.2 \text{ M}\Omega \times \text{cm}$ quality water. Potassium dihydrogen phosphate and ammonium sulfate were purchased from Merck (Overijse, Belgium), whereas dipotassium hydrogen phosphate was obtained from Panreac Chimica (Barcelona, Spain).

2.2. Peptide selection

Three different antimicrobial peptides, i.e., tyrothricin (cyclic polypeptide), vancomycin HCl (glycopeptide), bacitracin (cyclic polypeptide), and the chemotherapeutic peptide bleomycin sulfate (glycopeptide) were selected as test compounds, given their complex and multi-component drug nature and their closely structurally related impurity profiles. An overview of the selected peptides, as well as their major components and structural properties, is given in Table 1.




Fig. 1. Bleomycin sulfate (0.5 mg/mL) HPLC-UV (A) and UPLC-UV/MS (B) analysis with identification of major compounds and method comparison. The MS¹ spectra of the main compounds bleomycin A₂ and B₂, as well as impurities bleomycin A₅ and demethylbleomycin A₂, are given.

2.3. Sample preparation

Peptide samples were analyzed using HPLC-UV method derived from European Pharmacopoeia (Ph. Eur.), as well as by newly developed UPLC-UV/MS method [11–14]. As the latter also comprised MS detection, peptide sample solutions were prepared by replacing the buffer components with H₂O, whilst maintaining the same target concentrations as stipulated in the Ph. Eur. Details regarding peptide concentrations and solvents used for the samples are given in Table 1.

2.4. Chromatography

The HPLC apparatus consisted of a Waters Alliance 2695 separations module and a Waters 2996 Photodiode Array (PDA) detector with Empower 2 software for data acquisition (all Waters, Milford, MA, USA). Chromatographic analysis of the four selected peptides was based upon the methods as described in the Ph. Eur. [11–14]. Tyrothricin, vancomycin HCl and bleomycin sulfate were analyzed with a Lichrospher 100 RP-18 (250 mm × 4 mm, 5 µm) (Merck, Overijse, Belgium), whereas bacitracin was analyzed with a Hypersil ODS (250 mm × 4 mm, 5 µm) (Thermo Scientific,

Erembodegem-Aalst, Belgium). A brief overview of the chromatographic conditions is given in Table 2.

The UPLC equipment consisted of a Waters UPLC Quaternary Solvent Manager, a Waters ACQUITY Sample Manager, a Waters Ultra Performance LC (UPLC) PDA, a Waters ACQUITY Isocratic Solvent Manager and a Waters ACQUITY QDa detector (QDa) which is a compact single quad mass detector equipped with an electrospray ionization (ESI) interface. For data acquisition and instrument control, the Empower 3 FR 2 software was used. The ACQUITY UPLC CSH C₁₈ column (100 mm × 2.1 mm, 1.7 µm) (Waters, Zellik, Belgium) was maintained at 40 °C [15]. Mobile phase A consisted of a 95/5 H₂O/ACN+0.1% TFA, whilst mobile phase B consisted of a 5/95 composition of the same solvents. Different gradient methods were used for peptide analysis (Table 3). The flow rate was set at 0.5 mL/min and the injection volume was 2 µL. A post-column 10/1 PDA/QDa split ratio was employed, together with a post-column addition of 40/10/50 H₂O/propionic acid/2-propanol at a flow rate of 0.35 mL/min to the portion going to the QDa. This post-column addition neutralized the TFA ion suppression effect [16] and sustained a sufficient flow rate to the QDa detector. The QDa was operated in an electrospray positive ion mode by applying a voltage of 0.8 kV to the ESI capillary and the




Fig. 2. Tyrothricin (1 mg/mL) HPLC-UV (A) and UPLC-UV/MS (B) analysis with identification of major compounds and method comparison. Exemplary MS¹ spectra of main compounds tyrocidin C and gramicidin C₁, are given.

cone voltage was set at 15 V. The desolvation temperature was set at 600 °C. A full mass spectrum between *m/z* 100 and 1250 was acquired at a sampling rate of 2.0 points/sec.

The HPLC and UPLC methods were compared based upon (i) the system suitability tests (SST), i.e., resolution, peak to valley ratio (*p/v*), signal-to-noise ratio (*S/N*) and symmetry factor (*A_s*), listed in the Ph. Eur.; (ii) the quality limits of the active pharmaceutical ingredient (API); (iii) the limit of detection (LOD) based on *S/N* ratio; (iv) the number of peaks detected above reporting threshold (reported peaks); (v) the number of peaks (tentatively) identified; and (vi) the total run time.

3. Results

Chromatograms obtained from HPLC-UV and UPLC-UV/MS analysis of bleomycin sulfate are given in Fig. 1. Using the HPLC-UV method, bleomycin A₂ and B₂ were identified by their high normalized peak areas (62.79% and 34.05%) [11]. Impurity D (demethylbleomycin A₂) was identified by its relative retention time (RRT) to bleomycin A₂, i.e., 1.5 to 2.5. Although this HPLC-UV method

Criteria	Acceptance limits	Result	
		UPLC-UV/MS	HPLC-UV
System suitability test			
<i>p/v</i> Gramicidin A ₂ , A ₁	≥3	2.05	-
Quality requirements			
∑Gramicidin A ₁ , C ₁ , C ₂ , and A ₂	25%-50%	25.45%	-
∑Tyrocidins	50%-70%	66.28%	-
Gramicidins and Tyrocidins	Min. 85 %	91.73%	-
Sensitivity			
LOD	-	UV: 0.689 µg/mL (0.0698 % relative to 1mg/mL solution)* QDa [†] : 1.27 µg/mL (0.13 % relative to 1 mg/mL solution)	-
Extra Parameters			
Run time	-	16.8 min	60.0 min
Reported peaks	-	12	14 [°]
Identified peaks	-	7	0

* Calculated using the Tyrocidin C peak in the 0.02 mg/ml tyrothricin chromatogram

[°] no disregard threshold was used

[†] Selective ion monitoring using *m/z* value of 675.11

adheres to the individual monograph SST (resolution between bleomycin A₂ and B₂ is higher than 5; peak width at half height), significant tailing was observed, i.e., above the maximal value of 1.5 as given in the Ph. Eur. 2.2.46, and the amount of bleomycin B₂ surpassed the limit stipulated in the Ph. Eur. This is an example of how an averagely performing HPLC-UV QC method, combined with inadequate integration, can result in a false low quality conclusion. Using the UPLC-UV/MS method, which resulted in a better chromatographic separation of all compounds as evidenced by the reduction in tailing of the major bleomycin compounds and increased resolution with the present impurities, the amount of bleomycin B₂ was determined to be 29.0%, thus complying with the Ph. Eur. limit. Moreover, using MS detection, the identity of the bleomycin A₂ (*m/z* 708.08 and 472.44; *z*=2, 3) and B₂ (*m/z* 713.52 and 476.06; *z*=2, 3) was confirmed. Furthermore, two additional Ph. Eur.-listed impurities were identified by their mass spectra, i.e., bleomycinic acid (*m/z* 657.34; *z*=2) and bleomycin A₅ (*m/z* 720.42 and 480.73; *z*=2, 3). The identity of demethylbleomycin A₂ (impurity D) was also confirmed by MS (*m/z* 701.04; *z*=2). Finally, it was noted that the total run time of the UPLC-UV/MS method was reduced to 20 min, whereas the run time of HPLC-UV was 100 min.




Fig. 3. HPLC-UV (A) and UPLC-UV/MS (B) analysis of vancomycin HCl (2 mg/mL) with identification of major compounds and method comparison. The MS¹ spectra of the main compound, Vancomycin B, and two exemplary impurity peaks are given.

Chromatograms obtained from HPLC-UV and UPLC-UV/MS analysis of tyrothricin are given in Fig. 2. Due to the complex nature of the tyrothricin sample that contains different tyrocidins and gramicidins, formal identification was not possible with the HPLC-UV method without an additional analysis for retention time matching to reference standards. Moreover, in the absence of certified reference materials, calculation of the SST (peak to valley ratio between gramicidin A₁ and A₂) could not be performed to verify whether the method was Ph. Eur. compliant. Alternatively, using the UPLC-UV/MS method, tyrocidins A-D (*m/z*: 636.01, 655.67, 675.11 and 686.47, respectively; *z*=2) and gramicidin A₁, A₂ and C₁ (*m/z*: 941.65, 948.76 and 930.22, respectively; *z*=2) were identified based upon their MS spectra and the presence of K⁺ adducts.

Chromatograms obtained from HPLC-UV and UPLC-UV/MS analysis of vancomycin HCl are given in Fig. 3. Identification of vancomycin B using the traditional HPLC-UV method was done by its high relative peak area, whereas aglucovancomycin B (Ph. Eur. impurity C) was identified based upon its relative retention time to vancomycin B [17]. Analysis of the same sample using the UPLC-UV/MS method improved the resolution, resulting in the separation of more impurity peaks. Moreover, using the MS detection, the majority of these impurity peaks could be tentatively

Criteria	Acceptance limits	Result	
		UPLC-UV/MS	HPLC-UV
System suitability test			
S/N			
Vancomycin B (test sol c)	≥ 5.0	12.53	5.33
A _s			
Vancomycin B (test sol b)	0.8 - 1.6	1.083	0.888
Quality requirements			
Vancomycin B	Min. 93.0%	95.7%	93.3%
Total impurities	Max. 7.0%	4.3%	6.7%
Individual impurity	Max 4.0%	No impurity > 4.0 %	No impurity > 4.0 %
Sensitivity			
LOD	-	UV: 0.479 µg/mL (0.024 % relative to 2 mg/mL solution)*	1.13 µg/mL (0.057 % relative to 2 mg/mL solution)*
		QDa [†] : 0.220 µg/mL (0.011 % relative to 2 mg/mL solution)	
Extra Parameters			
Run time	-	30.0 min	35.0 min
Reported peaks	-	19	11
Identified peaks	-	13	2

* Calculated using the vancomycin B peak in the 0.002 mg/mL chromatogram

[†] Selective ion monitoring using *m/z* value of 725.44

identified based upon their MS¹ spectra and the work was performed by Diana et al. (Table 4) [17]. The mass spectra of two vancomycin-related impurities, which were tentatively identified as a demethylated impurity (peak 3) and aglucovancomycin B (peak 13), are given in Fig. 3. Even though their individual impurity amount, relative to the main compound vancomycin B, was well below 1%, i.e., 0.1% for peak 3 and 0.3% for peak 13, these mass spectra could be used for immediate identification purposes.

A similar observation was made for analysis of bacitracin by UPLC-UV/MS. Confirmation of the identity of the main compounds, as well as tentative identification of major impurity peaks was done with reference to the previous MS research performed by Govaerts et al. [18]. The mass spectra of the main compound bacitracin A, as well as of two small impurities, i.e., bacitracin E (1.5% relative to bacitracin A) and bacitracin H (1.7% relative to bacitracin A), are given in Fig. 4 and were positively used for identification purposes. The total run time of the UPLC-UV/MS method was also drastically reduced compared with the HPLC method, i.e., 16.8 min vs. 75 min.

A comparison between the performances of the traditional HPLC-UV and new UPLC-UV/MS method for the QC analysis of bleomycin sulfate, vancomycin HCl and bacitracin is given in Fig. 5.

Table 4
Tentatively identified vancomycin HCl related impurities.

#	RT (min)	m/z	Identification
1	1.134	661.24	Oxidized desvancosaminylvancosamin, Demethylleucylvancomycin B
2	5.028	734.78	Crystalline degradation product (demethylated), Deamidated vancomycin B (demethylated), N-demethylvancomycin B
3	5.095	718.40	Crystalline degradation product (demethylated), Deamidated vancomycin B (demethylated)
4	5.209	725.51	Oxidized vancomycin B (demethylated), Crystalline degradation product major, Crystalline degradation product minor
5	5.318	717.77	Crystalline degradation product (demethylated), Deamidated vancomycin B (demethylated), N-demethylvancomycin B
6	5.950	718.47	Crystalline degradation product (demethylated), Deamidated vancomycin B (demethylated)
7	6.009	724.95	Crystalline degradation product major, Crystalline degradation product minor
8	6.098	725.51	Oxidized vancomycin B (demethylated)
9	6.302	724.95	Crystalline degradation product major, Crystalline degradation product minor
10	6.369	725.44	Vancomycin B
11	6.518	725.92	Oxidized vancomycin B (demethylated)
12	6.707	725.58	Oxidized vancomycin B (demethylated)
13	12.694	1143.48	Aglucovancomycin B

As tyrothricin components could not be identified using the HPLC-UV method, performance comparison could not be made.

4. Discussion

The advantages of sub-2 μm chromatography (UPLC) over traditional HPLC chromatography for analysis of pharmaceutical compounds have already been extensively documented [19–23]. The use of smaller particle size with an optimized LC system gives improved efficiency, resulting in faster and higher resolution separations than those of traditional LC systems with column particle sizes in the 3–10 μm range. For the transfer of a method from traditional HPLC to UPLC (and vice versa), guidance about the adjustment of column size, injection volume, flow rate and time program for gradient elution are available [24]. However, this was not the main objective of this study. The goal of this study was to inspire pharmaceutical analysts to implement new technologies in method development. Currently, most LC-based quality control methods for APIs listed in the Ph. Eur. generally still utilize stationary phases with a particle size of 3–10 μm . Small variations in particle size are permitted for isocratic methods, i.e., reduction of 50% in particle size, whereas no alteration in particle size is allowed with gradient elutions [25]. As a result, the use of UPLC




Fig. 4. HPLC-UV (A) and UPLC-UV/MS (B) analysis of bacitracin (2 mg/mL) with identification of major compounds and method comparison. The MS¹ spectra of a main compound, bacitracin A, and two exemplary impurity peaks, are given.

Criteria	Acceptance limits	Result	
		UPLC-UV/MS	HPLC-UV
System suitability test			
p/v Bacitracin B ₁ , B ₂	≥ 1.2	21.801	2.83
Quality requirements			
∑Bacitracin A, B ₁ , B ₂ and B ₃	Min. 70%	72.14 %	67.70%
Bacitracin A	Min. 40.0%	44.07%	42.18%
Sensitivity			
LOD	-	UV: 0.149 $\mu\text{g/mL}$ (0.00745 % relative to 2 mg/mL solution) QDa [†] : 0.0207 $\mu\text{g/mL}$ (0.00103 % relative to 2 mg/mL solution)	4.00 $\mu\text{g/mL}$ (0.2 % relative to 2 mg/mL solution)
Extra Parameters			
Run time	-	16.8 min	75.0 min
Reported peaks	-	28	26
Identified peaks	-	14	4

* Calculated using the Bacitracin A peak from the 0.01 mg/ml bacitracin chromatogram

[†] Selective ion monitoring using m/z value of 712.12




Fig. 5. Relative performance comparison between UPLC-UV/MS and HPLC-UV method for bleomycin sulfate, vancomycin HCl and bacitracin QC analysis. HPLC-UV performance was set at 100%.

converting existing API quality control methods in Ph. Eur. into faster and more efficient sub-2 μm methods is currently not allowed without validation. Hence, simple transferring HPLC–UPLC methods in pharmaceutical QC method validation is currently not allowed. Moreover, it is emphasized that when a method is changed, the corresponding acceptance limits are also to be re-considered [26]. However, newly developed quality control methods have already used this sub-2 μm stationary phase as evidenced by the related substance methods for quetiapine fumarate and nevirapine hemihydrate [27,28]. Moreover, the revised Ph. Eur. 8.3 general chapter 2.2.29 liquid chromatography covers sub-2 μm particle technology [29].

Although MS is listed as an analytical technique in the Ph. Eur. [30], its use in the quality control of APIs is currently limited to a few cases, e.g. detection of impurity B of oseltamivir phosphate, detection of impurity F of imatinib mesilate, and determination of the interferon β -1a isoform distribution [31–33]. However, as the newly developed single quad MS detector is smaller, cheaper and easily used, it is expected that its application will increase exponentially. Coupled to UPLC separation module, it adequately answers the current shift in QC emphasis from API assay towards impurity profiling [34,35].

This study has demonstrated the applicability of identifying related impurities and/or multiple active components in complex API samples using a single quad MS detector. The current equipment set-up, using a post-column split, coupled with post-column addition of a propionic acid containing solvent, even permits the use of TFA to improve chromatographic peak shape [36], without compromising the MS detection too much. As liquid chromatographic methods with volatile buffers often show inferior chromatographic separation compared with analogous methods with non-volatile buffers [37], the use of the anionic ion-pairing reagent TFA for peptide separation is important to counterbalance the loss in chromatographic performance. As expected, resolution and LOD

were improved by the use of sub-2 μm column technology and the overall run time was significantly reduced, when compared with traditional 3–10 μm HPLC methods listed.

5. Conclusion

A newly developed, single quad MS detector was coupled to a UPLC separation module and used in routine quality control analysis of bleomycin sulfate, tyrothricin, vancomycin HCl and bacitracin peptide APIs. The results were compared with the results obtained by traditional HPLC-UV methods which were based upon the Ph. Eur. As expected, the UPLC separation resulted in a higher resolution and a lower limit of detection, as well as a significant reduction in run time. Furthermore, MS detector may enable to directly identify impurities or components even at low levels without the need of reference standards. Currently, MS is only sparsely applied in the Ph. Eur. However, recent technical advances deliver fit for using single quad MS detectors to the pharmaceutical field, and its applications in QC analysis are expected to rise in the near future.

Acknowledgments

This research was funded by PhD grants of 'Institute for the Promotion of Innovation through Science and Technology in Flanders (IWT-Vlaanderen)' (No. 101529 (MD) and 121512 (BG)). We would also like to thank Jonas Poppe for his experimental help in this study and Laurence Van Oudenhove, Stephen McDonald and Jan Claeriboudt (Waters) for the critical reading of this manuscript.

References

- [1] R. Aebersold, M. Mann, Mass spectrometry-based proteomics, *Nature* 422 (2003) 198–207.
- [2] A. Pandey, M. Mann, Proteomics to study genes and genomes, *Nature* 405 (2000) 837–846.
- [3] B. Domon, R. Aebersold, Mass spectrometry and protein analysis, *Science* 312 (2006) 212–217.
- [4] S.-E. Ong, M. Mann, Mass spectrometry-based proteomics turns quantitative, *Nat. Chem. Biol.* 1 (2005) 252–262.
- [5] X. Zhang, D. Wei, Y. Yap, et al., Mass spectrometry-based “omics” technologies in cancer diagnostics, *Mass Spectrom. Rev.* 26 (2007) 403–431.
- [6] M.S. Lee, E.H. Kerns, LC/MS applications in drug development, *Mass Spectrom. Rev.* 18 (1999) 187–279.
- [7] R. Kostianinen, T. Kotiaho, T. Kuuranne, et al., Liquid chromatography/atmospheric pressure ionization-mass spectrometry in drug metabolism studies, *J. Mass Spectrom.* 38 (2003) 357–372.
- [8] W. Xu, N.E. Manicke, G.R. Cooks, et al., Miniaturization of mass spectrometry analysis systems, *J. Lab. Autom.* 15 (2010) 433–439.
- [9] S.E. Hamilton, F. Mattrey, X. Bu, et al., Use of miniature mass spectrometer to support pharmaceutical process chemistry, *Org. Process Res. Dev.* 18 (2014) 103–108.
- [10] X. Bu, J. Yang, X. Gong, et al., Evaluation of compact mass spectrometer for routine support of pharmaceutical chemistry, *J. Pharm. Biomed. Anal.* 94 (2014) 139–144.
- [11] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Bleomycin sulfate, Strasbourg, France, 2015.
- [12] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Tyrothricin, Strasbourg, France, 2015.
- [13] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Vancomycin hydrochloride, Strasbourg, France, 2015.
- [14] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Bacitracin, Strasbourg, France, 2015.
- [15] M.A. Lauber, S.M. Koza, S.A. McCall, et al., High-resolution peptide mapping separations with MS-friendly mobile phases and charge-surface-modified C18, *Anal. Chem.* 85 (2013) 6936–6944.
- [16] A. Apffel, S. Fischer, G. Goldberg, et al., Enhanced sensitivity for peptide mapping with electrospray liquid chromatography–mass spectrometry in the presence of signal suppression due to trifluoroacetic acid-containing mobile phases, *J. Chromatogr. A* 712 (1995) 177–190.
- [17] J. Diana, D. Visky, J. Hoogmartens, et al., Investigation of vancomycin and related substances by liquid chromatography/ion trap mass spectrometry, *Rapid Commun. Mass Spectrom.* 20 (2006) 685–693.
- [18] C. Govaerts, C. Li, J. Orwa, et al., Sequencing of bacitracin A and related minor components by liquid chromatography/electrospray ionization ion trap tandem mass spectrometry, *Rapid Commun. Mass Spectrom.* 17 (2003) 1366–1379.
- [19] D. Guillaume, J. Ruta, S. Rudaz, New trends in fast and high-resolution liquid chromatography: a critical comparison of existing approaches, *Anal. Biochem. Chem.* 397 (2010) 1069–1082.
- [20] T. Wu, C. Wang, X. Wang, et al., Comparison of UPLC and HPLC for analysis of 12 phthalates, *Chromatographia* 68 (2008) 803–806.
- [21] S.A.C. Wren, P. Tchelitcheff, Use of ultra-performance liquid chromatography in pharmaceutical development, *J. Chromatogr. A* 1119 (2006) 140–146.
- [22] N. Wu, A.M. Clausen, Fundamental and practical aspects of ultrahigh pressure liquid chromatography for fast separations, *J. Sep. Sci.* 30 (2007) 1167–1182.
- [23] L. Nováková, L. Matysová, P. Solich, Advantages of application of UPLC in pharmaceutical analysis, *Talanta* 68 (2006) 908–918.
- [24] D. Guillaume, D.T.T. Nguyen, S. Rudaz, et al., Method transfer for fast liquid chromatography in pharmaceutical analysis: Application to short columns packed with small particle. Part II: gradient experiments, *Eur. J. Pharm. Biopharm.* 8 (2008) 430–440.
- [25] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, 2.2.46 Chromatographic separation techniques, Strasbourg, France, 2015.
- [26] J.O. De Beer, B.M.J. De Spiegeleer, J. Hoogmartens, et al., Relationship between content limits and assay methods: an interlaboratory statistical evaluation, *Analyst* 117 (1992) 933–940.
- [27] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Quetiapine fumarate, Strasbourg, France, 2015.
- [28] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Nevirapine hemihydrate, Strasbourg, France, 2015.
- [29] European Directorate for the Quality of Medicines and Healthcare, *Pharmeuropa* 24, Strasbourg, France, 4 October 2012.
- [30] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, 2.2.43 Mass spectrometry, Strasbourg, France, 2015.
- [31] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Oseltamivir phosphate, Strasbourg, France, 2015.
- [32] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.5, Imatinib mesilate, Strasbourg, France, 2015.
- [33] European Directorate for the Quality of Medicines and Healthcare, European Pharmacopoeia 8.2, Interferon β -1a concentrated solution, Strasbourg, France, 2015.
- [34] S. Görög, The importance and the challenges of impurity profiling in modern pharmaceutical analysis, *Trends Anal. Chem.* 25 (2006) 755–757.
- [35] S. Baertschi, Analytical methodologies for discovering and profiling degradation-related impurities, *Trends Anal. Chem.* 25 (2006) 758–767.
- [36] D.V. McCalley, Effect of buffer on peak shape of peptides in reversed-phase high performance liquid chromatography, *J. Chromatogr. A* 1038 (2004) 77–84.
- [37] L. Van den Bossche, A. Van Schepdael, S. Chopra, et al., Identification of impurities in polymyxin B and colistin bulk sample using liquid chromatography coupled to mass spectrometry, *Talanta* 83 (2011) 1521–1529.