High-Speed Imaging of Amoeboid Movements Using Light-Sheet Microscopy

Daisuke Takao^{1,9}, Atsushi Taniguchi^{2,9}, Takaaki Takeda³, Seiji Sonobe², Shigenori Nonaka¹*

1 Laboratory for Spatiotemporal Regulations, National Institute for Basic Biology, Okazaki, Aichi, Japan, 2 Department of Life Science, Graduate School of Life Science, University of Hyogo, Kamigori-cho, Hyogo, Japan, 3 National Astronomical Observatory of Japan, Mitaka, Tokyo, Japan

Abstract

Light-sheet microscopy has been developed as a powerful tool for live imaging in biological studies. The efficient illumination of specimens using light-sheet microscopy makes it highly amenable to high-speed imaging. We therefore applied this technology to the observation of amoeboid movements, which are too rapid to capture with conventional microscopy. To simplify the setup of the optical system, we utilized the illumination optics from a conventional confocal laser scanning microscope. Using this set-up we achieved high-speed imaging of amoeboid movements. Three-dimensional images were captured at the recording rate of 40 frames/s and clearly outlined the fine structures of fluorescent-labeled amoeboid cellular membranes. The quality of images obtained by our system was sufficient for subsequent quantitative analysis for dynamics of amoeboid movements. This study demonstrates the application of light-sheet microscopy for high-speed imaging of biological specimens.

Citation: Takao D, Taniguchi A, Takeda T, Sonobe S, Nonaka S (2012) High-Speed Imaging of Amoeboid Movements Using Light-Sheet Microscopy. PLoS ONE 7(12): e50846. doi:10.1371/journal.pone.0050846

Editor: Alexandre J. Kabla, University of Cambridge, United Kingdom

Received July 30, 2012; Accepted October 25, 2012; Published December 5, 2012

Copyright: © 2012 Takao et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Funding: This work was supported by the Human Frontier Science Program (HFSP, http://www.hfsp.org/), the Ministry of Education, Culture, Sports, Science and Technology (MEXT, http://www.jsps.go.jp/english/), the Japan Society for the Promotion of Science (JSPS, http://www.jsps.go.jp/english/index.html) and CREST (http://www.jst.go.jp/kisoken/crest/en/index.html). The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript.

Competing Interests: The authors have declared that no competing interests exist.

* E-mail: snonaka@nibb.ac.jp

• These authors contributed equally to this work.

Introduction

Selective plane illumination microscopy (SPIM) and its derivative, digital-scanned light-sheet microscopy (DSLM), were developed for fluorescence imaging of live biological specimens [1-6]. In SPIM, light-sheet illumination is achieved by selectively illuminate a single plane using a cylindrical lens [1]. The objective, which contains a CCD camera that acquires images, is placed with its optical axis perpendicular to path of illumination. The principle of DSLM is similar to that of SPIM except for the synchronization of the 'apparent' light-sheet by the scanning excitation beam and image acquisition [4]. The advantages of DSLM over SPIM include [4]: (1) a relatively constant intensity profile along the width of the 'apparent' light-sheet, which is important for quantitative imaging of large specimens; (2) a reduction in the optical aberrations; (3) an illumination efficiency of 95% for DSLM compared to approximately 3% for conventional SPIM. The major advantages of SPIM and DSLM are their efficiency of illumination, low phototoxicity, and high-speed image acquisition. In fact, the rate of image acquisition by SPIM is increased by 20 frames/s across a 360×360 pixel field of view [6].

Light-sheet microscopy is a very powerful tool for the biological studies. However, for many biologists who are not familiar with optics it is difficult to construct such an optical system on their own. In this regard, to simplify the system, we developed a DSLM that utilized the illumination optics from a conventional confocal laser-scanning microscope (CLSM). We referred to it as "ezDSLM" to reflect its ease of setup. The ezDSLM system makes light-sheet microscopy available to any laboratory because it is based on a conventional CLSM.

A shortcoming of conventional light-sheet microscopy is its limitations for high-speed imaging in some situations. For example, the specimen is typically embedded in agarose gel and moved through the light-sheet to obtain three-dimensional images. Since the position of light-sheet is fixed during imaging, the specimen must be moved using a sample holder. This can be problematic if the specimen needs to be observed without embedding in agarose gel, as is the case for the observation of the movements of amoeba on substrates or ciliary movements of cultured cells. Moreover, rapid movements or shaking of a specimen within medium may impact cell-substrate contacts or ciliary movements due to substrate vibrations or extracellular liquid flow. To overcome the limitations of standard light-sheet microscopy we modified the way that three-dimensional images are obtained by optically scanning the light-sheet instead of moving the specimen.

To demonstrate the ease of set-up for ezDSLM and its utility for biological imaging we observed the movements of *Amoeba proteus*. Although amoeboid movements can be observed by fluorescently labeling their plasmalemmas, live imaging has been limited to single optical planes. This was largely due to the fact that amoeba change their shapes too rapidly to obtain three-dimensional images by CLSM [7]. Despite being able to overcome this limitation with conventional DSLM, it would introduce the new problem of sample shaking as described above. Herein we describe the ezDSLM to capture clear images of rapidly moving amoeba with both high spatial and temporal resolutions that are sufficient for subsequent quantitative analyses.

Results and Discussion

Overview of the ezDSLM

We developed ezDSLM. Schematic drawings of image acquisition by ezDSLM are shown in Figure 1 and Movie S1. As in conventional DSLM [4], the line scan of a laser beam is used to generate a light-sheet that is synchronized with image acquisition by a CCD camera (Fig. 1A,B). In contrast to conventional DSLM, ezDSLM uses a two-dimensional scan of illumination to move the 'apparent' light-sheet along the axis of objective for emission detection. To keep the light-sheet on the focal plane of the objective for emission detection, the objective is mounted to the motor stage and synchronized to the 'apparent' light-sheet (Fig. 1C). Since it is unnecessary to move specimens during imaging, fast scans for three-dimensional image acquisition is achieved without shaking the specimens. Therefore, freely-moving specimens in medium can be observed without outside interference.

To simplify the ezDSLM system we used the illumination optics from a CLSM system (Fig. 2A). The FV1000 system (Olympus) was used as a master of signaling and electric devices were fully synchronized (Fig. S1). Laser scanning for illumination is controlled by FV1000, which outputs electric signals for synchronization. The signals are then used to trigger movement of the objective mounted on the motor stage and image acquisition by the CCD camera. The ezDSLM system is very simple because it utilizes the illumination optics of the CLSM, which transforms it into a light-sheet microscope. It is also highly advantageous to biologists because it makes light-sheet microscopy available to people without special skills or expertise.

Observation of Amoeboid Movements by ezDSLM

To evaluate the performance of ezDSLM, we applied it to highspeed imaging of amoeboid movements on substrates. A setup with standard conditions for the observation of the movements of *Amoeba proteus* [7] was constructed. Fluorescently labeled amoebae were placed on the glass surface of the sample holder (Fig. 2B,C). Amoebae have special membrane structures known as plasmalemma, which is composed of multiple layers including the plasma membrane, amorphous layer, and filamentous layer. The plasmalemma was labeled with lipophilic dye DiI. The chamber unit was filled with KCM medium, which is the standard medium for observation of amoebae. In conventional CLSM, it is difficult to observe the amoeboid movements with sufficient spatial and temporal resolutions for quantitative analyses. The images have low contrast and are blurred because amoebae change their shapes rapidly during the acquisition of images for a single XYZ stack. In ezDSLM, by contrast, we successfully obtained a time-series of XYZ-stacked images of whole amoebae (Fig. 3). The threedimensional images clearly outlined the cell shapes of moving amoebae.

Spatial and Temporal Resolutions Achieved by ezDSLM

The image quality, with regard to contrast or focusing, of ezDSLM was comparable to that in conventional DSLM, and thus ezDSLM worked well as a light-sheet microscope. The typical condition for image acquisition used here is as follows: Objectives for excitation, 5×, NA 0.10, air (Olympus), and for emission detection, $40 \times$, NA 0.80, water dipping (Olympus); 336×256 pixels (corresponding to $197.0 \times 228.2 \ \mu m$) in XY resolutions; 2.48 µm intervals for each optical slice in Z stacks; frame rate, 40 frames/s. Under these conditions, it takes 2.5 s to acquire a Z stack of 100 X-Y slices (~250 µm in Z dimension) that typically cover the whole body of a single amoeba (~200 µm in Z dimension; Fig. 3D). Compared to the frame rate of the CLSM (2.39 frames/s with the same resolutions), ezDSLM achieved a 17-times higher rate. Although the frame rate of CLSM can be improved by using spinning disc or resonant scanner, light-sheet microscopy has advantages such as greater penetration depth of illumination and reduced photobleaching and phototoxicity, as described [1-6]. In addition, our observation of amoeba with the same preparation using a spinning disk microscope often produced false signal that was produced by out-of-focus-light passed neighboring pinholes (data not shown). The limiting step in our high-speed imaging system is the recording rate of the CCD camera. Therefore, the imaging rate of ezDSLM would be remarkably improved with a higher performance camera.

Using CLSM, the outlines of amoebae were frequently blurred since they moved rapidly during image acquisition. In addition, the displacement and changes in cell shapes of amoebae between each time frame were too large to follow and analyze (Fig. 3E,F, Movie S2). In contrast, images obtained with ezDSLM clearly outlined the cell shapes in -XYZ dimensions (Fig. 3A–D, Fig. S2, Movie S3), as well as in the time series of the Z-stack (Fig. 3G,H, Movie S4, Movie S5). Changes in cell shape, such as shrinkages and extensions of protrusions, were captured in detail with sufficient spatial and temporal resolutions. We also made timelapse three-dimensional CG movies of membrane structures, using



Figure 1. Schematic of ezDSLM (A, B). Line scan (A) to make an 'apparent' light-sheet (B). This step is similar to that in the conventional DSLM. (C) Scan of the light-sheet to obtain three-dimensional image of specimen. The objective for emission detection is synchronized with the movement of the light-sheet scanning. Ex, objective for excitation; Em, objective for emission detection; Sp, specimen. doi:10.1371/journal.pone.0050846.g001



Figure 2. Setup of ezDSLM for observation of amoeboid movements. (A) Side view of the instrument (not to scale). Optical axes of objectives for excitation (Ex) and emission (Em) are shown in dashed lines. The illumination optics were derived from the confocal laser scanning microscope (CLSM). Specimens were placed on the sample holder (SH) in the chamber unit (CU) filled with medium. The CU containing the objective for emission detection is mounted on a movable motor stage. BF, barrier filter; TL, tube lens; CCD, CCD camera. (B) Photograph of the chamber unit. (C) Photograph of the sample holder. Amoebae were placed on the cover glass for imaging (arrow head). doi:10.1371/journal.pone.0050846.g002

4-dimensional digital universe (4D2U) system [8]. These movies depicted the motion of the membrane structures (Movie S6). Therefore, the image acquisition rate of ezDSLM was sufficiently high enough to follow three-dimensional shape changes of moving amoebae. Furthermore, clusters of accumulated dyes on the plasmalemma were traceable in the ezDSLM images (Fig. 3G,H). By using these clusters as landmarks, it is possible to analyze not only cell shapes but also membrane dynamics in detail.

Quantitative Analysis

To show that ezDSLM is applicable for subsequent quantitative analyses, we measured dynamics of protrusions in amoeboid movements. Since protrusions extend three-dimensionally and changes in their shapes are rapid as described above, for precise analyses, three-dimensional images with both spatially and temporally sufficient resolutions are needed. In this regard, ezDSLM would be a suitable tool.

Using data obtained by ezDSLM, mean velocities of extension and retraction of protrusions were measured and compared (Fig. 4A). The mean velocity of extension was significantly higher than that of retraction, which coincides with the observations. We then measured time-evolved velocity of extension (Fig. 4B). Extensions of protrusions were accelerated for the early phase and they reached max velocities within ~10 s. This may reflect the underlying mechanism in which actomyosin plays a role for extension of protrusion [9–12]. A small protrusion termed "bleb" is formed by contraction of the actomyosin cortex in the early phase, and this explosive force would accelerate the extension. Although further investigation is needed to understand biological significance, we showed that ezDSLM is applicable to such analyses.

Biological Significance of Application of ezDSLM

In this study, we showed the application of ezDSLM for high-speed imaging and quantitative analysis of movements of Amoeba proteus. Amoeboid movements and membrane dynamics have been well characterized in dendritic cells [13] and leucocytes [14]. For Amoeba proteus, on the other hand, observation of cell movements has been difficult due to its high motility: $\sim 100 \,\mu\text{m/min}$ compared to $< 10 \,\mu\text{m/min}$ for fibroblasts. Therefore, imaging of motile Amoeba proteus has been limited to be two dimensions [7,15]. Although some excellent studies of three-dimensional images of moving Dictyostelium discoideum were reconstructed from differential interference contrast (DIC) micrographs [16,17], the contrast was relatively low compared to fluorescence images. In addition, traceable clusters of dyes were observable in cells stained with DiI (Fig 3G,H) making fluorescence imaging beneficial for the analyses of membrane dynamics. Collectively, high-speed imaging achieved by ezDSLM would serve as a powerful tool



Figure 3. Images of freely moving Dil-labeled amoebae. Images were obtained using the ezDSLM (A–D, G, H) and CLSM (E, F). (A, B) X–Y optical sections from a Z-stack. Z35 and Z45 indicate slice numbers in Z-stacks with 2.48 μm intervals, i.e., actual distances are 86.8 and 111.6 μm, respectively. A complete series of Z-stack is shown in Fig. S2. (C, D) Front (C) and side (D) views of a reconstructed three-dimensional image. (E, F) Max projection images of Z-stacks at time-points T3 (E) and T4 (F) from time-lapse images obtained using CLSM. Intervals for each time-point were 0.54 s. It is difficult to follow the changes in shapes (arrowheads). (G, H) Max projection images from Z-stacks at time points T25 (G) and T30 (H) from time-lapse images obtained by ezDSLM. Intervals for each time-point were 0.04 s. Extensions of protrusions (arrowheads) and traceable clusters of dye (arrows) are depicted in detail. Note that the four more images (T26–T29) are present between these time points (not shown). Scale bars, 50 μm. doi:10.1371/journal.pone.0050846.g003

to investigate amoeboid movements as well as membrane dynamics. Further investigation using ezDSLM would provide a breakthrough in our understanding of the mechanisms of amoeboid movements. High-speed imaging with ezDSLM is also applicable to many other biological studies, such as the observation of ciliary movements. The most important features of ezDSLM include the observation of specimens without shaking and its ease of setup. These features of ezDSLM facilitate the application of light-sheet microscopy.

Conclusions

We developed ezDSLM, which we used to capture remarkably clear images of moving amoebae. This process did not require the specimen to be shaken because the light-sheet is optically scanned to obtain three-dimensional images. We demonstrated the clarity and speed of ezDSLM to outline the movements of *Amoeba proteus*. The ezDSLM substantially extended the application of light-sheetbased microscopy in biological field. Finally, since our system uses



Figure 4. Dynamics of protrusions during amoeboid movements. By using the images obtained by ezDSLM, velocities of extending or retracting protrusions were measured. (A) Mean velocities of protrusions during extension or retraction are shown in box plots. Blue and black dots indicate mean values and outliers, respectively. The mean velocities of extension and retraction were significantly different by t-test (p<0.01). Data from 25 protrusions for each category are shown. (B) Time-evolved velocities of extending protrusions. Data from 5 protrusions are shown as mean \pm SD.

doi:10.1371/journal.pone.0050846.g004

a conventional CLSM system for illumination optics, setup is very simple and easily adaptable to other laboratories.

Materials and Methods

Cell Culture and Fluorescent Labeling of Amoeba proteus

Amoeba proteus was cultured in KCM culture medium (7 mg of KCl, 8 mg of CaCl₂, 8 mg of MgSO₄-7H₂O per liter) at 25°C and fed on *Tetrahymena pyriformis* as previously described [18]. For live imaging of amoeboid movements the lyophilic dye DiI (Invitrogen) was used to label the plasmalemma. Amoebae were stained for 10 min in KCM medium containing 5.36 μ M DiI (2.68 mM stock in DMSO) and washed three times with KCM medium. For imaging, the amoebae were placed on the sample holder (Fig. 2C) and rested for 15 min prior to imaging.

Setup of ezDSLM

Setup of ezDSLM is schematically shown in Fig. 2, Fig. S1, and Movie S1. The Olympus FV1000 was used as the illumination optics as described above. The excitation wave length was 559 nm and the barrier filter BA570-625HQ (Olympus) was selected for detection of DiI. The motorized stage (M-111.1DG) and controller (C-862) for moving the objective were obtained from PI. The resolution of the stage is high (minimal incremental motion; 50 nm) and it has long travel range (15 mm). The movement of the stage was controlled by the original macros on the software, Mikromove (PI). The ORCA-R2 (Hamamatsu) was the CCD camera used for image acquisition. The camera was controlled by the HC Image software (Hamamatsu). The controllers were connected to a computer that is distinct from that in the FV1000 system. All electric devices were connected through their I/O interfaces and the FV1000 system was used as a master of synchronization signals. Parameters for laser scanning (laser power, scan rate, spatial interval, etc) were controlled by a computer in the FV1000 system. The parameters were adjusted for imaging conditions, e.g., the increments and range of movements for the light-sheet were typically 2.48 µm and 0.3 mm, respectively. When the FV1000 start scanning, it outputs synchronization signals that are active during each line scan ("line active" signal) or frame scan ("frame active" signal) as shown in Fig. S1. Exposures of CCD camera were synchronized to the line active signals so that each line scan makes apparent light-sheet. The motion of the motorized stage was triggered by the frame active signal and it moved at the constant rate that was synchronized to the movement of light-sheet. Trial-and-error processes may be required to adjust parameters at first, which can be achieved by imaging fluorescent microbeads or reflection of glass surface as in the case of conventional light-sheet microscopy. Images were processed using Image J, and in some cases, using the 4D2U system that is originally developed for astronomical studies [8].

Measurements of Velocities of Protrusions

Coordinates (X, Y, Z) of extending or retracting edges of protrusions on the acquired images were manually tracked using ImageJ. By calculating displacements in each time interval, the velocities were obtained.

Supporting Information

Figure S1 Flow chart and waveforms of the synchronization signals. The flow of the synchronization signals from the FV1000 is shown in upper panel. PC1 and PC2 are the computers that regulate laser scanning and other functions including image acquisition, respectively. PC1 is a master of the synchronization signal. CCD camera and motorized stage are regulated by PC2 through the controllers. Waveforms of the synchronization signals and responses of the devices are schematically shown in lower panel. "Line active" signal is active while the FV1000 scans lines to make the apparent light-sheet. "Frame active" signal is active during single movements of the apparent light-sheet. (TIF)

Figure S2 Complete series of Z-stack shown in Fig. 3A, **B.** Slice numbers are labeled. Actual interval for each slice is 2.48 μm. (TIF)

Movie S1 Schematic animation of three-dimensional image acquisition by ezDSLM. The objective for emission detection that is incorporated in the chamber unit is mounted on the motorized stage as shown in Fig. 2. The range of scanning is enlarged to clarify the detail of movement. Actual scanning range was typically small (~ 0.3 mm) and the movement of the chamber unit did not affect amoeboid movements. (MOV)

Movie S2 Images of freely moving DiI-labeled amoebae, obtained using the CLSM. Max projection images of Zstacks were shown.

(MOV)

Movie S3 Images of freely moving DiI-labeled amoebae, obtained using the ezDSLM.

(MOV)

Movie S4 Images of freely moving DiI-labeled amoebae, obtained using the ezDSLM. (MOV)

References

- Huisken J, Swoger J, Del Bene F, Wittbrodt J, Stelzer EH (2004) Optical sectioning deep inside live embryos by selective plane illumination microscopy. Science 305: 1007–1009.
- Huisken J, Stainier DYR (2007) Even fluorescence excitation by multidirectional selective plane illumination microscopy (mSPIM). Opt Letters 32: 2608–2610.
- Keller PJ, Pampaloni F, Stelzer EH (2007) Three-dimensional preparation and imaging reveal intrinsic microtubule properties. Nat Methods 4: 843–846.
- Keller PJ, Schmidt AD, Wittbrodt J, Stelzer EH (2008) Reconstruction of zebrafish early embryonic development by scanned light sheet microscopy. Science 322: 1065–1069.
- Planchon TA, Gao L, Milkie DE, Davidson MW, Galbraith JA, et al. (2011) Rapid three-dimensional isotropic imaging of living cells using Bessel beam plane illumination. Nat Methods 8: 417–423.
- Wu Y, Ghitani A, Christensen R, Santella A, Du Z, et al. (2011) Inverted selective plane illumination microscopy (iSPIM) enables coupled cell identity lineaging and neurodevelopmental imaging in Caenorhabditis elegans. Proc Natl Acad Sci U.S.A. 108: 17708–17713.
- Nishihara E, Shimmen T, Sonobe S (2007) New aspects of membrane dynamics of Amoeba proteus contractile vacuole revealed by vital staining with FM 4–64. Protoplasma 231: 25–30.
- Kokubo E, Hayashi M, Kato T, Takeda T, Kaifu N, et al. (2005) 4-dimensional digital universe project. J Korean Astron Soc 38: 153–155.
- 9. Boulbitch A, Simson R, Simson DA, Merkel R, Hackl W, et al. (2000) Shape instability of a biomembrane driven by a local softening of the underlying actin

Movie S5 Images of freely moving DiI-labeled amoebae, obtained using the ezDSLM.

Movie S6 Images of freely moving DiI-labeled amoebae, obtained using the ezDSLM. The images were processed using the 4D2U system. (MOV)

Acknowledgments

We thank T. Murata for technical assistance and discussion.

Author Contributions

Conceived and designed the experiments: DT AT SS SN. Performed the experiments: DT AT. Analyzed the data: DT AT TT. Contributed reagents/materials/analysis tools: DT AT TT. Wrote the paper: DT.

cortex. Phys Rev E Stat Phys Plasmas Fluids Relat Interdiscip Topics 62: 3974–3985.

- Charras GT, Yarrow JC, Horton MA, Mahadevan L, Mitchison TJ (2005) Nonequilibration of hydrostatic pressure in blebbing cells. Nature 435: 365–369.
- Langridge PD, Kay RR (2006) Blebbing of Dictyostelium cells in response to chemoattractant. Exp Cell Res 312: 2009–2017.
- Charras G, Paluch E (2008) Blebs lead the way: how to migrate without lamellipodia. Nat Rev Mol Cell Biol 9: 730–736.
- Renkawitz J, Schumann K, Weber M, Lämmermann T, Pflicke H, et al. (2009) Adaptive force transmission in amoeboid cell migration. Nat Cell Biol 11: 1438– 1443.
- Lee J, Gustafsson M, Magnusson KE, Jacobson K (1990) The direction of membrane lipid flow in locomoting polymorphonuclear leukocytes. Science 247: 1229–1233.
- Grebecki A (1986) Two-directional pattern of movements on the cell surface of amoeba proteus. J. Cell Sci. 83: 23–35.
- Wessels D, Voss E, Bergen NV, Burns R, Stites J, et al. (1998) A computerassisted system for reconstructing and interpreting the dynamic three-dimensional relationships of the outer surface, nucleus and pseudopods of crawling cells. Cell Motil. Cytoskeleton 41: 225–246.
- Soll DR, Wessels D, Kuhl S, Lusche DF (2009) How a cell crawls and the role of cortical myosin II. Eukaryot Cell 8: 1381–1396.
- Griffin JL (1960) An improved mass culture method for the large, freeliving amebae. Exp Cell Res 21: 170–178.