

CLINICAL RESEARCH

e-ISSN 1643-3750 © Med Sci Monit, 2022; 28: e935812 DOI: 10.12659/MSM.935812

Received:	2021.12.13
Accepted:	2022.01.31
Available online:	2022.02.17
Published:	2022.03.27

Mar

A Faster, Novel Technique to Detect COVID-19 Neutralizing Antibodies

Author Satatis Data Ir uscrip Liter Fun	s' Contribution: Study Design A ta Collection B tical Analysis C nterpretation D t Preparation E rature Search F ds Collection G	ABCDEF 1 C 1 DF 2-5 BF 6 BD 7 DF 6 B 8-10 AG 1	Wei-Huai Chiu D Wei-Yi Kong D Yu-Jang Su D Jyun-Wei Wen Ciao-Ming Tsai D Chitsung Hong Pang-Yen Chen D Cheng-Hao Ko D		 Graduate Institute of Automation and Control, National Taiwan University of Science and Technology, Taipei, Taiwan Poison Center, Department of Emergency Medicine, MacKay Memorial Hospita Taipei, Taiwan Yuanpei University of Medical Technology, HsinChu, Taiwan Department of Medicine, MacKay Medical College, Taipei, Taiwan MacKay Junior College of Medicine, Nursing and Management, Taipei, Taiwan SpectroChip, Inc., Hsinchu, Taiwan Department of Power Mechanical Engineering, National Tsing Hua University, Hsinchu, Taiwan Department of Emergency Medicine, MacKay Memorial Hospital, Taipei, Taiwa Department of Emergency Medicine, MacKay Memorial Hospital, Taipei, Taiwa Department of Environmental and Occupational Health Sciences, National Yang M Chiao Tung University College of Medicine, Taipei, Taiwan Institute of Public Health, National Yang Ming Chiao Tung University College Medicine, Taipei, Taiwan 	al, Aing of	
Corresponding Author: Financial support: Conflict of interest:			Cheng-Hao Ko, e-mail: Kevin.no.ntust.2@gmail.com, dr.ko.chenghao@gmail.com This research was funded by SpectroChip, MacKay Memorial Hospital, Kaohsiung Medical University Chung-Ho Memorial Hospital, and the Ministry of Science and Technology, Taiwan (MOST 110-2218-E-011-008 and MOST 109-2327-B-037-001) None declared				
-	Bacl	kground:	The COVID-19 pandemic has spread globally in a short period of time. It is known that antibody (nAb) level can effectively predict vaccine efficacy, which leads to the exploration of vaccine trials for efficacy assessment. Thus, the current study aimed to develop a platform to quantify nAb levels faster, at lower cost, and with better efficiency.				
	Material/Methods: A total of 69 sera samples were collected for the research, 28 of which were from unvaccinated participants. The other 27 samples and the remaining 14 samples were from the participants who had received the first and second dose, respectively, of AZ vaccine 1 month before. With cPass assays (Genscript cPass nAb ELISA as- say) used as a criterion standard and lateral flow immunoassay kit (Healgen Scientific – LFIA test kit) coupled with a spectrometer (LFIA+S) for checking each specimen, we aimed to detect the presence of neutralizing an- tibodies in sera and to confirm the relationship between the inhibition rate from cPass assays and the nAb in- dex from the LFIA+S						
	Results: Data analysis of the research were taken from the certified ELISA and LFIA+S, which indicated a high of tanget (Paarcan's r=0.864, ICC=0.00139) between the 2 methods.						
	Conclusions: The dataset demonstrated that LFIA+S was affordable, had a strong correlation with results of the cPass nAbs detection kit, and has potential clinical applications, with an exclusive feature that allows non-experts to use it with ease. It is believed that the proposed platform can be promoted in the near future.						
	Ke	eywords:	Antibodies, Neutralizing • COVID-19 • Self Efficacy • Spectrum Analysis				
	Full-	text PDF:	https://www.medscimonit.com/a	bstract/index/idArt/	ndex/idArt/935812		
			■ 2887 ■ 21	⊥ _{2 10} □	วิ 35		



e935812-1

Background

Since the initial outbreak was reported, COVID-19 has been identified as a public health issue and has caused millions of infections and deaths globally [1-4]. However, there is a broad agreement that eliminating the virus is no longer feasible [5-7]. Increasing the vaccination rate has been recognized as the most effective measure against the global pandemic. Therefore, critical issues such as rapid evaluation of the efficacy of the vaccine are of urgent concern [8-10]. Recently, it has been reported that levels of neutralizing antibodies (nAbs) are highly predictive of immune protection from symptomatic SARS-CoV-2 infection. This is because nAbs can inhibit the interaction between the receptor-binding domain of the spike-1 protein (S1-RBD) and human angiotensin-converting enzyme 2 (h-ACE2) [11,12]. Currently, there are 2 standard methods available for nAbs detection. One is the 50% plaque reduction neutralization test (PRNT₅₀) [13-15], while the other is the FDA-approved SARS-CoV-2 Neutralization Antibody Detection Kit cPassTM (Genscript cPass neutralization antibody assay) [16-18]. In addition, which is a semi-guantitative detection method and the first surrogate neutralization assay to be commercially available. However, both methods mentioned above have difficulties, such as maintaining virus viability, biosafety level 3 (BSL 3) laboratory requirement, and being time-consuming to reveal the result [19-22].

On the other hand, the paper-based SARS-CoV-2 Neutralizing Antibody Test Kit (Healgen Scientific - LFIA test kit, Healgen Scientific, LLC, Houston, USA) shows several advantages, such as low cost, short processing time, and ease of use (it can be used by a non-expert), which make LFIA an ideal commercial kit to detect the nAbs [23,24]. However, as the LFIA test provides only qualitative results, the determination of nAbs level relies on the experimental operator's subjective cognition. Therefore, a novel platform that integrates LFIA with a portable spectrometer (LFIA+S) was proposed in this research to increase sensitivity and accuracy. Furthermore, the detection method was upgraded from qualitative level to semi-quantitative level [25,26].

To determine if LFIA+S is a highly reliable method for nAbs detection, statistical analyses were performed to compare the effects of the proposed integrated platform and the normal criterion standard cPass neutralization antibody assay [27]. We found that LFIA+S results are strongly correlated with those of cPass and LFIA+S only costs around US \$5 per test, which is affordable for general clinics or areas with insufficient medical care. Our aim was to provide a highly reliable platform to make it easier to detect the level l of nAbs for COVID-19 vaccine efficacy.

Material and Methods

Ethics Approval

The protocol for this research was approved by the Institutional Review Board (IRB) of MacKay Memorial Hospital (IRB No: 21MMHIS141e) and Kaohsiung Medical University Chung-Ho Memorial Hospital (IRB No: KMUHIRB-F(I)-20200141). All of the participants understood the purpose of the study and agreed to sign the participant consent form prior to their participation.

Serum Samples

A total of 69 serum samples were studied in a retrospective cohort study from June 2021 to November 2021. Results of the observational study conducted with medical personnel of MacKay Memorial Hospital and Kaohsiung Medical University Chung-Ho Memorial Hospital were included in this study. Among them, 28 serum samples were collected from the people who had not received vaccination yet; the other 27 sera and 14 sera were obtained from people who had received the first and second dose, respectively, of the AZ vaccine 1 month earlier, before the experiment.

Spectral Analysis System

The spectral analysis platform is equipped with a cassette designed for the LFIA test kit to detect the absorbance spectrum. The spectral chip provides high spectral resolution (3-5 nm), the results of which were demonstrated with a certain wavelength range (300 to 1100 nm). The real-time spectrum is used as the experimental reference for the primary absorbance wavelengths, which was to detect the use of this platform from 500 nm to 600 nm and 680 nm.

SARS-CoV-2 Neutralizing Antibody Lateral Flow Immunoassay and Detection (LFIA Test Kit)

A commercially available SARS-CoV-2 Neutralizing Antibody Test Kit (Healgen Scientific – LFIA test kit, Healgen Scientific, LLC, Houston, USA) was used to detect the presence of neutralizing antibodies in sera. First, 25 μ L of serum sample was loaded and mixed with 2 drops of the diluted buffer in the well for 10 min. The loaded test solutions can be driven through a conjugated pad by capillary force. As a competitive assay, the presence of neutralizing antibodies can be detected based on the signal intensity at the test line (T-line) and the control line (C-line), as shown in (**Figure 1**). This is because the neutralizing antibodies can inhibit the binding between the hACE-2 proteins (coated at T-line) and S1-RBD proteins (conjugated to colloidal gold nanoparticles and stored in a conjugation pad). In other words, the increase in neutralizing antibodies within a serum sample could lead to less binding of hACE-2 and



Figure 1. Lateral flow immunoassay test strip with C-line and T-line.



Figure 2. Lateral flow immunoassay+portable spectrometer (LFIA+S).

S1-RBD/gold nanoparticles and make the T-line lighter. The testing results were considered invalid if the C-line did not appear during the 20-min waiting time.

Novel Technique Platform of Lateral Flow Immunoassay Coupled with Spectrometer

In this research, a patent-pending testing platform was designed to demonstrate/display the result of LFIA, as shown in **Figure 2**. A specifically designed optical system was installed on the platform to reveal the absorbance spectrum of the gold nanoparticles (**Figure 3**) displayed on the test kit (T-line). The absorbance intensity (ΔA) of the gold nanoparticles can be calculated by the following formula and was later used to quantify the amounts of neutralizing antibodies in sera:

The signal intensity of the T and C-line can be calculated as ΔA T-line and ΔA C-line, respectively. The value of nAb can be derived by comparing the ΔA value of both test and control line:

$$n_{Ab} = 1 - (\frac{\Delta A_{T-line}}{\Delta A_{C-line}})$$

The amount of neutralizing antibody is proportional to the nAb value. Then, nAb value was transferred to the cPass ELISA inhibition percentage with a regression fitting curve. The formula of the regression curve is:

Inhibition percentage (%)=56.2873×n_{Ab}+73.5053

We compared the inhibition percentage derived from nAb with the cPass ELISA results.



Figure 3. Absorbance spectrum of gold nanoparticles.

SARS-CoV-2 Neutralizing Antibody ELISA test

A commercial FDA-approved SARS-CoV-2 neutralizing antibody ELISA kit (Genscript cPass neutralization antibody assay) as a criterion standard was compared with the LFIA kit. Serum samples and positive and negative controls were first diluted 10 times and individually mixed with a volume ratio of 1: 1 HRP-RBD solution. The mixed samples were then incubated at 37°C for 30 min. We loaded 100 µL of each mixed sample to a microplate well pre-coated with hACE-2 proteins. The microplate was then incubated at 37°C for 30 min and washed with washing buffer 4 times. We added 100 µL of TMB (3, 3', 5, 5'-Tetramethylbenzidine) substrate solution to each well in the microplate and then incubated it in the dark at room temperature for 15 min to start the enzymatic reaction. We added 50 µL of stop solution to each well, and the microplate was then read by an ELISA reader (Tecan Sunrise™). Since the neutralizing antibodies inhibit the bindings between RBD and hACE-2 proteins, the lower the optical density (O.D.) at 450 nm, the more neutralizing antibody there is. The absorbance of positive control must be less than 0.3, and the absorbance of negative control must be more than 1. The inhibition percentage is obtained from the following formula:

Inhibition percentage (%)=
$$(1 - \frac{OD_{sample}}{OD_{negative control}}) \times 100\%$$

The serum samples with less than 30% inhibition percentage were identified as negative cases.

Statistical Analysis

t Test

A 2-tailed *t* test was used to evaluate whether there was any statistically significant difference between the pre-vaccinated and vaccinated groups with the LFIA+S method. In this study, results yielding a *P* value of 0.05 were identified as borderline

of statistical significance, *P* values under 0.01 were statistically significant, and *P* values below 0.001 were highly statistically significant.

Receiver Operating Characteristic Curve

The receiver operating characteristic (ROC) is an effective way to evaluate the overall accuracy of the proposed LFIA+S method by calculating the area under the curve (AUC). In general, the value of AUC ranges from 0 to 1, where a values of 0.9 to 1 reflect an outstanding result, values of 0.8 to 0.9 is considered excellent, and 0.8 to 0.7 is acceptable.

Pearson's r

Pearson's r was used to evaluate the correlation between cPass ELISA and LFIA+S. Pearson's r between 0 to 0.25 was the absence of correlation, Pearson's r between 0.25 and 0.5 was poor correlation, Pearson's r between 0.5 and 0.75 was moderate to good correlation, and Pearson's r between 0.75 and 1 was very good to excellent correlation.

Bland-Altman Analysis

Bland-Altman analysis was used to evaluate the agreement between cPass ELISA and LFIA+S. The upper limit of the 95% confidence interval was defined by adding 1.96 times standard deviations (SD) to the mean difference between these 2 methods; the lower limit of the 95% confidence interval was defined by subtracting the 1.96 times SDs from the mean difference. Data with 95% confidence interval were agreements between these 2 methods.

Interclass Correlation Analysis

We used intraclass correlation coefficient (ICC) analysis to estimate the reliability between cPass ELISA and the proposed LFIA+S. Coefficient of ICC below 0.5 was poor reliability, coefficient between 0.5 and 0.75 was moderate reliability, coefficient between 0.75 and 0.9 was good reliability, and coefficient above 0.9 was excellent reliability.

Results

The results are shown in **Figure 4**. The increasing concentration of neutralizing antibodies could replace the S-RBD and bind with h-ACE2, which leads to a lighter color of the T-line. However, no obvious relationship was found between the antibody concentration and C-line color. Therefore, 4 categories were defined based on the visual observation: T-line deeper than C-line (negative), T-line as deep as C-line (weak positive), T-line lighter than C-line (medium positive), and T-line fade away (high positive).



Figure 5. Spectral results of the competitive test kit from different concentration of neutralizing antibodies specimen. Spectrum of C-line shown with green lines. Spectrum of T-line shown with red lines.

In addition, absorbance was utilized to reveal the spectrum result of the test kits. The spectral peak intensity ranges from 520 nm to 540 nm and depended on the colloidal gold, which is a determining marker for S-RBD. **Figure 5** illustrates the results of the spectral observation, which once again confirms that the deeper the color of the T-line (HPOS), the higher the

absorbance of the colloidal gold, and the C-line (NEG) does not consistently change with concentration variation.

After scanning the test kits with a spectrometer, an nAb index was calculated by a spectral result formula. To evaluate the diagnostic ability of nAb, we first considered those specimens

e935812-5



Figure 6. The *t* test nAb results of negative and positive sera are significantly different (*P* value <0.001).



Figure 7. ROC curve for index nAb. The area under the curve is 0.9901.

with an inhibition rate greater than 30% in the cPass kit as positive samples, and we viewed the others as negative samples. We used the *t* test to compare the nAb value between positive and negative samples. As shown in **Figure 6**, the *P* value was less than 0.001, which shows that these 2 independent groups were significantly different. Then, we calculated the area under the receiver operating characteristic curve (ROC curve). This value was greater than 0.9, which means that nAb is a good predictor for distinguishing negative results from positive ones (**Figure 7**).

We performed both cPass assays and lateral flow test for each specimen. To confirm the relationship between the inhibition rate from cPass assays and the nAb index from lateral flow test, we performed linear regression analysis, showing a correlation coefficient of 0.864 (**Figure 8**).



Figure 8. The linear regression fitting between cPass ELISA inhibition percentage and LFIA+S inhibition percentage.



Figure 9. The Bland-Altman analysis between cPass ELISA inhibition percentage and LFIA+S inhibition percentage.

We converted the inhibition percentage of cPass ELISA and LFIA+S. The Bland-Altman analysis showed a 94.2% (65 out of 69) confidence interval between cPass ELISA and LFIA+S (**Figure 9**). ICC analysis showed excellent reliability between the cPass and LFIA+S, and the coefficient was 0.90138, greater than 0.9 (**Figure 10**). With the information above, we are more confident that LFIA+S is a reliable solution to predict nAbs.

Discussion

With the outbreak of the COVID-19 pandemic and the threat of SARS-CoV-2 virus variants, there have been more reinfection cases [1,5,28,29]. Even fully vaccinated people still have a high risk of getting infections. It is a necessity that the vaccine effectiveness (VE) be tested, which can be analyzed by neutralizing antibody level (nAb level). As a result, it is important to



Figure 10. The interclass correlation analysis between cPass ELISA inhibition percentage and LFIA+S inhibition percentage.

detect the nAb level in sera after being vaccinated [7,30,31]. In currently used detection methods, PRNT₅₀ is a criterion standard for verifying the nAb, which was first used to detect neutralizing antibodies in 2014. Although it has higher accuracy and sensitivity than other tests, such as hemagglutination and many commercial enzyme immunoassays, the main disadvantage is that it takes days to get results [10,15,18,22,32]. In addition, because live viruses are used while conducting PRNT₅₀, it requires extremely strict operating procedure and high requirements for biosafety level and personnel operations. Therefore, it cannot completely meet a large number of testing requirements. To address the shortcomings of PRNT₅₀, the first USFDA-EUA-approved commercial neutralizing antibody ELISA kit (Genscript cPass neutralization antibody assay) has turned out to be a better alternative way to quantify nAb since 2020 [13,14,16]. However, cPass still needs to be conducted by specialists and can only be used in the laboratory, and it takes at least 2 h to complete the detection. According to the discussions above, it is known that either PRNT₅₀ or

cPass requires professional specialists to operate and a rigorous experimental environment, which needs a huge budget for large-scale testing. In the light of the defects of existing methods, the purpose of the current study was to develop a method that has the same high accuracy as cPass but is easier to use, quickly provides results, and is more affordable [33-35].

The LFIA+S method proposed in the current study is believed to be a more convenient way to get the nAb result. It takes only 3 min to perform, and it is easy to carry out because of the small size of the kit. In the current research, the LFIA+S method could significantly differentiate between the positive and negative sera at P<0.001. The area of the ROC curve was 0.9901, which indicated excellent accuracy. Comparing the proposed method with the cPass ELISA kit, Pearson's r was 0.864, which meant that these 2 methods are highly correlated. The 69 clinical sera data were transferred from nAb value to percentage inhibition with the regression fitting formula. When comparing these serum data with Bland-Altman analysis, the result showed that 65 sera were in the critical interval, and only 4 entries were out of the interval. The agreement between these 2 methods was 94.2% (65 out of 69). ICC analysis showed the coefficient of ICC of the 2 methods was 0.90138, which demonstrated excellent reliability. In terms of the cost, LFIA+S does not require expensive reagents, and it only needs colloidal gold test paper for detecting neutralizing antibodies, at US \$5 per test paper.

Table 1 provides a simplified comparison and analysis. The results of the statistical analyses show that LFIA+S has an accuracy highly similar to that of cPass in the detection of nAb, and also has high specificity and sensitivity. Our results show that LFIA+S has great potential to be a new reliable method for quantifying nAb.

This current study has several limitations. The first limitation was the design of the optical platform; all the databases we established for comparing the concentration of gold nanoparticles to absorption spectrum information can only be conducted using

Current methods	PRNT ₅₀	cPass	LFIA+S
Since year	2014	2020	2021
Country developed	America	America	Taiwan
Cost/each test	>US\$100/each test	>US\$100/each test	US\$5/each test
Time to get report	3-5 days	3-5 days	3 minutes
Country used	Global	Global	Taiwan
Accuracy	Very high	High	High
Operational complexity	Very high	High	Low
Laboratory requirements	Biosafety 3 level	Biosafety 2 level	Anywhere

 Table 1. Comparison of current methods for detecting neutralizing antibodies.

the LFIA test kit produced by Healgen Scientific. In other words, using testing kits of other brands requires rebuilding the database. To assess all brands will require a significant storage space that makes real-time data and cloud computing impractical. Furthermore, it the number of collected samples should be larger; the insufficient sample size turned out to be the primary obstacle to fine-tuning the platform sensitivity. In the future, it will be essential that more vaccinated people be recruited so as to provide medical institutions with more valid data on vaccine efficacy.

Conclusions

As COVID-19 vaccine coverage increases, it is essential to investigate whether the vaccine efficacy is sufficient to prevent spread of the pandemic. In reality, the existing methods for detecting nAbs level, $PRNT_{so}$ and cPass, rely heavily on professional personnel, strict experimental laboratory settings, and expensive detection costs (more than US \$100/per test), yet the testing efficiency is low and it takes 3~5 days to get results. However, the results of the current study demonstrated that LFIA+S has a high correlation with cPass; furthermore, LFIA+S has several other advantages that the existing methods lack. LFIA+S features low testing cost (less than US \$5/per test), higher effectiveness and efficiency (results can be shown in 3 min), and is more portable and easier to use. Ultimately, we hope that this study can serve as a

References:

- 1. Al-Rohaimi AH, Al Otaibi F. Novel SARS-CoV-2 outbreak and COVID19 disease; A systemic review on the global pandemic. Genes Dis. 2020;7(4):491-501
- 2. Su YJ, Chang CW, Chen MJ, Lai YC. Impact of COVID-19 on liver. World J Clin Cases. 2021;9(27):7998-8007
- 3. Bolcato M, Aurilio MT, Aprile A, et al. Take-home messages from the COVID-19 pandemic: strengths and pitfalls of the Italian National Health Service from a medico-legal point of view. Healthcare. 2021;9(1): 17.
- 4. Perrone MA, Feola A, Pieri M, et al. The effects of reduced physical activity on the lipid profile in patients with high cardiovascular risk during COVID-19 lockdown. Int J Environ Res Public Health. 2021;18(16):8858
- GeurtsvanKessel CH, Okba NMA, Igloi Z, et al. An evaluation of COVID-19 serological assays informs future diagnostics and exposure assessment. Nat Commun. 2020;11(1):3436
- Hansen CH, Michlmayr D, Gubbels SM, et al. Assessment of protection against reinfection with SARS-CoV-2 among 4 million PCR-tested individuals in Denmark in 2020: A population-level observational study. Lancet. 2021;397(10280):1204-12
- Bolcato M, Rodriguez D, Feola A, et al. COVID-19 Pandemic and equal access to vaccines. Vaccines. 2021;9(6):538
- Thomas SJ, Moreira ED, Kitchin N, et al. Six month safety and efficacy of the BNT162b2 mRNA COVID-19 vaccine. medRxiv. 2021;2021:21261159
- 9. Baraniuk C. What do we know about China's COVID-19 vaccines? BMJ. 2021;373:n912
- Ju B, Zhang Q, Ge J, et al. Human neutralizing antibodies elicited by SARS-CoV-2 infection. Nature. 2020;584(7819):115-19
- Haas EJ, Angulo FJ, McLaughlin JM, et al. Impact and effectiveness of mRNA BNT162b2 vaccine against SARS-CoV-2 infections and COVID-19 cases, hospitalisations, and deaths following a nationwide vaccination campaign in Israel: An observational study using national surveillance data. Lancet. 2021;397(10287):1819-29

valuable source that provides a new method for nAbs detection. Our proposed method shows great promise as a testing tool for evaluating the level of nAbs. We expect that the clinical application of LFIA+S can provide adequate data support for CDC to control the pandemic and implement reasonable vaccination standards.

Acknowledgements

The authors thank Dr. Yuan-Hui Chueh from National Tsing Hua University and Dr. Chao-Hung Kuo and Dr. Chih-Hsing Hun from Kaohsiung Medical University Chung-Ho Memorial Hospital for valuable technical discussions.

Statement

The staff of SpectroChip provided experimental technical support in this research. In addition, SpectroChip provided a spectrometer (the ONE InstantCare chromogenic rapid screening analyzer) for experimental use. SpectroChip had no other direct or indirect financial involvement in the study

Declaration of Figures' Authenticity

All figures submitted have been created by the authors, who confirm that the images are original with no duplication and have not been previously published in whole or in part.

- 12. Khoury DS, Cromer D, Reynaldi A, et al. Neutralizing antibody levels are highly predictive of immune protection from symptomatic SARS-CoV-2 infection. Nat Med. 2021;27(7):1205-11
- Frenck RW, Jr., Klein NP, Kitchin N, et al. Safety, immunogenicity, and efficacy of the BNT162b2 COVID-19 vaccine in adolescents. N Engl J Med. 2021;385(3):239-50
- 14. Bewley KR, Coombes NS, Gagnon L, et al. Quantification of SARS-CoV-2 neutralizing antibody by wild-type plaque reduction neutralization, microneutralization and pseudotyped virus neutralization assays. Nat Protoc. 2021;16(6):3114-40
- 15. Xie X, Nielsen MC, Muruato AE, et al. Evaluation of a SARS-CoV-2 lateral flow assay using the plaque reduction neutralization test. Diagn Microbiol Infect Dis. 2021;99(2):115248
- Taylor SC, Hurst B, Charlton CL, et al. A new SARS-CoV-2 dual-purpose serology test: Highly accurate infection tracing and neutralizing antibody response detection. J Clin Microbiol. 2021;59(4):e02438-20
- 17. Jung K, Shin S, Nam M, et al. Performance evaluation of three automated quantitative immunoassays and their correlation with a surrogate virus neutralization test in coronavirus disease 19 patients and pre-pandemic controls. J Clin Lab Anal. 2021;35(9):e23921
- Jung J, Rajapakshe D, Julien C, Devaraj S. Analytical and clinical performance of cPass neutralizing antibodies assay. Clin Biochem. 2021;98:70-73
- 19. Lau EHY, Tsang OTY, Hui DSC, et al. Neutralizing antibody titres in SARS-CoV-2 infections. Nat Commun. 2021;12(1):63
- Papenburg J, Cheng MP, Corsini R, et al. Evaluation of a commercial culture-free neutralization antibody detection kit for severe acute respiratory syndrome-related coronavirus-2 and comparison with an antireceptorbinding domain enzyme-linked immunosorbent assay. Open Forum Infect Dis. 2021;8(6):ofab220

e935812-8

Indexed in: [Current Contents/Clinical Medicine] [SCI Expanded] [ISI Alerting System] [ISI Journals Master List] [Index Medicus/MEDLINE] [EMBASE/Excerpta Medica] [Chemical Abstracts/CAS]

- 21. Nandakumar V, Profaizer T, Lozier BK, et al. Evaluation of a surrogate enzyme-linked immunosorbent assay-based severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) cPass neutralization antibody detection assay and correlation with immunoglobulin G commercial serology assays. Arch Pathol Lab Med. 2021;145(10):1212-20
- 22. Jung BK, Yoon J, Bae J-Y, et al. Performance evaluation of the BZ COVID-19 neutralizing antibody test for the culture-free and rapid detection of SARS-CoV-2 neutralizing antibodies. Diagnostics. 2021;11(12):2193
- Bao H, Yuan M, Xiao C, et al. Development of a signal-enhanced LFIA based on tyramine-induced AuNPs aggregation for sensitive detection of danofloxacin. Food Chem. 2022;375:131875
- Hnasko RM, Jackson ES, Lin AV, et al. A rapid and sensitive lateral flow immunoassay (LFIA) for the detection of gluten in foods. Food Chem. 2021;355:129514
- Chen P-Y, Ko C-H, Wang CJ, et al. The early detection of immunoglobulins via optical-based lateral flow immunoassay platform in COVID-19 pandemic. PLoS One. 2021;16(7):e0254486
- Hung KF, Hung CH, Hong C, et al. Quantitative spectrochip-coupled lateral flow immunoassay demonstrates clinical potential for overcoming coronavirus disease 2019 pandemic screening challenges. Micromachines (Basel). 2021;12(3):321
- Zou KH, O'Malley AJ, Mauri L. Receiver-operating characteristic analysis for evaluating diagnostic tests and predictive models. Circulation. 2007;115(5):654-57

- Tan ST, Chen TH, Yang HW, Su YJ. Changes in poisoning during the COVID-19 pandemic worldwide. Am J Emerg Med. 2021 [Online ahead of print]
- Su YJ, Lai YC. Comparison of clinical characteristics of coronavirus disease (COVID-19) and severe acute respiratory syndrome (SARS) as experienced in Taiwan. Travel Med Infect Dis. 2020;36:101625
- Bartsch SM, O'Shea KJ, Ferguson MC, et al. Vaccine efficacy needed for a COVID-19 coronavirus vaccine to prevent or stop an epidemic as the sole intervention. Am J Prev Med. 2020;59(4):493-503
- Edara VV, Hudson WH, Xie X, et al. Neutralizing antibodies against SARS-CoV-2 variants after infection and vaccination. JAMA. 2021;325(18):1896-98
- Girardin RC, Dupuis AP, Payne AF, et al. Temporal analysis of serial donations reveals decrease in neutralizing capacity and justifies revised qualifying criteria for coronavirus disease 2019 convalescent plasma. J Infect Dis. 2021;223(5):743-51
- Xu L, Li D, Ramadan S, Li Y, Klein N. Facile biosensors for rapid detection of COVID-19. Biosens Bioelectron. 2020;170:112673
- Miyakawa K, Jeremiah SS, Kato H, et al. Rapid detection of neutralizing antibodies to SARS-CoV-2 variants in post-vaccination sera. J Mol Cell Biol. 2022;13(12):918-20
- Moriyama S, Adachi Y, Sato T, et al. Temporal maturation of neutralizing antibodies in COVID-19 convalescent individuals improves potency and breadth to circulating SARS-CoV-2 variants. Immunity. 2021;54(8):1841-1852.e4