



Photoconversion of Alloreactive T Cells in Murine Peyer's Patches During Acute Graft-Versus-Host Disease: Tracking the Homing Route of Highly Proliferative Cells *In Vivo*

Katja J. Jarick^{1,2}, Zeinab Mokhtari¹, Lukas Scheller^{1,2}, Julia Hartweg^{1,2}, Sina Thusek¹, Duc-Dung Le¹, Maria Ranecky^{1,2}, Haroon Shaikh^{1,2}, Musga Qureischi^{1,2}, Katrin G. Heinze³ and Andreas Beilhack^{1,2,4*}

¹ Interdisciplinary Center for Clinical Research (IZKF) Laboratory for Experimental Stem Cell Transplantation, Department of Internal Medicine II, University Hospital, Würzburg, Germany, ²Graduate School of Life Sciences, University of Würzburg, Würzburg, Germany, ³Rudolf Virchow Center, University of Würzburg, Würzburg, Germany, ⁴Department of Pediatrics, University Hospital, Würzburg, Germany

OPEN ACCESS

Edited by:

Ulrike Koehl, Hannover Medical School, Germany

Reviewed by:

Mette D. Hazenberg, University of Amsterdam, Netherlands Benedetto. Bruno, Università degli Studi di Torino, Italy

> *Correspondence: Andreas Beilhack beilhack_a@ukw.de

Specialty section:

This article was submitted to Alloimmunity and Transplantation, a section of the journal Frontiers in Immunology

> Received: 22 January 2018 Accepted: 13 June 2018 Published: 27 June 2018

Citation:

Jarick KJ, Mokhtari Z, Scheller L, Hartweg J, Thusek S, Le D-D, Ranecky M, Shaikh H, Qureischi M, Heinze KG and Beilhack A (2018) Photoconversion of Alloreactive T Cells in Murine Peyer's Patches During Acute Graft-Versus-Host Disease: Tracking the Homing Route of Highly Proliferative Cells In Vivo. Front. Immunol. 9:1468. doi: 10.3389/fimmu.2018.01468 The regulation of immune cell migration throughout the body is essential to warrant immunosurveillance and to maintain immune homeostasis. Marking and tracking of these cells has proven important to study mechanisms of immune cell trafficking and cell interaction in vivo. Photoconversion is a well-suited technique for intravital application because it enables contactless time- and location-specific marking of cells in the tissue without surgically manipulating the microenvironment of the cells in question. However, in dividing cells the converted fluorescent protein may decline quickly. Here, we provide a detailed description of the photoconversion technique and its applicability to tracking highly proliferating T cells from the priming site of T cell activation to peripheral target organs of effector function in a preclinical model. Dendra2⁺ T cells were photoconverted in the Peyer's patches during the initiation phase of acute graft-versus-host disease (GvHD) and tracked through the mesenteric lymph nodes and the peripheral blood to the small intestine with flow cytometry and intravital two-photon microscopy. Photoconverted alloreactive T cells preserved the full proliferative capacity, homing, and migration of alloreactive T cells in the intestinal lamina propria. We conclusively proved that photoconversion of highly proliferative alloreactive T cells in the Peyer's patches is an effective tool to study trafficking of alloreactive T cells under physiologic conditions and to GvHD target tissues. This technique can also be applied to the study of immune cell tracking under inflammatory and non-inflammatory conditions.

Keywords: T cell migration, acute graft-versus-host disease, mouse models, photoconversion, Dendra2, Peyer's patch, *in vivo* cell tracking, lymphocyte homing

INTRODUCTION

Immune cells can migrate to distant locations within the body to warrant systemic immunosurveillance. This efficient immune cell homing is a prerequisite to protect against intruders or to regulate immune responses at any given location throughout the body. Specific immune surveillance becomes particularly impressive when comparing the volume of a T cell (35–95 fL) to the dimensions of the entire body and considering the small population size of a particular T cell clone with specific recognition capability (1). T cell homing plays an important role in many immunological reactions—therefore, delineating this process is central for understanding T-cell-mediated immunity. From a therapeutic perspective, it is essential to better understand the underlying mechanism of lymphocyte trafficking. This would allow specific blocking or fostering of homing routes in organ-specific inflammatory conditions such as autoimmune diseases, graft rejection after solid organ transplantation, or graft-versus-host disease (GvHD) following allogeneic hematopoietic cell transplantation.

Pioneering work in tracking of cell dynamics *in vivo* utilized transfer of radioactively labeled lymphocytes into rats, sheep, and other animals (2–6). Detecting radioactivity in different bodily fluids and organs indirectly proved the labeled cells' presence. Later, the discovery of congenic markers simplified transfer studies because it enabled the detection of adoptively transferred cells *in vivo* without needing to label the cells before transfer (7, 8). The introduction of biocompatible fluorescent dyes permitted tracking of fluorescently labeled, adoptively transferred cells and their division cycle *in vivo* (9). When using cells expressing fluorescent proteins, they require no labeling before transfer or before detection analysis at all (10, 11), and the fluorescence can even be induced at a specific time (12). Still, labeling dyes are widely used as a flexible means to label cells of interest before adoptive transfer (13).

However, none of these labeling methods are time- and site-specific at the same time, therefore it was not possible to mark specific populations of transferred cells after the transfer. Photoconversion is an excellent technique to mark cells *in vivo* in a specific location, because it enables contactless labeling without surgical manipulation of the organ of interest itself. This facilitates studying T cell homing without changing the homing properties of neither the studied cells nor the tissue surrounding the cells of interest. It is pivotal to minimize tissue perturbance and to avoid an experimental bias in the T cell homing process. By contrast, local introduction of a dye may lead to undesirable tissue perturbance and inflammation. Also, this would non-specifically label all cells present. Therefore, photoconversion is an expediant technique to introduce a time- and site-specific label specifically to transferred cells.

To date, several studies have employed photoconversion to monitor immune cell trafficking. The first (14) and most of the following studies employed the photoconvertible protein Kaede. These studies ranged from monitoring subcellular trafficking of single molecules (15) over organelles (16) to whole body trafficking of cell populations. Tracked cells comprised different precursor cells during embryogenesis (17, 18), as well as immune cell populations (19) and pathogens (20). However, there are two limitations when transferring the findings of the mentioned studies to tracking T cells *in vivo*.

First, due to the development of new photoconvertible proteins, there are now even better choices for *in vivo* imaging. Due to its higher pK, the Dendra2 protein is photoconverted 20 times more efficiently than other convertible fluorescent proteins like Kaede and mEOS (21). This renders Dendra2 more suitable for low-illumination and deep-tissue photoconversion. Gurskaya and colleagues derived the photoconvertible fluorescent protein Dendra from the octocoral *Dendronephthya sp* (22). A subsequent A224V substitution in this protein resulted in Dendra2 with improved maturation and brighter fluorescence both before and after photoswitching (23). Photoconversion irreversibly switches Dendra2 green fluorescence from green to red. However, cell metabolism and cell division dilute the converted protein over time, especially in highly proliferating cells.

Second, immune cell trafficking was up to now studied mostly in non-proliferating cells or under steady state conditions (24–26). More limited are studies providing data from infection models (27, 28) or cancer (29), where only a few specific T cell clones divide. Yet, at present, there are neither detailed reports on the photoconversion procedure nor systematic evaluation of possible side effects on T cell functions *in vivo*. Taken together, it remained unclear whether photoconversion would still be applicable to track an immune response with highly proliferative immune cells. Therefore, we tested whether we could employ photoconversion to track alloreactive T cells from the priming site to the peripheral target organ.

Here, we provide the first in-depth description of the intravital photoconversion procedure. We performed *in vivo* mouse experiments to study the homing from the Peyer's patches as priming sites to the mucosa of the small intestine as target site. We transplanted T cells expressing the photoconvertible protein Dendra2 in mitochondria (16). Subsequently, we photoconverted anatomical sites of initial T cell activation (30, 31) just before trafficking starts, then tracked the cells for the subsequent 24 h to the target site.

We show that photoconverted alloreactive T cells can be tracked from the Peyer's patches through the mesenteric lymph nodes and the peripheral blood to the small intestines as target site, using flow cytometry. Furthermore, we demonstrate that photoconversion is well suited to observe T cell migration in target tissues with intravital microscopy.

Photoconversion does not impede cell division, homing, and migration in alloreactive T cells. Hence it is a suitable technique to study trafficking of highly proliferative cells *in vivo*. This protocol provides insight into the procedures and caveats of photoconverting alloreactive T cells in the Peyer's patches during acute GvHD. The photoconversion technique opens new opportunities for effective and exact tracking of cell distribution *in vivo*. It enables researchers to prove directly the time- and location-specific origin of a marked cell in any given tissue analyzed.

MATERIALS AND EQUIPMENT

Materials

Materials are listed in **Table 1** corresponding to the experimental procedures.

Preparation of Reagents

 Anesthetic: add 2 ml of Ursotamin[®] (100 mg/ml, Serumwerk) and 2 ml of Xylavet[®] (20 mg/ml, CP-pharma) to 21 ml DPBS. Inject 10 μl/g body weight to reach desired concentrations (Ketamin 80 mg/kg, Xylazin 16 mg/kg).

TABLE 1 | Materials used for the photoconversion and associated techniques.

Dependables Section PDL 99500 PDL 995000 PDL 9950000 PDL 9950000 PDL 9950000 PDL 99500000 PDL 995000000 PDL 99500000000000000000000000000000000000	Material	Catalog #	Company
Transportation PD4 36200 Pan Botech Old Southout Control 2403 Harmmetein Call conting humber (Nuture Mouse Toells K() Type Hub solution Harmmetein Transport Busing Control Second Southout Southo	Dispensables		
DFRS windle Ca ¹¹ Mg1 ¹⁶ P04 3000P04 isolochCal contrig drawtor Nuclausel2403HaransteinT cell encharter kit (2011). CO11622, C0248, and114190Theme FelterTrain 19, publicade Linkouchel Mouse T cells (3)Theme FelterTheme FelterTrain 19, publicade Linkouchel Mouse T cells (3)Theme FelterTheme FelterTrain 19, publicade Linkouchel Mouse T cells (3)Theme FelterTheme FelterTrain 19, publicade Linkouchel Mouse T cells (3)Statistica (3)Theme FelterAmethetics SolutionThe surgest 2000 (3)The surgest 2000 (3)Theme FelterT rain surgest 2000 (3)Statistica (3)Statistica (3)The surgest 2000 (3)Theme FelterAndgefer Investignin2000 (1)NeolaBarnelTheme FelterAndgefer Investignin2000 (1)NeolaBarnelTheme FelterStatistica (3)Statistica (3)Statistica (3)Statistica (3)Statistica (3)Theme FelterSt	Transplantation		
Cal outsine Ki (C111b, C115k2, C145R, and 11413D Thems Fisher Tar-118, prabada Lincanch Mause Toells Kit) Tipon blue Studion Fachards Lincanch Mause Toells Kit Tipon blue Studion Fachards Lincanch Mater Fachards Lincanch Mause Toells Kit Tipon blue Studion Fachards Lincanch Mater Fachards Lincanch Mater Fachar	DPBS without Ca ²⁺ /Mg ²⁺	P04-36500	Pan Biotech
T adl archithment kil (2011b., CD1602-, CD48-, and 11413D Therm Pikher Trails - Dynoback binoched Mouse T cells K0 Filter 15, Dynoback binoched Mouse T cells K0 Train Structure Evidence Eviden	Cell counting chamber (Neubauer)	ZK03	Hartenstein
Jer 19. J. Jandenski Linkoufed Kouse I dels Kitj Tipan bla soliton Partonic Latter Acadebiai e soliton Pater and partalon Pater and partalon	T cell enrichment kit (CD11b ⁻ , CD16/32 ⁻ , CD45R ⁻ , and	11413D	Thermo Fisher
Ingan Due Souldon Experiment Subsection Series Areasite souldon Experiment (Suparition*) Bayer Experiment (Suparition*) Bootn Dickinson Bootn Dickinson In Instruct (Suparition*) Bootn Dickinson Bootn Dickinson In Instruct (Suparition*) Series Bootn Dickinson Series (Suparition*) Series (Suparition*) Bootn Dickinson Series	Ter-119 ⁻ , Dynabeads Untouched Mouse T cells Kit)		
Encomment Control Protection Specific	Trypan blue solution		
Protectioners Based on Diskingson Based on Diskingson Protectioners Based on Diskingson Based on Diskingson Protectioners Stand Diskingson Brain Protectioners Stand Diskingson Brain Protectioners 2500 (11 1) Noclab Stand Protectioners 704028 Koramik Schmidt Stand Protectioners 704028 Koramik Schmidt Stand Protectioners 704028 Koramik Schmidt Stand Protectioners Reserving Schmidt Koramik Schmidt Stand Protectioners Reserving Schmidt Koramik Schmidt Stand Protectioners Reserving Schmidt Reserving Schmidt Protectioners Reserving Schmidt Reserving Schmidt Stand Protectioners Reserving Schmidt Reserving Schmidt Protectioners Reserving Schmidt Reserving Schmidt Stand Protectioners Reserving Schmidt Reserving Schmidt Stand Protectioners Reserving Schmidt Reserving Schmidt Stand Protectioners Reserving Schmidt Reserving Schmidt P	Anosthatia aduitian		
For United Start Spartner (*) Bayer 1 nd springs 28GA x 29F (*) Add mm x 10 mm), BD Plasipak ^M 300015 Becton Dokinoon Areguint (*) Novelger) Samuel Samuel Pasteur Spects 2001111 Novelger) Satur Spects 704028 Samuel Satur Spects 704028 Bestpringer) Satur Spects 704028 Bestpringer) Satur Spects 704028 Califier Control Satur Spects 70304 Califier Control Satur Spects 70304 Califier Control Satur Spects 70304 Satur Spects Satur Spects Satur Spects 704040 Satur Spects Satur Spects <td>Photoconversion encretion</td> <td></td> <td></td>	Photoconversion encretion		
n m gering 2004, x (39 ° (x) (34 ° m x 10 °m), BD Plasipak [™] 200115 Barban 12 ° m haulin gering 2005, x (34 ° m x 12 ° m), Omnican* 100 9151141 Barban 12 ° m of the stand 12 ° m of the	Eve ointment (Benanthen®)		Baver
Initial spring 005x x12* (0.3 mm x12 mm), Omnean* 100 9151141 Beau Analgatis (Novalgim) 2001111 NeoLab Stanie roachs 2001111 NeoLab Stanie roachs 704028 Body products, Palax Phama i Neometh GmbH Stanie roachs 704028 Commeth GmbH Stanie roachs 704028 Commeth GmbH Stanie roachs Commeth GmbH Editional All Reacher Stanie Stanie roachs Stanie roachs Commeth GmbH Stanie roachs Editional All Reacher Stanie Commeth GmbH Stanie roachs Editional All Reacher Stanie Coll Reacher Stanie Stanie roachs Editional All Reacher Stanie Editional Coll Reacher Stanie Stanie roachs Stanie roachs Editional Coll Reacher Stanie Stanie roachs Stanie roachs Editional Reacher Stanie Stanie roachs Stanie roach Editional Reacher St	1 ml svringe 26GA x $3/8''$ (0.45 mm x 10 mm) BD Plastinak TM	300015	Becton Dickinson
Analgetic (Nondig)** Sanoti Sanoti Destaur pierter Postur pierter Nonda ba Sterie razors 704028 Body products. Palax Pharma. Sterie guaze sval, Bettack* ster. 2) 12780 Lohman & Bauscher Sterie guaze sval, Bodtabor) BH21 Destacher Bauscher Sterie guaze sval, Bodtabor) BARD Balaxcher Balaxcher Of to val Dodda00 BARPEN, Molitabor healthcare Balaxcher Of to val Dodda00 BARPEN, Molitabor healthcare Balaxcher Onichaber 70% expapated T331.3 Carl Both Gel Balax Doublaber 10% with beveled needle 3001.6 Bern Bestacher Bioscabit Bardoga patentis 760180, 607180, 607180, 606180 Grien Bio-one Acou API pro Discabit Bardoga patentis 04270000, 642752433, WR WR Discabit Bardoga patentis 04270000, 642752433, WR Discabit Bardoga patentis 1042714080, 94271102 Griene Bio-one Discabit Bardoga patentis 1042701 Griene Bio-one Dacot Dibordit Bardon <td>1 ml insulin svringe 30Gx x $1/2''$ (0.3 mm x 12 mm). Omnican[®] 100</td> <td>9151141</td> <td>Braun</td>	1 ml insulin svringe 30 Gx x $1/2''$ (0.3 mm x 12 mm). Omnican [®] 100	9151141	Braun
Parefar properties2600111NecLabStarle nazors70028800 products. Rate PramaraStarle nazors70028800 products. Rate PramaraStarle dissociality swab (Satpack*ize 2)127201014001Starle dissociality swab (Satpack*ize 2)127201014001Starle dissociality swab (Rotlabc')1421201401Of towal80043080043084AFIER, Moinicke healthcareOf towal80043080043084AFIER, Moinicke healthcareOr towal201509700110161Or towal2016087001Starle scape bladks, fasther #108851087401Disposible serological potters6000400Griener Elo-oneAcourter Tro60000600180Griener Elo-oneCall starler, 70 µm EASYstrainer^M542070Griener Elo-oneCall starler, 70 µm EASYstrainer^M542070Griener Elo-oneCall starler, 70 µm EASYstrainer^M222261Griener Elo-oneCall starler, 70 µm EASYstrainer^M542070Griener Elo-oneCall starler, 70 µm EASYstrainer^M542070Griener Elo-oneCall starler, 70 µm EASYstrainer542070Griener Elo-one	Analgetic (Novalgin®)		Sanofi
Sherie racons70428Body modules, Fielar Pharman, PSherie disserting swab (Stopck" size 2)12780Lohman A RauscherSherie disserting swab (Stopck" size 2)12780Lohman A RauscherSherie paus swab (Rollako T)800430EARPIER, Mohiloke healthcareO'P towal800430EARPIER, Mohiloke healthcareO'R towal70% 2 (propand)EARPIERCulckpad" 70% 2 (propand)3884085EARIANCulckpad" 70% 2 (propand)3884085EARIANCulckpad" 70% 2 (propand)3884085EARIANCulckpad" 70% 2 (propand)Status 6.0 (MB Augurt 70% 2 (propand)EhlorCulckpad" 70% 2 (propand)3884085EhlorStatus 6.0 (WB Augurt 70% 2 (propand)Status 6.0 (WB Augurt 70% 2 (propand)Disposable stockpad blacks, feather #10BBS10B.FrainDisposable stockpad blacks, feather #10200300, 640725433WRDisposable stockpad blacks, feather #10242700080, 64275433WRDisposable stockpad blacks, feather #10242700080, 6427543754,WRCall statuser, 70 µm EASYstrainer**242700Graine Elio-oneMortubes 1.5 ml27201Graine Elio-oneStatusTobe 50 ml27271Graine Elio-oneStatusFyrthorycon protoche yas buffer130-110-917Mittery ElootenCall strainer, 70 µm EASYstrainer**130-110-917Mittery ElootenCall strainer, 100 µm (MACS* Smart Strainers)15140-122Gibco ⁴ Call strainer21875-034Gibco ⁴ Call strainer<	Pasteur pipettes	2600111	NeoLab
Selied dissecting web (Selpack*ize 2)Kasmink KombH Shelie aguze swab (Gapack*ize 2)13985Lohman A Rauscher Shelie aguze swab (Gapack*ize 2)13985Lohman A Rauscher Shelie aguze swab (Gapack*ize 2)H12.1Cal Roth OP towal800430B4458ER, Moliikke healingare OP towal1391.3Cal Roth Op dispack*70% 2 prognardHelisch Medizinprodukte Gmbil Povitare iodine (Brauch*7.5% solution)3840.065BFAUN Stune, 6-0 with bewied needleV3016Ethicon Cal IsociatorB8510B. Braun Shelie scapite baddes, feather #10B8510Gener Elio-one Disposable serological pipettesQ270080, 047180, 606180Gener Elio-one Cal Istainer, 70 µm EASYstrainer™542070Genire Elio-one Cal Istainer, 70 µm EASYstrainer™542070Genire Elio-one Tube Somi272261Genire Elio-one Cal Istainer, 70 µm EASYstrainer™542070Genire Elio-one Tube Somi272261Genire Elio-one Cal Istainer, 70 µm EASYstrainer™542070Genire Elio-one Tube Somi272261Genire Elio-one Cal Istainer, 70 µm EASYstrainerM542070Genire Elio-one Tube Somi27200Genire Elio-one Tube Somi27201Genire Elio-one Cal Istainer, 70 µm EASYstrainerM542070Genire Elio-one Tube Somi27201Genire Elio-one Tube Somi27201Genire Elio-one Cal Istainer, 70 µm EASY	Sterile razors	704028	Body products, Relax Pharma u.
Skelie dissekting svab (kapta/k size 2)12780Lohman & RaukcherSkelie guate svab (kaptir 5 or x 5 orn)39696Lohman & RaukcherOf PowelBM12.1Carl RohOf PowelBM12.1Carl RohOdikapad "70% 2 propanol1931.3Carl RohOdikapad "70% 2 propanolStatute, or S ornBM201ER, Monitocia healthoursOdikapad "70% 2 propanol100.000BM201EROdikapad "70% 2 propanolStatute, or S ornBM201EOdikapad "70% 2 propanol100.00180, 600180Greiner Bio-oneStatute, or S orn26300, 607180, 600180Greiner Bio-oneDisposable service(carl properties)26300, 607180, 600180Greiner Bio-oneAccuriet "a norm26300, 642752433,WVRAccuriet "a norm242704050, 94271708, 342733754,WVRMicropipettas242070Greiner Bio-oneMicropipettas242704164, 942711302Greiner Bio-oneMicropipettas24270164, 942711302Greiner Bio-oneMicropipettas25272464, 942711302Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100.000180Greiner Bio-oneMicropipettas100			Kosmetik GmbH
Sherie gaure swab (Gaite '5 cm' s 5 cm)13895Lohman & BauscherSherie oton swab (Foliabo')EH12.1Carl RothO'P towal800430BARHIER, Molnicke healthcareO'R ethanol1931.3Carl RothQuickpad' 70% 2 propenolHoltsch Medizinprodukte GmbilQuickpad' 70% 2 propenolHoltsch Medizinprodukte GmbilQuickpad' 70% 2 propenolHoltsch Medizinprodukte GmbilQuickpad' 70% 2 propenolHoltsch Medizinprodukte GmbilCall BoaltonBB510B BraunCall BoaltonBB510Griener Bio-oneSterle scälpel blades, feather #10BB510Griener Bio-oneDisposable serological pipettes40270030Greiner Bio-oneAcou-jetr 'Bro20300BrantMicropipettes404720454, 402711302Cell strainer, 70 µm EAS'strainer ^{AM} 42000Greiner Bio-oneTabes 50 ml227261Greiner Bio-oneTabes 50 ml22706Greiner Bio-oneTabes 50 ml227061Greiner Bio-oneTabes 50 ml2175-034Greiner Bio-oneCell strainer, 100 µm (MACS' Smart Strainers)130-110-917Profiler 100 mos 50 ml21675-034<	Sterile dissecting swab (Setpack® size 2)	12780	Lohmann & Rauscher
Skrile cotton swab (Robibio') EH12.1 Carl Roth 0P lowel 800430 BARNER, Molinicke neshtaren 70% ethand 1931.3 Carl Roth Ouickpad" 70% 2-propanol B31.3 Carl Roth Povidone lother (Braund" 7.5% solution) 3864065 BRAIN Status, e.0 with bevelde needie 2016 Encon Cell isolation B5510 B. Braun Bisposable serological pipettes 760180, 607180, 605180 Greiner Bio-one Accurgit* pro 042760280, 642752433, WNR Microppettes 0427741788, 342733754, 042741788, 342733754, Microppettes 0427741788, 342733754, 042741788, 342733754, Microppettes 0427741788, 342733754, 042741788, 342733754, Cell istainer, 70 µm EASYstrainer™ 542070 Greiner Bio-one Tube 50 ml 227261 Greiner Bio-one Cell istainer, 70 µm (MACS* Smart Strainers) 130-110-917 Mittery Biotec Erythoropic lysis buffer Erythoropic lysis buffer Better Division Dissociation buffer Erythoropic lysis Coll face of Coll face of Coll face of Coll fa	Sterile gauze swab (Gazin [®] 5 cm × 5 cm)	13695	Lohmann & Rauscher
OP towal 600430 EARIET, Mönlick nealhtraar Op Sig ethanol 1931.3 Carl Roh Oukkpadr 70% 2-prognal Hotsch Madkingrodukte Gmbl Povidone loadne (Braunol 7.5% solution) 3846065 ERAUN Sature, 6-0 with beveled neede V301G Erbicon Cell isolatom Erbicon Erbicon Sature, 6-0 with beveled neede 700180, 607180, 606180 Gener Bo-one Disposable serological pipetes 700180, 607180, 606180 Gener Bo-one Call straine, 70 µm EASVstrainer™ 642700 Greiner Bo-one Micropipetes 642720454, 942711302 Greiner Bo-one Tube 50 ml 227206 Greiner Bo-one Tube 50 ml 22706 Greiner Bo-one Tube 50 ml 22706 Greiner Bo-one Tube for mages	Sterile cotton swab (Rotilabo®)	EH12.1	Carl Roth
70% stand T931.3 Call Roth Povidenci Odine (Graunoff 7,5% solution) 3884066. BRAUN Sture, 6.0 with bowled needle V301G Ethicon Cell isolation BB510 B. Braun Disposable scroppical pipettes 760180, 607180, 607180, 606180 Greiner Bio-one Accur-jeft pro 26300 Brandt Micropipettes 042760330, 642752433, WR 9427417488, 342733754, 9427417408, 342733754, 0427804780, 642752433, 042720454, 942711302 Greiner Bio-one 10455 0ml 222761 Greiner Bio-one 10456 0ml 2272044, 942711302 Cell strainer, 70 µm EASYstrainerM 1045 0ml 227261 Greiner Bio-one 10456 0ml 227261 Greiner Bio-one 10456 0ml 227261 Greiner Bio-one 10456 0ml 22781 Greiner Bio-one 10450 0ml 21875-034 Greiner Bio-one 1050010ml MACS® Smart Strainers) 130-110-917 Millenyi Biotec Profileration assay E F Periodianin, straptomycin (PenStrep) 15140-122 Glaco [®] 1041000 21875-034 Glaco [®] Politeration tracking dye (call trace violet) 33922.500 Sarstedt 1040	OP towel	800430	BARRIER, Mölnlicke healthcare
Quickparl 70% 2-proganol Hottsch Medizingrodukte Gmbil Povidene iodine (Raunoll [®] 7.5% solution) 3864065 BRALNN Stature, 60 with beveled needle V301G Bhioon Cell isolation BE510 B. Braun Disposable serological pipettes 760180, 607180, 606180 Greiner Bio-one Acu-jief" pro 26300 Brandt Outzroposality and the service of the servi	70% ethanol	T931.3	Carl Roth
Povidenci iddine (Braunoff 7.5% solution) 3884066 BRAUN Cell isolation Ethics (name in the version is the intervention of the	Quickpad® 70% 2-propanol		Holtsch Medizinprodukte GmbH
Suture, 6-0 with beveled needle V301G Ethicon Source scalped blades, feather #10 B8510 B. Braun Disposable servelogical pipettes 760180, 607180, 606180 Greiner Bio-one Accur-jet* pro 26300 Brandt Micropipettes 042705243, WR 042720454, 942711302 Greiner Bio-one Cell strainer, 70 µm EASYstrainer ^M 542070 Greiner Bio-one Micropipettes 042720454, 942711302 Greiner Bio-one Micropipettes 042720454, 942711302 Greiner Bio-one Micropipettes 042720454, 942711302 Greiner Bio-one Micropipettes 227281 Greiner Bio-one Strainer, 100 µm (MACS* Straat Strainers) 130-110-917 Mitterry Biotec Strainer, 100 µm (MACS* Straat Strainers) 130-110-917 Mitterry Biotec Periforation transbudde 21875-034 Gibco* Poliferation transbudde 21875-034 Gibco* Pendelling, straptomycin (PenStrep) 15140-122 Gibco* Fold works gredy (edi trace violet) 34557 Lif Echinhologies Anti	Povidone iodine (Braunol® 7.5% solution)	3864065	BRAUN
Cell isolation Strain Scala Bistop Bistop <th< td=""><td>Suture, 6-0 with beveled needle</td><td>V301G</td><td>Ethicon</td></th<>	Suture, 6-0 with beveled needle	V301G	Ethicon
Stein Scape Databs, Rainfer in 10 besite Disposable service Scape Databs, Rainfer in 10 besite Disposable service Disposable service Scape Disposable service Disposable dispo	Cell isolation		
Displase Profes Profes Orted Out Poil Out Poil Out Poil Micropipettes 042770030, 0427520430, WNP 042721054, 942711302 042711002 Out Poil	Sterile scalpel blades, feather #10	BB510	B. Braun
Accorge pro 20000 Linkin Micropipetties 0427200300, 642752433, WWR 942741768, 342733754, 042720454, 942711302 Call strainer, 70 µm EASYstrainer™ 542070 Greiner Bio-one Microtubes 1.5 ml 227261 Greiner Bio-one Microtubes 1.5 ml 227261 Greiner Bio-one S ml syringe BD Discardt™ II Becton Dickinson Call strainer, 100 µm (MACS® Smart Strainers) 130-110-917 Mittenyi Biotec Fythroxyte lysis buffer E E Dissociation buffer Proliferation assay I Proliferation assay E Gibco® Proliferation assay E Gibco® Proliferation assay E E Proliferation assay E Gibco® P	Accuriate pro	26200	Brandt
Mulcipipetities 04270306, 0427036, 042703764, 042720454, 942711302 042720454, 942711302 Cell strainer, 70 µm EASYstrainer™ 542070 Greiner Bio-one Tube 50 ml 227261 Greiner Bio-one Microtubes 1.5 ml 72.706 Sarstedt 5 ml syringe Boliscardit™ II Becton Dickinson Cell strainer, 100 µm (MACS® Smart Strainers) 130-110-917 Mittenyi Biotec Erythrocyte lysis buffer Dissociation buffer Profferation assay Profiferation assay 64257-034 Gibco® Pencillin, streptomycin (PenStrep) 15140-122 Gibco® Poliferation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2011) 553058 BD Pharmingen™ Recombinant human IL-2 589106 BioLegend Flow cytometry U U BioLegend VabitOby (Si-6.7) coupled to APC-Cy7 100722 BioLegend Anti-CD4 antibody (Si-6.7) coupled to APC-Cy7 100722 BioLegend Solure pruly asaya L220 BioLegend	Micropipettes	20300	
Od2720054, 942711302 Cell stainer, 70 µm EASYstrainer™ 542070 Greiner Bio-one Microtubes 1.5 ml 72.706 Sarstedt S ml syringe BD Discardt™ II Becton Dickinson Betton Dickinson Cell strainer, 100 µm (MACS* Smart Strainers) 130-110-917 Miltenyi Biotec Profileration assay Profileration assay Becton Dickinson RPMI Medium1640 21875-034 Gibco® Penicillin, streptomycin (PenStrep) 15140-122 Gibco® L-Gutamine 25030081 Gibco® Foliferation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human IL-2 639106 BioLegend Cytotoxidy saya 23101 BioLegend V-botom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (154-50 coupled to APC-Cy7 100722 BioLegend Vabitity dye (Zombie aqua) 423101 BioLegend Cytotoxidy sasay E E E Solum pyruvata, 100 nM, 100x 1380-039 BioLegend Cytotoxidy sasay E	Micropipeites	942741768 342733754	V V V 1
Cell strainer, 70 µm EASYstrainer™ 542070 Greiner Bio-one Tube 50 ml 227261 Greiner Bio-one Microtubes 1.5 ml 72.706 Sarstedt 5 ml syringe BD Discardit™ II Becton Dickinson Cell strainer, 100 µm (MACS* Smart Strainers) 130-110-917 Mitteryi Biotec Cell strainer, 100 µm (MACS* Smart Strainers) 130-110-917 Mitteryi Biotec Dissociation buffer T T Profiferation assay E E Profiferation sass Gibco* Gibco* Penicillin, streptomycin (PenStrep) 15140-122 Gibco* 1-Giutamine 25030081 Gibco* Fold greation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human IL-2 S9106 BioLegend Flow cytometry U U U U-bottom 96-well plate S022,500 Sarstedt Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Vability dye (Zombie aqua) 423101 BioLegend Ottototity assay E E D-Lucferin, firefly, potassium salt E220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco*		042720454 942711302	
Tube 50 ml 227261 Greiner Bio-one Microtubes 1.5 ml 72.706 Sarstedt 5 ml syinge BD Discardit™ II Becton Dickinson Cell strainer, 100 µm (MACS® Smart Strainers) 130-110-917 Mitteryi Biotec Erythrocyte lysis buffer Dissociation buffer F Dissociation buffer 72.706 Gibco® Proliferation assay E E Proliferation assay Gibco® Gibco® Penicillin, streptomycin (PenStrep) 15140-122 Gibco® 1-Gultarnine 25030081 Gibco® Fetal bovine serum 10270-106 Gibco® Proliferation tracking dye (cell trace violet) 34557 Life Technologies Arti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human L-2 Sey106 BioLegend Plow cytometr U U Differation differ U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD3 antibody (M4-5) coupled to APC-Cy7 100526 BioLegend Vability dye (Zombie aqua) 423101 BioLegend Vability dye (Zombie aqua) Life Technologies Gibco® Vability dye (Zombie aqua) BioLegend SolLegend Sodium pyruvate, 100 mM, 100x 1360-039	Cell strainer. 70 um EASYstrainer™	542070	Greiner Bio-one
Microtubes 1.5 ml 72.706 Sarstedt 5 ml syringe BD Discrit™ II Becton Dickinson Cell strainer, 100 µm (MACS® Smart Strainers) 130-110-917 Brythrocyte lysis buffer Iteravitation assay Proliferation assay Iteravitation assay RPMI Medium1640 21875-034 Gibco® Penicilini, streptomycin (PenStrep) 15140-122 Gibco® - Calutamine 25030081 Gibco® - Fetal bovine serum 10270-106 Gibco® Proliferation tracking dye (cell trace violet) 34557 Life Technologies Arti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human IL-2 589106 BioLegend Poilor (Combie aqua) 33322.500 Sarstedt Arti-CD3 antibody (63-6.7) coupled to APC-Cy7 100526 BioLegend Yability dye (Zombie aqua) 83392.500 Sarstedt Arti-CD3 antibody (63-6.7) coupled to APC-Cy7 100526 BioLegend Yability dye (Zombie aqua) 8220 BioLegend Vability dye (Zombie aqua) 8220 BioLegend Yability dye (Zombie aqua) 1360-039 Gibco® Yability dye (Zombie aqua) 1360-039 Gibco® Yability dye (Zombie aqua) 6005182 PerkInEImer <	Tube 50 ml	227261	Greiner Bio-one
5 ml syringe BD Discardit™ II Becton Dickinson Cell strainer, 100 µm (MACS® Smart Strainers) 130-110-917 Miltenyi Biotec Erythrocycle tysis buffer Dissociation buffer Profileration assay RPMI Medium1640 Penicilin, steptomycin (PenStrep) 15140-122 Gibco® Glutamine 25030081 Gibco® Fetal bovine serum 0270-106 Fold penicilin, steptomycin (VenStrep) 34557 Dislocation tracking dye (cell trace violet) 34557 Arti-CD3 antbody (145-2C11) 553058 Becombinant human IL-2 Bio2egend Flow cytometry U U-bottom 96-well plate 83.3922.500 Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend BioLegend Vability dye (Zombie aqua) 423101 Bick 96-well plate BioLegend D-Luciferin, firefly, potassium sait L8220 Bick 96-well plate BioSynth Sodium pryruate, 100 mM, 100x 11360-039 Gibco® Gibco® Sorgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 fexible needle holder, 1 needle holder Karl Hammacher GmbH and m 1 fexible needle holder, 1 needle h	Microtubes 1.5 ml	72.706	Sarstedt
Cell strainer, 100 µm (MACS® Smart Strainers) 130-110-917 Mittenyi Biotec Erythrocyte lysis buffer <	5 ml syringe BD Discardit™ II		Becton Dickinson
Erythrocyte lysis buffer Dissociation buffer Proliferation assay RPMI Medium1640 21875-034 Gibco® Penicillin, streptomycin (PenStrep) 15140-122 Gibco® 1-Glutamine 25030081 Gibco® Fetal bovine serum 10270-106 Gibco® Proliferation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2C11) 53058 BD Pharmingen™ Recombinant human IL-2 589106 BioLegend Flow cytometry U U Sarstedt Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD4 antibody (S3-6.7) coupled to PE-Cy7 100722 BioLegend Vabitry dye (Zombie aqua) 423101 BioLegend Vabitry dye (Zombie aqua) 423101 BioLegend Vabitry dye (Zombie aqua) BioLegend BioLegend Vabitry dye (Zombie aqua) CP-160 BioLegend Vabitry dye (Zombie aqua) Gibco® Sarstedt D-Luciferin, firefly, potassium sait Eg220 BioLegend	Cell strainer, 100 µm (MACS [®] Smart Strainers)	130-110-917	Miltenyi Biotec
Dissociation buffer Proliferation assay FMMI Medium1640 21875-034 Gibco® Penicillin, streptomycin (PenStrep) 15140-122 Gibco® Penicillin, streptomycin (PenStrep) 2503081 Gibco® Petal bovine serum 2503081 Gibco® Proliferation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human L-2 S99106 BioLegend Flow cytometry U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD4 antibody (G3-6.7) coupled to APC-Cy7 100526 BioLegend Anti-CD8 antibody (G3-6.7) coupled to APC-Cy7 100722 BioLegend Anti-CD8 antibody (G3-6.7) coupled to PE-Cy7 100722 BioLegend Sodium pyruvate, 100 mM, 100x Kingend	Erythrocyte lysis buffer		
Proliferation assayCillatorCillatorRPMI Medium 164021875-034Gibco®RPMI Medium 164021875-034Gibco®Penicillin, streptomycin (PenStrep)15140-122Gibco®1-Glutamine25030081Gibco®Fetal bovine serum0270-106Gibco®Proliferation tracking dye (cell trace violet)34557Life TechnologiesAnti-CD3 antibody (145-2C11)553058BD Pharmingen™Recombinant human IL-2589106BioLegendFlow cytometryUUSarstedtV-bottom 96-well plate83.3922.500SarstedtAnti-CD4 antibody (S3-6.7) coupled to APC-Cy7100526BioLegendAnti-CD4 antibody (S3-6.7) coupled to APC-Cy7100722BioLegendVability dye (Zombia equa)023101BioLegendVability dye (Zombia equa)6005182PerkinElmerD-Luciferin, firefly, potassium saltL3200BioSynthSodium pyruvate, 100 mM, 100x1580-039Gibco®Sodium pyruvate, 100 mM, 100xCP-160FaxitronSurgery tools: 2 fine forceps, 1 pair of small scissors, 1 fiscible needle holderKarl Hammacher GmbH and m1 fiscible needle holder11801001Laservision1 fiscible needle holderFisPLL051001LaservisionHigh-power UV LED lamp, fiber-coupledSilver LED-405 nmFizervision	Dissociation buffer		
RPMI Medium1640 21875-034 Gibco® Penicillin, streptomycin (PenStrep) 15140-122 Gibco® - Giltamine 25030081 Gibco® real bovine serum 10270-106 Gibco® Proliferation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2C1) 589106 BioLegend Recombinant human IL-2 589106 BioLegend Fbox cytometry U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (184-50) coupled to APC-Cy7 100526 BioLegend BioLegend Anti-CD4 antibody (53-6.7) coupled to APC-Cy7 100722 BioLegend Silegend Viability dye (Zombia aqua) 423101 BioLegend Silegend Viability dye (Zombia aqua) 6005182 PerkinEImer Silegend D-Luciferin, firefly, potassium salt L8200 BioLsgend Silegend Sodium pruvate, 100 mM, 100x 11360-039 Gibco® Silegend Yara irradiation source CP-160 Faxitron Silegend Surgery tools: 2 fine forceps, 1 pair of small scissors, 1 Karl Hammacher GmbH and m 118xibe needle holder	Proliferation assay		
Pencillin, streptomycin (PenStrep)15140-122Gibco®L-Glutarmine25030081Gibco®Proliferation tracking dye (cell trace violet)34557Life TechnologiesProliferation tracking dye (cell trace violet)34557Life TechnologiesAnti-CD3 antibody (145-2C11)553058BD Pharmingen™Recombinant human IL-2589106BioLegendFlow cytometryUUU-bottom 96-well plate83.3922.500SarstedtAnti-CD4 antibody (RM4-5) coupled to APC-Cy7100526BioLegendAnti-CD8 antibody (53-6.7) coupled to PE-Cy7100722BioLegendViability dye (Zombie aqua)423101BioLegendCytotxcity assayBioLegendSioSynthBlack 96-well plate6005182PerkinElmerD-Luciferin, firefly, potassium saltL8220BioSynthSodium pyruvate, 100 mM, 100×11360-039Gibco®Equipment1Kart Hammacher GmbH and m1 flexible needle holder, 1 needle holder1Faxitron1 flexible needle holder, 1 needle holder1Fi8P1L051001LaservisionViprotection glassesF18P1L051001LaservisionHigh-power UV LED lamp, fiber-coupledSilver LED-405 mmFizentixViprotection glassesF18P1L051001LaservisionHigh-power UV LED lamp, fiber-coupledSilver LED-405 mmFizentixHigh-power UV LED lamp, fiber-coupledSilver LED-405 mmFizentixHigh-power UV LED lamp, fiber-coupledSilver LED-405 mmFizentix	RPMI Medium1640	21875-034	Gibco®
L-Glutamine 25030081 Gibco [®] Fetal bovine serum 0200081 Fetal bovine serum 10270-106 Gibco [®] Proliferation tracking dye (cell trace violet) 34557 Life Technologies Anti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human IL-2 589106 BioLegend Flow cytometry U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD4 antibody (S3-6.7) coupled to APC-Cy7 100526 BioLegend Anti-CD4 antibody (S3-6.7) coupled to PE-Cy7 100722 BioLegend Yiability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay BioLegend Cytotoxicity assay D-Luciferin, firefly, potassium salt 6005182 PerkinElmer D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco [®] Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, 1 1 flexible needle holder, 1 1 flexible needle holder, 1 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH and m 1 flexible needle holder, 1 Uv protection glasses FisP1L051001 Laservision High-power UV LED lamp, fiber-coupled CP	Penicillin, streptomycin (PenStrep)	15140-122	Gibco®
Fetal bovine serum102/0-106Gibco®Proliferation tracking dye (cell trace violet)34557Life TechnologiesAnti-CD3 antibody (145-2C11)553058BD Pharmingen™Recombinant human IL-2599106BioLegendFlow cytometryUUU-bottom 96-well plate83.3922.500SarstedtAnti-CD4 antibody (RM4-5) coupled to APC-Cy7100526BioLegendAnti-CD8 antibody (53-6.7) coupled to PE-Cy7100722BioLegendVability dye (Zombie aqua)423101BioLegendCytotoxicity assayBioLegendCytotoxicity assayBlack 96-well plate6005182PerkinElmerD-Luciferin, firefly, potassium saltL8220BioSynthSodium pyruvate, 100 mM, 100x11360-039Gibco®Karly irradiation sourceSurgery tools: 2 fine forceps, 1 pair of small scissors, 1 flexible needle holder, 1 needle holderCP-1602 heating mats (20 cm x 30 cm)76085Trixie Heimtierbedarf GmbHUV protection glassesF18P1L051001LaservisionHigh-power UV LED lamp, fiber-coupledSilver LED-405 nmPrizmatix		25030081	Gibco®
Proliferation tracking dye (cell trace violet) 34557 Life lechnologies Anti-CD3 antibody (145-2C11) 553058 BD Pharmingen™ Recombinant human IL-2 589106 BioLegend Flow cytometry U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (M4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD8 antibody (53-6.7) coupled to APC-Cy7 100722 BioLegend Viability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay BioLegend Sologend Black 96-well plate 6005182 PerkinElmer D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, 1 fiex/sologend Yes Yes 1 fiex/ible needle holder, 1 needle holder Yes Yes 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED l	Fetal bovine serum	10270-106	
Anti-CDS antibody (145-2011) 553058 BD Pharmingen ^{1,44} Recombinant human IL-2 589106 BioLegend Flow cytometry U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD8 antibody (53-6.7) coupled to APC-Cy7 100722 BioLegend Viability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay Black 96-well plate 6005182 PerkinElmer Black 96-well plate 8005182 PerkinElmer D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, CP-160 Faxitron 1 flexible needle holder 1 flexible needle holder 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Proliferation tracking dye (cell trace violet)	34557	Life Technologies
Recombinant interact Discrete Flow cytometry BioLegend U-bottom 96-well plate 83.3922.500 Sarstedt Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD8 antibody (53-6.7) coupled to PE-Cy7 100722 BioLegend Viability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay BioLegend Cytotoxicity assay Black 96-well plate 6005182 PerkinElmer D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco* Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, 1 fine forceps, 1 pair of small scissors, 1 Fine Heimtierbedarf GmbH and m 1 flexible needle holder, 1 needle holder 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Anti-GD3 antibody (145-2G11)	530106	BD Pharmingen M
Provide of the field of the		369100	BioLegena
Orbitol if 30 weig plate 00.3322,300 Galated Anti-CD4 antibody (RM4-5) coupled to APC-Cy7 100526 BioLegend Anti-CD8 antibody (53-6.7) coupled to PE-Cy7 100722 BioLegend Viability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay 6005182 PerkinElmer Black 96-well plate 6005182 BioSynth D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder	Libottom 96-well plate	83 3022 500	Saretadt
Anti-CD8 antibody (53-6.7) coupled to PE-Cy7 100722 BioLegend Viability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay 6005182 PerkinElmer Black 96-well plate 6005182 BioSynth D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment C P-160 Faxitron X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Anti-CD4 antibody (BM4-5) coupled to APC-Cv7	100526	Biol egend
Viability dye (Zombie aqua) 423101 BioLegend Cytotoxicity assay 6005182 PerkinElmer Black 96-well plate 6005182 BioSynth D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Anti-CD8 antibody (53-6.7) coupled to PE-Cv7	100722	BioLegend
Cytotoxicity assay 6005182 PerkinElmer Black 96-well plate 6005182 PerkinElmer D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder	Viability dye (Zombie aqua)	423101	BioLegend
Black 96-well plate 6005182 PerkinElmer D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m Mathematical Science 1 flexible needle holder, 1 needle holder 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Cytotoxicity assay		Ŭ
D-Luciferin, firefly, potassium salt L8220 BioSynth Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder Faxitron 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Black 96-well plate	6005182	PerkinElmer
Sodium pyruvate, 100 mM, 100x 11360-039 Gibco® Equipment CP-160 Faxitron X-ray irradiation source Karl Hammacher GmbH and m 1 1 flexible needle holder, 1 needle holder Karl Hammacher GmbH and m 1 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	D-Luciferin, firefly, potassium salt	L8220	BioSynth
Equipment X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder 76085 Trixie Heimtierbedarf GmbH 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Sodium pyruvate, 100 mM, 100×	11360-039	Gibco®
X-ray irradiation source CP-160 Faxitron Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder 76085 Trixie Heimtierbedarf GmbH 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Equipment		
Surgery tools: 2 fine forceps, 1 pair of small scissors, Karl Hammacher GmbH and m 1 flexible needle holder, 1 needle holder Trixie Heimtierbedarf GmbH 2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	X-ray irradiation source	CP-160	Faxitron
1 flexible needle holder, 1 needle holder 2 heating mats (20 cm × 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	Surgery tools: 2 fine forceps, 1 pair of small scissors,		Karl Hammacher GmbH and mearo
2 heating mats (20 cm x 30 cm) 76085 Trixie Heimtierbedarf GmbH UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	1 flexible needle holder, 1 needle holder		
UV protection glasses F18P1L051001 Laservision High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	2 heating mats (20 cm × 30 cm)	76085	Trixie Heimtierbedarf GmbH
High-power UV LED lamp, fiber-coupled Silver LED-405 nm Prizmatix	UV protection glasses	F18P1L051001	Laservision
	High-power UV LED lamp, fiber-coupled	Silver LED-405 nm	Prizmatix
Powermeter PM1100 discontinued I horlabs	Powermeter PM100	discontinued	Thorlabs

TABLE 1 | Continued

	Photoconversion	n of F	P	During	GvHD
--	-----------------	--------	---	--------	------

Material	Catalog #	Company		
Equipment				
Glass fiber, 1,500 μ m, NA = 0.5		Prizmatix		
Collimator, 405 nm, $f = 4.02$ mm, NA = 0.6 SMA fiber collimation	F671SMA-405	THORLABS		
Laboratory stand with bosshead and ring clamp		VWR		
Infrared lamp	BF 27	Beurer		
Thermometer (dual thermo max/min)	E609790	Amarell Electronic		
Centrifuge (Megafuge 40R)		Thermo Scientific		
Laminar flow hood Hera safe	KS 18	Thermo Scientific		
CO ₂ incubators	150i	Thermo Scientific		
Water bath	WNB 14	Memmert		
FACS Canto II equipped with 405, 488, and 633 nm lasers and		Becton Dickinson		
a high-throughput sampler (HTS)				
Filter sets for FACS 488 channel: 735 LP + 780/60, 655LP + 760LP,				
610LP, 556LP + 585/42, 520LP + 530/30, 488/10; 405 channel:				
502 LP + 510/50, 450/50; 633 laser: 735LP + 780/60, 685LP, 660/20				
TrimScope II equipped with a titanium sapphire laser (Chameleon Ultra II,		Lavision Biotec		
Coherent), beam splitters at 500, 570, and 655 nm, bandpass				
filters 420/50, 535/50, 605/70, and three photomultipliers				
IVIS Spectrum	124262	PerkinElmer		
Software				
FlowJo	Version X	TreeStar		
Imaris	Versions 7.7.2 and 8	Bitplane		
Living image®	Version 4.0	PerkinElmer		
Matlab	Version R2016a	Version R2016a Mathworks		
Mice				
B6;129S-Gt(ROSA)26Sor ^{tm1.1(CAG-COX8A/Dendra2)Dcc} /J	018397	Jackson Laboratories		
BALB/c (BALB/cAnNCrl)	Strain Code 028	Charles River		
C57BL/6 (C57BL/6NCrl)	Strain Code 027	Charles River		

Materials listed once may be required also for other techniques explained.

- Analgetics: metamizol (Novalgin[®]): add 266 μl of Novalgin[®] to 100 ml of drinking water (1.33 mg/ml). In addition, mice may receive 200 mg/kg of Novalgin injection solution. Inject 10 μl of Novalgin injection solution per g body weight, prepared from 8 μl of Novalgin mixed with 192 μl of NaCl, subcutaneously at recovery from anesthesia.
- Erythrocyte lysis buffer: dissolve 89.9 g NH₄Cl, 10 g KHCO₃, and 0.37 g EDTA in 1 l of deionized water.
- Trypan blue solution: dissolve 1 g of trypan blue in 100 ml of PBS. Dilute 1:10 in PBS to receive working solution to mix at equal volume with cell suspension.
- Enrichment buffer: dissolve 0.5 g BSA and 0.375 g EDTA in 500 ml DPBS, sterile-filter.
- Dissociation buffer: per mouse, prepare 70 ml. Dissolve 55.836 mg of EDTA in 66.5 ml HBSS without Ca²⁺/Mg²⁺. Add 3.5 ml FCS to receive HBSS with 5% FCS, 2 mM EDTA.
- Proliferation medium: RPMI medium1640 with 10% FCS, 1% PenStrep solution, 30 ng/ml anti-CD3 antibody, and 50 IU/ml IL-2.

STEPWISE PROCEDURES

This protocol focusses on the *in situ* photoconversion technique and subsequent flow cytometry analysis of converted donor T cells in distant sites. Therefore, these sections are in greater detail, whereas additional procedures are outlined in brief.

This study was carried out in accordance with the recommendations and regulations of the Regierung von Unterfranken. The protocol was approved under the protocol number 55.2 2532-2-48.

Transplantation

Donor T Cell Enrichment

T cells were isolated from the spleen of Dendra²⁺ donor mice (16) at the age of 8–12 weeks. Splenocyte suspensions were enriched for T cells according to the manufacturer's protocol, counted by trypan blue exclusion, and adjusted to a density of 12×10^6 cells/ml in PBS. Typical T cell yields lay between 15 and 30% of the splenocyte input with a final T cell purity of 85–95%.

Donor Bone Marrow Isolation

Bone marrow cells were isolated from hind legs (femura and tibiae) of 8- to 12-week-old C57BL/6 mice. Cell numbers were determined by trypan blue exclusion, and the cell concentration was adjusted to 50×10^6 cells/ml. Typical bone marrow cell yields were $1-1.5 \times 10^8$ cells per mouse. Both the T cells and the bone marrow preparations can be stored at 4°C for up to 2 h until transplantation.

Transplantation

At the time point of transplantation, recipient mice should ideally have more than 20 g body weight to be sufficiently resilient to partial body weight loss. BALB/c mice at the age of 10–14 weeks were myeloablatively irradiated (8 Gy), and 1.2×10^6 Dendra2⁺ donor T cells were intravenously injected *via* the retro-orbital venous plexus together with 5×10^6 bone marrow cells in a total volume of 200 μ l PBS (**Figure 1A**). The drinking water was supplied with Baytril (Enrofloxacin, 0.05%) for 7 days after transplantation to avoid infections, and mice were monitored twice daily within this period. GvHD was scored clinically and body weight was measured daily.

Photoconversion

We photoconverted the cells 4 days after transplantation. This is the time point when the majority of donor cells are still in the secondary lymphoid organs and the first donor T cells leave the secondary lymphoid organs to enter the target tissues. This time point is optimal because short intervals between marking the cells and detecting them in the target tissue avoid dilution of the photoconverted protein before detection.

Preparation

Two sterile working spaces were prepared: one for animal preparation and one for the operation and photoconversion. For each of them, a heating mat was pre-warmed and covered with a sterile OP towel. All needed instruments and reagents were sterilized and laid out.

An operation cover was prepared by cutting an OP towel to a size of 15 cm \times 15 cm, and an oval window was cut into the center with a diameter of 1.5 cm \times 2.5 cm to be placed on the operation site.

A gauze swab was prepared by folding it to 7.5 cm \times 5 cm and incising a 2 cm long central slit.

The collimator was connected to the UV lamp *via* the glass fiber and fixed on the laboratory stand.

A clean recovery cage was placed under an infrared lamp, and the sensor of a thermometer was placed on the bedding.

Operation and Photoconversion

Mice were anesthetized by injecting 10 μ l anesthetic solution (Ketamine and Xylazine) per gram body weight intraperitoneally, and placed on the sterile, heated OP towel for preparation. Eyes were protected from dehydration by applying eye ointment. The ventral fur was wetted with 70% ethanol using a dissecting swab and shaved; the skin was disinfected with povidone iodine in an outward circling motion.

After 10 min, anesthetic depth was ensured to be stadium III.2 (surgical tolerance) by pinching the hind paw. If anesthesia was insufficient, further anesthetic was provided retro-orbitally. Not more than $20-50 \mu l$ was injected at a time to avoid over-dosage.

The prepared animal was transferred onto a fresh sterile and heated OP towel and the prepared OP-cover was applied. A 1–1.5 cm midline incision was made in the skin and in the peritoneum on the linea alba, no bleeding should be visible. The wound was covered with the prepared gauze swab, the slit overlaid on the incision. The swab was soaked with PBS to prevent drying of externalized intestinal tissue. Next, two cotton swabs were soaked with PBS and used to handle the intestinal contents. The swabs were inserted into the peritoneum to localize and gently exteriorize the cecum, which is located left craniolateral to the incision. The intestinal tissues were handled with great care to prevent a postoperative ileus. First, the cecum was placed on the gauze swab left lateral to the incision, and the intestinal tissue was kept moist at all times by dripping PBS on it using a Pasteur pipette. The small intestine was then gently pulled out, starting from the cecum, until the first Peyer's patch appeared (**Figure 1B**). This was placed on the gauze swab facing upwards and covered with sterile tin film. A central hole in the foil with a diameter of 2 mm exposed the Peyer's patch but covered the surrounding intestinal tissue.

The collimator of the UV lamp was placed at a distance of 5 mm over the Peyer's patch, and the tissue was illuminated for 2 min at maximum power (**Figure 1C**). An additional tinfoil shield protected the experimenters from scattered UV light, and UV protection glasses were used whenever the illumination was turned on.

The tin foil was removed, the intestine was gently pulled out until the next Peyer's patch appeared, and the procedure was repeated until six to seven patches were converted. In order to later be able to identify the converted patches during isolation, the distances in centimeter to the cecum and between the individual converted spots were documented for all converted patches.

After photoconversion, the small intestine and the cecum were gently placed back into the peritoneum as depicted in **Figure 1D**. Correct organization of intestinal loops was important to prevent a postoperative ileus. The peritoneum was closed with three non-consecutive stitches. Subsequently, the skin was closed with four stitches. The overall operation procedure took 30–45 min per animal. The animal was placed into a clean cage and held warm with an infrared lamp until recovered from anesthesia. The local temperature around the animal was monitored not to rise above 30°C to avoid overheating. When conscious, the animals were transferred back into their group cage. Analgesic was supplied in the drinking water.

Isolation of T Cells

At different time points after photoconversion (0, 4, 12, and 24 h), mice were sacrificed, and the blood, mesenteric lymph nodes, Peyer's patches, and small intestine were harvested. One spleen was isolated for control stainings.

Peripheral Blood

Mice were anesthetized with Ketamine and Xylazin. After 10 min, the anesthetic depth was ensured to be stadium III.2 (surgical tolerance) by pinching the hind paw. The skin of the animals was disinfected, the animals were fixed on a pad and ventrally opened. The intestine was laid to the animal's left to expose the left renal vein and the vena cava inferior. An insulin syringe containing 200 µl of lysis buffer was inserted into the vena cava just below the fusion with the renal vein (Figure S1 in Supplementary Material), with the opening facing cranial and ventral. Then, the needle opening was turned to the animal's left and gently pushed to the ventral side of the vena cava. Blood was collected very slowly to prevent collapsing of the vessel. Up to 700 µl were collected per animal. The erythrocytes were lysed in 30 ml of lysis buffer for 10 min at room temperature. The lysis was stopped by adding 20 ml of DPBS. The cells were pelleted for 5 min at $330 \times g$ at 4°C and resuspended in 5 ml of PBS. All processed cell suspensions can be stored for up to 4 h at 4°C in DPBS until further analysis.

Afterward, the mouse was killed by cervical dislocation. The intestinal tract was removed from duodenum to rectum and placed on PBS-moistened paper towels. The mesenteric lymph nodes were excised, and all converted Peyer's patches were



FIGURE 1 | Procedure of Peyer's patch photoconversion. (A) Time line of transplantation and photoconversion procedure. BALB/c recipients were maeloablatively irradiated (8 Gy) and transplanted with wild type bone marrow and Dendra2⁺ T cells. Donor T cells were photoconverted in the Peyer's patches 4 days after transplantation. (B) Photograph of an externalized Peyer's patch placed on a PBS-soaked sterile gauze. Cecum, colon, peritoneal opening, and small intestine are circled with black lines. (C) Application of a tinfoil stencil to expose a Peyer's patch to UV light illumination using a high-power LED coupled to a glass fiber. (D) The duodenum (a) and the stomach-proximal part of the jejunum (b) were not externalized when photoconverting the six distal Peyer's patches closest to the cecum. After photoconversion, (1) the proximal jejunum has to be placed cranial to the cecum. Subsequently, (2) the ileum (d) is placed opposite of the cecum, after which (3) the distal jejunum (c) is placed ventral of the ileum. Finally, the cecum is placed back at the location where it laid on the gauze.

excised generously and trimmed afterward. All non-converted Peyer's patches were removed from the small intestine.

Mesenteric Lymph Nodes

All surrounding tissue (fat and vessels) was carefully removed from the mesenteric lymph nodes. The lymph nodes were kept

in a 1.5 ml tube containing 200 μl of PBS for up to 30 min on ice until further processing.

Peyer's Patches

Peyer's patches surrounded by intestinal tissue were spread out on PBS-moistened paper towels, and all surrounding tissue was precisely cut off using a scalpel. Exact trimming was fundamental for a good T cell recovery. The patches were kept in a 1.5 ml tube containing 200 µl of PBS for up to 30 min on ice until further processing.

Mesenteric lymph nodes and Peyer's patches were thoroughly minced through a PBS-pre-wet 70 μ m strainer using the sterile plunger of a 5 ml syringe, and rinsed through with 5 ml of icecold PBS. Cells were pelleted at 330 × *g* for 5 min at 4°C, and the supernatant was removed. Lymph node cells were resuspended in the residual volume and directly transferred completely onto a 96-well plate for staining. Cells from the Peyer's patch were washed once with 5 ml of ice-cold PBS, centrifuged, and then transferred onto the 96-well plate in the residual volume.

Intestine

Cells were isolated from the intestine using a protocol modified from Geem and colleagues (32). The small intestinal tube was opened longitudinally and vigorously dragged through 300 ml of PBS to remove fecal contents. It was cut into 2 cm segments and transferred into a 50 ml tube containing 30 ml of dissociation buffer. The tube was horizontally fixed on an orbital shaker set to 250 rpm for 20 min at 37°C. The solution was passed through a 100-µm strainer into a 50 ml tube. The incubation and filter step were repeated with another 30 ml of dissociation buffer. The first dissociation fraction was centrifuged at $330 \times g$ for 5 min at 4°C and the pellet was washed in 30 ml of PBS. The two filtered cell suspensions were united and the PBS washing step was repeated. In initial experiments, the lamina propria fraction was additionally processed according to Geem et al. (32), but all donor T cells were found in the epithelial wash fraction (Figure S2 in Supplementary Material).

Spleen

A 70-µm cell strainer on a 50-ml Falcon was pre-wet with 2 ml of erythrocyte lysis buffer. The organ was placed on the strainer, incised crosswise a few times, minced with the plunger of a 5-ml syringe, and passed through the strainer by rinsing with 8 ml of erythrocyte lysis buffer. Erythrocytes lysis was stopped after 2 min by adding 10 ml of DPBS through the strainer. Cells were pelleted for 5 min at $330 \times g$ at 4°C and resuspended in 5 ml of PBS.

Flow Cytometry

Up to 1×10^6 cells were stained per well. Cells were resuspended in 100 µl blocking buffer and incubated for 5 min at 4°C. 100 µl of antibody mix was added, and cells were stained for 30 min at 4°C in the dark. Cells were pelleted at 330 × g for 5 min at 4°C and resuspended in 200 µl of PBS. The stained cells can be measured up to 6 h after staining when stored in the dark at 4°C. Longer pauses may require preceding fixation of the stained cells.

Fluorescence signal was acquired at 1,000–3,000 events/s, and washing wells were measured between different sample groups to exclude cross-contamination. The Dendra2 signal is extremely bright, especially in the FITC/GPF channel, and care was taken not to saturate any PMT. Table S1 in Supplementary Material contains a representative compensation matrix.

Donor T cells were identified by gating on lymphocytes, live cells, singlets, Dendra2-green, and/or Dendra2-red positive cells (Figure S3 in Supplementary Material).

Proliferation Assay

Splenocytes were isolated from Dendra2⁺ mice as described in Section Spleen. Half of the cells were photoconverted in a black 24-well dish for 2 min with a rotating motion of the glass fiber and without the collimator. 5×10^6 cells/ml were stained with 1 µl of cell trace violet for 6 min at room temperature. 500 µl of FCS and 5 ml of RPMI were added successively, and samples were incubated at 37°C in a water bath for 5 min. The samples were pelleted at 330 × g for 5 min at room temperature and resuspended in 5 ml of PBS. Cells were seeded on a 96-well round bottom plate at 200,000 cells/well in 200 µl of proliferation medium and left to proliferate for 3–4 days.

In Vivo Cell Migration

Mice were anesthetized and an intestinal loop was exposed as described in the photoconversion operation procedure (see Operation and Photoconversion). Instead of wetted gauze and OP towels, cling foil was used to avoid lint sources in the microscope. The mouse was positioned on a heating pad and the intestinal loop was positioned under a glass cover slip using two custommade holders. The tissue was kept moist using sterile 0.9% NaCl solution. Fluorophores were excited at a wavelength of 840 nm. The light intensity was increased as the square of penetration depth along the *z*-axis, between 5 and 30%. Images were acquired every 30 s in a field of view of 495 μ m \times 495 μ m with a resolution of 256 \times 256 pixels and a distance of 3 μ m between two z-planes.

Cytotoxicity Assay

Lymphocytes were isolated from the axillary, brachial, cervical, inguinal, and mesenteric LN and the spleen of naïve Dendra2+ mice. The containing T cells were stimulated in a mixed lymphocyte reaction (MLR) with allogenic BALB/c splenocytes that had been irradiated with 30 Gy or that had been T-cell-depleted. 8×10^{6} B6; Dendra²⁺ responder cells were incubated with 4×10^{6} BALB/c stimulator cells in 2 ml/well of RPMI with 10% FCS, 1% Pen/Strep, and 1 mM sodium pyruvate in 24-well plates for 4-5 days at 37°C, 5% CO₂. Cytotoxicity was assessed in round bottom 96-well plates with 50,000 Luciferase-expressing MOPC-315.BMP2.FUGLW target cells (33) and 50,000 to 5×10^{6} MLR effector cells containing cytotoxic T lymphocytes. Target cells only were used as a negative control, 0.1% Triton was used as a positive control. Survival of the target cells was quantified from the bioluminescence signal. The cells were pelleted at $330 \times g$ for 5 min and resuspended in 110 μ l of 2.75 mg/ml Luciferin in PBS. The cells were incubated for 3 min at room temperature and the signal was measured in a black 96-well plate. Percent of lysis was calculated with the formula (experimental signal - background signal)/(maximal signal – background signal) \times 100.

RESULTS

T Cells Are Photoconverted in Peyer's Patches

First, we determined the efficiency of *in situ* photoconversion of Dendra2⁺ T cells in the Peyer's patches. The LED UV lamp transferred 73 mW of UV light onto the Peyer's patch through a

stencil with a diameter of 2–3 mm (**Figure 2A**). The *in situ* photoconversion efficiently switches the majority of Dendra2⁺ donor cells in the Peyer's patch from green to red. About 90% (range: 75–98%) of donor cells isolated immediately after photoconversion contained photoconverted Dendra2 protein (**Figure 2B**). We detected photoconverted cells only in the Peyer's patches, and no photoconverted cells leaked into the blood or the small intestine at this time point (**Figure 2C**). Only single cells resided in the mesenteric lymph nodes immediately after the photoconversion procedure, which took up to 30 min per mouse. In conclusion, photoconversion of the Peyer's patches reliably marks the majority of currently residing Dendra2⁺ donor T cells in the converted tissue area.

T Cells Retain Their Proliferative Capacity

To address whether photoconversion of T cells would affect key biological functions such as T cell proliferation, we compared the proliferative capacity of photoconverted T cells versus untreated T cells *in vitro*. Both populations proliferated to the same extent (Figures 3A,B). To find out how long the photoconverted Dendra2-red remains detectable in rapidly dividing cells, we measured the decay of red fluorescent signal during proliferation. Stimulated splenocytes lost the photoconverted Dendra2-red signal at the same rate as the proliferation dye, with half of the red signal being retained after each division cycle. The photoconverted population separated from the untreated control for at least four divisions before the two populations merged (Figure 3C). In this setting, CD8 T cells had completed two to five or more cycles and proliferated more vigorously as compared to the CD4 T cells, which had divided one to four times after 4 days proliferation (Figure 3D). We measured the activation marker CD44 before and after stimulation on converted versus unconverted T cells. Photoconversion did not affect the expression of CD44 with or without stimulation after 4 days in culture (Figure 3E). Overall, photoconversion does not affect the cells' ability to proliferate. Notably, non-dividing photoconverted Dendra2-red⁺ T cells retained their fluorescence intensity at least 4 days in cell culture.

T Cells Traffic From Peyer's Patches to the Mesenteric Lymph Nodes, Blood, and Intestine

To evaluate whether photoconversion would affect T cell motility or efficient homing, we analyzed the homing behavior of photoconverted Dendra2⁺ donor T cells at different time points after photoconversion (**Figure 4**). After 4 h, first T cells had left the Peyer's patches and resided in the mesenteric lymph nodes, where 0.5% of donor lymphocytes originated from the photoconverted Peyer's patches. In the peripheral blood, only 0.04% of donor lymphocytes originated from the Peyer's patches.

12 h after photoconversion, cell numbers had increased in both sites, mesenteric lymph nodes and peripheral blood, as compared to 4 h after photoconversion, and reached even higher numbers after 24 h. Notably, the vast majority of donor T cells in the Peyer's patch still remained positive for the red fluorescent signal at this time point. We readily detected T cells in the small intestine after 24 h at a frequency of 4.3% (mean).



FIGURE 2 | Efficiency of Dendra2⁺ donor cell photoconversion in the Peyer's patches. (A) Output power of LED lamp as measured by power meter. The plot shows total output (without stencil) and applied power on Peyer's patch (with stencil). (B) Photoconversion efficiency of donor T cells in the Peyer's patch. Before photoconversion, 0% of the cells contained red protein, after photoconversion, 89% (mean) of the cells were red. (C) Representative FACS plots of donor cells isolated from the Peyer's patches, mesenteric lymph nodes, peripheral blood, and the small intestine before and right after photoconversion. Peyer's patches contain a clear photoconverted population after photoconversion, whereas all other organs contain unconverted cells at this time point.



proliferate equally well (non-parametric unpaired test, Mann–Whitney). (C) Flow-cytometric signal of the photoconverted protein plotted against the proliferation dye. Each cell division leads to a 50% reduction of fluorescence intensity of the photoconverted fluorescent protein Dendra2-red. (D) Comparison of CD4 and CD8 cell proliferation upon polyclonal activation of photoconverted Dendra2⁺ T cells. More rapidly dividing CD8 T cells lose faster their fluorescent cell trace violet (CTV) signal because they proliferate more than the CD4 T cells under this *in vitro* condition. (E) CD44 expression as a marker of T cell activation before and after stimulation. CD44 levels are comparable in photoconverted cells (red) versus not photoconverted cells (green). All exp: (n = 3 mice).

Photoconversion Does Not Impair T Cell Migration Within the Target Organ or Cytotoxicity

After entering the target organs by extravasation, alloreactive T cells migrate in the tissue to find their target cells and mediate the potentially lethal graft-versus-host response. We tested whether the photoconverted T cells performed equally well in migrating in the target organ, and found that the photoconverted cells were able to migrate as well as the green cells in the lamina propria of the small intestine (**Figure 5A**). Migration speed as a basic measure of the ability to migrate, and the distribution of turning angles as a measure of directionality were not affected by photoconversion

(**Figures 5B,C**). Furthermore, we assessed whether photoconversion affects the cells' cytotoxicity. The photoconverted cells performed similarly in lysing the target cells *in vitro* (**Figure 5D**). These results indicate that homing, migration, and effector functions of alloreactive T cells are preserved after photoconversion. We can track the cells on their way from the priming to the target site without perturbing the course of the immune reaction.

In summary, photoconverted T cells fully retained the ability to proliferate, to traffic to distant sites, and to migrate within tissues at a time point between the initiation and effector phase of acute GvHD. Therefore, we conclude that photoconversion is suitable to study adaptive immune responses such as the



FIGURE 4 | Time course of T cell homing from Peyer's patches through the mesenteric lymph nodes and the blood to the small intestine. (A) Peyer's patches contain high percentages of photoconverted cells within the 24 h period observed. (B) Single cells enter the mesenteric lymph nodes during the photoconversion procedure, and low numbers populate the mesenteric lymph nodes as early as 4 h after photoconversion. 24 h after photoconversion, on average 5.7% of donor cells are from the photoconverted Peyer's patches. (C) First photoconverted cells appear in the blood after 4 h and increase to an average of 0.78% (range: 0.04-1.41%) 24 h after photoconversion. (D) Cells were readily detected in the intestine 24 h after photoconversion. In some mice, single red cells show up in the intestinal sample, presumably due to contamination. One data point represents one mouse. Statistics are according to an unpaired non-parametric Kruskal–Wallis test, *p < 0.001, ***p < 0.001, and ****p < 0.0001.

homing and migration of heavily proliferative alloreactive T cells in GvHD.

DISCUSSION

Our findings demonstrate that Dendra2-expressing T cells can be efficiently photoconverted in the Peyer's patches during the induction phase of acute GvHD. Photoconversion yields proliferating cells that traffic *via* their regular homing route to the intestine. There, they migrate within the tissue to find their target antigen. Hence, photoconversion proved to be applicable to this highly proliferative cell population for at least four cell divisions. The main advantages of this technique are

- Contactless and efficient labeling of cells.
- Effective and exact tracking of cells in vivo.
- Direct evidence about time and location of origin in any given organ analyzed.

We photoconverted the cells toward the end of the initiation phase. This is the time when the first cells start entering the small intestinal lamina propria, and therefore, it is a suitable time point to track the early effector phase in the target organ. At the time point of photoconversion, more than 90% of cells in the Peyer's patches express CD44 and the majority have divided (34). We readily detected the first infiltrating cells in the lamina propria 24 h after photoconversion, which corresponds to day +5 after transplantation. At this time point, the donor cells had a larger size than when transplanted as measured by flow cytometry (data not shown), indicating a T cell blast phenotype. Since at comparable time points, almost all donor T cells proliferated (35), one can expect that more than 90% of the cells detected have undergone multiple divisions. This means that low numbers of proliferating cells can be detected with high specificity. Care has to be taken to exclude cell debris and dead cells in order to receive results that are biologically relevant (see also Figure S3 in Supplementary Material). Although there is a contamination of other cell types in the transplanted enriched T cell population (5-15%) which will also be photolabeled if present in the Peyer's patches, the vast expansion of the donor T cells in an allogenic setting outgrows the contamination by far. Furthermore, contaminating non-immune cells will most likely not home to the host lymphoid organs after transplantation. Therefore, we expect that all donor cells detected 5 days after transplantation are indeed donor T lymphocytes. The majority of these cells is still positive for the photoconverted



fluorophore at the time point of detection 24 h after photoconversion. Therefore, later time points during the effector phase are likely to yield similarly positive results, since more cells will leave the Peyer's patches during the course of disease.

It is not possible to determine the maximum time after conversion that converted cells remain detectable in the intestine. This is because the mice die from GvHD about 4 days after photoconversion (about 8 days after transplantation) in this model. Nevertheless cells proliferate greatly in GvHD. It is thus likely that the photoconversion signal would be diluted out after five division cycles in these cells (Figure 3). This is expected to be 24-48 h after photoconversion assuming a division rate of 6 h (36). This is distinctly different from the lack of dilution of photoconverted Dendra2 protein in a non-proliferative setting. In such steady state circumstances, it is possible to follow the fate of photoconverted cells for several months (28). A further problem arises if the cells to be converted are in large or highly pigmented organs. The converting light may then not penetrate deeply enough. Hence, it is necessary to ensure efficient penetration by the converting light before performing homing analyses to grant complete conversion of the residing cells. Furthermore, highly pigmented organs heat up more than normal tissues under the UV illumination. In such organs, careful measures are needed

to avoid tissue overheating and drying. Extended troubleshooting guidelines are displayed in **Table 2**. More factors may have to be adjusted when using other photoconvertible proteins which photoconvert less efficiently. Finally, there were no visible alterations of the microenvironment, but we cannot exclude minor phototoxicity to cell populations that were not relevant to this study. In brief, limitations of this technique are

- Limited accessibility of very thick or highly pigmented organs.
- Potential phototoxicity affecting cell populations that were not analyzed in this study.

Notably, we found that even immediately after photoconversion, single converted cells can be detected in the mesenteric lymph nodes. This was surprising given that the two organs are far apart. We postulate that these cells may have resided in the efferent lymphatics of one of the first converted Peyer's patches, and had traveled to the mesenteric lymph nodes within the 30-40 min of photoconversion and blood collection. Such rapid migration is entirely realistic assuming a lymph flow rate of $100 \ \mu m/s$ (38).

Our photoconversion protocol generated a high percentage of converted cells (90%). This is considerably higher than when the intestinal lymphocytes were converted endoscopically (30%) or

TABLE 2 | Troubleshooting guide.

Step	Problem	Possible reason	Solution
Anesthesia	Mice react to footpad pinching	Anesthesia underdosed	Slowly add small doses of anesthetic intravenously until desired anesthetic depth is achieved
	Mice take longer than 1 h to recover from anesthesia	Anesthesia overdosed	Do not use more anesthesia than needed to reach desired anesthetic depth
		Low body temperature	Monitor the body temperature with a rectal probe. Place the animal on a heating pad during anesthesia recover
Lethality up to 24 h Posta after operation ileus despite recovery from bowe anesthesia throu mous	Postoperative ileus (a postoperative ileus can be detected by widened bowel loops becoming apparent through the ventral skin of the mouse)	Rough handling of the bowel	Handle bowels gently and with care. Use pre-wetted cotton swabs to handle the bowels. Refrain from pinching the bowel with forceps. Do not apply pressure or tension to the bowel. Keep bowels lubricated at all times
	_	Bowels were not placed back into the peritoneal cavity in the right organization	Follow Figure 1D for placing back the intestinal loops. The intestinal loop organization may vary to some extend individually
	Peritonitis	Too intensive light exposure leading to necrosis which in turn leads to intestinal perforation	Reduce UV illumination intensity or exposure time to UV illumination. Carefully lubricate the tissue to avoid heat development
		Introduction of infectious agents during the operation procedure	Work aseptically according to (37)
		Introduction of infectious agents after wound closure	Suture both the peritoneum and the skin sufficiently
Cell detection	Only dead cells are detected in flow cytometry	Contaminating digestive enzymes or intestinal contents in the cell suspension after isolation	Wash the pellets of isolated cells with PBS one or two times after obtaining single-cell suspension
		Too intensive light exposure leading to apoptosis and necrosis	Reduce UV illumination intensity or exposure time of UV illumination. Carefully lubricate the tissue to avoid heat development
	No converted protein signal can be detected	Insufficient UV light exposure for photoconversion	Increase duration and/or intensity of UV illumination
		Incorrect wavelength of converting light	The Dendra2 protein converts optimally with 405 nm wavelength illumination
		Insufficient excitation of red protein form	Use a flow cytometer equipped with a 561 nm laser
			Increase the sensitivity of the PMT detecting the red fluorescent signal
		Photoconverted protein diluted out by proliferation	Measure at earlier time points after photoconversion
	Converted cells are detected in the intestine 0 h after photoconversion	Photoconversion is not confined to the Peyer's patches	Reduce the size of the tin foil stencil
		Peyer's patch tissue is included in small intestinal sample	Note down the position of converted patches during the conversion procedure. Carefully remove all Peyer's patches from the intestine before tissue digestion
			Discard the tissue closely surrounding the Peyer's patches

where cells in cervical lymph nodes were converted through the skin (46%) (24). The numbers are comparable to another study converting Dendra2 in the Peyer's patches (99%). In contrast to our study, their protocol included much shorter but more intense ultraviolet light exposure (27).

We conclude that photoconversion of highly proliferative alloreactive T cells in the Peyer's patches as presented in this protocol is a suitable and effective tool for studying physiological trafficking of T cells during inflammatory disease conditions such as acute GvHD. This protocol can be directly applied to the study of T cells under physiological steady state settings, in cancer models or under infectious conditions, to investigate cell migration. With some adjustments, it can also be applied to other cell types or to other animal models.

ETHICS STATEMENT

This study was carried out in accordance with the recommendations and regulations of the Regierung von Unterfranken. The protocol was approved under the protocol number 55.2 2532-2-48.

AUTHOR CONTRIBUTIONS

KJ designed, planned, and carried out the experiments, analyzed and interpreted the data, and wrote the manuscript. ZM analyzed and interpreted the cell migration tracks. LS, JH, and ST helped in the photoconversion and cell isolation. D-DL, HS, MR, and MQ helped in cell isolation. KH helped with the establishment of ultraviolet illumination and gave feedback on the project. AB designed the study, interpreted the data, and revised the manuscript. All authors read and discussed the manuscript.

ACKNOWLEDGMENTS

We thank Tim Steinfatt and Dr. Julia Delgado-Tascón for advice with the cell isolation and stimulation, and Mike Friedrich and Jürgen Pinnecker for help with the photometer and the two-photon microscopy setup. We thank Dr. Niklas Beyersdorf for advice on the mixed lymphocyte reactions. We would like to express our gratitude to Drs. Andreas Brandl, Julia Delgado-Tascón, Michael Sixt, Manfred Lutz, Eleni Koutsilieri, Andrew Davis, and the whole Beilhack lab for valuable advice and discussion.

FUNDING

This project was funded by unrestricted research grants to AB; KJ, ZM, and ST were funded by the Interdisciplinary Center for Clinical Research (IZKF) University of Würzburg, Project B233; LS, MR, MQ, JH, and KJ receive funds from the Graduate

REFERENCES

- Beilhack A, Rockson SG. Immune traffic: a functional overview. *Lymphat Res Biol* (2003) 1(3):219–34.
- 2. Gowans JL. The recirculation of lymphocytes from blood to lymph in the rat. *J Physiol* (1959) 146(1):54–69. doi:10.1113/jphysiol.1959.sp006177
- 3. Heath TJ, Lascelles AK, Morris B. The cells of sheep lymph. *J Anat* (1962) 96(Pt 3):397–408.392.
- Hall JG, Morris B. The lymph-borne cells of the immune response. Q J Exp Physiol Cogn Med Sci (1963) 48(3):235–47. doi:10.1113/expphysiol.1963. sp001660
- Hall JG, Smith ME. Homing of lymph-borne immunoblasts to the gut. *Nature* (1970) 226(5242):262–3. doi:10.1038/226262a0
- Halstead TE, Hall JG. The homing of lymph-borne immunoblasts to the small gut of neonatal rats. *Transplantation* (1972) 14(3):339–46. doi:10.1097/ 00007890-197209000-00009
- Shen FW, Saga Y, Litman G, Freeman G, Tung JS, Cantor H, et al. Cloning of Ly-5 cDNA. Proc Natl Acad Sci U S A (1985) 82(21):7360–3. doi:10.1073/ pnas.82.21.7360
- Sykes M, Chester CH, Sundt TM, Romick ML, Hoyles KA, Sachs DH. Effects of T cell depletion in radiation bone marrow chimeras. III. Characterization of allogeneic bone marrow cell populations that increase allogeneic chimerism independently of graft-vs-host disease in mixed marrow recipients. *J Immunol* (1989) 143(11):3503–11.
- Lyons AB. Analysing cell division in vivo and in vitro using flow cytometric measurement of CFSE dye dilution. *J Immunol Methods* (2000) 243(1): 147–54. doi:10.1016/S0022-1759(00)00231-3
- Chalfie M, Tu Y, Euskirchen G, Ward WW, Prasher DC. Green fluorescent protein as a marker for gene expression. *Science* (1994) 263(5148):802–5. doi:10.1126/science.8303295
- Manfra DJ, Chen S-C, Yang T-Y, Sullivan L, Wiekowski MT, Abbondanzo S, et al. Leukocytes expressing green fluorescent protein as novel reagents for adoptive cell transfer and bone marrow transplantation studies. *Am J Pathol* (2001) 158(1):41–7. doi:10.1016/S0002-9440(10)63942-9

School of Life Sciences, University of Würzburg. KJ and ST received funding from the EFRE Program of the European Union (Center for Personalized Molecular Immunotherapy), D-DL was funded by the DFG Transregio 221,B09 and ST by Z01. KH is funded by the Rudolf Virchow Center of the University of Würzburg.

SUPPLEMENTARY MATERIAL

The Supplementary Material for this article can be found online at https://www.frontiersin.org/articles/10.3389/fimmu.2018.01468/full#supplementary-material.

FIGURE S1 | Scheme displaying the collection of blood from the vena cava. After ventral opening of the mouse, the intestines were laid to the right and the vena cava was exposed. A 30 G syringe containing 200 µl of lysis buffer was inserted at the junction of the left renal vein with the inferior vena cava and blood is slowly collected from the vessel.

FIGURE S2 | Comparison of donor T cell yield from EDTA wash fraction and subsequent enzymatic digestion of the tissue. (A) Representative flow cytometry plots showing the percentage of donor cells within live single lymphocyte population. All donor cells in the intestine are washed out by the extensive washing steps (left), and no T cells are found after digesting the remaining tissue (right). (B) Quantification of donor T cells in the two fractions.

FIGURE S3 | Exemplary flow cytometry gating scheme to identify the percentage of photoconverted cells among live donor lymphocytes. From left to right: lymphocytes are gated by size and granularity, live cells are gated negative for dead stain, and single cells were selected by plotting height and area of the forward scatter. Donor cells were identified by being positive for the green and/or the red fluorescent Dendra2 protein, and percentages of cells with red Dendra2 fluorescence were quantified within this population.

- Foudi A, Hochedlinger K, Van Buren D, Schindler JW, Jaenisch R, Carey V, et al. Analysis of histone 2B-GFP retention reveals slowly cycling hematopoietic stem cells. *Nat Biotechnol* (2008) 27:84. doi:10.1038/nbt.1517
- Progatzky F, Dallman MJ, Lo Celso C. From seeing to believing: labelling strategies for in vivo cell-tracking experiments. *Interface Focus* (2013) 3(3):20130001. doi:10.1098/rsfs.2013.0001
- Tomura M, Yoshida N, Tanaka J, Karasawa S, Miwa Y, Miyawaki A, et al. Monitoring cellular movement in vivo with photoconvertible fluorescence protein "Kaede" transgenic mice. *Proc Natl Acad Sci U S A* (2008) 105(31): 10871–6. doi:10.1073/pnas.0802278105
- Woods E, Courtney J, Scholz D, Hall WW, Gautier VW. Tracking protein dynamics with photoconvertible Dendra2 on spinning disk confocal systems. *J Microsc* (2014) 256(3):197–207. doi:10.1111/jmi.12172
- Pham AH, McCaffery JM, Chan DC. Mouse lines with photo-activatable mitochondria to study mitochondrial dynamics. *Genesis* (2012) 50(11):833–43. doi:10.1002/dvg.22050
- Hatta K, Tsujii H, Omura T. Cell tracking using a photoconvertible fluorescent protein. Nat Protoc (2006) 1(2):960–7. doi:10.1038/nprot.2006.96
- Nowotschin S, Hadjantonakis A-K. Use of KikGR a photoconvertible green-to-red fluorescent protein for cell labeling and lineage analysis in ES cells and mouse embryos. *BMC Dev Biol* (2009) 9(1):49. doi:10.1186/1471-213X-9-49
- Nemet I, Ropelewski P, Imanishi Y. Applications of phototransformable fluorescent proteins for tracking dynamics of cellular components. *Photochem Photobiol Sci* (2015) 14(10):1787–806. doi:10.1039/C5PP00174A
- Müller AJ, Aeschlimann S, Olekhnovitch R, Dacher M, Späth, Gerald F, et al. Photoconvertible pathogen labeling reveals nitric oxide control of *Leishmania major* infection in vivo via dampening of parasite metabolism. *Cell Host Microbe* (2013) 14(4):460–7. doi:10.1016/j.chom.2013.09.008
- Adam V, Nienhaus K, Bourgeois D, Nienhaus GU. Structural basis of enhanced photoconversion yield in green fluorescent protein-like protein Dendra2. *Biochemistry* (2009) 48(22):4905–15. doi:10.1021/bi900383a
- 22. Gurskaya NG, Verkhusha VV, Shcheglov AS, Staroverov DB, Chepurnykh TV, Fradkov AF, et al. Engineering of a monomeric green-to-red photoactivatable

fluorescent protein induced by blue light. Nat Biotechnol (2006) 24:461. doi:10.1038/nbt1191

- Chudakov DM, Lukyanov S, Lukyanov KA. Using photoactivatable fluorescent protein Dendra2 to track protein movement. *Biotechniques* (2007) 42(5): 553, 555, 557 passim. doi:10.2144/000112470
- Morton AM, Sefik E, Upadhyay R, Weissleder R, Benoist C, Mathis D. Endoscopic photoconversion reveals unexpectedly broad leukocyte trafficking to and from the gut. *Proc Natl Acad Sci U S A* (2014) 111(18):6696–701. doi:10.1073/pnas.1405634111
- Hampton HR, Bailey J, Tomura M, Brink R, Chtanova T. Microbe-dependent lymphatic migration of neutrophils modulates lymphocyte proliferation in lymph nodes. *Nat Commun* (2015) 6:7139. doi:10.1038/ncomms8139
- Kitano M, Yamazaki C, Takumi A, Ikeno T, Hemmi H, Takahashi N, et al. Imaging of the cross-presenting dendritic cell subsets in the skin-draining lymph node. *Proc Natl Acad Sci U S A* (2016) 113(4):1044–9. doi:10.1073/ pnas.1513607113
- Schulz O, Ugur M, Friedrichsen M, Radulovic K, Niess JH, Jalkanen S, et al. Hypertrophy of infected Peyer's patches arises from global, interferonreceptor, and CD69-independent shutdown of lymphocyte egress. *Mucosal Immunol* (2014) 7(4):892–904. doi:10.1038/mi.2013.105
- Marriott CL, Dutton EE, Tomura M, Withers DR. Retention of Ag-specific memory CD4+ T cells in the draining lymph node indicates lymphoid tissue resident memory populations. *Eur J Immunol* (2017) 47(5):860–71. doi:10.1002/eji.201646681
- Torcellan T, Hampton HR, Bailey J, Tomura M, Brink R, Chtanova T. In vivo photolabeling of tumor-infiltrating cells reveals highly regulated egress of T-cell subsets from tumors. *Proc Natl Acad Sci U S A* (2017) 114(22):5677–82. doi:10.1073/pnas.1618446114
- Beilhack A, Schulz S, Baker J, Beilhack GF, Wieland CB, Herman EI, et al. In vivo analyses of early events in acute graft-versus-host disease reveal sequential infiltration of T-cell subsets. *Blood* (2005) 106(3):1113–22. doi:10.1182/blood-2005-02-0509
- Bäuerlein CA, Riedel SS, Baker J, Brede C, Garrote A-LJ, Chopra M, et al. A diagnostic window for the treatment of acute graft-versus-host disease prior to visible clinical symptoms in a murine model. *BMC Med* (2013) 11(1):134. doi:10.1186/1741-7015-11-134

- Geem D, Medina-Contreras O, Kim W, Huang CS, Denning TL. Isolation and characterization of dendritic cells and macrophages from the mouse intestine. *J Vis Exp* (2012) (63):4040. doi:10.3791/4040
- Trebing J, Lang I, Chopra M, Salzmann S, Moshir M, Silence K, et al. A novel llama antibody targeting Fn14 exhibits anti-metastatic activity in vivo. *MAbs* (2014) 6(1):297–308. doi:10.4161/mabs.26709
- Beilhack A, Schulz S, Baker J, Beilhack GF, Nishimura R, Baker EM, et al. Prevention of acute graft-versus-host disease by blocking T-cell entry to secondary lymphoid organs. *Blood* (2008) 111(5):2919–28. doi:10.1182/ blood-2007-09-112789
- Chung J, Ebens CL, Perkey E, Radojcic V, Koch U, Scarpellino L, et al. Fibroblastic niches prime T cell alloimmunity through Delta-like Notch ligands. J Clin Invest (2017) 127(4):1574–88. doi:10.1172/JCI89535
- Maury S, Salomon BT, Klatzmann D, Cohen JL. Division rate and phenotypic differences discriminate alloreactive and nonalloreactive T cells transferred in lethally irradiated mice. *Blood* (2001) 98(10):3156–8. doi:10.1182/blood. V98.10.3156
- Hoogstraten-Miller SL, Brown PA. Techniques in aseptic rodent surgery. *Curr Protoc Immunol* (2008) Chapter::12–.14. doi:10.1002/0471142735. im0112s82
- Blatter C, Meijer EF, Nam AS, Jones D, Bouma BE, Padera TP, et al. In vivo label-free measurement of lymph flow velocity and volumetric flow rates using Doppler optical coherence tomography. *Sci Rep* (2016) 6:29035. doi:10.1038/ srep29035

Conflict of Interest Statement: The authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Copyright © 2018 Jarick, Mokhtari, Scheller, Hartweg, Thusek, Le, Ranecky, Shaikh, Qureischi, Heinze and Beilhack. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.