


## **Protocol**

# Multiplexed proteomic profiling of cysteine reactivity and ligandability in human T cells



Differential amino acid reactivity with chemical probes can provide valuable information on the functionality and ligandability of proteins in native biological systems. Here, we present a quantitative, multiplexed chemical proteomic protocol for in-depth reactivity and ligandability profiling of cysteines in proteins in quiescent and stimulated T cells. This protocol illuminates dynamic immune state-dependent alterations in cysteine reactivity, revealing chemoselective and stereoselective small-molecule interactions with cysteines in structurally and functionally diverse proteins that lack chemical probes.

Ekaterina V. Vinogradova, Benjamin F. Cravatt

vinograd@rockefeller.edu (E.V.V.) cravatt@scripps.edu (B.F.C.)

#### Highlights

Cysteine reactivity profiling to compare biochemical changes across T cell states

Cysteine ligandability profiling to discover small-moleculeprotein interactions

Stereoisomeric covalent probes facilitate mechanismof-action studies

Vinogradova & Cravatt, STAR Protocols 2, 100458 June 18, 2021 © 2021 The Authors.

https://doi.org/10.1016/ j.xpro.2021.100458





#### Protocol

## Multiplexed proteomic profiling of cysteine reactivity and ligandability in human T cells

Ekaterina V. Vinogradova<sup>1,2,3,\*</sup> and Benjamin F. Cravatt<sup>1,4,\*</sup>

<sup>1</sup>Department of Chemistry, The Scripps Research Institute, La Jolla, CA 92037, USA

<sup>2</sup>Laboratory of Chemical Immunology and Proteomics, The Rockefeller University, 1230 York Avenue, New York, NY 10065, USA

<sup>3</sup>Technical contact

<sup>4</sup>Lead contact

\*Correspondence: vinograd@rockefeller.edu (E.V.V.), cravatt@scripps.edu (B.F.C.) https://doi.org/10.1016/j.xpro.2021.100458

#### **SUMMARY**

Differential amino acid reactivity with chemical probes can provide valuable information on the functionality and ligandability of proteins in native biological systems. Here, we present a quantitative, multiplexed chemical proteomic protocol for in-depth reactivity and ligandability profiling of cysteines in proteins in quiescent and stimulated T cells. This protocol illuminates dynamic immune state-dependent alterations in cysteine reactivity, revealing chemoselective and stereoselective small-molecule interactions with cysteines in structurally and functionally diverse proteins that lack chemical probes.

For complete details on the use and execution of this protocol, please refer to Vinogradova et al. (2020).

#### **BEFORE YOU BEGIN**

#### Naïve T cell isolation

© Timing: 2 h

- 1. Isolate peripheral blood mononuclear cells (PBMCs) over standard Lymphoprep (STEMCELL Technologies) gradient using slightly modified manufacturer's instructions.
  - a. Layer 25 mL of freshly isolated blood on top of 12.5 mL of Lymphoprep in a 50 mL Falcon tube avoiding mixing of blood with Lymphoprep.
  - b. Centrifuge the tubes at room temperature (931 g, 20 min, 23°C) with brake off.
  - c. Transfer the plasma and Lymphoprep layers containing PBMCs to new 50 mL Falcon tubes with a 2:1 dilution with PBS.
  - d. Pellet the cells (524 g, 8 min, 4°C) and combine the cells into 2 Falcon tubes by washing the rest of the tubes with PBS (10 tubes combined into 1 using 40 mL of PBS).
  - e. Combine pellets from 2 Falcon tubes into 1 Falcon tube with PBS (20 mL and 20 mL wash), count the cell number, and spin down the cells (524 g, 8 min, 4°C).

*Optional:* In case red blood cell contamination is significant (the pellet is completely red), use RBC lysis buffer (Thermo Fisher) before proceeding to the next step.

Isolate T cells from fresh PBMCs using EasySep Human T Cell Isolation Kit (STEMCELL Technologies, negative selection) according to manufacturer's instructions.







#### T cell activation for mass-spectrometry analysis (activated T cells)

#### © Timing: 3-4 days

- 3. Pre-coat non-tissue culture treated 6-well plates using a cocktail of  $\alpha$ CD3 (5  $\mu$ g/mL, BioXCell) and  $\alpha$ CD28 antibodies (2  $\mu$ g/mL, BioXCell) in PBS (2 mL/well) and keep at 4°C overnight (or 37°C for 2 h on the day of the experiment).
- 4. The next day, transfer the plates to a  $37^{\circ}$ C incubator for 30 min 1 h followed by a wash with PBS (2 × 5 mL/well).
- 5. Re-suspend freshly isolated T cells obtained in step 2 in RPMI media supplemented with 10% FBS, L-glutamine (2 mM), penicillin (100 U/mL), and streptomycin (100  $\mu$ g/mL) at 1 × 10<sup>6</sup> cells/mL, plate into the pre-coated 6-well plates (3–7 mL/well) and keep at 37°C in a 5% CO<sub>2</sub> incubator for 3 days.
- 6. Following this incubation period, combine the cells in 50 mL Falcon tubes, pellet (524 g, 5 min, 4°C), and wash with PBS (10 mL). Pellet the cells after the wash (524 g, 5 min, 4°C).
- 7. Transfer the cells into an Eppendorf tube in 1 mL of PBS, pellet and either re-suspend in fresh RPMI media for *in situ* treatments or flash-freeze and keep at -80°C until further analysis (*in vitro* treatments). The cell count typically increases 1.5–2 fold compared to the originally plated cells.

#### T cell expansion for mass-spectrometry analysis (control T cells)

#### © Timing: 13-15 days

- 8. Pre-coat a non-tissue culture treated 6-well plate with  $\alpha$ CD3 (1.5  $\mu$ g/mL) antibody in PBS (3 mL/well) and keep at 4°C overnight.
- 9. The following day transfer the plates to a  $37^{\circ}$ C incubator for 30 min 1 h followed by a wash with PBS (2 × 5 mL/well).
- 10. Re-suspend freshly isolated T cells isolated in step 2 in RPMI media (10% FBS, L-glutamine (2 mM), penicillin (100 U/mL), streptomycin (100  $\mu$ g/mL)), containing  $\alpha$ CD28 antibody (1  $\mu$ g/mL) at 1 × 10<sup>6</sup> cells/mL, plate into the pre-coated 6-well plate (3–7 mL/well) and keep at 37°C in a 5% CO<sub>2</sub> incubator for 3 days. Typically, 30–60 × 10<sup>6</sup> of freshly isolated T cells were selected for an expansion from each donor to afford 800–1600 × 10<sup>6</sup> expanded T cells on average.
- 11. Combine the cells in 50 mL Falcon tubes, pellet (524 g, 5 min, 4°C), and wash with PBS (10 mL).
- 12. Re-suspend the cells in RPMI media containing recombinant IL-2 (10 U/mL) at  $2.5 \times 10^5$  cells/mL and keep at  $37^{\circ}$ C in a 5% CO<sub>2</sub> incubator for 10–12 days, splitting the cells every 3–4 days with IL-2 containing media to keep cell density below  $2 \times 10^6$  cells/mL.
- 13. After 13–15 days of cell expansion (3 day pre-activation on the plate + 10–12 days of expansion), pellet the cells (524 g, 5 min, 4°C), and combine the cells into 2 Falcon tubes by washing the rest of the tubes with PBS.
- 14. Combine pellets from 2 Falcon tubes into 1 Falcon tube with PBS (20 mL and 20 mL wash), count the cell number, and spin down the cells (524 g, 5 min, 4°C).
- 15. Re-suspend the cells in fresh RPMI media for *in situ* treatments or flash-freeze and keep at  $-80^{\circ}$ C until further analysis (*in vitro* treatments).

Optional: For in situ treatments, re-suspend the cells in fresh RPMI media (3 × 10<sup>6</sup> cells/mL) and treat with stereoisomeric probes (5  $\mu$ M or 20  $\mu$ M, Figure 1) at 37°C in a 5% CO<sub>2</sub> incubator for 3 h. Following the treatment, pellet the cells (524 g, 5 min, 4°C), wash with PBS (10 mL), transfer into Eppendorf tubes (1 mL PBS), spin down (524 g, 5 min, 4°C), flash-freeze and keep at  $-80^{\circ}$ C until further analysis.

**Note:** T cell activation can be evaluated and monitored using flow cytometry-based analysis of activation markers CD25 and CD69 or ELISA analysis of secreted proinflammatory cytokines (e.g., IL-2, IFN- $\gamma$ ). Viability of cells can be measured using Trypan blue or flow cytometry-based approach using near-IR LIVE-DEAD stain (Invitrogen).



## Stereoisomeric acrylamide compounds

Figure 1. Stereoisomeric acrylamide electrophiles

Structures of four stereochemically defined acrylamide compounds, where only one compound (EV-96) stereoselectively inhibits T cell activation. The stereoisomeric relationships are shown in blue (diastereomers) and red (enantiomers). Red color in chemical structures indicates the acrylamide reactive group. (Adapted from Vinogradova et al., 2020))

Note: T cell viability at the end of expansion should be >80% as measured by Trypan blue.

#### **KEY RESOURCES TABLE**

Media should be stored at  $4^{\circ}$ C and warmed to  $37^{\circ}$ C prior to use. IL-2 should be added to the media right before use. Lymphoprep should be stored at  $4^{\circ}$ C and used cold for better separation during the blood layering step.

REAGENT or RESOURCE	SOURCE	IDENTIFIER
Antibodies		
PB anti-human CD4	BD Pharmingen	Cat#: 558116
PerCP/Cyanine5.5 anti-human CD8a	BioLegend	Cat#: 301032
FITC anti-human CD69	BioLegend	Cat#: 310904
PE anti-human CD25	BioLegend	Cat#: 302606
Monoclonal mouse anti-human CD3 (clone OKT3)	Bio X Cell	Cat#: BE0001-2
Monoclonal mouse anti-human CD28 (clone 9.3)	Bio X Cell	Cat#: BE0248
Biological samples		
Human blood	Scripps Research Normal Blood Donor Services	https://nbds.scripps.edu

(Continued on next page)



Continued		
REAGENT or RESOURCE	SOURCE	IDENTIFIER
Chemicals, peptides, and recombina	nt proteins	
Recombinant IL-2 (TECIN™ Teceleukin)	Hoffmann-La Roche, Inc.	Bulk Ro 23-6019
DMSO	Corning	Cat#: 25-950-CQC
Fetal Bovine Serum, Australian	Omega Scientific	Cat#: FB-21
Lymphoprep	STEMCELL	Cat#: 07861
EasySep <sup>TM</sup> Human T cell isolation kit	STEMCELL	Cat#: 17951
EDTA (0.5M, pH 8.0)	Invitrogen	Cat#: AM9260G
Desthiobiotin polyethyleneoxide iodoacetamide	Santa Cruz Biotechnology	Cat#: sc-300424
Urea	Fisher Scientific	Cat#: M1084871000
lodoacetamide	Sigma-Aldrich	Cat#: I1149-25G
Dithiothreitol (DTT)	Fisher Bioreagents	Cat#: BP172-25
Sequencing Grade Modified Trypsin	Promega	Cat#: V5111
Streptavidin Agarose Resin	Fisher Scientific	Cat#: 20349
Micro Bio-Spin column	Bio-Rad	Cat#: 7326204
Nonidet <sup>TM</sup> P40 substitute (Igepal <sup>TM</sup> CA-630)	USB Corporation	Cat#: 19628
Triethylammonium bicarbonate buffe	r Sigma-Aldrich	Cat#: T7408-500ML
TMT10plex™ Isobaric Label Reagent Set	Thermo Scientific	Cat#: 90406
Acetonitrile	VWR Intl	Cat#: BJAS017-0100
Hydroxylamine solution	Sigma-Aldrich	Cat#: 467804-10ML
Formic acid, ~98%, for mass spectrometry	Honeywell Fluka	Cat#: 94318-250ML-F
Acetonitrile LC-MS Grade	Fisher Scientific	Cat#: A955
Water LC-MS Grade	Fisher Scientific	Cat#: W6-4
Ammonium bicarbonate	Sigma-Aldrich	Cat#: 09830-500G
Protease Inhibitor Cocktail	Sigma-Aldrich	Cat#: 05892791001
Scout fragment KB02	Synthesized in-house	N/A
Scout fragment KB05	Synthesized in-house	N/A
Critical commercial assays		
Micro BCA™ Protein Assay Kit	Thermo Scientific	Cat#: 23235
DC Protein Assay Kit I	Bio-Rad	Cat#: 5000111
LIVE/DEAD™ Fixable Near-IR Dead Cell Stain Kit	Invitrogen	Cat#: L34976
Software and algorithms		
RAW Converter	Version 1.1.0.22; 2004 release	http://fields.scripps.edu/rawconv/
Integrated Proteomics Pipeline (IP2) and ProLuCID	Integrated Proteomics Applications, Inc	http://goldfish.scripps.edu/
FlowJo (v10.0.7)	TreeStar Inc.	https://www.flowjo.com/
Other		
ELISA microplates	R&D Systems	Cat#: 991427
'	3M Empore	Cat#: 143863
	Fisher Scientific	Cat#: E925000090
Osonica Probe Sonicator	Fisher Scientific	Cat#: 15338281

#### **MATERIALS AND EQUIPMENT**

IA-DTB stock solution (store in the dark at $-80^{\circ}$ C)		
Reagent	Final concentration	Amount
IA-DTB	10 mM	5 mg
DMSO	n/a	921.8 μL

## **Protocol**



DTT solution (make fresh before use)		
Reagent	Final concentration	Amount
DTT	10 mM	6.2 mg
Urea	9 M	2.2 g
TEAB buffer (1.0M)	50 mM	200 μL
LC-MS grade H <sub>2</sub> O	n/a	2 mL

lodoacetamide solution (make fresh before use)		
Reagent	Final concentration	Amount
lodoacetamide	500 mM	23.1 mg
LC-MS grade H <sub>2</sub> O	n/a	0.25 mL

TEAB buffer (Buffer A, store at 4°C)		
Reagent	Final concentration	Amount
TEAB (1.0 M, pH 8.5)	50 mM	1 mL
LC-MS grade H₂O	n/a	19 mL

EPPS buffer (Buffer B, store at 20°C-23°C)		
Reagent	Final concentration	Amount
EPPS	200 mM	2.52 g
LC-MS grade H <sub>2</sub> O	n/a	50 mL

#### $\triangle$ CRITICAL: EPPS pH needs to be adjusted after re-suspension.

Ammonium bicarbonate buffer (Buffer C, store at 20°C–23°C)		
Reagent	Final concentration	Amount
NH <sub>4</sub> HCO <sub>3</sub>	10 mM	1.58 g
LC-MS grade H <sub>2</sub> O	n/a	1000 mL

## $\triangle$ CRITICAL: Ammonium bicarbonate causes eye, skin, and respiratory tract irritation. Avoid inhalation and contact with eyes, skin, and clothing.

Trypsin solution (store at $-80^{\circ}$ C)		
Reagent	Final concentration	Amount
Trypsin	0.25 μg/μL	20 µg
CaCl <sub>2</sub> (100 mM)	25 mM	20 μL
Trypsin buffer	n/a	60 μL

2× wash buffer (store at 20°C-23°C)		
Reagent	Final concentration	Amount
TEAB (1.0 M, pH 8.5)	50 mM	2 mL
NaCl (3 M)	150 mM	2 mL
NP-40	0.2% (v/v)	80 μL
LC-MS grade H₂O	n/a	36 mL





Solvent A (store at 20°C-23°C)		
Reagent	Final concentration	Amount
LC-MS grade H <sub>2</sub> O	95% (v/v)	950 mL
LC-MS grade CH <sub>3</sub> CN	5% (v/v)	50 mL
Formic Acid	0.1% (v/v)	1 mL

Solvent B (store at 20°C–23°C)		
Reagent	Final concentration	Amount
LC-MS grade H <sub>2</sub> O	20% (v/v)	200 mL
LC-MS grade CH <sub>3</sub> CN	80% (v/v)	800 mL
Formic Acid	0.1% (v/v)	1 mL

Formic acid (store at 4°C)		
Reagent	Final concentration	Amount
LC-MS grade H <sub>2</sub> O	80% (v/v)	30 mL
Formic acid	20% (v/v)	7.5 mL

△ CRITICAL: Formic acid is highly irritating and corrosive. It can cause severe burns to eyes, skin, and the respiratory system. Wear appropriate protective gear (lab coat, gloves, eye goggles) and handle large volumes in a well-ventilated chemistry hood.

TMT tag stock solution (store in the dark at -80°C)			
Reagent	Final concentration	Amount	
TMT tag	20 μg/μL	5 mg	
Dry acetonitrile	n/a	250 μL	

 $\triangle$  CRITICAL: CH<sub>3</sub>CN is a flammable liquid that is harmful upon inhalation or in contact with skin and eyes. Wear appropriate protective gear (lab coat, gloves, eye goggles) when working with it.

#### **STEP-BY-STEP METHOD DETAILS**

The timing for each step is calculated for the preparation of one 10-plex mass-spectrometry sample (10 separate treatment channels, Figure 2).

#### Sample preparation and treatment

#### <sup>©</sup> Timing: 4 h

This step involves preparation of the proteome, treatment with the broad cysteine-reactive desthio-biotin-containing chemical probe (Figure 3), reduction and alkylation of free unreacted thiols.

- 1. Thaw the cell pellets on ice (calculate the amount of cells that will yield >1 mg of proteome), resuspend in cold PBS (550–600  $\mu$ L) and lyse by sonication (probe sonicator, 2 × 8 pulses, 40% duty cycle, output setting = 4).
- 2. Normalize protein concentration to 1.5-2.0 mg/mL using a standard DC protein assay (Bio-Rad).

#### **Protocol**






Figure 2. Ligandability profiling using scout fragments or stereoisomeric electrophiles

(A) Extended experimental workflow for chemical proteomic experiments measuring covalent electrophile-cysteine interactions in primary human T cells using isobaric tandem mass tags for mass differentiation and MS3-based quantification (TMT-ABPP).

(B) Schematic representation of several major steps in the protocol: (1) Protein precipitation after alkylation step; (2) TMT-tags used for peptide labeling and mass-differentiation; (3) 96-well plate used for high pH fractionation highlighting the wells that get combined into separate fractions; (4) IP2 analysis parameters for static and differential modifications.

*Optional:* Separate soluble and particulate proteomic fractions by ultracentrifugation  $(100,000 \, g, 45 \, min)$ .

Optional: Treat total proteome, soluble or particulate fractions with scout fragments KB02 or KB05 (Backus et al., 2016) (500  $\mu$ M final concentration) for 1 h at ambient temperature for broad ligandability analysis.

- 3. Treat the resulting proteomes (500  $\mu$ L, 1.5–2.0 mg/mL) with iodoacetamide polyethyleneoxide desthiobiotin (IA-DTB (Santa Cruz), 5  $\mu$ L of 10 mM stock in DMSO, final concentration: 100  $\mu$ M) for 1 h at ambient temperature.
  - $\triangle$  CRITICAL: IA-DTB is light sensitive, but can be stored for several months as DMSO aliquots at  $-80^{\circ}$ C in the absence of light in amber 1.5 mL Eppendorf tubes.
- 4. Add ice-cold MeOH (600  $\mu$ L), CHCl<sub>3</sub> (100  $\mu$ L), and H<sub>2</sub>O (100  $\mu$ L), vortex the mixture and centrifuge (10,000 g, 10 min, 4°C) to afford a protein disc at the interface of CHCl<sub>3</sub> and aqueous layers (Figure 2B).
- 5. Aspirate both layers without perturbing the disk.

△ CRITICAL: if the disk does not form properly, aspirate just the top layer.

6. Add cold methanol (1 mL) to the disk and vortex. Pellet the proteins (10,000 g, 10 min, 4°C), and re-suspend the resulting pellets in 90  $\mu$ L of buffer containing 9M urea, 10 mM DTT and 50 mM triethylammonium bicarbonate (1/20 dilution of 1.0 M stock solution, pH 8.5) by thorough pipetting up and down or sonication.






Figure 3. Chemical probes used in TMT-ABPP experiments to evaluate cysteine reactivity and ligandability Structures of electrophilic scout fragments (A) and broad reactivity probe (B) targeting cysteine. IA-DTB - Desthiobiotin polyethyleneoxide iodoacetamide.

 $\triangle$  CRITICAL: Avoid having excessive formation of bubbles, a cloudy solution is acceptable at this stage.

- 7. Heat the resulting mixture at  $65^{\circ}$ C for 20 min, then cool the sample to room temperature.
- 8. Add iodoacetamide (10  $\mu$ L, 500 mM solution in H<sub>2</sub>O), and incubate the samples in the dark at 37°C for 30 min with shaking.

 $\triangle$  CRITICAL: DTT and IA should be stored at 4°C and the stock solutions in water should be made fresh on the day of the experiment.

#### Trypsin digestion and streptavidin enrichment

© Timing: 5 h

This step describes trypsin digestion of the proteins, streptavidin enrichment of desthiobiotin-containing peptides, and elution of the enriched peptides from streptavidin resin.

- 9. Following the labeling with iodoacetamide, dilute the samples with 305  $\mu$ L of Buffer A (TEAB buffer, 50 mM, pH 8.5) to afford a 2.0 M final urea concentration, add trypsin (4  $\mu$ L of 0.25  $\mu$ g/ $\mu$ L trypsin in trypsin buffer, containing 25 mM CaCl<sub>2</sub>), and digest the proteins at 37 °C for 2 h.
  - △ CRITICAL: Dilution of urea to 2.0 M concentration is critical to ensure proper enzymatic activity of trypsin.

III Pause point: The tryptic digest step can be performed overnight in a 37°C incubator with agitation.

- 10. The following day, dilute the samples with  $2\times$  wash buffer, 400  $\mu$ L, 50 mM TEAB, pH 8.5, 150 mM NaCl, 0.2% NP-40), add streptavidin-agarose beads (50% slurry in  $2\times$  wash buffer) to each sample (50  $\mu$ L/sample) and rotate the bead mixture for 2 h at rt.
- 11. For a 10-plex sample, wash streptavidin-agarose bead slurry in  $2\times$  wash buffer (550  $\mu$ L, 50% slurry;  $2\times1$  mL, 50 mM TEAB, pH 8.5, 150 mM NaCl, 0.1% NP-40) and bring up to the initial volume in the  $2\times$  wash buffer prior to the addition to the sample.
- 12. After incubation, pellet the beads by centrifugation (2,000 g, 1 min), transfer the beads to BioSpin columns and wash (3  $\times$  1 mL wash buffer, 3  $\times$  1 mL PBS, 3  $\times$  1 mL H<sub>2</sub>O).
- 13. Elute the peptides by addition of 300  $\mu L$  of 50% aqueous CH<sub>3</sub>CN containing 0.1% FA.
- 14. Evaporate the eluate to dryness using SpeedVac vacuum concentrator.

△ CRITICAL: Excessive washes of the beads after the enrichment step are critical for removal of non-specific binding peptides, which detrimentally affect the coverage of labeled peptides.

III Pause point: The samples can be stored at  $-20^{\circ}$ C or  $-80^{\circ}$ C for several months.

#### Protocol



#### TMT tag labeling

#### © Timing: 2 h

- 15. Re-suspend the peptides in 100  $\mu$ L Buffer B (200 mM EPPS, pH 8.0) with 30% CH<sub>3</sub>CN, vortex and spin down (2,000 g, 1 min).
- 16. Add TMT tags (3  $\mu$ L/channel in dry CH<sub>3</sub>CN, 20  $\mu$ g/ $\mu$ L; Figure 2B) to the corresponding tubes and allow the reaction to proceed for 60 min.
- 17. Quench the reaction by the addition of 5% hydroxylamine (3  $\mu$ L per sample), vortex, spin down, and leave at room temperature for 15 min.
- 18. Add FA (5  $\mu$ L, Final concentration: 5% v/v) to each tube, vortex the tubes, spin down and combine all of the labeled channels in a single low binding 1.5 mL Eppendorf tube.
- 19. Dry the final combined sample in a SpeedVac vacuum concentrator and keep at  $-80^{\circ}$ C until the desalting and high pH fractionation step.

**Note:** The dry sample will look like viscous oil due to the presence of salts that will be removed in the following step.

△ CRITICAL: TMT tags can decompose over time. Re-suspend the solid tags in dry acetonitrile and avoid storing them at room temperature for prolonged time, as well as extensive number of freeze/thaw cycles. Avoid using Tris buffer, as TMT tags are unstable in the presence of primary amines.

#### Desalting and high pH fractionation

#### <sup>©</sup> Timing: 3 h

In this step the sample is desalted and fractionated to afford better LC/MS separation at the stage of the mass-spectrometry analysis.

- 20. Re-suspend the samples in solvent A (1 mL, 95%  $H_2O$ , 5%  $CH_3CN$ , 0.1% FA), acidify with 20  $\mu$ L of 20% FA in  $H_2O$ , and desalt by passing through Sep-Pak C18 cartridges (55–105  $\mu$ m, Waters) to avoid loading excessive salts on the HPLC column (300Extend-C18, 3.5  $\mu$ m, 4.6 × 250 mm).
  - a. Condition the cartridge with  $CH_3CN$  (3 × 1 mL);
  - b. Equilibrate the cartridge with solvent A (3  $\times$  1 mL);
  - c. Load the sample by dropwise addition and wash the cartridge with the flow-through from the sample (1  $\times$  0.5 mL);
  - d. Desalt by passing solvent A (3  $\times$  1 mL);
  - e. Elute the peptides by adding 1 mL of solvent B (80% CH<sub>3</sub>CN, 20% H<sub>2</sub>O, 0.1% FA);
- 21. Evaporate the eluent to dryness using SpeedVac vacuum concentrator.

Note: At this stage there should be no visible salts in the Eppendorf tube.

III Pause point: The samples can be stored at  $-20^{\circ}$ C or  $-80^{\circ}$ C for several months.

- 22. Re-suspend the sample in solvent A (500  $\mu$ L) with sonication and fractionate into a 96 deep-well plate using HPLC (Agilent).
  - a. Add 20  $\mu$ L of 20% FA to each well in a 96-well collection plate to acidify the eluting peptides.
  - b. Elute the peptides onto a capillary column (ZORBAX 300Extend-C18, 3.5 μm) and separate at a flow rate of 0.5 mL/min using the following gradient: 100% Ammonium bicarbonate (Buffer C; 10 mM) from 0–2 min, 0%–13% acetonitrile from 2–3 min, 13%–42% acetonitrile from 3–60 min, 42%–100% acetonitrile from 60–61 min, 100% acetonitrile from 61–65 min, 100%–0% acetonitrile from 65–66 min, 100% buffer C from 66–75 min, 0%–13% acetonitrile



from 75–78 min, 13%–80% acetonitrile from 78–80 min, 80% acetonitrile from 80–85 min, 100% buffer C from 86–91 min, 0%–13% acetonitrile from 91–94 min, 13%–80% acetonitrile from 94–96 min, 80% acetonitrile from 96–101 min, and 80%–0% acetonitrile from 101–102 min. Use the following time-based fraction collection: 5–75 min, 96 fractions (collection mode: A1-A12, B1-B12, C1-C12, D1-D12, E1-E12, F1-F12, G1-G12, H1-H12).

- c. Evaporate the eluent to dryness in the plate using SpeedVac vacuum concentrator.
- 23. Re-suspend the peptides in solvent B (80%  $CH_3CN$ , 20%  $H_2O$ , 0.1% FA; 75  $\mu$ L/well) and combine all the wells from each column in the 96-well plate (Columns 1–12 become fractions f1–f12) into a new clean Eppendorf tube (1.5 mL, 12 fractions total, Figure 2B), following an additional wash with the same solvent (200  $\mu$ L/column, washing wells in each row).
- 24. Remove the solvent using SpeedVac vacuum concentrator.

III Pause point: The samples are ready for mass-spectrometry analysis and can be stored at  $-80^{\circ}\text{C}$  for several days to months.

#### TMT-ABPP liquid chromatography mass-spectrometry (LC-MS/MS/MS) analysis

This step describes LC-MS/MS/MS analysis conditions for each proteomic fraction and the Orbitrap Fusion mass-spectrometer settings for the analysis.

- 25. Re-suspend the resulting 12 combined fractions (Figure 2B) in solvent A (10 or 14  $\mu$ L, depending on the injection loop size) and analyze on the Orbitrap Fusion mass-spectrometer (5  $\mu$ L or 10  $\mu$ L injection volume).
- 26. Analyze the samples by liquid chromatography tandem mass-spectrometry using an Orbitrap Fusion mass spectrometer (Thermo Scientific) coupled to an UltiMate 3000 Series Rapid Separation LC system and autosampler (Thermo Scientific Dionex).
  - a. Elute the peptides onto a capillary column (75  $\mu$ m inner diameter fused silica, packed with C18 (Waters, Acquity BEH C18, 1.7  $\mu$ m, 25 cm) and separate at a flow rate of 0.25  $\mu$ L/min using the following gradient: 5% CH<sub>3</sub>CN (0.1% FA) in H<sub>2</sub>O (0.1% FA) from 0–15 min, 5%–35% CH<sub>3</sub>CN (0.1% FA) from 15–155 min, 35%–95% CH<sub>3</sub>CN (0.1% FA) from 155–160 min, 95% CH<sub>3</sub>CN (0.1% FA) from 160–169 min, 95%–5% CH<sub>3</sub>CN (0.1% FA) from 169–170 min, and 5% CH<sub>3</sub>CN (0.1% FA) from 170–200 min.
  - b. The voltage applied to the nano-LC electrospray ionization source was 1.9 kV.
  - c. Data was acquired using an MS3-based TMT method adapted from Wang, Y. et al.(Wang et al., 2019)
    - Begin the scan sequence with an MS1 master scan (Orbitrap analysis, resolution 120,000, 400–1700 m/z, RF lens 60%, automatic gain control [AGC] target 2E5, maximum injection time 50 ms, centroid mode) with dynamic exclusion enabled (repeat count 1, duration 15 s).
    - ii. Select the top ten precursors for MS2/MS3 analysis.
    - iii. MS2 analysis consisted of: quadrupole isolation (isolation window 0.7) of precursor ion followed by collision-induced dissociation (CID) in the ion trap (AGC 1.8E4, normalized collision energy 35%, maximum injection time 120 ms).
    - iv. Following the acquisition of each MS2 spectrum, synchronous precursor selection (SPS) enabled the selection of up to 10 MS2 fragment ions for MS3 analysis.
    - v. MS3 precursors were fragmented by HCD and analyzed using the Orbitrap (collision energy 55%, AGC 1.5E5, maximum injection time 120 ms, resolution was 50,000).
    - vi. For MS3 analysis, use charge state–dependent isolation windows. For charge state z = 2, the MS isolation window was set at 1.2; for z = 3-6, the MS isolation window was set at 0.7.
- 27. Extract the MS2 and MS3 files from the raw files using RAW Converter (version 1.1.0.22; available at <a href="http://fields.scripps.edu/rawconv/">http://fields.scripps.edu/rawconv/</a>), upload to Integrated Proteomics Pipeline (IP2), and search using the ProLuCID algorithm (Figure 2B, publicly available at <a href="http://fields.scripps.edu/downloads.php">http://fields.scripps.edu/downloads.php</a>) using a reverse concatenated, non-redundant variant of the

#### Protocol



Human UniProt database (release-2012\_11; SwissProt and Trembl entries were combined and filtered to choose a single protein isoform for each Ensembl gene identifier such that there's one protein per gene. Reviewed entries (SwissProt) were chosen if available, otherwise the longest Unreviewed (Trembl) entry was used).

Note: For future applications, newer versions of the database should be generated and used.

- a. Search cysteine residues with a static modification for carboxyamidomethylation (+57.02146
  Da) and up to one differential modification for the desthiobiotin (DTB) tag (+398.2529 Da) (Figure 2B).
- b. Search the N-termini and lysine residues with a static modification corresponding to the TMT tag (+229.1629 Da).
- c. Trypsin was selected as the enzyme of choice. Peptides were required to be at least 6 amino acids long, to have at least one tryptic terminus, and to contain the DTB modification.
- d. Filter ProLuCID data through DTASelect (version 2.0) to achieve a peptide false-positive rate below 1%
- e. Perform the MS3-based peptide quantification with reporter ion mass tolerance set to 20 ppm with Integrated Proteomics Pipeline (IP2).

#### **EXPECTED OUTCOMES**

It is expected that this protocol will, on average, lead to enrichment and identification of >10000 cysteine-containing peptides (from >5000 proteins) per each 10-plex TMT-based experiment. This level of proteomic coverage allows for broad characterization of cysteine reactivity and/or small molecule-interactions across different proteomes. Using lower amounts of proteome is possible, but might lead to decreases in proteomic coverage. When combined with an unenriched quantitative proteomics workflow, this protocol allows for the proteome-wide discovery of biochemical changes in proteins that reflect dynamic alterations in cysteine reactivity.

The use of a 10-plex multiplexing strategy, in principle, allows for profiling up to 9 different electrophilic compounds in a single TMT-ABPP experiment (along with a DMSO control sample); however, in practice, we recommend performing each compound and DMSO control treatment in duplicate to facilitate interpretation of data, resulting in 4 compounds assayed per TMT-ABPP experiment. This approach informs on the extent of cysteine engagement by electrophilic compounds and proteomewide selectivity of these compounds, as well as structure-activity relationships for each quantified cysteine in a single experiment. This information facilitates mechanism-of-action studies related to phenotypic screening, where identifying structurally similar inactive compounds and their target landscapes is a powerful strategy for streamlined target identification. In particular, the deployment of compounds that can be stereochemically differentiated in their phenotypic effects and proteomic reactivity provides excellent active and inactive controls for mapping biologically relevant targets, identifying actionable sites of ligandability on diverse proteins in cells, and for performing follow-up functional studies on these proteins. A representative cysteine reactivity landscape for four stereoisomeric probes EV-96, EV-97, EV-98, and EV-99 is shown in Figure 4.

#### **QUANTIFICATION AND STATISTICAL ANALYSIS**

#### TMT-ABPP R value calculation for cysteine state-dependent reactivity dataset

This section of the protocol describes the filters applied at the stage of data compilation and processing.

- 1. At the individual TMT experiment level, the following filters should be applied to remove lowquality peptides:
  - a. Remove non-unique peptides;
  - b. Remove half-tryptic peptides;
  - c. Remove peptides with more than one internal missed cleavage site;




Figure 4. Target landscape of stereoisomeric acrylamide electrophiles

Heat map showing cysteines engaged >50% by stereoisomeric acrylamide compounds (5  $\mu$ M, 3 h). For inclusion in the map, cysteines were also required to show equivalent or greater engagement by the relevant stereoisomeric electrophile at 20  $\mu$ M. Stereoselective target engagement for EV-96 and EV-97 is highlighted in red and green, respectively. (Adapted from Vinogradova et al., (Vinogradova et al., 2020))

- d. Remove peptides with low (<10,000) sum of reporter ion intensities in both expanded or activated channels;
- e. Remove peptides with high variation (coefficient of variance >0.5) between the replicate expanded or activated channels if their sum of reporter ion intensities is greater than 5,000.
- f. Calculate R values (activated vs. expanded) for each peptide using the average reporter ion intensities of activated and expanded TMT channels.
- g. Once the R values are calculated, two types of grouping should be performed to aggregate peptide quantification data:
  - i. Group together overlapping peptides with the same modified cysteine (e.g., different charge states, high pH fractionation fractions, or tryptic termini), then average their R values, and report the shortest unique tryptic peptide;
  - ii. Group together multiple modified cysteines on a tryptic peptide, then report the averaged R values for further data processing.
- 2. Report the median value derived from at least two biological replicates as the final R value for each peptide with a maximum ratio cap of 20.

#### Protocol



#### Data processing and analysis for IA-DTB reactivity dataset

3. Proteins must have at least three unique quantified peptides in either particulate or soluble fraction in the TMT-ABPP experiments within the state-dependent dataset to be analyzed.

**Note:** Proteins with lower amount of quantified peptides in the reactivity profiling experiment can also be analyzed for reactivity changes, if 1) no subcellular fractionation into soluble and particulate fractions is performed in reactivity profiling and 2) a fully matched (donor, treatment conditions) unenriched proteomics data is available.

*Optional:* If the sample was fractionated into soluble and particulate fractions, the fraction with the most quantified unique peptides was selected for analysis for each protein. If a protein had an equal number of unique quantified peptides in both fractions, the peptide R ratios (activated vs. expanded) from both fractions were averaged.

- 4. To account for potential donor variations in protein expression level, proteins are required to have at least one peptide R ratio within 1.5-fold of the protein expression level measured in unenriched proteomics (TMT-exp) experiments (if available) and excluded from the analysis if all peptide R ratios are greater than 2.0 or less than 0.5.
- 5. For proteins with 5 or more quantified peptides, a cysteine can be considered for potential change in reactivity if its peptide R value differs more than two-fold from both the median R value of all quantified cysteines on the same protein and from the protein expression level measured in TMT-exp experiments (if available). For proteins with three or four quantified peptides, a cysteine was considered for potential change in reactivity if its peptide R value differed more than two-fold from the protein expression level measured by TMT-exp data, with an additional requirement that the maximum peptide R ratio differed more than 2-fold from the minimum peptide R ratio. All the cysteines that passed the initial filters described above were manually curated to remove low quality profiles.

#### **LIMITATIONS**

Like many untargeted mass spectrometry-based proteomic platforms, the cysteine reactivity profiling method described herein requires a substantial amount of protein (~0.5–1.0 mg of protein per cell state) for optimal sensitivity. This quantity of protein may be challenging to acquire for rare cell states, leading to lower sensitivity and cysteine coverage. Some dynamic cysteine reactivity changes may be lost upon cell lysis and adapted protocols where the IA-DTB probe is added immediately upon cell lysis might address this limitation. Low stoichiometry liganding events that produce < 50% engagement of a cysteine (i.e., < 50% reduction in IA-DTB reactivity for that cysteine) may be challenging to quantify with accuracy. Averaging the cysteine reactivity changes across several replicate experiments can increase confidence in lower stoichiometry liganding events, and this confidence can be further strengthened if such events show additional features indicative of specific binding (e.g., stereoselectivity).

#### **TROUBLESHOOTING**

#### **Problem 1**

Low intensity of the mass-spectrometry signal. (steps 25 and 26)

#### Potential solutions

Run a BSA digest standard to test the performance of the mass spectrometer; Increase the amount of proteome used;

Troubleshoot the enrichment step by performing quality control of all essential reagents (IA-DTB, streptavidin beads);





Sonicate the Eppendorf tube with the sample and solvent A before transferring the sample into the mass-spectrometry vial;

#### **Problem 2**

High levels of unlabeled peptides. (steps 12, 25, and 26)

#### **Potential solutions**

Add additional wash steps after the enrichment step; Perform quality control of essential reagents (IA-DTB, TMT tags);

#### Problem 3

High levels of channel variability between TMT channels. (steps 15-19)

#### **Potential solutions**

Make sure to fully close the Eppendorf tubes and carefully inspect them at each step to ensure absence of leaks;

Perform quality control of TMT tags;

#### **Problem 4**

High pressure and/or clogging of the mass-spectrometry column. (steps 25 and 26)

#### **Potential solutions**

Make sure you are not overloading the column;

When possible, use a trap column to avoid clogging of the main column;

If you see visible gooey residue after high pH fractionation step and removal of solvent using SpeedVac – pass each sample through a C18 tip to remove any salts that could clog the instrument. Re-make your buffers used for high pH fractionation of future samples;

If you see visible precipitate at the final step of the protocol, spin down the sample (10,000 g; 5 min) and avoid getting the precipitate into the sample loading vial;

#### Problem 5

Low peptide numbers from the IP2 search results. (step 27)

#### Potential solutions

Check if you are using the right FASTA database for the search (e.g., the correct organism); Double-check the search parameters, including the expected masses of static and differential modifications.

#### **RESOURCE AVAILABILITY**

#### Lead contact

Further information and requests for resources and reagents should be directed to and will be fulfilled by the lead contact, Benjamin F. Cravatt (cravatt@scripps.edu).

#### **Materials availability**

All chemical probes generated in this study are available from the lead contact with a completed Materials Transfer Agreement.

#### Protocol



#### Data and code availability

Raw proteomic data is available from the lead contact upon request. Processed proteomic data is provided in Data S1 of Vinogradova et al., (Vinogradova et al., 2020)

#### **ACKNOWLEDGMENTS**

This work was supported by NIH (CA231991 to B.F.C.), the Life Sciences Research Foundation (to E.V.V.), and Vividion Therapeutics. We thank the following individuals for contributing to the original manuscript, where this method was reported: Zhang, X.; Remillard, D.; Lazar, D.C.; Suciu, R.M.; Wang, Y.; Bianco, G.; Yamashita, Y.; Crowley, V.M.; Schafroth, M.A.; Yokoyama, M.; Konrad, D.B.; Lum, K.M.; Simon, G.M.; Kemper, E.K.; Lazear, M.R.; Yin, S.; Blewett, M.M.; Dix, M.M.; Nguyen, N.; Shokhirev, M.N.; Chin, E.; Lairson, L.; Melillo, B.; Schreiber, S.L.; Forli, S.; and Teijaro, J.R.

#### **AUTHOR CONTRIBUTIONS**

E.V.V. and B.F.C. conceived the project and drafted the manuscript.

#### **DECLARATION OF INTERESTS**

B.F.C. is a founder and scientific advisor to Vividion Therapeutics. B.F.C. and E.V.V. are co-inventors on a patent application related to this work.

#### **REFERENCES**

Backus, K.M., Correia, B.E., Lum, K.M., Forli, S., Horning, B.D., Gonzalez-Paez, G.E., Chatterjee, S., Lanning, B.R., Teijaro, J.R., Olson, A.J., et al. (2016). Proteome-wide covalent ligand discovery in native biological systems. Nature *534*, 570–574.

Vinogradova, E.V., Zhang, X., Remillard, D., Lazar, D.C., Suciu, R.M., Wang, Y., Bianco, G., Yamashita, Y., Crowley, V.M., Schafroth, M.A., et al. (2020). An activity-guided map of electrophile-cysteine interactions in primary human T cells. Cell *182*, 1009–1026.e9.

Wang, Y., Dix, M.M., Bianco, G., Remsberg, J.R., Lee, H.Y., Kalocsay, M., Gygi, S.P., Forli, S., Vite, G., Lawrence, R.M., et al. (2019). Expedited mapping of the ligandable proteome using fully functionalized enantiomeric probe pairs. Nat. Chem. 11, 1113–1123.