RESEARCH ARTICLE

WILEY

Analysis of current status of quantitative detection of biomarkers for liver fibrosis in Clinical labs in China

Chao Zhang 💿 🕴 Chuanbao Zhang

National Center for Clinical Laboratories, National Center of Gerontology, Institute of Geriatric Medicine, Beijing Hospital, Chinese Academy of Medical Sciences, Beijing, P. R. China

Correspondence

Chao Zhang and Chuanbao Zhang, National Center for Clinical Laboratories, Beijing Hospital, Beijing, China. Emails: zc_mdy@163.com (C. Z.), cbzhang@nccl.org.cn (C. Z.)

Funding information

This work was supported by the National Natural Science Foundation of China (no. 82003809) and Beijing Hospital Nova Project (BJ-2020-087).

Abstract

Aim: To explore the quality control and implementation of the quantitative detection of liver fibrosis biomarkers, laminin (LN), collagen IV (Col IV), procollagen III aminoterminal propeptide (PIIINP), hyaluronic acid (HA), and cholyglycine (CG), in China.

Methods: Two quality control products were measured in different laboratories using different measurement methods and reagents, and the acquired results were subjected to analysis. The quantitative detection technique was based on the conventional assessment criteria, with a target value $\pm 30\%$ being employed.

Results: Hundred labs were involved in the External Quality Assessment with 88 laboratories completing the assessment, and the pass rates were 84%, 80.2%, 67.5%, 77.3%, and 58.3% for HA, LN, PIIINP, Col IV, and CG, respectively. Chemiluminescence immunoassay was used most for HA (90.1%), LN (90.1%), PIIINP (87.9%), and Col IV (82.9%) determination, whereas the chemiluminescence immunoassay (31.6%), latex-enhanced immunoturbidimetry (36.7%), and homogeneous enzyme immunoassay (26.7%) were used for CG determination. The coefficients of variation for HA, LN, PIIINP, Col IV, and CG in different laboratories were 3.3%–19.49%, 1.74%–38.81%, 1.97%–41.29%, 2.85%–41.69%, and 2.71%–41.8%, respectively.

Conclusion: The clinical quantitative detection of liver fibrosis biomarkers is highly performed in China. The existing problems are that there are many manufacturers producing reagents and instruments, the quality of reagents is uneven, the specificity and sensitivity of reagents are greatly different, the comparability of results of various systems is poor, and the accuracy and consistency between different systems are lacking. All above underscores the critical importance of EQA in improving and monitoring the identification of biomarkers for liver fibrosis.

KEYWORDS

external quality assessment, internal quality control, liver fibrosis biomarkers, quality control products, quantitative detection and analysis

This is an open access article under the terms of the Creative Commons Attribution-NonCommercial-NoDerivs License, which permits use and distribution in any medium, provided the original work is properly cited, the use is non-commercial and no modifications or adaptations are made. © 2022 The Authors. *Journal of Clinical Laboratory Analysis* published by Wiley Periodicals LLC.

1 | INTRODUCTION

The characteristics of chronic liver disease include increasing fibrosis, which may result in portal hypertension, cirrhosis, liver synthesis damage, and hepatocellular carcinoma.^{1,2} Cholestatic and hepatotoxic chronic liver injuries are the two most prevalent forms of chronic liver injury.³ Hepatotoxic injury is induced by persistent liver cell damage, such as HBV and/or HCV infection, and non-alcoholic or alcoholic steatohepatitis.^{4,5} More than 240 million individuals worldwide are infected with chronic HBV, with China being one of the countries with a high incidence of HBV.^{6,7} Therefore, a large number of people with chronic HBV infection develop liver fibrosis in China.

Liver fibrosis is defined as the production of fibrous scarring caused by the accumulation of extracellular matrix (ECM) proteins (mainly type I and III crosslinked collagen) that substitute injured normal tissue.^{1,3} To be frank, the stage of liver fibrosis could only be determined via the biopsy of the liver at present. Nevertheless, a biopsy of the liver can cause serious problems to patients in up to 3%, as well as mortality occurring in 0.03% of the cases.⁵ For monitoring chronic liver disease progression, assessing fibrosis staging is critical, which determines prognosis and optimal treatment time, monitors treatment response, and assesses disease progression to reduce morbidity and mortality associated with cirrhosis sequelae. Thus, accurate non-invasive tests for liver fibrosis diagnosis and staging are required. Non-invasive fibrosis testing of fibrosis is gaining traction in clinical practice to diagnose liver fibrosis. For fibrosis assessment criteria, blood-based indicators provide considerable benefits than standard biopsy of the liver, such as accessibility, cost, and safety. Serum markers estimate fibrosis using markers that measure ECM degradation/fibrogenesis. An extensive spectrum of chronic liver disorders, such as non-alcoholic and alcoholic fatty liver disease, as well as chronic viral hepatitis, has been studied and verified with the aid of this technique.

Currently, five major liver fibrosis indicators are utilized in clinical laboratories, namely laminin (LN), hyaluronic acid (HA), procollagen III's N-terminal peptide (PIIINP), cholyglycine (CG), and collagen IV (Col IV).⁸⁻¹³ However, none of these biomarkers was included in the formal External Quality Assessment (EQA) plan. This means it is impossible to assess whether the measurement is accurate. The reliability and utility of disease diagnosis and treatment are determined by the accuracy of the clinical outcomes. The comparability and accuracy of findings across various methods employed in different laboratories, as a consequence, are critically essential considerations. EQA is a method in which the same samples are examined by various labs, and the findings are obtained, analyzed, and subjected to a comparison made by an external, independent organization. The detection and calibration capabilities of clinical laboratories are determined by the National Center for Clinical Laboratories (NCCL) via interlaboratory comparisons, and the laboratory progress is monitored. As a consequence, EQA contributes to the guarantee of comparability and accuracy of outcomes. Now that clinical laboratories are providing detection

services for liver fibrosis biomarkers, it is necessary for the NCCL to conduct EQA services for them to ensure the accuracy and comparability of test results.

2 | MATERIAL AND METHODS

2.1 | Respondents

In the EQA survey for the quantitative identification of liver fibrosis biomarkers, the responders comprised clinical labs that had consented to be involved in it.

2.2 | Questionnaire survey

In response to the query on "Investigations on the Quantitative Detection of Biomarkers of Liver Fibrosis," the clinical laboratories provided responses using the EQA system established by NCCL. In addition to the essential data, the questionnaire comprised of eight queries, including the quantitative identification of liver fibrosis biomarkers, the type of the sample, method of detection, the brand of the reagent, the brand of the instrument, the unit of the concentration used by the laboratory, and their readiness to be involved in the EQA survey (serum matrix) for the quantitative identification of liver fibrosis biomarkers (Table S1). After the data were obtained, a summary of the current state of the quantitative detection of liver fibrosis indicators in clinical labs was prepared.

2.3 | EQA survey

The NCCL handled the notification of the EQA investigation and the submission of the application. The quality control items were provided to the labs that had submitted their applications. 202011 and 202012 were the two samples (comprising HA, LN, PIIINP, Col IV, and CG biomarkers) with different concentrations (cut-off values and high values).

2.4 | Criteria of EQA evaluation

For the purpose of conducting the present research, the statistically robust mean \pm 30% was used as the evaluation criterion. The detection result was considered as "qualified" if the deviation of the two samples was within the set range.

2.5 | Statistical analysis

Based on the response and submitted findings, the accumulated number and the proportion of the questionnaire responses were calculated. The quantitative identification outcomes of the liver fibrosis biomarkers were classified in accordance with the reagents and evaluation made. Moreover, every reagent manufacturing company was named as reagents 1, 2, 3, and so on (Table S2). The coefficient of variation (CV), mean value, and bias between the target and the mean values of the reagent groups were derived when no less than two laboratories were in a single group.

3 | RESULTS

3.1 | The proportion of detection projects of liver fibrosis biomarkers

In the questionnaire, the proportion of the various detection items was evaluated. As shown in Table 1, 49.6% (112/226) of the laboratories quantitatively detected four liver fibrosis biomarkers (HA, LN, PIIINP, and Col IV); 39.8% (90/226) of the laboratories detected five liver fibrosis biomarkers (HA, LN, PIIINP, Col IV, and CG): 10.6% (24/226) of the laboratories detected only the CG biomarker. We assessed the proportion of distinct detection techniques with the aid of a questionnaire. The results revealed that the proportion of chemiluminescence immunoassay (CLIA) for HA, LN, PIIINP, and Col IV detection was 86.6%. Other detection methods included latex-enhanced immunoturbidimetry, up-conversion luminescence immunoassay, (ULIA), electrochemiluminescence, (ECLIA), and radioimmunoassay (RIA), accounting for approximately 13.4% of the total detection methods. The common methods for CG detection were CLIA, latex-enhanced immunoturbidimetry, and homogeneous enzyme immunoassay (HEI), accounting for 28.1%, 37.7%, and 28.9%, respectively. The other methods including ECLIA and RIA represented around 5.3% of the overall estimate.

TABLE 1 Summary of questionnaire information

3.2 | Group statistics on the basis of the methods of detection

For the quality control materials that were distributed, two concentrations (cut-off value and high value) were established. For the purpose of collecting more data, the quality control materials were distributed to the labs that quantitatively detected four or five biomarkers. Moreover, the clinical labs that were involved in the present research employed a variety of methods, procedures, and reagents. "Other category" was used to classify laboratories that did not describe their detection techniques in their application. As illustrated in Table 2, robust CV, standard uncertainty and the robust standard deviation are all in accordance with the ISO13528 standards. The results revealed significant differences in the robust mean values among the various identification methods of the laboratories, indicating the necessity for grouped statistics. For HA detection using RIA, the quality control substance concentration increased with an increase in robust CV (24.1% vs 70.8%), which can be attributed to identification procedures defects and small trial size. With increased concentration, robust CV of chemiluminescence immunoassay-HRP for HA, LN, PIIINP, and Col IV showed a downward trend, while the chemiluminescence immunoassay-Acridinium Ester showed an opposite trend. We tried to explain that maybe CLIA-HRP is prone to be stable in HA, LN, Col IV, and PIIINP detection when the sample concentration is increased, while CLIA-AE is more suitable for detection at lower concentrations. Conversely, for Col IV detection, the robust CV of CLIA-AE and microplate decreased with an increase in concentration. However, the detection of PIIINP, HA, LN, and Col IV via ULIA was relatively stable. For CG detection, the robust CV of CLIA decreased with increased sample concentration. The robust CV of

Project		Amount		Proportion
HA, LN, PIIINP, Col IV		112		49.6%
HA, LN, PIIINP, Col IV, and CG		90		39.8%
CG		24		10.6%
	Method		Amount	Proportion (%)
HA, LN, PIIINP, Col IV	Chemiluminescence immunoassay,	CLIA	175	86.6
	Latex-enhanced immunoturbidimet	ry	17	8.4
	Up-conversion luminescence immu	noassay, ULIA	5	2.5
	Electrochemiluminescence, ECLIA		3	1.5
	Radioimmunoassay, RIA		2	1.0
CG	Chemiluminescence immunoassay,	CLIA	32	28.1
	Electrochemiluminescence, ECLIA		3	2.6
	Radioimmunoassay, RIA		1	0.9
	Latex-enhanced immunoturbidimet	ry	43	37.7
	Homogeneous enzyme immunoass	ay, HEI	33	28.9
	Others		2	1.8

Note: The proportion of detection projects of liver fibrosis biomarkers and composition ratios of detection methods.

4 of 12 | WILEY-

TABLE 2 The result of grouped statistics according to detection methods

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV (%)
HA	202011	All	81	119.54	23.3	3.236	19.49
		Chemiluminescence immunoassay-AMPPD	8	97.2	26.66	12.594	27.42
		Chemiluminescence immunoassay-HRP	22	128.39	9.97	2.657	7.77
		Chemiluminescence immunoassay-Acridinium Ester	13	118.36	21.7	7.524	18.34
		Chemiluminescence immunoassay-ABEI	25	121.61	14.58	3.645	11.99
		Chemiluminescence immunoassay-Microplate	5	54.55	4.1	2.958	7.51
		Up-conversion luminescence immunoassay, ULIA	3	102.74	5.24	3.782	5.1
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	2	137.1	33.04	29.201	24.1
		Latex-enhanced immunoturbidimetry	1	-	-	-	-
		Others	1	-	-	-	-
	202012	All	81	303.35	66.8	9.278	22.02
		Chemiluminescence immunoassay-AMPPD	8	339.91	44.37	20.965	13.06
		Chemiluminescence immunoassay-HRP	22	343.84	20.89	5.568	6.08
		Chemiluminescence immunoassay-Acridinium Ester	13	360.37	113.51	39.353	31.5
		Chemiluminescence immunoassay-ABEI	25	288.06	23.93	5.982	8.31
		Chemiluminescence immunoassay-Microplate	5	87.69	13.33	9.618	15.2
		Up-conversion luminescence immunoassay, ULIA	3	126.98	1.44	1.04	1.14
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	2	341.9	239.6	211.774	70.08
		Latex-enhanced immunoturbidimetry	1	-	-	-	-
		Others	1	-	-	-	-
LN	202011	All	81	146.94	57.03	7.921	38.81
		Chemiluminescence immunoassay-AMPPD	8	120.43	40.49	19.127	33.62
		Chemiluminescence immunoassay-HRP	22	124.14	8.98	2.394	7.24
		Chemiluminescence immunoassay-Acridinium Ester	13	115.95	57.95	20.091	49.98
		Chemiluminescence immunoassay-ABEI	25	189.48	7.04	1.761	3.72
		Chemiluminescence immunoassay-Microplate	5	131.99	17.45	12.59	13.22
		Up-conversion luminescence immunoassay, ULIA	3	81.91	1.98	1.432	2.42
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	2	326	49.72	43.942	15.25
		Latex-enhanced immunoturbidimetry	1	-	-	-	-
		Others	1	-	-	-	-
	202012	All	81	457.63	178.68	24.816	39.04
		Chemiluminescence immunoassay-AMPPD	8	404.71	54.62	25.804	13.5
		Chemiluminescence immunoassay-HRP	22	395.11	23.49	6.261	5.95
		Chemiluminescence immunoassay-Acridinium Ester	13	363.45	266.71	92.464	73.38
		Chemiluminescence immunoassay-ABEI	25	607.1	34.2	8.55	5.63
		Chemiluminescence immunoassay-Microplate	5	334.05	19.9	14.359	5.96
		Up-conversion luminescence immunoassay, ULIA	3	73.31	2.81	2.026	3.83
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	2	1213.85	39.53	34.942	3.26
		Latex-enhanced immunoturbidimetry	1	-	-	-	-
		Others	1	-	-	-	-

ZHANG AND ZHANG

TABLE 2 (Continued)

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV (%)
PIIINP	202011	All	83	12.36	4.98	0.683	40.28
		Chemiluminescence immunoassay-AMPPD	8	10.32	3.92	1.852	37.99
		Chemiluminescence immunoassay-HRP	22	10.86	1.49	0.397	13.73
		Chemiluminescence immunoassay-Acridinium Ester	13	4.3	1.29	0.448	30.06
		Chemiluminescence immunoassay-ABEI	25	14.57	1.08	0.269	7.38
		Chemiluminescence immunoassay-Microplate	5	156.64	31.6	22.802	20.17
		Up-conversion luminescence immunoassay, ULIA	3	61.96	4.87	3.512	7.86
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	1	-	-	-	-
		Latex-enhanced immunoturbidimetry	3	4.19	3.01	2.175	71.99
		Others	1	-	-	-	-
	202012	All	83	30.85	12.74	1.747	41.29
		Chemiluminescence immunoassay-AMPPD	8	23.96	9.27	4.378	38.67
		Chemiluminescence immunoassay-HRP	22	35.22	3.87	1.032	11
		Chemiluminescence immunoassay-Acridinium Ester	13	10.23	4.86	1.686	47.55
		Chemiluminescence immunoassay-ABEI	25	29.37	2.68	0.671	9.13
		Chemiluminescence immunoassay-Microplate	5	475.43	355.83	256.801	74.84
		Up-conversion luminescence immunoassay, ULIA	3	59.72	4.49	3.242	7.52
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	1	-	-	-	-
		Latex-enhanced immunoturbidimetry	3	10.23	0.69	0.5	6.77
		Others	1	-	-	-	-
Col IV	202011	All	88	78.34	32.66	4.352	41.69
		Chemiluminescence immunoassay-AMPPD	8	144.65	10.95	5.173	7.57
		Chemiluminescence immunoassay-HRP	22	96.01	5.4	1.438	5.62
		Chemiluminescence immunoassay-Acridinium Ester	13	74.29	26	9.014	25
		Chemiluminescence immunoassay-ABEI	25	48.77	3.51	0.895	7.19
		Chemiluminescence immunoassay-Microplate	5	98.79	23.82	17.188	24.11
		Up-conversion luminescence immunoassay, ULIA	3	129.65	6.18	4.458	4.76
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	2	57.15	7.78	6.875	13.61
		Latex-enhanced immunoturbidimetry	8	61.82	25.69	10.704	41.55
		Others	1	-	-	-	-
	202012	All	88	306.87	110.22	14.686	35.92
		Chemiluminescence immunoassay-AMPPD	8	454.64	189.92	89.731	41.78
		Chemiluminescence immunoassay-HRP	22	387.81	17.78	4.739	4.59
		Chemiluminescence immunoassay-Acridinium Ester	13	282.14	71.14	24.665	25.22
		Chemiluminescence immunoassay-ABEI	25	258.77	28.29	7.219	10.93
		Chemiluminescence immunoassay-Microplate	5	426.51	40.5	29.231	9.5
		Up-conversion luminescence immunoassay, ULIA	3	327.88	9.73	7.024	2.97
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	2	100.1	13.95	12.332	13.94
		Latex-enhanced immunoturbidimetry	8	135.41	73.86	30.776	54.55
		Others	1	-	-	-	_

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV (%)
CG	202011	All	60	2.47	0.65	0.106	26.46
		Chemiluminescence immunoassay-AMPPD	2	2.02	1.14	1.007	56.39
		Chemiluminescence immunoassay-HRP	2	2.15	0.15	0.134	7.05
		Chemiluminescence immunoassay-ABEI	15	2.7	0.2	0.067	7.46
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	1	-	-	-	-
		Immunoturbidimetry	7	2.28	0.5	0.238	22.1
		Latex-enhanced immunoturbidimetry	15	2.03	0.23	0.074	11.3
		Homogeneous enzyme immunoassay, HEI	16	2.46	0.53	0.167	21.72
		Others	1	-	-	-	-
	202012	All	60	14.12	5.9	0.952	41.8
		Chemiluminescence immunoassay-AMPPD	2	10.2	3.37	2.977	33.02
		Chemiluminescence immunoassay-HRP	2	17.16	1.24	1.099	7.25
		Chemiluminescence immunoassay-ABEI	15	12.89	0.31	0.103	2.38
		Electrochemiluminescence, ECLIA	1	-	-	-	-
		Radioimmunoassay, RIA	1	-	-	-	-
		Immunoturbidimetry	7	14.84	8.22	3.885	55.42
		Latex-enhanced immunoturbidimetry	15	9.69	2.3	0.743	23.76
		Homogeneous enzyme immunoassay, HEI	16	17.76	7.32	2.287	41.22
		Others	1	-	_	_	_

If there is only one laboratory, robust standard deviation, standard uncertainty, and robust CV cannot be calculated.

immunoturbidimetry increased significantly. We suspect that CLIA may have high reliability in detecting high concentrations.

3.3 | Group statistics premised on detection reagents

Besides the detection techniques, the responses were classified according to the detection reagents used in view of the fact that the same detection technique can use different reagents; therefore, the classification by reagents increases accuracy. As shown in Table 3, reagents 1 and 13 represented the majority of the market share, indicating that they are frequently employed in most Chinese clinical labs for the quantitative identification of liver fibrosis biomarkers. Furthermore, the pass rates of the reagents 1 and 13 groups (100%) were higher than those of the other groups, suggesting that the larger sample size could increase the statistical accuracy. Although the reagent 1, 7, 10, 13, and 14 groups used CLIA for the detections of Col IV, LN, HA, and PIIINP, the results varied, thereby suggesting the need for grouping based on the use of reagents and improving the commutability of quality control materials. For CG detection, reagents 1, 13, and 26 were commonly used in most Chinese clinical laboratories. Reagents 1 and 13 were used in the CLIA method, and reagent 26 was used in the HEI method. With increased sample concentration, CLIA's robust CV decreased, whereas that of HEI was elevated, which

was consistent with the findings of grouping by detection methods (Table 2).

3.4 | Group statistics premised on the principle of grouping

According to the statistical data shown in Tables 2 and 3, it is reasonable to hypothesize that the detection reagents and techniques were employed in their entirety. The reagent produced by a particular manufacturer is typically matched to the detection technique used by that manufacturer. Regardless of the data for grouping based on equipment, procedures, or reagents, grouping was conducted according to the number of labs involved: greater than or equal to 18 or 12 were in one group under ISO 13528. The "other" category was formed if they did not fit into any of the other groupings. According to ISO 13528, the findings were examined using the robust mean, robust standard deviation, and robust CV of the involved labs. However, owing to the limited number of participating labs, the grouping criterion was unsuitable. Therefore, the participating laboratories were classified into separate groups based on the number of laboratories using the same reagent \geq 5 (Table 4). For the detection of PIIINP, the lowest values were obtained using reagent 7 and the highest using reagent 14. However, for the detection of LN, the results obtained on using reagent 7 and reagent 14 were opposite. Furthermore, reagents 7 and 14 were both used in the CLIA

6 of 12

TABLE 2 (Continued)

 TABLE 3
 The result of grouped statistics according to detection reagents

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV (%)
HA	202011	All	81	119.54	23.3	3.236	19.49
		Reagent 1	22	131.94	7.81	2.131	5.92
		Reagent 2	2	137.1	33.04	29.201	24.1
		Reagent 5	3	102.74	5.24	3.782	5.1
		Reagent 7	5	42.17	1.39	0.778	3.3
		Reagent 8	3	63.91	0.8	0.709	1.25
		Reagent 9	3	106.87	9.92	8.767	9.28
		Reagent 10	5	101.61	5.36	2.995	5.27
		Reagent 11	3	99.39	4.9	3.533	4.93
		Reagent 13	25	122.45	13.49	3.372	11.01
		Reagent 14	5	265.2	11.58	6.473	4.37
	202012	All	81	303.35	66.8	9.278	22.02
		Reagent 1	22	349.89	16.54	4.511	4.73
		Reagent 2	2	341.9	239.6	211.774	70.08
		Reagent 5	3	126.98	1.44	1.04	1.14
		Reagent 7	5	63.93	9.96	5.567	15.58
		Reagent 8	3	140.03	1.6	1.418	1.15
		Reagent 9	3	230.33	33.53	29.64	14.56
		Reagent 10	5	311.77	16.08	8.988	5.16
		Reagent 11	3	375.08	6.61	4.773	1.76
		Reagent 13	25	289.44	22.98	5.744	7.94
		Reagent 14	5	705.18	42.11	23.541	5.97
LN	202011	All	81	146.94	57.03	7.921	38.81
		Reagent 1	22	123.44	7.72	2.107	6.26
		Reagent 2	2	326	49.72	43.942	15.25
		Reagent 5	3	81.91	1.98	1.432	2.42
		Reagent 7	5	127.6	2.23	1.244	1.74
		Reagent 8	3	251.22	130.25	115.122	51.85
		Reagent 9	3	91.58	3.39	2.998	3.7
		Reagent 10	5	63.27	3.26	1.822	5.15
		Reagent 11	3	87.38	1.76	1.269	2.01
		Reagent 13	25	189.18	7.18	1.794	3.79
		Reagent 14	5	235.67	15.85	8.86	6.73
	202012	All	81	457.63	178.68	24.816	39.04
		Reagent 1	22	394.83	18.76	5.118	4.75
		Reagent 2	2	1213.85	39.53	34.942	3.26
		Reagent 5	3	73.31	2.81	2.026	3.83
		Reagent 7	5	360.99	35.81	20.016	9.92
		Reagent 8	3	707.88	372.38	329.136	52.61
		Reagent 9	3	266.35	3.79	3.345	1.42
		Reagent 10	5	117.53	6.56	3.668	5.58
		Reagent 11	3	392.16	17.59	12.697	4.49
		Reagent 13	25	608.39	32.44	8.111	5.33
		Reagent 14	5	624.26	25.49	14.252	4.08

8 of 12 | WILEY

TABLE 3 (Continued)

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV (%)
PIIINP	202011	All	83	12.36	4.98	0.683	40.28
		Reagent 1	22	10.54	1.18	0.321	11.17
		Reagent 2	2	233.95	26.38	23.318	11.28
		Reagent 5	3	61.96	4.87	3.512	7.86
		Reagent 7	5	182.71	31.28	17.488	17.12
		Reagent 8	3	16.42	9.85	8.704	59.97
		Reagent 9	3	10.07	0.1	0.085	0.95
		Reagent 10	5	3.21	0.16	0.089	4.96
		Reagent 11	3	7.36	0.4	0.287	5.39
		Reagent 13	25	14.78	0.99	0.247	6.69
		Reagent 14	5	4.2	0.24	0.132	5.65
		Reagent 16	2	2.69	0.92	0.66	34.02
	202012	All	83	30.85	12.74	1.747	41.29
		Reagent 1	22	34.76	3.13	0.854	9
		Reagent 2	2	1858.1	1460.67	1291.059	78.61
		Reagent 5	3	59.72	4.49	3.242	7.52
		Reagent 7	5	618.64	108.46	60.633	17.53
		Reagent 8	3	29.43	14.28	12.623	48.53
		Reagent 9	3	17.24	0.64	0.567	3.72
		Reagent 10	5	6.21	0.25	0.139	3.99
		Reagent 11	3	16.8	0.94	0.675	5.57
		Reagent 13	25	29.38	2.14	0.536	7.29
		Reagent 14	5	9.38	0.19	0.103	1.97
		Reagent 16	2	9.56	0.69	0.5	7.25
Col IV	202011	All	88	78.34	32.66	4.352	41.69
		Reagent 1	22	94.82	4.01	1.094	4.23
		Reagent 2	2	57.15	7.78	6.875	13.61
		Reagent 5	3	129.65	6.18	4.458	4.76
		Reagent 7	5	100	5.45	3.048	5.45
		Reagent 8	3	161.97	9.82	8.682	6.07
		Reagent 9	3	70.49	23.37	20.66	33.16
		Reagent 10	5	107.02	12.58	7.035	11.76
		Reagent 11	3	151.01	2.14	1.546	1.42
		Reagent 13	25	47.65	4.12	1.051	8.65
		Reagent 14	5	56.4	4.83	2.702	8.57
		Reagent 17	7	64	17.15	8.75	26.79
	202012	All	88	306.87	110.22	14.686	35.92
		Reagent 1	22	385.22	17.66	4.818	4.59
		Reagent 2	2	100.1	13.95	12.332	13.94
		Reagent 5	3	327.88	9.73	7.024	2.97
		Reagent 7	5	401.93	52.77	29.498	13.13
		Reagent 8	3	404.9	38.82	34.311	9.59
		Reagent 9	3	200.4	44.48	39.314	22.2
		Reagent 10	5	408.7	30.17	16.864	7.38
		Reagent 11	3	637.21	6.36	4.593	1
		Reagent 13	25	254.99	23.07	5.886	9.05
		Reagent 14	5	235.62	6.72	3.758	2.85
		Reagent 17	6	146.89	38.78	19.791	26.4

TABLE 3 (Continued)

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV (%)
CG	202011	All	60	2.47	0.65	0.106	26.46
		Reagent 1	5	3.44	0.71	0.395	20.58
		Reagent 7	4	1.3	0.02	0.021	1.85
		Reagent 13	15	2.74	0.18	0.056	6.4
		Reagent 17	3	2.2	1.06	0.935	48.09
		Reagent 19	2	1.88	0.29	0.255	15.37
		Reagent 24	4	1.94	0.07	0.046	3.77
		Reagent 25	4	4.06	0.1	0.085	2.37
		Reagent 26	10	2.46	0.19	0.074	7.61
		Reagent 29	3	2.1	0.13	0.095	6.24
	202012	All	60	14.12	5.9	0.952	41.8
		Reagent 1	5	24.02	4.17	2.333	17.37
		Reagent 7	4	7.82	0.45	0.397	5.74
		Reagent 13	15	12.93	0.35	0.113	2.71
		Reagent 17	3	29.4	25.68	22.695	87.33
		Reagent 19	2	11.96	7.1	6.272	59.36
		Reagent 24	4	8.21	1.14	0.714	13.91
		Reagent 25	4	27.78	2.92	2.58	10.51
		Reagent 26	10	17.28	3.04	1.2	17.58
		Reagent 29	3	8.71	0.54	0.392	6.24

There was only one case of reagent 3, 4, 6, 12, 15 for detection of HA, LN, PIIINP, and Col IV, so statistical calculation could not be carried out. There was only one case of reagent 2, 4, 12, 18, 20, 21, 22, 23, 27, 28 for detection of CG, so statistical calculation could not be carried out.

method but labeled with different substrates. This can be attributed to the different sensitivity of the various reagents for the detection items. When the EQA of the quantitative detection of liver fibrosis biomarkers is properly performed, the variability of the results may be minimized by expanding the number of labs that were included in the present research.

All included labs were divided into groups and subjected to analyses according to the reagents they employed, with the findings being reviewed in relation to the target value $\pm 30\%$. The desired outcome is a robust average value, with a dispersion of $\pm 30\%$ being the maximum permitted for the degree of dispersion of the data. A satisfactory result was judged to be achieved if the detection value fell within a certain range (Table 4). In the 202011 sample detection, the pass rates of HA, LN, PIIINP, Col IV, and CG in the "other group" were 70%, 30%, 31.8%, 36.45%, and 63.3%, whereas those of reagents 1, 7, 10, 13, and 14 were nearly 100%. Therefore, if laboratories using the same reagent can be individually grouped, the pass rates are relatively high. This situation can also be seen in the detection of high-concentration samples. The pass rate of the separate group was remarkably higher than that of the mixed group. The coefficient of variation for HA, LN, PIIINP, Col IV, and CG in different laboratories were 3.3%-19.49%, 1.74%-38.81%, 1.97%-41.29%, 2.85%-41.69%, and 2.71%-41.8%, respectively. Generally, the lower CV appeared in the separate group, while the higher CV appeared in the mixed group. Therefore, these results signify the importance of grouping. With the increase in concentration, the pass

rates of HA, LN, PIIINP, Col IV, and CG, apart from the other group, showed no significant difference. The pass rates of HA, PIIINP, and CG decreased with increasing concentrations, whereas that of LN and Col IV increased but not by much.

3.5 | Total pass rates

After calculating the pass rate of each detection project, the total pass rates of all projects were calculated. The total pass rates of HA, LN, PIIINP, Col IV, and CG were 84.0%, 80.2%, 67.5%, 77.3%, and 58.3%, respectively (ie pass for both low and high concentrations). Therefore, the pass rates of PIIINP, Col IV, and CG need improvement (Table 5).

4 | DISCUSSION

EQA (referred to as "proficiency testing") is employed for the purpose of assessing the capacity of a laboratory to perform tests. It is an essential tool for identifying problems in a clinical laboratory and designing the appropriate solutions. EQA is an essential external monitoring technique for quality assurance, especially with the lack of a reference method and reference material. Several novel protein detection markers have recently been introduced into the clinical testing environment. This contributes to China's fast expansion of the medical

10 of 12 | WILEY-

 TABLE 4
 The result and passing rates of grouped statistics according to the grouping principle

HA202011AlAl191470Reagent122131.947.8122.1315.922221100Reagent1512.171.3010.733.355100Reagent105101.415.5382.9755.2755100Reagent1450.62.211.886.4734.375100202012Al8303.256.6799.2782.20191050Reagent1563.939.5955.56715.585210050Reagent1563.979.2782.202191050100Reagent1563.979.2782.202191050Reagent15131.7716.078.6865.1655100Reagent15131.7716.078.6865.1655100Reagent15123.447.7322.1076.262.5100Reagent15123.447.7322.1076.262.5100Reagent15123.447.7311.7443.755100Reagent15123.45178.672.48161.7455100Reagent1523.64178.672.48161.7455100Reagent153.6093.5052.0144.7455 </th <th></th> <th>Batch number</th> <th>Group</th> <th>Total number</th> <th>Robust mean</th> <th>Robust SD</th> <th>Standard uncertainty</th> <th>Robust CV</th> <th>Total number</th> <th>Number of passing labs</th> <th>Pass rate (%)</th>		Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV	Total number	Number of passing labs	Pass rate (%)
Reagent 2 12 741 <th74< th=""> 741<td>HA</td><td>202011</td><td>All</td><td>81</td><td>119.54</td><td>23.3</td><td>3.236</td><td>19.49</td><td>19</td><td>14</td><td>70</td></th74<>	HA	202011	All	81	119.54	23.3	3.236	19.49	19	14	70
Reagent542.171.3910.7783.355100Reagent105101.615.3582.9955.2755100Reagent145252.211.886.4734.3755100202012All301303.566.7999.2782.202191050Reagent125434.886.1534.5114.732221100Reagent105311.7716.0798.9885.165480Reagent132525.275.56715.585480Reagent145705.1842.11123.5415.9755100Reagent1325705.1842.11123.5415.9755100Reagent145123.447.0327.92138.8119630Reagent1325123.457.0327.92138.8119630Reagent145123.457.1751.7747.79255100Reagent10563.273.261.4225.155100Reagent1325136.358.866.7355100Reagent1325136.973.5652.541001010Reagent1325160.9725251001010Reagent145117.536.5623.6685.88 </td <td></td> <td></td> <td>Reagent 1</td> <td>22</td> <td>131.94</td> <td>7.812</td> <td>2.131</td> <td>5.92</td> <td>22</td> <td>21</td> <td>100</td>			Reagent 1	22	131.94	7.812	2.131	5.92	22	21	100
Reagent10510.15.3582.9955.2755100Reagent132512.24513.4863.37211.01252496Reagent14526.211.586.4734.3755100202012All0103.036.6799.27822.02191050Reagent10531.0716.5764.5144.732221100Reagent10531.0716.0798.98851.655100Reagent10521.1716.0798.98851.655100Reagent105015.1342.11123.415.9755100Reagent105127.627.213.84119630Reagent105127.622.211006.262221100Reagent105127.6127.6127.41.741.7455100Reagent105127.6127.61.25.71.558.846.735100Reagent105192.617.551.558.846.735100Reagent105192.617.558.5655100Reagent10517.558.565.585100Reagent10517.538.562.5414.254.0855100Reagent10517.			Reagent 7	5	42.17	1.391	0.778	3.3	5	5	100
Reagent 13 25 12.45 13.48 3.372 11.01 25 24 96 Reagent 14 5 262.2 11.58 6.473 4.37 5 5 100 Reagent 1 22 349.88 16.536 4.511 4.73 22 21 100 Reagent 0 5 63.93 9.295 5.567 15.58 5 4 80 Reagent 10 5 11.77 16.079 8.988 5.567 5 5 100 Reagent 10 5 11.77 16.079 8.988 5.577 5 5 100 Reagent 10 5 11.74 2.52 12.4 1.749 25 25 100 Reagent 13 25 182.4 7.723 2.107 6.26 2.2 100 Reagent 10 5 127.6 2.226 1.244 1.745 1.74 1.74 5 5 100 Reagent 10 5			Reagent 10	5	101.61	5.358	2.995	5.27	5	5	100
Reagent 14 5 265.2 11.58 6.473 4.37 5 5 100 202012 All 81 303.3 66.79 9.278 22.02 19 10 50 Reagent 1 25 349.88 16.536 4.511 4.73 5 4 80 Reagent 10 5 31.77 16.07 8.988 5.16 5 5 100 Reagent 13 25 289.4 22.05 5.744 7.94 5 5 100 Reagent 14 5 705.18 42.111 23.541 7.97 5 5 100 Reagent 1 22 12.64 7.032 7.921 38.81 19 6 30 Reagent 1 5 12.76 2.226 1.822 51.5 5 100 Reagent 13 25 18.18 7.75 1.753 3.66 7.97 2 2 100 Reagent 13 25 18.			Reagent 13	25	122.45	13.486	3.372	11.01	25	24	96
202012 All 81 30.3.5 6.6.799 9.278 22.02 19 10 50 Reagent 1 22 34.98 16.536 4.511 4.73 22 21 100 Reagent 10 5 311.77 16.079 8.988 5.16 5 50 100 Reagent 13 25 28.04 22.971 5.744 7.74 2.5 2.5 1000 LN Reagent 1 21 12.344 7.723 2.107 6.26 2.2 2.1 100 Reagent 10 5 0.53.7 3.26 1.244 1.74 5 5 100 Reagent 10 5 0.52.7 1.75.7 7.75 5 5 100 Reagent 10 5 0.53.7 3.26 1.24 1.74 5 5 100 Reagent 10 5 0.52.7 1.75.8 8.66 6.73 5 100 Reagent 10 5 1.75.3			Reagent 14	5	265.2	11.58	6.473	4.37	5	5	100
Reagent1 22 349.88 16.33 4.511 4.73 22 21 100 Reagent7 5 63.03 9.95 5.567 15.88 5 4 80 Reagent13 25 289.44 22.975 5.744 7.94 25 25 100 Reagent13 25 289.44 22.975 5.744 7.94 25 25 100 Reagent14 5 705.18 42.111 23.541 5.7 5 100 Reagent14 51 127.6 2.226 1.244 1.74 5 5 100 Reagent13 25 199.18 7.175 1.774 3.7 2.5 100 Reagent13 22 199.18 7.175 1.744 5 5 100 Reagent13 22 199.18 7.175 2.666 5.58 5 5 100 Reagent13 22 199.13 5.562 3.666 3.668		202012	All	81	303.35	66.799	9.278	22.02	19	10	50
Reagent 7563.9399595.56715.885480Reagent 10531.17716.0798.9885.1655100Reagent 132529.9444.201123.5415.9755100Reagent 145705.1842.11123.5415.9755100INReagent 12212.3447.7232.1076.262221100Reagent 10563.273.261.8225.1555100Reagent 1325189.187.1751.7943.792525100Reagent 14523.56715.858.866.7355100Reagent 1325189.187.1751.7943.792525100Reagent 14536.09935.80520.0169.9255100Reagent 1325608.393.24423.6835.5855100Reagent 145642.6225.44414.2524.0855100Reagent 1325668.393.2423.6834.0272173.88PIINP20211All8312.364.7780.6834.0272173.88Reagent 1325668.390.2476.692525100100100100Reagent 13251.6241.7471.74			Reagent 1	22	349.88	16.536	4.511	4.73	22	21	100
Reagent 105311.7716.0798.9885.1655100Reagent 1325297.422.9755.747.942525100LNAll81146.9457.0327.92138.8119630Reagent 7562221100100100100100Reagent 7562221100100100100100Reagent 85632.73.261.82251555100Reagent 1325189.187.1751.7943.792525100Reagent 145235.6715.858.866.7355100Reagent 105137.672.8169.7072.515100Reagent 105137.531.7943.792525100Reagent 105137.536.5623.6685.5855100Reagent 105117.536.5623.6685.5855100Reagent 10256.63.93.2428.1184.752.22209.52100Reagent 1056.622.54.9411.224.085100 </td <td></td> <td></td> <td>Reagent 7</td> <td>5</td> <td>63.93</td> <td>9.959</td> <td>5.567</td> <td>15.58</td> <td>5</td> <td>4</td> <td>80</td>			Reagent 7	5	63.93	9.959	5.567	15.58	5	4	80
Reagent 1325289.4422.9755.7447.942525100Reagent 145705.1842.11123.5415.97555100LNAllAll81146.94570.327.2138.6119630Reagent 1563.273.261.2441.7455100Reagent 10563.273.261.8225.1555100Reagent 1125189.187.7151.7943.702525100Reagent 1225158.568.866.7355100Reagent 122394.8318.7645.1184.752221100Reagent 122394.8318.7645.1184.752221100Reagent 125608.3932.4428.1115.332525100Reagent 1325618.273.68020.0169.2955100Reagent 145624.2625.44414.2524.0855100Reagent 1325251001.733.645.185.55100Reagent 14532.411.1481.7125.551001.781.78Reagent 132532.11.1590.0974.9655100Reagent 14532.10.1590.0974.965 <td></td> <td></td> <td>Reagent 10</td> <td>5</td> <td>311.77</td> <td>16.079</td> <td>8.988</td> <td>5.16</td> <td>5</td> <td>5</td> <td>100</td>			Reagent 10	5	311.77	16.079	8.988	5.16	5	5	100
Reagent 1575.184.2.112.3.5415.97555100LNAll81146.9457.0327.92138.8119630Reagent 722123.447.7232.1076.262221100Reagent 105127.62.261.2441.7455100Reagent 12563.273.261.8225.1555100Reagent 135235.671.5.858.866.7355100Reagent 145135.761.7.872.8465.7355100Reagent 105137.511.7.943.792.52.5100Reagent 105137.518.6623.6685.5855100Reagent 105630.9135.8052.0169.9255100Reagent 105630.9135.8052.0169.9255100Reagent 105630.913.5802.0169.9255100Reagent 105630.913.5802.0169.9255100Reagent 105640.215.411.1780.32111.172.22.09.52PIINP81.821.2741.7481.71255100Reagent 1051.821.7370.1821.7474.1255			Reagent 13	25	289.44	22.975	5.744	7.94	25	25	100
LN202011All81146.9457.0327.92138.8119630Reagert 122123.447.7232.1076.262221100Reagert 1563.273.261.8225.1555100Reagert 1325189.187.1751.7943.792525100Reagert 14525.6715.858.866.7355100Reagert 122394.8317.867724.81639.0419840Reagert 122394.8317.867724.81639.0520.01100100Reagert 1530.0935.00520.0169.9255100Reagert 105117.536.5623.6685.585100Reagert 1325608.3932.4428.1115.332525100Reagert 1455100Reagert 13256.9332.4428.1115.33255100Reagert 1325608.3932.4428.1115.33255100 <t< td=""><td></td><td></td><td>Reagent 14</td><td>5</td><td>705.18</td><td>42.111</td><td>23.541</td><td>5.97</td><td>5</td><td>5</td><td>100</td></t<>			Reagent 14	5	705.18	42.111	23.541	5.97	5	5	100
Reagent 1 22 123.44 7.723 2.107 6.26 22 21 100 Reagent 7 5 127.6 2.226 1.244 1.74 5 5 100 Reagent 10 5 63.27 3.26 1.822 5.15 5 5 100 Reagent 10 5 235.67 15.85 8.86 6.73 5 5 100 Reagent 1 52 394.83 18.764 5.118 4.75 22 21 100 Reagent 1 52 394.83 18.764 5.118 4.75 22 21 100 Reagent 1 52 360.99 35.805 20.016 9.92 5 5 100 Reagent 1 51 662.26 2.668 5.88 5 100 Reagent 1 51 10.54 1.178 0.321 11.17 22 20 95.2 Reagent 1 5 14.78 0.989 0.247	LN	202011	All	81	146.94	57.032	7.921	38.81	19	6	30
Pill NP 202012Reagent 75127.61.2261.2441.7455100Reagent 10563.273.261.8225.1555100Reagent 1325189.187.1751.7943.792525100Reagent 145235.6715.858.866.7355100Reagent 175360.9715.858.866.732221100Reagent 105177.577.8.164.752221100Reagent 132566.02935.80520.0169.9255100Reagent 132563.6935.4228.115.3355100Reagent 132564.2625.49414.2524.0855100Reagent 145624.2625.49414.2524.0855100Reagent 132110.541.1780.32111.17222095.2Reagent 132518.27131.28417.48817.1255100Reagent 1453.210.1590.08317.1255100Reagent 132514.780.9890.2476.692525100Reagent 1453.120.1590.98592221100Reagent 15510.011.74741.2921418.2 <t< td=""><td></td><td></td><td>Reagent 1</td><td>22</td><td>123.44</td><td>7.723</td><td>2.107</td><td>6.26</td><td>22</td><td>21</td><td>100</td></t<>			Reagent 1	22	123.44	7.723	2.107	6.26	22	21	100
Reagent 10563.273.261.8225.1555100Reagent 1325189.187.1751.7943.792525100Reagent 145235.6715.858.866.7355100202012All81457.63178.67724.81639.0419840Reagent 151360.9935.80520.0169.9255100Reagent 105117.536.5623.6685.5855100Reagent 1325608.3932.4428.1115.332525100Reagent 145624.2625.49414.2524.0855100PIIINP20211All8312.364.9780.683117.17222095.2Reagent 132516.5417.170.2111.17222095.2100Reagent 1453.210.1590.0894.0272525100Reagent 132514.780.9890.2476.692525100Reagent 1453.210.1590.0894.9655100Reagent 132514.780.9890.2476.692525100Reagent 1454.20.2370.1325.655100Reagent 1524.1420.5361.74741.2920			Reagent 7	5	127.6	2.226	1.244	1.74	5	5	100
Reagent 1325199.187.1751.7943.792525100Reagent 45235.6715.858.866.7355100202012All81457.63178.67724.81639.0419840Reagent 122394.8318.7645.1184.752221100Reagent 0536.0935.6520.0169.9255100Reagent 10517.536.523.6485.585100Reagent 105608.3932.4428.1115.332525100Reagent 1210.5425.49414.2524.0855100PIINP202011All8312.364.9780.68340.2721731.8Reagent 12110.541.1780.32111.17222095.2100Reagent 353.210.590.2426.69255100Reagent 453.210.1790.890.2476.692525100Reagent 32514.763.1290.85492221100Reagent 454.20.2370.1325.655480202012All8330.8512.7361.74741.2921418.2Reagent 559.80.1290.85492			Reagent 10	5	63.27	3.26	1.822	5.15	5	5	100
Reagent 145235.6715.858.866.7355100202012All81457.63178.67724.81639.0419840Reagent 122394.8318.7645.1184.752221100Reagent 05360.9735.80520.0169.9255100Reagent 1325608.3932.4428.1115.332525100Reagent 145624.2625.49414.2524.0855100Reagent 132110.541.1780.32111.17222095.2Reagent 1053.210.1590.08340.2721731.8Reagent 1053.210.1590.08340.2721731.8Reagent 1053.210.1590.0894.9655100Reagent 132514.780.9290.2476.692525100Reagent 132514.780.9290.2476.692525100Reagent 132514.780.9290.2476.69255100Reagent 132529.3812.7361.74741.292148.2202012All8330.8512.7361.74741.292148.2Reagent 13259.380.1280.1331.7535			Reagent 13	25	189.18	7.175	1.794	3.79	25	25	100
202012All81457.63178.67724.81639.0419840Reagent 122394.8318.7645.1184.752221100Reagent 75360.9935.80520.0169.9255100Reagent 105117.536.5623.6685.5855100Reagent 1325608.3932.4428.1115.332525100Reagent 145624.2625.49414.2524.0855100Reagent 1210.541.1780.62111.17222095.2Reagent 1210.541.1780.21111.172225100Reagent 1053.210.1590.0894.9655100Reagent 132514.780.9890.2476.692525100Reagent 132514.780.9890.2476.692525100Reagent 1454.20.2370.1325.655100Reagent 1524.120.2370.1325.655100Reagent 165618.64108.4660.63317.535100Reagent 165618.64108.4660.63317.535100Reagent 1659.380.1031.775480Reagent 1759.380.103 <td></td> <td></td> <td>Reagent 14</td> <td>5</td> <td>235.67</td> <td>15.85</td> <td>8.86</td> <td>6.73</td> <td>5</td> <td>5</td> <td>100</td>			Reagent 14	5	235.67	15.85	8.86	6.73	5	5	100
Reagent 122394.8318.7645.1184.752221100Reagent 75360.9935.80520.0169.9255100Reagent 105117.536.5623.6685.5855100Reagent 1325608.3992.4428.1115.332525100Reagent 145624.2625.49414.2524.0855100PIINP202011All8312.364.9780.68340.2721731.8Reagent 75182.7131.28417.180.12111.7222095.7Reagent 1053.210.1590.0894.9655100Reagent 132514.780.9890.2476.692525100Reagent 1453.210.1590.0894.96554080Reagent 132514.780.9890.2476.692525100Reagent 1454.20.2370.1325.655480Reagent 132529.381.27361.74741.2921418.2Reagent 132529.380.1250.5655100Reagent 132529.382.1420.5367.29252080Reagent 14550.1031.75355100<		202012	All	81	457.63	178.677	24.816	39.04	19	8	40
Reagent 75360.9935.80520.0169.9255100Reagent 105117.536.5623.6685.5855100Reagent 1325608.3932.4428.1115.332525100PIINP202011All8312.364.9780.68340.2721731.8Reagent 12110.541.1780.32111.17222095.2Reagent 153.210.1590.0894.9655100Reagent 1053.210.1590.0894.9655100Reagent 112514.780.9890.2476.692525100Reagent 1254.20.2370.1325.65100Reagent 14534.763.1290.85492221100Reagent 12134.763.1290.85492221100Reagent 12134.763.1290.85492221100Reagent 1259.381.0846.0331.75355100Reagent 13259.380.1650.1031.975480Col IV202011All8878.3332.664.35241.6922836.4Reagent 13259.380.1631.97548036.4 </td <td></td> <td></td> <td>Reagent 1</td> <td>22</td> <td>394.83</td> <td>18.764</td> <td>5.118</td> <td>4.75</td> <td>22</td> <td>21</td> <td>100</td>			Reagent 1	22	394.83	18.764	5.118	4.75	22	21	100
Reagent 105117.536.5623.6685.5855100Reagent 1325608.3932.4428.1115.332525100Reagent 145624.2625.49414.2524.0855100PIINP202011All8312.364.9780.68340.2721731.8Reagent 12110.541.1780.32111.17222095.2Reagent 75182.7131.28417.48817.1255100Reagent 1053.210.1590.0894.9655100Reagent 132514.780.9890.2476.692525100Reagent 1454.20.2370.1325.655480202012All8330.8512.73611.7741.2921482Reagent 12134.763.1290.85492221100Reagent 1450.2382.1420.5367.29252080Reagent 1459.380.1850.1031.9755100Reagent 1459.380.1850.1031.975480Col IV202014All8878.3332.664.35241.692221100Reagent 1459.380.1850.1031.97			Reagent 7	5	360.99	35.805	20.016	9.92	5	5	100
Reagent 13 25 608.39 32.442 8.111 5.33 25 25 100 PIIINP 202011 All 83 12.36 4.978 0.683 40.27 21 7 31.8 Reagent 1 21 10.54 1.178 0.321 11.17 22 20 95.2 Reagent 7 5 182.71 31.284 17.488 17.12 5 5 100 Reagent 10 5 3.21 0.159 0.089 4.96 5 5 100 Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 0.85 12.736 1.747 41.29 21 4.0 100 <t< td=""><td></td><td></td><td>Reagent 10</td><td>5</td><td>117.53</td><td>6.562</td><td>3.668</td><td>5.58</td><td>5</td><td>5</td><td>100</td></t<>			Reagent 10	5	117.53	6.562	3.668	5.58	5	5	100
Reagent 145624.2625.49414.2524.0855100PIIINP202011All8312.364.9780.68340.2721731.8Reagent 12110.541.1780.32111.17222095.2Reagent 75182.7131.28417.48817.1255100Reagent 1053.210.1590.0894.9655100Reagent 132514.780.9890.2476.692525100Reagent 1454.20.2370.1325.655480202012All8330.8512.7361.74741.2921418.2Reagent 12134.763.1290.85492221100Reagent 1256.210.2480.1393.9955100Reagent 13259.380.1650.1331.75352080Reagent 13259.380.1650.1331.975480Col IV202014All8878.3332.664.35241.69228136.4Reagent 159.380.1650.1031.975480Col IV202014All8878.3332.664.35241.692221100Reagent 159.4824.011			Reagent 13	25	608.39	32.442	8.111	5.33	25	25	100
PHINP 202011 All 83 12.36 4.978 0.683 40.27 21 7 31.8 Reagent 1 21 10.54 1.178 0.321 11.17 22 20 95.2 Reagent 7 5 182.71 31.284 17.488 17.12 5 5 100 Reagent 10 5 3.21 0.159 0.089 4.96 5 5 100 Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagen			Reagent 14	5	624.26	25.494	14.252	4.08	5	5	100
Reagent 1 21 10.54 1.178 0.321 11.17 22 20 95.2 Reagent 7 5 182.71 31.284 17.488 17.12 5 5 100 Reagent 10 5 3.21 0.159 0.089 4.96 5 5 100 Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 31.29 0.854 9 22 21 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.48	PIIINP	202011	All	83	12.36	4.978	0.683	40.27	21	7	31.8
Reagent 7 5 182.71 31.284 17.488 17.12 5 5 100 Reagent 10 5 3.21 0.159 0.089 4.96 5 5 100 Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All			Reagent 1	21	10.54	1.178	0.321	11.17	22	20	95.2
Reagent 10 5 3.21 0.159 0.089 4.96 5 5 100 Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 7 5 618.64 108.464 60.633 17.53 5 5 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All			Reagent 7	5	182.71	31.284	17.488	17.12	5	5	100
Reagent 13 25 14.78 0.989 0.247 6.69 25 25 100 Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 7 5 618.64 108.464 60.633 17.53 5 5 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1<			Reagent 10	5	3.21	0.159	0.089	4.96	5	5	100
Reagent 14 5 4.2 0.237 0.132 5.65 5 4 80 202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 7 5 618.64 108.464 60.633 17.53 5 5 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 10<			Reagent 13	25	14.78	0.989	0.247	6.69	25	25	100
202012 All 83 30.85 12.736 1.747 41.29 21 4 18.2 Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 7 5 618.64 108.464 60.633 17.53 5 5 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.45 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 <td></td> <td></td> <td>Reagent 14</td> <td>5</td> <td>4.2</td> <td>0.237</td> <td>0.132</td> <td>5.65</td> <td>5</td> <td>4</td> <td>80</td>			Reagent 14	5	4.2	0.237	0.132	5.65	5	4	80
Reagent 1 21 34.76 3.129 0.854 9 22 21 100 Reagent 7 5 618.64 108.464 60.633 17.53 5 5 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.45 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 14 5		202012	All	83	30.85	12.736	1.747	41.29	21	4	18.2
Reagent 7 5 618.64 108.464 60.633 17.53 5 5 100 Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.45 6 5 83.3 Reagent 7 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14			Reagent 1	21	34.76	3.129	0.854	9	22	21	100
Reagent 10 5 6.21 0.248 0.139 3.99 5 5 100 Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.45 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 14 5 56.4 4.834 2.702 9.57 25 22 22 100			Reagent 7	5	618.64	108.464	60.633	17.53	5	5	100
Reagent 13 25 29.38 2.142 0.536 7.29 25 20 80 Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.45 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14 5 56.4 4.834 2.702 9.57 25 22 22 22 24 22 24 100			Reagent 10	5	6.21	0.248	0.139	3.99	5	5	100
Reagent 14 5 9.38 0.185 0.103 1.97 5 4 80 Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.455 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100			Reagent 13	25	29.38	2.142	0.536	7.29	25	20	80
Col IV 202011 All 88 78.33 32.66 4.352 41.69 22 8 36.4 Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.455 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14 5 56.4 4.834 2.702 9.57 25 22 0.5 9			Reagent 14	5	9.38	0.185	0.103	1.97	5	4	80
Reagent 1 22 94.82 4.011 1.094 4.23 22 21 100 Reagent 7 5 100 5.452 3.048 5.455 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14 5 56.4 4.834 2.702 9.57 25 22 05.9	Col IV	202011	All	88	78.33	32.66	4.352	41.69	22	8	36.4
Reagent 7 5 100 5.452 3.048 5.45 6 5 83.3 Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14 5 56.4 4.834 2.702 8.57 25 22 05.9			Reagent 1	22	94.82	4.011	1.094	4.23	22	21	100
Reagent 10 5 107.02 12.584 7.035 11.76 5 5 100 Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14 5 56.4 4.834 2.702 8.57 25 22 05.9			Reagent 7	5	100	5.452	3.048	5.45	6	5	83.3
Reagent 13 25 47.65 4.121 1.051 8.65 5 5 100 Reagent 14 5 56.4 4.834 2.702 8.57 25 22 05.9			Reagent 10	5	107.02	12.584	7.035	11.76	5	5	100
Reagent 1/ 5 5/ / / 83/ 2702 9.57 25 22 05.0			Reagent 13	25	47.65	4.121	1.051	8.65	5	5	100
			Reagent 14	5	56.4	4.834	2.702	8.57	25	23	95.8
Reagent 17 6 64 17.146 8.75 26.79 5 5 100			Reagent 17	6	64	17.146	8.75	26.79	5	5	100
202012 All 88 306.87 110.217 14.686 35.92 22 10 45.5		202012	All	88	306.87	110.217	14.686	35.92	22	10	45.5
Reagent 1 22 385.22 17.664 4.818 4.59 22 21 100			Reagent 1	22	385.22	17.664	4.818	4.59	22	21	100
Reagent 7 5 401.93 52.768 29.498 13.13 6 4 66.7			Reagent 7	5	401.93	52,768	29.498	13.13	6	4	66.7
Reagent 10 5 408.7 30.168 16.864 7.38 5 5 100			Reagent 10	5	408 7	30,168	16.864	7.38	5	5	100
Reagent 13 25 254.99 23.069 5.886 9.05 5 5 100			Reagent 13	25	254.99	23.069	5.886	9.05	5	5	100
Reagent 14 5 235.62 6.722 3.758 2.85 25 24 100			Reagent 14	5	235.62	6,722	3.758	2.85	25	24	100
Reagent 17 6 146.89 38.783 19.791 26.4 5 5 100			Reagent 17	6	146.89	38.783	19.791	26.4	5	5	100

TABLE 4 (Continued)

	Batch number	Group	Total number	Robust mean	Robust SD	Standard uncertainty	Robust CV	Total number	Number of passing labs	Pass rate (%)
CG	202011	All	60	2.47	0.654	0.106	26.46	30	19	63.3
		Reagent 1	5	3.44	0.707	0.395	20.58	5	4	80
		Reagent 13	15	2.74	0.175	0.056	6.4	15	15	100
		Reagent 26	10	2.46	0.187	0.074	7.61	10	10	100
	202012	All	60	14.12	5.9	0.952	41.8	30	8	26.7
		Reagent 1	5	24.02	4.173	2.333	17.37	5	5	100
		Reagent 13	15	12.93	0.35	0.113	2.71	15	15	100
		Reagent 26	10	17.28	3.037	1.2	17.58	10	9	90

TABLE 5 The total passing rates

Project	Total number	Number of passing laboratories	Pass rate (%)
НА	81	68	84.0
LN	81	65	80.2
PIIINP	83	56	67.5
Col IV	88	68	77.3
CG	60	35	58.3

diagnostic industry; however, it also highlights the significance of conducting the necessary interlaboratory quality assessment method. Therefore, the NCCL carried out an EQA survey for the quantitative identification of liver fibrosis biomarkers, which contributes to the preparation for a formal EQA program for the identification of liver fibrosis indicators.

There are three accepted methods for determining the evaluation criteria of EQA. The first method is based on clinical use-value, which is based on the experience of clinicians (eg glycosylated hemoglobin). The second method is based on biological variation comprising a few indices that have been studied and can be found on the Westgard (https://www.westgard.com/minimum-biodatabas e1.htm) and EFLM (https://biologicalvariation.eu/) websites (eg cholesterol). The third method is based on existing detection technology.¹⁴⁻¹⁸ The detection of liver fibrosis biomarkers uses the third method. For the data collected, the pass rates according to the target values of $\pm 10\%$, 15%, 20%, 25%, and 30% were calculated. These results help to determine the appropriate criteria for existing detection capabilities.

According to the results of the EQA study conducted in China, the CLIA technique accounted for the greatest proportion of the quantitative detection of HA, LN, PIIINP, and Col IV, representing around 90%, whereas, for the identification of CG, the CLIA, HEI, and latex-enhanced immunoturbidimetry methods accounted for approximately one-third each (Table 1). This suggested that despite the fact that there are multiple detection techniques, they are considerably concentrated. On analyzing the reported data, the results revealed significant differences based on the detection techniques employed, even in different orders of magnitude

(Table 2). Even with the similar detection techniques but varied reagents, the detection values varied 3-100-fold (Tables 3 and 4; reagents 1, 7, 10, 13, and 14). Although using a similar reagent and procedure, the numerical variance across various labs differed by more than twofold. These differences could be attributed to the following reasons. The reagents used were different, which adopt different methodologies (some manufacturers use the competition method, whereas others use the sandwich method). Additionally, the binding sites and affinity for quality control products along with sensitivity differ as different reagent manufacturers select different antibodies. The nature of the raw materials of quality control products (natural antigen or recombinant antigen) also contributes to the varied results. The quality control products employed in the present research were obtained from human serum, which might include a number of different antibody proteins in addition to the indicators for liver fibrosis. Hence, it is reasonable to hypothesize that the significant disparities in test findings across labs are due to changes in reagents, procedures, and characteristics, such as anti-interference effects and specificity. In accordance with the existing regulations for medical device registration, the detection system that is formed of the reagents and their calibration products provided by the reagent manufacturers, as well as the "suitable equipment" listed on the kit's specifications, can serve as a supporting system, which may be defined as the "open" supporting system. The reasons for the substantial differences in outcomes between the "closed" and "open" matching methods are complicated. The development of various analytical systems, such as the data reading techniques, configuration of absorbance wavelengths, and built-in calibrations, could be a possibility. Besides, plenty of invalid data were removed (such as wrong filling, missing file, and unit error) when analyzing the results, suggesting that the data reporting behavior of laboratory staff needs to be standardized. Therefore, the human factor in the data submission process cannot be ignored. According to the present EQA data, it is impossible to isolate a specific factor.

A large number of labs engaged in the quantitative identification of liver fibrosis biomarkers in China. Although the pass rate for HA and LN was relatively high in this survey, it does not suggest that the quality of this marker satisfied the clinical standards (Tables 4 and 5) in view of the fact that the evaluation criterion of target value set

^{12 of 12} WILEY

at \pm 30% is not very strict. At the same time, we found that the pass rate of the individual group (grouped according to ISO 13528) was substantially greater than the one of the other groups. This not only shows the necessity for grouping but also suggests that reagents with a small market share are at a disadvantage in grouping (Table 4). In addition to objective factors, many subjective factors also affect the results including data input error, reagent or method selection error, unit error, incorrect sample sequence use, internal quality control, and personnel operation error. All these factors restrict the accuracy of EQA results.

5 | CONCLUSION

The present research is the first report to reveal the current state of the quantitative identification of liver fibrosis biomarkers in China, as far as we know. Despite the small number of labs involved in the present research, our findings indicate that there are a few drawbacks in the quantitative detection of liver fibrosis biomarkers. Commercially accessible quantitative detection kits need methodological study, and quality control procedures must be improved for the purpose of better understanding the reasons for differences in detection findings. Overall, the present research establishes the basis for the development of a formal EQA for the quantitative identification of liver fibrosis biomarkers. At the same time, it also provides guidance for the selection of instruments, methods, and reagents for laboratory liver fibrosis marker detection.

AUTHOR CONTRIBUTIONS

Chao Zhang was involved in acquisition and analysis of data. Chao Zhang and Chuanbao Zhang participated in management of data. Chao Zhang designed the study and drafted the manuscript.

ACKNOWLEDGMENTS

We thank Autobio for providing quality control products.

CONFLICT OF INTEREST

The authors declare no competing interests.

DATA AVAILABILITY STATEMENT

Not applicable.

ORCID

Chao Zhang () https://orcid.org/0000-0002-3836-5107

REFERENCES

- Adams LA. Biomarkers of liver fibrosis. J Gastroenterol Hepatol. 2011;26:802-809.
- Caballeria L, Toran P, Caballeria J. Markers of hepatic fibrosis. Med Clin (Barc). 2018;150:310-316.

- Kisseleva T, Brenner D. Molecular and cellular mechanisms of liver fibrosis and its regression. Nat Rev Gastroenterol Hepatol. 2021;18:151-166.
- Dhar D, Baglieri J, Kisseleva T, Brenner DA. Mechanisms of liver fibrosis and its role in liver cancer. *Exp Biol Med (Maywood)*. 2020;245:96-108.
- Masuzaki R, Kanda T, Sasaki R, et al. Noninvasive assessment of liver fibrosis: current and future clinical and molecular perspectives. *Int J Mol Sci.* 2020;21(14):4906.
- Tong S, Revill P. Overview of hepatitis B viral replication and genetic variability. J Hepatol. 2016;64:S4-S16.
- Tang LSY, Covert E, Wilson E, Kottilil S. Chronic hepatitis B infection: a review. JAMA. 2018;319:1802-1813.
- El-Mezayen HA, Habib S, Marzok HF, Saad MH. Diagnostic performance of collagen IV and laminin for the prediction of fibrosis and cirrhosis in chronic hepatitis C patients: a multicenter study. *Eur J Gastroenterol Hepatol.* 2015;27:378-385.
- Enomoto H, Bando Y, Nakamura H, Nishiguchi S, Koga M. Liver fibrosis markers of nonalcoholic steatohepatitis. World J Gastroenterol. 2015;21:7427-7435.
- Neuman MG, Cohen LB, Nanau RM. Hyaluronic acid as a noninvasive biomarker of liver fibrosis. *Clin Biochem*. 2016;49:302-315.
- Mak KM, Mei R. Basement membrane type IV collagen and laminin: an overview of their biology and value as fibrosis biomarkers of liver disease. Anat Rec (Hoboken). 2017;300:1371-1390.
- 12. Adams LA, Wang Z, Liddle C, et al. Bile acids associate with specific gut microbiota, low-level alcohol consumption and liver fibrosis in patients with non-alcoholic fatty liver disease. *Liver Int*. 2020;40:1356-1365.
- Loomba R, Adams LA. Advances in non-invasive assessment of hepatic fibrosis. Gut. 2020;69:1343-1352.
- Kenny D, Fraser CG, Hyltoft Petersen P, Kallner A. Consensus agreement. Scand J Clin Lab Invest. 1999;59(7):585.
- Kallner A, McQueen M, Heuck C. The Stockholm Consensus Conference on quality specifications in laboratory medicine, 25-26 April 1999. Scand J Clin Lab Invest. 1999;59(7):475.
- Ceriotti F, Fernandez-Calle P, Klee GG, et al. Criteria for assigning laboratory measurands to models for analytical performance specifications defined in the 1st EFLM Strategic Conference. *Clin Chem Lab Med.* 2017;55:189-194.
- 17. Kallner A. Quality specifications based on the uncertainty of measurement. *Scand J Clin Lab Invest*. 1999;59:513-516.
- Sandberg S, Thue G. Quality specifications derived from objective analyses based upon clinical needs. Scand J Clin Lab Invest. 1999;59:531-534.

SUPPORTING INFORMATION

Additional supporting information may be found in the online version of the article at the publisher's website.

How to cite this article: Zhang C, Zhang C. Analysis of current status of quantitative detection of biomarkers for liver fibrosis in Clinical labs in China. *J Clin Lab Anal*. 2022;36:e24490. doi:10.1002/jcla.24490