# Crosstalk between AMPK activation and angiotensin II-induced hypertrophy in cardiomyocytes: the role of mitochondria

Jessica Soto Hernández <sup>a</sup>, Giselle Barreto-Torres <sup>a</sup>, Andrey V. Kuznetsov <sup>b</sup>, Zaza Khuchua <sup>c</sup>, Sabzali Javadov <sup>a, \*</sup>

<sup>a</sup> Department of Physiology, School of Medicine, University of Puerto Rico, San Juan, PR, USA <sup>b</sup> Department of Cardiac Surgery, Cardiac Surgery Research Laboratory, Innsbruck Medical University, Innsbruck, Austria <sup>c</sup> Cincinnati Children's Research Foundation, University of Cincinnati, Cincinnati, OH, USA

Received: August 1, 2013; Accepted: November 28, 2013

# Abstract

AMP-kinase (AMPK) activation reduces cardiac hypertrophy, although underlying molecular mechanisms remain unclear. In this study, we elucidated the anti-hypertrophic action of metformin, specifically, the role of the AMPK/eNOS/p53 pathway. H9c2 rat cardiomyocytes were treated with angiotensin II (AngII) for 24 hrs in the presence or absence of metformin (AMPK agonist), losartan [AngII type 1 receptor (AT1R) blocker], N $_{
m O}$ -nitro-L-arginine methyl ester (L-NAME, pan-NOS inhibitor), splitomicin (SIRT1 inhibitor) or pifithrin- $\alpha$  (p53 inhibitor). Results showed that treatment with metformin significantly attenuated AngII-induced cell hypertrophy and death. Metformin attenuated AngII-induced activation (cleavage) of caspase 3, Bcl-2 down-regulation and p53 up-regulation. It also reduced AngII-induced AT1R up-regulation by 30% (P < 0.05) and enhanced AMPK phosphorylation by 99% (P < 0.01) and P-eNOS levels by 3.3-fold (P < 0.01). Likewise, losartan reduced AT1R up-regulation, indicating that metformin mediated its effects *via* AMPK activation. Beneficial effects of metformin and losartan converged on mitochondria that demonstrated high membrane potential ( $\Delta \psi_m$ ) and low permeability transition pore opening. Thus, this study demonstrates that the anti-hypertrophic effects of metformin are associated with AMPK-induced AT1R down-regulation and prevention of mitochondrial dysfunction through the SIRT1/ eNOS/p53 pathway.

Keywords: AMP kinase • metformin • mitochondria • angiotensin II receptors • hypertrophy

# Introduction

AMP-kinase (AMPK), an important metabolic sensor, regulates energy metabolism in organs and tissues including the liver, adipose tissue, and skeletal and cardiac muscles [1, 2]. In response to various intra- and extracellular stimuli, it enhances fatty acid uptake, and oxidative metabolism, glucose transport, glycolysis and mitochondrial biogenesis. Metformin, a biguanide drug with pleiotropic actions widely used for the treatment of patients with type 2 diabetes [3, 4], is a powerful activator of AMPK. Apart from its anti-hyperglycaemic effects, metformin significantly improves cardiac function in patients and various animal models of heart failure [5], although underlying

\*Correspondence to: Sabzali JAVADOV, M.D., Ph.D., Department of Physiology, School of Medicine, University of Puerto Rico, San Juan 00936-5067, PR, USA. Tel.: 787-758-2525; ext. 2909 Fax: 787-753-0120 mechanisms of action remain unclear. For example, in mice with heart failure induced by left coronary artery occlusion and subsequent reperfusion, metformin improved LV function and attenuated cardiac hypertrophy [6]. Similarly, metformin attenuated oxidative stress-induced cardiomyocyte apoptosis and prevented the progression of heart failure in dogs [7]. However, metformin did not exert anti-hyper-trophic effects on mice lacking functional AMPK or eNOS, which suggests that its effects are mediated *via* activation of the AMPK/eNOS pathway [6].

Cardiomyocyte hypertrophy, a major consequence of pressure or volume overload, plays a central role in the progression of heart failure. Various paracrine and autocrine factors are involved in pathogenesis and regulation of cardiomyocyte hypertrophy, including angiotensin II (AngII), a key player of the renin–angiotensin system that induces cell hypertrophy, differentiation and apoptosis through activation of various intracellular signalling molecules including Gq protein, calcineurin, mitogen-activated protein kinases (MAPK) and several transcription factors [8]. AngII type 1 (AT1R)

E-mail: sabzali.iavadov@upr.edu

Journal of Cellular and Molecular Medicine published by John Wiley & Sons Ltd and Foundation for Cellular and Molecular Medicine. This is an open access article under the terms of the Creative Commons Attribution License, which permits use, distribution and reproduction in any medium, provided the original work is properly cited.

doi: 10.1111/jcmm.12220

<sup>© 2014</sup> The Authors.

and type 2 (AT2R) G protein-coupled receptors have been shown to mediate the effects of circulating and local (intracellular) AngII [9, 10]. AT1R mediates pro-hypertrophic effects of AngII, but AT2R, in contrast, attenuates AT1R activation-induced cell growth. Recent studies demonstrated that pro-hypertrophic effects of AngII can be mediated through mitochondria and induce cell death [11]. Although the exact underlying mechanisms remain unclear, AngIIinduced depolarization of the mitochondrial membrane and increased production of mitochondrial reactive oxygen species (ROS) are associated with cardiomyocyte autophagy and hypertrophy [12].

Metformin, a promising pharmacological agent, may be used for the treatment of heart failure. Most previous studies were conducted *in vivo* on intact hearts. However, the heart contains many non-cardiomyocyte cell types, including fibroblasts, vascular endothelial cells, smooth muscle cells and immune cells, with cardiomyocytes accounting for only 30–40% of total heart cells. Accordingly, new studies using cultured cardiomyocytes are required to establish whether metformin exerts a direct anti-hypertrophic effect on these cells.

In this study, we elucidated the role of AMPK in AngII-induced hypertrophy in cultured H9c2 cardiomyocytes, providing additional evidence that AMPK activation improves mitochondrial function. The results demonstrate that metformin exerted anti-hypertrophic effects on cardiomyocytes and prevented AngII-induced cell death. We observed a negative reciprocal interaction between AMPK activation and AT1R levels: metformin activated AMPK, annulling AngII-induced up-regulation of AT1R, whereas losartan (AT1R antagonist) enhanced AMPK activation in AngII-treated H9c2 cardiomyocytes. Furthermore, metformin attenuated mitochondrial dysfunction and hypertrophy through the eNOS/SIRT1/p53 pathway.

## Material and methods

#### **Cell culture**

H9c2 embryonic rat cardiomyocytes were purchased from the American Type Culture Collection (Manassas, VA, USA) and cultured according to the manufacturer's protocol. In short, the cells were cultured in medium containing DMEM/Ham's F-12 (Invitrogen, Carlsbad, CA, USA), supplemented with 10% foetal bovine serum, 10  $\mu$ g/ml transferrin, 10  $\mu$ g/ml insulin, 10 ng/ml selenium, 1% penicillin and streptomycin, 2 mg/ml bovine serum albumin, 5  $\mu$ g/ml linoleic acid, 3 mM pyruvic acid, 0.1 mM minimum essential medium non-essential amino acids, 10% MEM vitamin, 0.1 mM 5-bromo-2'-deoxyuridine, 100  $\mu$ M L-ascorbic acid and 30 mM HEPES, pH 7.1, and maintained in 95% air and 5% CO<sub>2</sub> at 37°C. Prior to all experiments, cells were serum-starved for 24 hrs. Cells with 85–90% confluence from passages 3–20 were used for experiments.

#### **Experimental protocol**

Cells were treated with 200 nM AngII for 24 hrs in the presence or absence of 2 mM metformin (AMPK activator), 10  $\mu$ M losartan (AngII

type 1 receptor antagonist), 5  $\mu M$  compound C (AMPK inhibitor), 300  $\mu M$  N $_{0}$ -nitro-L-arginine methyl ester (L-NAME, pan-NOS inhibitor), 10  $\mu M$  splitomicin (SIRT1 inhibitor) or 10  $\mu M$  pifithrin- $\alpha$  (p53 inhibitor). All drugs were purchased from Sigma-Aldrich (St. Louis, MO, USA) and added to the culture medium 45 min. before AngII administration.

#### **Cell viability**

Cell viability was determined using the trypan blue exclusion method. The cells were cultured at a density of 5  $\times$  10<sup>5</sup> cells on 100-mm dishes and exposed to AngII (0 nM, 200 nM, 500 nM, 1  $\mu$ M or 5  $\mu$ M) for 24 hrs in the presence or absence of 2 mM metformin. After treatment, cells were rinsed with PBS, detached with Hyclone Trypsin (Thermo Fisher Scientific, Waltham, MA, USA) and counted using the TC20 Automated Cell Counter (Bio-Rad, Hercules, CA, USA). The per cent of viable and dead cells was calculated from a total number of counted cells.

#### Cell area measurement

Plated cells ( $3 \times 10^4$  cells per 6-cm dish) were treated in accordance with the experimental protocol. At the end of the experiment, the cells were viewed using a Leica DMIL LED inverted microscope (Wetzlar, Germany), equipped with a 3.2 MP scientific grade digital microscopy camera. Six to eight random photographs were taken from each dish, and at least, 60 individual cell size measurements were made from each group using Accu-Scope Micrometrics SE Premium software.

#### **SDS-PAGE and Western blotting**

Following the treatments, cells were washed twice with cold PBS and lysed in buffer containing 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 10 mM tetrasodium phosphate, 40 mM β-glycerophosphate, 2 mM ethylenediaminetetraacetate (EDTA), 2 mM ethylene glycol-bis(2-aminoethylether)-N,N, N',N'-tetraacetate (EGTA), 1% Triton X, 10% glycerol, 50 mM sodium fluoride, 0.2 mM sodium orthovanadate and protease inhibitors. The cell lysates were transferred to 1.5-ml Eppendorf tubes, homogenized and centrifuged at 10,000  $\times$  g for 5 min. at 4°C. The supernatant was transferred to a fresh tube, and protein concentration was determined using the Bio-Rad Protein Assay Reagent (Bio-Rad). Equal amounts of protein samples were resolved on SDS-PAGE and transferred onto Amersham Hvbond ECL 0.45 µM nitrocellulose membranes (GE Healthcare Bio-Sciences, Piscataway, NJ, USA). The membranes were immunoblotted with actin, AMPK, P-AMPK $\alpha^{Thr172}$ , Bcl-2, caspase 3 full length, eNOS, PeNOS (Cell Signaling, Danvers, MA, USA), AT1R, AT2R or p53 (Santa Cruz, Dallas, TX, USA) antibodies, followed by secondary antibodies. The signals were visualized using Thermo Scientific Pierce ECL Western Blotting Detection reagents (Thermo Fisher Scientific) at VersaDoc 3000 Gel Imaging System (Bio-Rad).

#### Mitochondrial membrane potential

To monitor  $\Delta \Psi_{m}$ , cells (4  $\times$  10<sup>5</sup> cells/well) plated in a 24-well culture plate were incubated for 30 min. with the membrane potential-sensitive dye JC-1 (10 µg/ml, 5,5',6,6'-tetraethyl-benzimidazolylcarbocyanine

iodide, Molecular Probes, Eugene, OR, USA). The intensity of fluorescence was immediately measured using a Spectramax M3 microplate reader (Molecular Devices) at 527 and 590 nm for emission and 488 nm for excitation. For confocal microscopy, the cells, plated on glass-bottom dishes, were loaded with JC-1 for 30 min. at 37°C, and images were captured using a Zeiss LSM510 META (Carl Zeiss, Oberkochen, Germany) microscope.

#### Mitochondrial permeability transition pore (PTP)

To quantify mitochondrial PTP opening, cardiomyocytes were loaded with 5  $\mu$ M calcein-acetoxymethylester (calcein-AM, Molecular Probes). Cobalt chloride (5 mM) was added to the culture medium to quench cytosolic and nuclear calcein. Cells were visualized with an Olympus IX73 inverted fluorescence microscope (Center Valley, PA, USA), and fluorescence was captured using an Olympus DP73 high-performance Peltier cooled digital colour camera. Olympus CellSens Dimension Imaging software was used to analyse calcein fluorescence in cells.

#### **Complex I activity**

Complex I activity was determined by measuring the decrease in the concentration of NADH at 340 nm and 30°C. Cells were resuspended in hypotonic phosphate buffer containing 5 mM MgCl<sub>2</sub> and 0.55 mg/ml saponin and, then, freeze-thawed three times to give the substrate access to the inner mitochondrial membrane. The assay was performed in phosphate buffer containing 1 mM KCN, 5 mM MgCl<sub>2</sub>, 2.5 mg/ml BSA, 2  $\mu$ M antimycin, 100  $\mu$ M decylubiquinone and 100  $\mu$ M NADH, pH 7.4. Specific activities were determined by calculating the slope of the reaction in the linear range and normalizing per mg of protein.

#### **Total ROS levels**

Cells were trypsinized with a 0.25% trypsin–EDTA solution (Thermo Scientific HyClone, Logan, UT, USA) and centrifuged at  $500 \times g$  for 5 min. at room temperature. The pellet was resuspended in culture medium, and cells were incubated with 20  $\mu$ M 2'7 dichlorofluorescein diacetate (DCF-DA; Alexis Biochemicals) for 30 min. at 37°C. The cells were again centrifuged at  $500 \times g$  for 5 min. to remove medium with excess dye, and the pellets were resuspended in PBS and added to a 96-well plate. Fluorescence intensity was measured using a Spectramax M3 plate reader (Molecular Devices, Sunnyvale, CA, USA) at an excitation of 485 nm and emission of 530 nm.

#### **Mitochondrial ROS levels**

To quantify mitochondrial ROS, cells (2 × 10<sup>5</sup> cells/well) in the 96-well culture plate were incubated with the mitochondrial superoxide-sensitive fluorescent dye MitoSOX Red (1  $\mu$ M, Invitrogen) for 20 min. at 37°C. The intensity of fluorescence was measured using a Spectramax M3 microplate reader (Molecular Devices) at 510 and 580 nm for excitation and emission respectively.

#### Statistical analysis

Data are presented as means  $\pm$  SEM. Statistical significance was evaluated using one-way ANOVA, followed by Tukey's multiple comparison *post hoc* test or an unpaired two-tailed Student's *t*-test. Differences were considered to be statistically significant when P < 0.05.



Fig. 1 The effects of metformin (Met) on cell size. Cell size was measured in H9c2 cells 24 hrs after treatment with AngII (200 nM) in the presence or absence of losartan (Los), Met, N $_{00}$ -nitro-L-arginine methyl ester (L-NAME) and splitomicin (SP) or pifithrin- $\alpha$  (Pf). (A) Representative images of individual cells, (B) cell size ( $\mu$ m), (C) cell size presented as percentage of control. \**P* < 0.05, \*\**P* < 0.01 *versus* control; \**P* < 0.05 *versus* AngII. *n* = 3–4 per group.

 $\odot$  2014 The Authors.

Journal of Cellular and Molecular Medicine published by John Wiley & Sons Ltd and Foundation for Cellular and Molecular Medicine.



**Fig. 2** The effects of metformin (Met) on AnglI-induced cell death. (**A**) The effect of AnglI (24 hrs), at concentrations ranging from 0 to 5  $\mu$ M, on cell death, in the presence or absence of 2 mM Met was measured by Trypan blue. (**B** and **C**) The effects of Met on caspase 3 (full size) and Bcl-2 protein levels in cells treated with 200 nM or 1  $\mu$ M AnglI. \**P* < 0.05, \*\**P* < 0.01 *versus* control; \**P* < 0.05 *versus* AnglI. *n* = 3–4 per group (for apoptotic proteins) and 6–8 (for cell death) per group.



**Fig. 3** The effect of losartan (Los), metformin (Met), N $_{\odot}$ -nitro-L-arginine methyl ester (L-NAME) or splitomicin (SP) on AngII type 1 (AT1R, **A**) and type 2 (AT2R, **B**) receptors in cell hypertrophy. Western blot data were normalized to actin and expressed as fold change relative to control groups. \*P < 0.05 versus control; \*P < 0.05 versus Met. n = 6-8 per group.



Fig. 5 The effects of AnglI on phospho-eNOS (P-eNOS, **A**) and total-eNOS (T-eNOS, **B**) in cardiomyocytes in the presence or absence of losartan (Los), metformin (Met), N $\infty$ -nitro-L-arginine methyl ester (L-NAME), and/or splitomicin (SP). Western blot data were normalized to actin and expressed as fold change relative to control groups. \*P < 0.05, \*\*P < 0.01 versus control; \*P < 0.05 versus AnglI;  $^{\&}P < 0.05$  versus AnglI+Met. n = 8-12 per group.

### Results

# Metformin attenuates AnglI-induced hypertrophy and cell death

In the first set of experiments, we examined whether metformin exerts a direct anti-hypertrophic effect on AngII-treated cardiomyocytes. As

previous studies demonstrate that metformin at high concentrations ( $\geq 5$  mM) reduces the energy state of intact hearts [4], but at 2 mM, it exerts no direct effect on complex I activity [13] and high-energy phosphates [14], we used metformin at a concentration of 2 mM. Likewise, this study revealed no detrimental effects of 2 or 5 mM metformin on complex I activity in H9c2 cardiomyocytes (see below).

As shown in Figure 1, cells treated for 24 hrs with 200 nM Angll demonstrated significant hypertrophy and a 22% (P < 0.01) increase

![](_page_5_Figure_0.jpeg)

**Fig. 6** The effects of AngII on the level of p53 in cardiomyocytes in the presence or absence of losartan (Los), metformin (Met), N $\omega$ -nitro-L-arginine methyl ester (L-NAME) and/or splitomicin (SP). Western blot data were normalized to actin and expressed as fold change relative to control groups. \**P* < 0.05 *versus* control; \**P* < 0.05 *versus* AngII. *n* = 10–12 per group.

in cell size. Both metformin and losartan prevented cell hypertrophy significantly, reducing cell size in the presence of AngII. Cell size values were not significantly affected in cells pre-treated with L-NAME, splitomicin or pifithrin- $\alpha$ . To assess the effect of metformin on cell death, the cells were pre-treated with 2 mM metformin in the presence of AngII at different concentrations (0 nM, 200 nM, 500 nM, 1  $\mu$ M or 5  $\mu$ M). AngII-induced cell death was dose-dependent and, when treated with 5  $\mu$ M AngII, reached a peak rate of 13.3% (Fig. 2A). Metformin prevented cell death at all concentrations of AngII and, interestingly, demonstrated anti-apoptotic effects on cells treated with 200 nM AngII, preventing caspase 3 cleavage and BcI-2 down-regulation. However, when the cells were treated with 1  $\mu$ M AngII, the presence or absence of metformin did not alter caspase 3 or

Bcl-2 protein levels (Fig. 2B and C), suggesting that cell death occurs through necrosis, rather than apoptosis, at high concentrations of Angll. Thus, these data demonstrate that metformin exerts a direct, antihypertrophic effect on cardiomyocyte hypertrophy and prevents Angll-induced cell death.

#### Metformin decreases AnglI-induced up-regulation of AT1R

AT1R is a predominant AnglI receptor in cardiac cells that is typically coupled to  $G\alpha_{q/11}$  and activates PKC, Ras, MAPKs, JAK-STAT and other pro-hypertrophy signalling molecules [15]. In the first set of experiments, we examined whether metformin affected AT1R levels in AnglI-treated H9c2 cells. As shown in Figure 3A, AnglI up-regulated AT1R and AT2R by 38% (P < 0.05) and 40% (P < 0.05) respectively. The presence of metformin or losartan completely prevented AT1R up-regulation. Also, losartan and metformin further increased AT2R levels by 94% (P = N.S.) and 124% (P < 0.05), respectively, in Angll-treated cells, compared with the control group (Fig. 3B). Neither the presence of compound C, splitomicin, nor L-NAME significantly affected AT2R levels. Thus, our data demonstrate that metformin reduces AT1R and increases AT2R in AngII-treated cells. Compound C (AMPK inhibitor) abolished these effects in the presence of metformin, suggesting that anti-hypertrophic effects of metformin require AMPK activation.

# Inhibition of AT1R is associated with AMPK activation

In the next group of experiments, we examined whether AT1R inhibition induces AMPK activation in cardiomyocytes. Phosphorylation of AMPK was accepted as an indicator of enzymatic activity, based on previous data showing that phosphorylation at Thr172 correlates with AMPK activation [16]. AngII had no significant effect on P-AMPK $\alpha^{Thr172}$ , indicating no activation (Fig. 4A). Pretreatment with losartan and metformin increased AMPK phosphorylation by 54% (P < 0.05) and 99% (P < 0.05), respectively, compared with the control group. Neither treatment with splitomicin

**Fig. 7** Effects of losartan (Los), metformin (Met), N<sub> $\omega$ </sub>-nitro-L-arginine methyl ester (L-NAME), splitomicin (SP) or pifithrin- $\alpha$  (Pf) on mitochondrial membrane potential ( $\Delta \Psi_m$ ) and permeability transition pores (PTP) in AngII-treated cardiomyocytes. (**A**) Representative images of JC-1 fluorescence. The images were obtained using a Zeiss LSM510 META (Carl Zeiss) microscope from cells after incubation with JC-1 (10 µg/ml). The cells were treated with Los, Met, SP, L-NAME or Pf in the presence or absence of 200 nM AngII. In addition, to demonstrate specificity of the fluorescence signal, the cells were treated with FCCP (carbonylcyanide-p-trifluoromethoxyphenylhydrazone), a chemical uncoupler of electron transport and oxidative phosphorylation that induces depolarization of the mitochondrial inner membrane. Red fluorescent images of dye aggregates ('J-aggregates') indicate high  $\Delta \Psi_m$ , whereas green images of monomeric dye ('JC-1 monomers') show low  $\Delta \Psi_m$ . (**B**) Quantitative results of JC-1 fluorescence intensity, measured using a Spectramax M3 microplate reader (Molecular Devices). Data were normalized to cells and expressed as a percentage of the control group. \**P* < 0.05 *versus* control. *n* = 6–8 per group. (**C**) Representative images of calcein fluorescence. Cells treated with Los, Met, SP, L-NAME or Pf in the presence or absence of 200 nM AngII for 24 hrs were co-loaded with cobalt chloride (5 mM) and calcein-AM (5 µM) and, then, imaged using an Olympus IX73 (Center Valley, PA, USA) inverted fluorescence microscope. (**D**) Quantitative results of calcein fluorescence, normalized to individual cells and expressed as a percentage of the control group. \**P* < 0.05 *versus* control; *n* = 4 per group.

![](_page_6_Figure_1.jpeg)

nor L-NAME affected AMPK phosphorylation (Fig. 4A and B). It should be noted that AngII, in the presence or absence of the pharmacological agents, did not affect total AMPK levels (Fig. 4C and D). Thus, these findings demonstrate that losartan may reduce AngII-induced cardiomyocyte hypertrophy, at least partially, *via* AMPK activation.

#### Metformin and losartan increase eNOS activation and reduce p53 in AnglI-induced hypertrophy

Previous studies demonstrate that AMPK phosphorylates eNOS at Ser1177 and increases nitric oxide production in rat hearts [7]. Nitric oxide has opposite effects on cell metabolism and mitochondrial function; it protects cells against oxidative stress at low concentrations but exerts detrimental effects at high concentrations [17, 18]. Recent studies demonstrate that eNOS activation (P-eNOS<sup>Ser1177</sup>) plays a central role in H<sub>2</sub>O<sub>2</sub>- and AngII-induced signalling in cardiac cells [19]. In our experiments, incubation of cells with AngII enhanced eNOS phosphorylation by 2.1 fold (P < 0.05), compared with the control group (Fig. 5A). Surprisingly, both losartan and metformin further increased eNOS activation by 3.3 (P < 0.05) and 3.6 fold (P < 0.05) respectively. Pre-treatment with splitomicin or L-NAME completely prevented AngII-induced phosphorylation. AngII did not affect the total level of eNOS in the presence or absence of the pharmacological agents (Fig. 5B).

Recent studies identified p53, a well-known tumour suppressor protein, as a key player in heart failure. It is activated in response to cellular stress, promoting apoptosis and development of heart failure [20]. AnglI-induced hypertrophy up-regulated p53 and pro-apoptotic signalling pathways in cardiomyocytes [11, 21]. To the best of our knowledge, there are no studies about the role of AMPK activation in AnglI-induced p53 up-regulation. Our studies demonstrated that AnglI increased the p53 levels by 58% (P < 0.05), after 24 hrs of treatment (Fig. 6). Both metformin and losartan prevented p53 up-regulation in AnglI-treated cells. SIRT1 regulates p53 activity through acetylation/deacetylation [22, 23]. Pre-treatment of cells with splitomicin or L-NAME did not significantly change p53 levels (Fig. 6). Thus, these findings demonstrate that metformin reduces AnglI-induced hypertrophy in cardiomyocytes through the regulation of the AMPK/SIRT1/eNOS/p53 pathway.

#### Metformin prevents mitochondrial membrane depolarization and increases the electron transfer chain (ETC) complex I activity in AnglI-treated cardiomyocytes

Hypertrophic signalling pathways mostly converge on mitochondria to regulate energy metabolism in response to various cellular stresses. AMPK plays a central role in mitochondrial metabolism and ATP production by stimulating glucose uptake, fatty acid oxidation and glycolysis. We examined whether metformin prevents AngII-induced mitochondrial dysfunction in H9c2 cells. As shown in Figure 7A and B. Angll reduced the membrane potential of mitochondria by 26% (P < 0.05). All drugs, except splitomicin, prevented AnglI-induced membrane depolarization. Particularly, pifithrin- $\alpha$  prevented membrane depolarization, indicating that AngII-induced mitochondrial dysfunction may be associated with p53 activation. Indeed, p53 translocation to mitochondria may be associated with mitochondrial dysfunction and cell death, according to numerous studies performed by various groups [24, 25]. Furthermore, physical interaction between p53 and cyclophilin D (CyP-D), a regulatory protein of the mitochondrial PTP, was found in the mitochondria of cancer cells [24] and in brains subjected to ischaemia/reperfusion [25]. However, our experiments with immunoprecipitation of p53 followed by immunoblotting of CyP-D revealed no physical interaction between these two proteins or other possible PTP components, such as adenine nucleotide translocase and voltage-dependent anion channels (data not shown). In addition, metformin and other pharmacological agents, except L-NAME, prevented AnglI-induced PTP opening (Fig. 7C and D). Interestingly, however, AnglI did not induce a significant increase in either total or mitochondrial ROS production in the presence or absence of the tested pharmacological agents, except L-NAME, which significantly increased total ROS levels (Fig. 8A and B). Metformin may mediate its beneficial effects through direct inhibition of the ETC complex I activity. To elucidate this hypothesis, we assessed a direct effect of metformin on complex I activity in untreated cells (not pretreated with Angll or pharmacological agents), at concentrations ranging from 0 to 20 mM. Low concentrations (2 and 5 mM) of metformin did not affect complex I activity, but 10 mM metformin inhibited it (Fig. 9A). Interestingly, AnglI alone did not significantly change complex I activity. However, pre-treatment of AnglI-treated cells with metformin increased the activity of the complex (Fig. 9B). Thus, metformin prevents depolarization of mitochondria, inhibits PTP opening and improves ETC activity in AnglI-treated cells. The beneficial effects of metformin on mitochondria are not associated with its direct interaction with complex I of the ETC.

## Discussion

The results of this study demonstrate the following: (i) In response to AngII, there is a negative reciprocal relationship between AMPK and AT1R. Metformin activates AMPK, inducing down-regulation of AT1R, whereas losartan blocks AT1R and stimulates AMPK activation. (ii) Both AMPK activation and AT1R inhibition increase eNOS phosphorylation and down-regulate p53. (iii) The beneficial effects of AMPK activation and AT1R inhibition may be associated with the prevention of mitochondrial dysfunction in cardiomyocytes. In Figure 10, we propose a mechanism for the anti-hypertrophic effect of metformin on AnglI-induced cardiomyocyte hypertrophy. Previous studies demonstrating anti-hypertrophic and anti-remodelling effects of metformin were mostly performed in vivo on intact hearts [7, 26]. However, along with cardiomyocytes, intact hearts contain other cell types, including endothelial and smooth muscle cells, fibroblasts and other connective tissue cells, mast cells and immune system-related cells. This study

![](_page_8_Figure_1.jpeg)

**Fig. 8** Effects of losartan (Los), metformin (Met), N<sub>D</sub>-nitro-L-arginine methyl ester (L-NAME) or splitomicin (SP) on total and mitochondrial reactive oxygen species (ROS) levels in AngII-treated H9c2 cells. (**A**) Total ROS levels were determined in the cells incubated with the ROS-specific probe, DCF-DA, after treatment with 200 nM AngII in the presence or absence of Los, Met, L-NAME or SP. In additional groups, the cells were treated with 50  $\mu$ M H<sub>2</sub>O<sub>2</sub> (positive control) or 1  $\mu$ M AngII for 24 hrs. (**B**) Mitochondrial ROS levels were measured under the same conditions using the mitochondrial-sensitive fluorescence dye, MitoSOX Red. In additional groups, the cells were treated with 50  $\mu$ M H<sub>2</sub>O<sub>2</sub> or 10  $\mu$ M antimycin (positive control) for 24 hrs. Results were normalized per 10<sup>6</sup> cells (for total ROS) and 10<sup>3</sup> cells (for mitochondrial ROS) and expressed as a percentage of the control group. \**P* < 0.05 *versus* control; *n* = 6–8 per group.

![](_page_8_Figure_3.jpeg)

**Fig. 9** Effects of losartan (Los), metformin (Met),  $N_{\omega}$ -nitro-L-arginine methyl ester (L-NAME) or splitomicin (SP) on the activity of ETC. complex I, in the presence or absence of AngII. (**A**) Met, at concentrations ranging from 0 to 20 mM, was added directly to cell lysates isolated from control H9c2 cardiomyocytes, and the activity of complex I was measured. (**B**) The complex I activity was measured in cardiomyocytes treated with 200 nM AngII, in the presence or absence of Los, Met, L-NAME or SP. Complex I activity was measured by monitoring the decrease in NADH concentration, normalized to mg of protein and expressed as a percentage of the control group.

attempts to clarify the role of cardiomyocytes in the anti-hypertrophic action of metformin.

Metformin is a well-known AMPK activator, and accumulated data demonstrate that AMPK activation plays a central role not only

in cardioprotection against ischaemia/reperfusion injury [27] but also in induction of anti-hypertrophic/anti-remodelling effects [28]. Currently, there are no studies elucidating the relationship between AMPK and AT1R in cardiomyocyte hypertrophy. Our data demon-

![](_page_9_Figure_0.jpeg)

Fig. 10 Proposed mechanism of the anti-hypertrophic effect of metformin on AngII-induced hypertrophy.

strate that metformin reduced the amount of AT1R in cultured cardiomyocytes; thus, it blocked the hypertrophic effects of Angll. The effect of metformin was attributed to AMPK activation, as the presence of compound C down-regulated AT1R. AMPK may reduce AT1R levels through down-regulation of genes encoding the receptor, followed by P-AMPK translocation to the nucleus. In response to stress, AMPK subunits accumulate in nuclei, independent of kinase activation [29], where the nuclear AMPK $\alpha$ 2 complex regulates gene transcription [30]. In addition, active AMPK may abrogate the effect of AnglI by inhibiting NAD(P)H oxidase, a downstream target of AT1R, as shown in endothelial cells [31]. In our experiments, losartan stimulated AMPK phosphorylation, suggesting that a negative reciprocal relationship may exist between AMPK and AT1R. It should be noted that beneficial effects of metformin may also be mediated through AT2R stimulation, which would counteract AT1R and abrogate AT1R-mediated growth responses in cardiac hypertrophy [32]. Previous studies demonstrate an increase [33], decrease [34] or no change [35] in Angliinduced activation of AMPK in cardiac cells. We observed that AnglI did not remarkably affect AMPK phosphorylation; however, losartan significantly stimulated AMPK phosphorylation by a yet unknown mechanism. Our data are consistent with previous studies, which show that administration of the AT1R blockers, telmisartan and candesartan stimulated AMPK activity in cultured myotubes [36] and rat hypothalamuses [37]. It remains unclear how AT1R inhibition stimulates AMPK. We are tempted to speculate that AT1Rinduced AMPK activation may be mediated through regulation of the upstream molecules, LKB1 and/or CamKK, and internalization and interaction of AT1R with AMPK.

The beneficial effects of AMPK activation and AT1B inhibition on mitochondria may be mediated through the same mechanism(s) or independent pathways. Angll binds to AT1R, activating the transcription factor, p53, and apoptosis. p53 binding to AT1R gene promoter regions is primarily responsible for up-regulating the receptor [38]. We found that losartan prevention of AT1R up-regulation was associated with decreased p53 levels, indicating that the latter mediates the AT1R-induced apoptotic pathway. On the other hand, metformin- or losartan-induced AMPK activation can increase NAD\*/NADH, which, in turn, stimulates the NAD<sup>+</sup>-dependent deacetylases sirtuins, including SIRT1, a regulator of p53 activity. According to previous studies. when metformin induces activation of AMPK and SIRT1, an inverse relationship exists between metformin-triggered AMPK-SIRT1 signalling and p53 protein levels in high glucose-exposed HepG2 cells [39]. Accordingly, in the presence of metformin, AMPK-induced p53 downregulation may prevent AnglI-induced AT1R up-regulation [38].

Phosphorylation of eNOS at Ser1177 may be a downstream target of the AMPK/SIRT1 pathway [40]. Remarkably, both metformin and losartan increased P-eNOS<sup>Ser1177</sup> levels in our studies. Splitomicin, however, prevented eNOS activation, indicating that SIRT1 is an upstream molecule of eNOS. Although eNOS inhibition enhanced AMPK activation in endothelial cells [41], our data revealed that treatment with L-NAME did not appreciably alter P-AMPK levels in cardiomyocytes. This difference can be due to cell specificity, as different cell types were used in these studies. Activation of the AMPK/eNOS pathway and subsequent metformin- or losartan-enhanced nitric oxide production may have various, further cardioprotective effects. In addition to coronary vasodilatation, nitric oxide regulates cell metabolism, including mitochondrial function in cardiomyocytes [42].

A number of studies [6, 43], including ours, [13] demonstrate that metformin therapy attenuates mitochondrial dysfunction during ischaemia/reperfusion and heart failure via yet unknown mechanisms. Consistent with these reports, we found that metformin prevented Angll-induced depolarization of the mitochondrial membrane, an index of mitochondrial function in cardiomyocytes. Likewise, both losartan and pifithrin- $\alpha$  may improve mitochondrial function, indicating that the effects of AMPK activation and AT1R inhibition converge on mitochondria. Various mechanisms can mediate the beneficial effects of AMPK activation on mitochondria. For example, metformin protects the heart through PI3K- [43] and PPAR $\alpha$ - [13] mediated inhibition of the mitochondrial PTP. On the other hand, in response to oxidative stress, p53 translocates to mitochondria and forms the CyP-D-p53 complex, which is associated with mitochondrial PTP opening [24, 25]. In the presence or absence of metformin, we found no difference between control and Angll-treated groups with regard to the interaction of CyP-D with p53, although pifithrin- $\alpha$  abrogated AnglI-induced membrane depolarization. These results indicate that metformininduced activation of the AMPK/SIRT1/eNOS pathway prevents Angll-induced, p53-dependent membrane depolarization and PTP opening in mitochondria, which is not associated with protein-protein interaction between p53 and CyP-D. It appears that further studies are necessary to clarify the role of p53 in hypertrophy-induced mitochondrial dysfunction. At high concentrations (>10 mM), metformin may improve mitochondrial function, independent of AMPK activation via direct inhibition of complex I [44]. However, this study excluded

that possibility, because 2 mM metformin did not inhibit complex I activity (Fig. 9A).

In conclusion, metformin activates AMPK, preventing AT1R up-regulation, and, *vice versa*, losartan down-regulates AT1R, stimulating AMPK. As such, this study demonstrates a negative reciprocal relationship between AMPK activation and AT1R levels in AngII-induced cardiomyocyte hypertrophy. In response to AngII administration, the AMPK/SIRT1/eNOS pathway abrogation of mitochondrial dysfunction mediates, at least partially, the beneficial effects of metformin and losartan. Our studies provide further evidence that AMPK can be targeted to promote cardiomyocyte survival *via* activation of several signalling pathways in hypertrophy and heart failure.

# Acknowledgements

The authors thank Bryan Agostini and Alan Preston for critical reading of the manuscript, and Gabriela Bergollo for her assistance with cell size experiments. This study was supported by the National Heart, Lung, And Blood Institute of the National Institutes of Health through Research Grant SC1HL118669 (to SJ), and in part, by the RCMI National Center for Research Resources NIH Grant G12RR-03051.

# **Conflicts of interest**

The authors confirm that there are no conflicts of interest.

# References

- Hardie DG, Hawley SA, Scott JW. AMP-activated protein kinase–development of the energy sensor concept. J Physiol. 2006; 574: 7–15.
- Arad M, Seidman CE, Seidman JG. AMPactivated protein kinase in the heart: role during health and disease. *Circ Res.* 2007; 100: 474–88.
- Bertrand L, Ginion A, Beauloye C, et al. AMPK activation restores the stimulation of glucose uptake in an *in vitro* model of insulin-resistant cardiomyocytes *via* the activation of protein kinase B. Am J Physiol Heart Circ Physiol. 2006; 291: H239– 50.
- Zhang L, He H, Balschi JA. Metformin and phenformin activate AMP-activated protein kinase in the heart by increasing cytosolic AMP concentration. Am J Physiol Heart Circ Physiol. 2007; 293: H457–66.
- Kirpichnikov D, McFarlane SI, Sowers JR. Metformin: an update. Ann Intern Med. 2002; 137: 25–33.
- Gundewar S, Calvert JW, Jha S, et al. Activation of AMP-activated protein kinase by metformin improves left ventricular function and survival in heart failure. *Circ Res.* 2009; 104: 403–11.
- Sasaki H, Asanuma H, Fujita M, et al. Metformin prevents progression of heart failure in dogs: role of AMP-activated protein kinase. *Circulation*. 2009; 119: 2568– 77.
- Molkentin JD, Dorn GW 2nd. Cytoplasmic signaling pathways that regulate cardiac hypertrophy. *Annu Rev Physiol.* 2001; 63: 391–426.
- Reid IA. Interactions between ANG II, sympathetic nervous system, and baroreceptor reflexes in regulation of blood pressure. *Am J Physiol.* 1992; 262: E763–78.

- Siragy HM. AT(1) and AT(2) receptors in the kidney: role in disease and treatment. Am J Kidney Dis. 2000; 36: S4–9.
- Liu Q, Wang G, Zhou G, et al. Angiotensin II-induced p53-dependent cardiac apoptotic cell death: its prevention by metallothionein. *Toxicol Lett.* 2009; 191: 314–20.
- Dai DF, Rabinovitch P. Mitochondrial oxidative stress mediates induction of autophagy and hypertrophy in angiotensin-II treated mouse hearts. *Autophagy*. 2011; 7: 917–8.
- Barreto-Torres G, Parodi-Rullan R, Javadov S. The role of PPARalpha in metformininduced attenuation of mitochondrial dysfunction in acute cardiac ischemia/reperfusion in rats. *Int J Mol Sci.* 2012; 13: 7694– 709.
- Hawley SA, Gadalla AE, Olsen GS, et al. The antidiabetic drug metformin activates the AMP-activated protein kinase cascade via an adenine nucleotide-independent mechanism. Diabetes. 2002; 51: 2420–5.
- Cook JL, Re RN. Lessons from *in vitro* studies and a related intracellular angiotensin II transgenic mouse model. *Am J Physiol Regul Integr Comp Physiol*. 2012; 302: R482– 93.
- Lei B, Matsuo K, Labinskyy V, et al. Exogenous nitric oxide reduces glucose transporters translocation and lactate production in ischemic myocardium in vivo. Proc Natl Acad Sci USA. 2005; 102: 6966–71.
- Bolli R. Cardioprotective function of inducible nitric oxide synthase and role of nitric oxide in myocardial ischemia and preconditioning: an overview of a decade of research. *J Mol Cell Cardiol.* 2001; 33: 1897–918.
- Jugdutt BI. Nitric oxide in heart failure: friend or foe. *Heart Fail Rev.* 2002; 7: 385–9.
- 19. Sartoretto JL, Kalwa H, Pluth MD, *et al.* Hydrogen peroxide differentially modulates

cardiac myocyte nitric oxide synthesis. *Proc Natl Acad Sci USA*. 2011; 108: 15792–7.

- Fujita T, Ishikawa Y. Apoptosis in heart failure. -The role of the beta-adrenergic receptor-mediated signaling pathway and p53-mediated signaling pathway in the apoptosis of cardiomyocytes. *Circ J.* 2011; 75: 1811–8.
- Grishko V, Pastukh V, Solodushko V, et al. Apoptotic cascade initiated by angiotensin II in neonatal cardiomyocytes: role of DNA damage. Am J Physiol Heart Circ Physiol. 2003; 285: H2364–72.
- Langley E, Pearson M, Faretta M, et al. Human SIR2 deacetylates p53 and antagonizes PML/p53-induced cellular senescence. EMBO J. 2002; 21: 2383–96.
- Vaziri H, Dessain SK, Ng Eaton E, et al. hSIR2(SIRT1) functions as an NAD-dependent p53 deacetylase. *Cell.* 2001; 107: 149–59.
- Wolff S, Erster S, Palacios G, et al. p53's mitochondrial translocation and MOMP action is independent of Puma and Bax and severely disrupts mitochondrial membrane integrity. Cell Res. 2008; 18: 733–44.
- Vaseva AV, Marchenko ND, Ji K, et al. p53 opens the mitochondrial permeability transition pore to trigger necrosis. *Cell.* 2012; 149: 1536–48.
- Yin M, van der Horst IC, van Melle JP, et al. Metformin improves cardiac function in a nondiabetic rat model of post-MI heart failure. Am J Physiol Heart Circ Physiol. 2011; 301: H459–68.
- Russell RR 3rd, Li J, Coven DL, et al. AMPactivated protein kinase mediates ischemic glucose uptake and prevents postischemic cardiac dysfunction, apoptosis, and injury. J Clin Invest. 2004; 114: 495–503.
- Fu YN, Xiao H, Ma XW, et al. Metformin attenuates pressure overload-induced car-

diac hypertrophy *via* AMPK activation. *Acta Pharmacol Sin.* 2011; 32: 879–87.

- Kodiha M, Rassi JG, Brown CM, et al. Localization of AMP kinase is regulated by stress, cell density, and signaling through the MEK-ERK1/2 pathway. Am J Physiol Cell Physiol. 2007; 293: C1427–36.
- Salt I, Celler JW, Hawley SA, et al. AMPactivated protein kinase: greater AMP dependence, and preferential nuclear localization, of complexes containing the alpha2 isoform. Biochem J. 1998; 334: 177–87.
- Schuhmacher S, Foretz M, Knorr M, et al. alpha1AMP-activated protein kinase preserves endothelial function during chronic angiotensin II treatment by limiting Nox2 upregulation. Arterioscler Thromb Vasc Biol. 2011; 31: 560–6.
- van Kesteren CA, van Heugten HA, Lamers JM, et al. Angiotensin II-mediated growth and antigrowth effects in cultured neonatal rat cardiac myocytes and fibroblasts. J Mol Cell Cardiol. 1997; 29: 2147–57.
- Nagata D, Takeda R, Sata M, et al. AMPactivated protein kinase inhibits angiotensin II-stimulated vascular smooth muscle cell proliferation. *Circulation*. 2004; 110: 444– 51.

- Stuck BJ, Lenski M, Bohm M, et al. Metabolic switch and hypertrophy of cardiomyocytes following treatment with angiotensin II are prevented by AMP-activated protein kinase. J Biol Chem. 2008; 283: 32562–9.
- Jiang SY, Xu M, Ma XW, et al. A distinct AMP-activated protein kinase phosphorylation site characterizes cardiac hypertrophy induced by L-thyroxine and angiotensin II. *Clin Exp Pharm Physiol.* 2010; 37: 919–25.
- Feng X, Luo Z, Ma L, et al. Angiotensin II receptor blocker telmisartan enhances running endurance of skeletal muscle through activation of the PPAR-delta/AMPK pathway. J Cell Mol Med. 2011; 15: 1572–81.
- Marques MB, Ribeiro-Oliveira A Jr, Guimaraes J, et al. Modifications in basal and stress-induced hypothalamic AMP-activated protein kinase (AMPK) activity in rats chronically treated with an angiotensin II receptor blocker. Stress. 2012; 15: 554–61.
- Pierzchalski P, Reiss K, Cheng W, et al. p53 Induces myocyte apoptosis via the activation of the renin-angiotensin system. Exp Cell Res. 1997; 234: 57–65.
- Nelson LE, Valentine RJ, Cacicedo JM, et al. A novel inverse relationship between metformin-triggered AMPK-SIRT1 signaling

and p53 protein abundance in high glucoseexposed HepG2 cells. *Am J Physiol Cell Physiol*. 2012; 303: C4–13.

- Ota H, Eto M, Ogawa S, et al. SIRT1/eNOS axis as a potential target against vascular senescence, dysfunction and atherosclerosis. J Atheroscler Thromb. 2010; 17: 431–5.
- Jin BY, Sartoretto JL, Gladyshev VN, et al. Endothelial nitric oxide synthase negatively regulates hydrogen peroxide-stimulated AMP-activated protein kinase in endothelial cells. Proc Natl Acad Sci USA. 2009; 106: 17343–8.
- Rakhit RD, Marber MS. Nitric oxide: an emerging role in cardioprotection? *Heart*. 2001; 86: 368–72.
- Bhamra GS, Hausenloy DJ, Davidson SM, et al. Metformin protects the ischemic heart by the Akt-mediated inhibition of mitochondrial permeability transition pore opening. Basic Res Cardiol. 2008; 103: 274–84.
- El-Mir MY, Nogueira V, Fontaine E, et al. Dimethylbiguanide inhibits cell respiration via an indirect effect targeted on the respiratory chain complex I. J Biol Chem. 2000; 275: 223–8.