

Immunoblotting protocol is included under Supplementary Information section.

Please note: SDS PAGE analyses of lysates used for immunoblotting utilized prestained MW markers (Precision Plus Protein Dual Color Standards, #1610374, from BIO-RAD).

Following transfer, nitrocellulose membranes were cropped according to migration of prestained MW markers.

Fig 2A



Fig 2E




Fig 2F




Fig 3A




Fig4A






Fig4F




Fig5F




Fig5G




Fig5H







Fig8A




FigS5A




FigS5E



FigS6B



FigS6C



FigS7C

