How to read biparametric MRI in men with a clinical suspicious of prostate cancer: Pictorial review for beginners with public access to imaging, clinical and histopathological database

Acta Radiologica Open 10(11) 1–22 © The Author(s) 2021 Article reuse guidelines: sagepub.com/journals-permissions DOI: 10.1177/20584601211060707 journals.sagepub.com/home/arr (\$SAGE

Ivan Jambor^{1,2}, Alberto Martini³, Ugo G Falagario⁴, Otto Ettala⁵, Pekka Taimen⁶, Juha Knaapila⁵, Kari T Syvänen⁵, Aida Steiner^{1,7}, Janne Verho^{1,7}, Ileana M Perez^{1,8}, Harri Merisaari^{1,8}, Paula Vainio⁶, Tarja Lamminen⁴, Jani Saunavaara^{1,9}, Giuseppe Carrieri⁴, Peter J Boström⁵ and Hannu J Aronen^{1,3}

Abstract

Prostate Magnetic Resonance Imaging (MRI) is increasingly being used in men with a clinical suspicion of prostate cancer (PCa). Performing prostate MRI without the use of an intravenous contrast (IV) agent in men with a clinical suspicion of PCa can lead to reduced MRI scan time. Enabling a large array of different medical providers (from mid-level to specialized radiologists) to evaluate and potentially report prostate MRI in men with a clinical suspicion of PCa. The aim of this pictorial review is to provide an insight into acquisition, quality control and reporting of prostate MRI performed without IV contrast agent in men with a clinical suspicion of PCa, aimed specifically at radiologists starting reporting prostate MRI, urologists, urology/radiology residents and mid-level medical providers without experience in reporting prostate MRI. Free public access (http://petiv.utu.fi/improd/and http://petiv.utu.fi/multiimprod/) to complete datasets of 161 and 338 men is provided. The imaging datasets are accompanied by clinical, laboratory and histopathological findings. Several topics are simplified in order to provide a solid base for the development of skills needed for an unsupervised review and potential reporting of prostate MRI in men with a clinical suspicion of PCa. The aire of PCa. The current review represents the first step towards enabling a large array of different medical providers to review and report accurately prostate MRI performed without IV contrast agent in men with a clinical suspicion of PCa.

Corresponding author:

Ivan Jambor, Department of Diagnostic Radiology, University of Turku, Kiinamyllynkatu 4-8, Turku 20521, Finland. Email: ivjamb@utu.fi



Creative Commons Non Commercial CC BY-NC: This article is distributed under the terms of the Creative Commons Attribution-NonCommercial 4.0 License (https://creativecommons.org/licenses/by-nc/4.0/) which permits non-commercial use,

reproduction and distribution of the work without further permission provided the original work is attributed as specified on the SAGE and Open Access pages (https://us.sagepub.com/en-us/nam/open-access-at-sage).

¹Department of Diagnostic Radiology, University of Turku, Turku, Finland

²Department of Radiology, Icahn School of Medicine at Mount Sinai, New York, NY, USA

³Department of Oncology/Unit of Urology, Urological Research Institute, IRCCS Ospedale San Raffaele, Milan, Italy

⁴Department of Urology and Organ Transplantation, University of Foggia, Foggia, Italy

⁵Department of Urology, University of Turku and Turku University Hospital, Turku, Finland

⁶Institute of Biomedicine, University of Turku and Department of Pathology, Turku University Hospital, Turku, Finland

⁷Medical Imaging Centre of Southwest Finland, Turku University Hospital, Turku, Finland

⁸Turku Brain and Mind Center, University of Turku, Turku, Finland

⁹Department of Medical Physics, Turku University Hospital, Turku, Finland

Keywords

IMPROD biparametric magnetic resonance imaging, diffusion weighted imaging, prostate cancer, prostate specific antigen, multi-institutional trial, prostate cancer screening, review

Received 23 May 2021; Accepted 1 November 2021

Introduction

Prostate cancer (PCa) is characterized by a spectrum of aggressiveness that ranges from minimal clinical significance to highly aggressive disease.¹ Biopsy procedures as well as the treatments of prostate cancer are not devoid of short- and long-term complications and adverse effects.^{2,3} Thus, optimal risk stratification of men with a clinical suspicion of PCa is of utmost importance to limit unnecessary biopsies and treatment-related complications while appropriately detecting clinically significant PCa (csPCa), for example a commonly used definition of Gleason score >3+3/Gleason Grade Group >1. After decades of insufficient ability of prostate specific antigen (PSA) and systematic transrectal ultrasound (TRUS) biopsy for detecting csPCa, prostate Magnetic Resonance Imaging (MRI) and prostate MRI targeted biopsy are currently increasingly being used in men with a clinical suspicion of prostate cancer (PCa). Both AUA and EAU guidelines recommend the use of prostate MRI and targeted biopsy in biopsy-naive men as well as in men with prior negative prostate biopsy with a clinical suspicion of prostate cancer (PCa) (https://uroweb.org/guidelines/; https://www.auanet. org/guidelines).

Multiple trials⁴ (PROMIS, 05/2012–10/2015⁵; PRECI-SION, 02/2016-08/2017⁶ and MRI-FIRST, 07/2015-08/ 2016⁷) demonstrated that multiparametric MRI performed using intravenous contrast agents as a triage test could reduce unnecessary biopsies, over-detection of insignificant PCa and increase detection of csPCa compared to the standard systematic biopsy pathway. However, adopting multiparametric MRI as triage pathway for all men with a clinical suspicion of PCa will require a lot of imaging time and radiology resources. The associated costs and ensuring MRI availability to all men with a clinical suspicion of PCa are major issues of concern. Several single (IMPROD, 03/2013-02/2015⁸; BIDOC, 10/2015-07/2017^{9,10}) and multicenter (MULTI-IMPROD, 02/2015-05/2017¹¹; STHLM3-MR trial 05/2016-05/2017,¹²) trials as well retrospective studies¹³ demonstrated similar findings for prostate MRI performed without an intravenous contrast agent.

Currently, prostate MRI in routine clinical practice is reported by radiologists. However, the use of an optimized prostate MRI acquisition and reconstruction protocol at 1.5 Tesla (T) or 3 T with continuous quality control and performance measures could hypothetically lead to a wider array of prostate MRI readers and ultimately lead to reduction of cost associated with wide-spread adoption of prostate MRI in men with a clinical suspicion of PCa. Enabling urologists and mid-level medical providers such as dedicated prostate MRI trained physician assistants to report prostate MRI could potentially lead to a cost reduction and population-based application of prostate MRI in men with a clinical suspicion of PCa. However, this is purely hypothetically and a robust evidence to support cost effectiveness of prostate MRI reported by and MRI targeted biopsy performed by dedicated prostate MRI trained physician assistants remains to be explored.

The aim of this review assay is to provide an insight into the reporting of prostate MRI performed without an intravenous contrast agent in men with a clinical suspicion of PCa. Specifically, this pictorial review article with free public access to full prostate MRI datasets is aimed radiologist starting to report prostate MRI, radiology residents, urologists, urology residents and mid-level medical providers such as physician assistants without experience in reporting prostate MRI. Several topics will be simplified and geared towards this audience in order to provide a solid base for development of skills needed for an evaluation and reporting of prostate MRI in men with a clinical suspicion of PCa. After reading this pictorial assay, the readers should gain knowledge in acquisition, quality control and reporting of prostate MRI performed without an intravenous contrast agent in men with a clinical suspicion of PCa.

Where to start: access to prostate MRI acquisition protocol and imaging datasets

Application of prostate MRI in the diagnostic pipeline for csPCa starts with high quality prostate MRI acquisition and reconstruction protocol. Using high quality prostate MRI is of utmost importance especially for less experienced readers who are the target audience of this pictorial review. Current clinical MRI scanners provide very homogeneous magnetic field with high gradient performance. Thus, end-users can obtain high quality prostate MRI with a wide range of different MRI scanners from different MRI vendors. The current review article provides access to IMPROD bpMRI protocol which was an independent and a parallel development to the European Society of Urogenital Radiology and PI-RADs guidelines since the development put strong emphasis on MRI physics.14-17 The IMPROD bpMRI (Figure 1) protocol was developed in 2012-2013 before the start of the IMPROD trial (NCT01864135) and free public access to the imaging datasets and protocol is provided at http://petiv.utu.fi/improd/and http://petiv.utu.fi/ multiimprod/. These study servers provide centres with little experience with prostate MRI a unique opportunity to import with minimal modifications a standardized fast acquisition protocol of prostate MRI performed without the use of intravenous contrast agents,¹⁸ developed and validated in prospective clinical trial settings. Additionally, IMPROD and MULTI-IMPROD servers (http://petiv.utu.fi/improd/and http://petiv. utu.fi/multiimprod/) are freely publicly accessible and offer the unique opportunity to function as a learning place with further standardization of MRI reporting systems. By the use of complete MRI datasets with corresponding structured MRI reports stored in the server and corresponding histopathology findings (including scanned whole mount prostatectomy section), prostate MRI readers of various experience could learn prostate MRI reporting, thus enabling easier adaptation for everyday clinical use. Acquisition protocols for prostate MRI of all three major MRI vendors are provided in supporting material enabling readers to implement those protocols at home institutions.

Basics of MRI physics with emphasis on prostate MRI

Magnetic resonance imaging uses radiofrequency energy (pulse) emitted by the RF coil to generate a magnetic field with subsequent excitation and relaxation processes. Details are well beyond the scope of the current review. Briefly, 'most common' (hydrogen proton) MRI works mainly because of hydrogen protons, which are placed in an external static magnetic field, are subjected to a turning force (torque), causing them to generate their own small magnetic field, or magnetic moment. Multiple different tissue contrasts can be generated. Most commonly used for prostate images are: T2-weighted images and T1-weighted images (Figures 2 and 3).

T2-weighted images and T1-weighted images: Anatomic magnetic resonance imaging of the prostate

T2-weighted images enable to depict prostate zonal anatomy (Figure 2). The normal peripheral zone has homogeneous high signal intensity on T2-weighted images in older men at least partly due to its relatively high-water content. However, normal peripheral zone in younger men can have diffuse intermediate to low signal intensity on T2-weighted images so it is important not to confuse it with a prostate cancer. The central portion of the prostate is called the central gland and consists of the central zone and transition zone; both have a lower signal intensity compared to the peripheral zone and a heterogeneous signal intensity at least partly due to its stroma-rich composition (Figure 2). Central gland, consisting of central zone and transition zone, is separated from the peripheral zone by a hypointense 'pseudocapsule'. The neurovascular bundles are predominantly located posterolaterally to the 'true' outer prostate capsule at about the 5- and 7-o'clock positions on an axial images. The seminal vesicles are hyperintense on T2weighted images and appear as tubular structures visible posteriorly to the base of the prostate, degree of distension varies based on amount of stored fluid. However, seminal vesicles can be hypointense on T2-weighted if they contain, for example, blood products.

No zonal anatomy can usually be seen on T1-weighted images. T1-weighted images are useful for the detection of




Figure 1. Teaching point - Example of all sequences acquired using a routine mid-to lower end 1.5 Tesla MRI scanner in scan time under 15 min (without shimming and calibration): IMPROD bpMRI consists of localizer, sagittal T2-weighted imaging (T2W Sag), axial T2-weighted imaging (T2W Axial) and three separate Diffusion Weighted Imaging (DWI5b500, DWI2b1500, DWI2b2000) acquisitions. Shown images represent imaging data sets of IMPROD study number 99, all images available at http://petiv.utu.fi/improd/.




Figure 2. Teaching point - Normal zonal anatomy: Axial T2-weighted, sagittal T2-weighted and transversal T1-weighted images of a 67year-old man with a PSA value of 4.1 ng/mL. In this cross-section the normal PZ has a high signal intensity and the central gland has (central and transition zone) an intermixed signal intensity with presence of benign prostatic hyperplasia. No zonal anatomy is visible on the T1-weighted image (C). Axial images were acquired with no image tilt.




Figure 3. *Teaching point* - Prostate cancer in peripheral zone: Transversal and sagittal T2-weighted images from the same 67-year-old man as in Figure 1. In this cross-section a low-signal-intensity lesion, indicated by the arrow, is present in the left peripheral zone at the level of the apex. The lesion was confirmed to be Gleason score 3+4 prostate cancer at targeted biopsy.

blood products which can appear hyperintense in T1-weighted images depending on the type of blood products. T1-weighted images should be obtained if a patient has history of recent prostate biopsy. In general, prostate MRI should not be performed within 8–12 weeks following prostate biopsy to allow sufficient time for resolution of blood products related to prostate biopsy, although large variations in presence and resolution of post-biopsy blood products in individual patients have been reported.^{19–21} In the presence of blood products, prostate cancer can appear as hypointense area what some authors call as 'haemorrhage exclusion sign' which is defined as presence of a well-defined area of low signal intensity surrounded by areas of high signal intensity on T1-weighted images.

Diffusion weighted imaging

Diffusion behaviour, at the tissue level, can be measured using MRI and special gradient acquisition schemes; this method is called Diffusion weighted imaging (DWI). Essentially all clinical MRI scanners can acquire DWI. The strength of diffusion weighting is controlled by the so-called b-value (Figure 4). In simple terms, DWI consists of two parts (Figure 5): A. gradient scheme which controls bvalues and B. read-out method which measures signal following radiofrequency excitation. Prostate DWI is most commonly acquired using a so-called echo planar (epi) read-out. This is a very fast approach but epi read-out is susceptible to artifacts caused for example by A. main magnetic field inhomogeneities which can be caused by air in the rectum (Figure 6) or hip prosthesis, B. insufficient fat suppression (Figure 7). These limitations can be overcome by different methods which are outside of scope of this review.²² It is important to stress that, DWI acquisition methods which are more robust for main magnetic field inhomogeneities (e.g. caused by air in rectum), takes significantly more time to be acquired (Figure 8) and some




Figure 4. *Teaching point* - Each individual b-values in prostate DWI can provide different tissue contrast and signal-to-noise ratio, even images with the same b-values can look totally different, thus, attention to details of prostate MRI acquisition protocol is needed: Multiple different Diffusion Weighted Imaging b-values with corresponding whole mount prostatectomy section, axial T2-weighted images and axial T2-relaxation values maps of a patients with Gleason score 3+4 prostate cancer in the left transition zone.

of these methods are not available at older scanner hardware/software versions. Thus, the simplest approach is to ensure patients do not have rectal gas. A sample instruction on patient preparation before prostate MRI is provided at IMPROD and MULTI-IMPROD servers (http:// petiv.utu.fi/improd/and http://petiv.utu.fi/multiimprod/). Thus, it is important to limit rectal gas if prostate DWI is performed using epi read-out such as DWIs in IMPROD bpMRI protocol. Diffusion weighted imaging is cornerstone of prostate MRI and obtaining high quality DWI is of utmost importance. If a lot of rectal gas is present, DWI performed using epi read-out will be non-diagnostic (Figure 9). Please note, b-value is just one of many parameters affecting tissue contrast between PCa and normal tissues. Same b-values could be acquired with many different sets of parameters and lead to very different contrast between PCa and normal tissues (Figure 10). The general idea of higher the b-value, higher the contrast between PCa and normal tissues is not always valid since trace DWI can reach level of background noise and contrast between PCa and normal tissues is affected by many parameters other than b-value (Figures 10 and 11). We urge readers not to fall in a trap of oversimplification of prostate DWI and believe that prostate DWI is a measure of prostate cancer cellularity. In fact, signal decay of prostate DWI is affected by a large array of parameters and interested readers are referred to simplified overview of prostate DWI by Roger Bourne and Eleftheria Panagiotaki.²

Each b-value of prostate DWI is collected individually and called trace DWI image (Figures 1 and 4). These trace DWI images are designated by their b-value which is commonly expressed in units of s/mm². There are many different ways how these individual b-values can be postprocessed. Details of those methods are we well beyond the scope of this manuscript and interested readers are referred elsewhere.^{16,23} The simplest and most commonly used method to post-process trace DWI images is monoexponential function which leads to generation of so-called Apparent Diffusion Confidents (ADC or ADCm) map (Figures 1 and 12).²⁴ Some authors prefer to use ADCm to specifically denote the Apparent Diffusion Coefficient maps calculated using monoexponential function since Apparent Diffusion Coefficient maps can be calculated using different mathematical function, for example, kurtosis (ADCk).

How to view prostate MRI images

Prostate MRI images are commonly collected in Digital Imaging and Communications in Medicine (DICOM) file format and stored in a Picture Archiving and Communication System (PACS). There are a lot of different clinically used DICOM viewers as well as a large arrays of free DICOM viewers. Weasis (used at IMPROD and MULTI-IMPROD server) is one of many options for Windows users while Osirix for Mac OS users.

How to quality control prostate biparametric MRI: minimal requirements

Quality control of prostate MRI images is of pivotal importance and any radiologist, radiology resident, urologist, urology resident, mid-level medical providers should be able to understand differences between good quality MRI and bad quality MRI especially when performed within and outside




Figure 5. Teaching point - Quality of DWI images is impacting by both gradient scheme as well as read-out: Representative Diffusion Weighted Imaging acquisition scheme.

her/his own institution. Unfortunately, PIRADs recommendation provides only a very limited general suggestion regarding MRI acquisition protocols and image quality can vary across different institutions. Quality control of prostate biparametric MRI could be performed in multiple different ways, a sample example is detailed below:

A. Make sure all MRI sequences were collected (Figure 1). A complete study should have the following:

1. T2-weighted images (axial and sagittal), 2. DWI, if IMPROD bpMRI protocol was used, then the following three DWI acquisitions: 2a. DWI with 5 b-values (0–500 s/mm²), 2b. DWI with 2 b-values (0–1500 s/mm²), 2c. DWI with 2 b-values (0–2000 s/mm²). If patient had prior history of prostate biopsy, a fast T1-weighted imaging should be performed as well. Fast T1-weighted imaging performed using volume interpolated gradient echo can be added to




Figure 6. *Teaching point* - Degree of prostate distortion due to susceptibility depends on the type of DWI acquisition protocol: Degree of susceptive artifacts depends on acquisition parameters: In this example, Diffusion Weighted Imaging of prostate was performed using echo planar (epi) read-out and a turbo spin echo (TSE) read-out. Due to rectal gas and resulting main magnetic field inhomogeneities, images acquired using epi read-out are severely distorted in contrast to images collected using TSE. Yellow arrow points to sharply delineated peripheral zone with a small müllerian duct cyst which cannot be seen on images acquired using epi read-out.




Figure 7. Teaching point - Degree of imaging artefacts, such as fat shift, depends on the type of DWI acquisition protocol: Example of different fat suppression methods for prostate Diffusion Weighted Imaging: In the middle row (Sh SE-EPI STIR) bright line in the middle of images is due to insufficient fat suppression.

routine IMPROD bpMRI protocol (Figure 1) to enable detection of blood products. T1-weighted sequences can be acquired in under 10 s covering whole prostate.²⁵ Thus, this is well-invested time in men with prior history of prostate biopsy. B. Review all individual trace DWI images. Routinely clinical MRI scanners calculate ADC/ADCm map automatically but a person reading prostate MRI cannot blindly rely on ADC/ADCm maps. ADC/ADCm maps are as good




Figure 8. *Teaching point* - Acquisition time for different DWI acquisition protocols can be very different: Comparison of Diffusion Weighted Imaging performed using a commonly used single-shot spin echo planar read-out (Single-shot SE-EPI) versus multi-shot spin echo planar read-out (Multi-shot SE-EPI) which is more robust for B0 (main magnetic field) variation, leading to less susceptibility artefacts, but can tends have longer acquisition time.




Figure 9. *Teaching point* - Impact of rectal air on prostate DWI performed using echo planar read-out, rectal air related susceptibility artefacts can limit prostate cancer detection: Due to rectal air (part F), trace DWI of b-value 2000 s/mm² acquired using echo planar read-out is distorted and these distortions are obscuring Gleason score 4+3 prostate cancer in the left peripheral zone (encircled in green on whole mount prostatectomy section – part I). Sagittal T2-weighted image collected during the same acquisition as the part F trace DWI of b-value 2000 s/mm² demonstrate rectal air in contrast to sagittal T2-weighted image (part B) collected as part of IMPROD bpMRI protocol (parts B, C, D, H). Images in parts A, F, G and E were collected using on a different day and using a different MRI scanner than images in parts B, C, D, and H. Following collection of images in part F, patient was asked to go on toilette and new DWI images were collected (part G) which did clearly show PCa in the left peripheral zone without presence of distortions.

as individual trace DWI images. Thus, the person reading prostate MRI needs to ensure the following: 1. Each individual trace DWI image has sufficient signal-to-noise (e.g. part C in Figure 10 has too low signal-to-noise), 2. Compare prostate shape on each axial individual trace DWI image with axial T2-weighted images – if axial trace DWI images are distorted it suggests the presence of susceptibility artifacts if DWI was performed using an echo planar read-out, for example, part F of Figure 9 demonstrates rectal distension with distortions of prostate shape. If individual trace DWI images are of low quality, calculated ADC/ADCm maps will be of low quality as well. Since ADC/ADCm maps calculation commonly takes into account all acquired trace DWI data, it is important to review each individual trace DWI series to ensure good data quality.

How to read prostate biparametric MRI

Each readers of prostate MRI can create her/his own systems and search pattern. Below is presented a standard way to view prostate MRI images:

- 1. Create 2×2 image window matrix
- 2. Perform quality control of data as stated in 'How to quality control prostate biparametric MRI: minimal requirements'.
- If T1-weighted images were collected, review those images first to ensure there are no haemorrhage related products which, depending on timing, would appear as T1 hyperintense areas.
- Once T1-weighted images have been reviewed, place axial T2-weighted images in the upper left-hand corner, ADC/ADCm of first IMPROD bpMRI DWI,




Figure 10. *Teaching point* - Same b-values can provide different tissue contrast and signal to noise, thus, end user should be familiar with NMR physics and/or use protocol specifically developed for prostate MRI: Multiple different trace DWI images of two prostate cancer lesions (right peripheral zone Gleason score 3+4 and left peripheral zone 3+3 - whole mount prostatectomy section, part E). Trace DWI image in part C had the same b-value as image in part D, however, tissues contrast is different. The right peripheral zone Gleason score 3+4 lesion can be seen on part B, F and G while cannot be seen in part C (due low signal-to-noise ratio), part D (due to low prostate cancer to normal tissue contrast) and part H (due to low prostate cancer to normal tissue contrast). The small left peripheral zone Gleason score 3 + 3 lesion was not detected on any DWI trace images.




Figure 11. Teaching point - Different appearance of prostate cancer depending on MRI acquisition protocol: Images in parts A, B, C, D and E were collected using a different MRI protocol and scanner than images in parts F, G, H and I. As can be seen, tissues contrast between prostate cancer (Gleason score 3+4 prostate cancer encircled on whole mount images, part J) is different even though the same b-values were collected. However, all images were collected at 3 T. Similar relative image windowing was applied between corresponding DWI trace images. The shown images in parts F, G, H and I are imaging data sets of IMPROD study number 90 and all MRI images available at http://petiv.utu.fi/improd/.

DWI5b500, in the upper right-hand corner and in the lower corners place the trace DWI images of b-values 1500 and 2000 s/mm², respectively (Figure 13). Some DICOM Viewer use image sync between series, when reviewing images first, disable this feature since prostate can move in cranio-caudal direction and the use of image sync between series could be misleading.

 Review DWI images of b-values 1500 and 2000 s/ mm² first. All high signal areas in these images need to be explained – there are four basic options: 1. PCa, 2. Artifact, 3. Benign Prostatic Hyperplasia, 4. Other rare benign conditions such as granulomatous prostatitis especially following *Bacillus* Calmette-Guérin (BCG) treatment.

 Differentiate peripheral zone from central gland (transition zone and central zone). A simplified scheme is shown in Figure 14. Normal peripheral zone (PZ) has high signal intensity on T2-weighed images (Figures 2, 3 and 15) while central gland (transition zone and




Figure 12. *Teaching point* - Typical imaging appearance of prostate cancer in both peripheral zone and anterior fibromuscular stroma: Axial T2-weighted image (part A), trace DWI images (parts B, C, D) and corresponding Apparent Diffusion Confidents (ADC/ADCm) maps (parts F, G, H). The shown images represent imaging data sets of IMPROD study number 99. All MRI images are available at http:// petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 99 are presented in supporting material. Please note, an example of viewing hanging protocol for IMPROD bpMRI is shown in Figures 15–20 and 23–28.




Figure 13. *Teaching point* - Typical imaging appearance of prostate cancer in peripheral zone: IMPROD study number 3: Focal area of low signal on T2-weighted (T2W Axial), low signal intensity of apparent diffusion coefficient maps (ADCm5b500) and high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm² was present in the right peripheral zone corresponding to Gleason score 3+4 prostate cancer. Prostatectomy was performed with wide margins on the right side and negative surgical margins were found on whole mount prostatectomy sections. All MRI images available at: http://petiv.utu.fi/ improd/ Additional images and videos of IMPROD study number 3 are presented in supporting material.

central zone) can have heterogeneous appearance related to BPH nodules (Figure 15).

7. Use/develop own and/or standardized reporting systems. It is important to communicate with all stakeholders in the diagnostic pipeline of csPCa to ensure that reporting system is well understood and provides the best diagnostic performance. One option is to use PIRADs version 2.1 systems in addition to other systems. In the current review, all findings are reported using PIRADs version 2.1 system as well as IMPROD bpMRI Likert scoring systems.²⁶

It is crucial to develop a standardized way to communicate imaging findings in a reporting system which can be easily interpreted by all stake-holders, including men undergoing prostate MRI. This is especially important if a person reading MRI scan is different from a person performing biopsy. Moreover, creating and maintaining a database of imaging, clinical and histopathological findings is an important part of quality-control/improvement process. As an example of a publicly available database, readers are referred to a PRO-MOD database.²⁷

Reporting findings in the peripheral zone

Details on peripheral zone findings using PIRADs version 2.1 are presented elsewhere. Two recent studies compared different prostate MRI reporting systems and did not find any major differences between prostate MRI studies reported using different reporting systems (PIRADsv2.0 versus Likert reported using clinical information, PI-RADsv2.1 vs IMPROD bpMRI Likert) when applied by the same readers.^{26,28} The major difference in using IMPROD bpMRI Likert scoring systems is integrating all available information rather than relying on 'dominant sequences' as per PIRADs version 2.1. In studies acquired using IM-PROD bpMRI protocol, there was rather small difference between PIRADs version 2.1 and IMPROD bpMRI Likert scoring system.²⁶ Specifically, in 161 patients who took part in the IMPROD trial, the decision to perform systematic biopsy or systematic targeted biopsy would be different depending on reporting system only in 4 men (2%, 4/161). All 4 cases are presented in this 15-tweet series: https://twitter. com/jambor ivan/status/1185272940671180800

Either extremes of reporting spectrum are rather easy to detect – score 1: no suspicion areas (Figure 15); score 4 or 5: Prostate cancer appears as focal area of decreased signal on T2-weighed images with decreased signal on ADC/ADCm images and increased signal on trace DWI b-values acquired with strong diffusion weighting, such as 1500 and 2000 s/ mm² (Figures 13 and 16). However, abscess and blood products (which could be hyperintense on T1-weighted images depending on the timing of biopsy/hemorrhage)




Figure 14. *Teaching point* - It is important to communicate imaging finding with referring physicians and a simplified prostate imaging reporting scheme can be of help: Simplified prostate zonal anatomy. Peripheral zone is marked as hashed light grey, central zone as light grey and transition zone of dark grey.

could also appears as decreased signal on ADC/ADCm images and increased signal on trace DWI b-values acquired with strong diffusion weighting, such as 1500 and 2000 s/mm². Please note, the b-value is just one of many parameters affecting contrast between tissues, thus, images with the same b-values can have very different tissue contrast.

The presence of prostatitis or prior episodes of prostatitis and/or biopsy can lead to diffuse decrease signal on T2-weighted images without a focal decreased signal on ADC/ADCm images and no focal increased signal bvalues acquired with strong diffusion weighting, in IM-PROD bpMRI protocol those trace DWI images are bvalues of 1500 and 2000 s/mm² (Figure 17). These areas could be reported as PIRADs version 2.1 and/or IMPROD bpMRI Likert scoring systems score of one or 2, implying low probability of SPCa and routinely does not require biopsy. When in doubts regarding the nature of the lesion, follow up MRI scans in 6-12-24 months, could be considered taking into account lesion appearance and location on MRI




Figure 15. *Teaching point* - Typical BPH: IMPROD trial number 167: Normal peripheral zone and normal central gland with multiple nodules of benign prostatic hyperplasia were present on T2-weighted images (T2W Axial), apparent diffusion coefficient maps (ADCm5b500) and trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm², in the central zone. All central gland (transitional and central zone) lesions have a well-defined capsule. No prostate cancer was found in any cores of systematic biopsy. All MRI images are available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 167 are presented in the supporting material.




Figure 16. Teaching point - Typical prostate cancer in peripheral zone: IMPROD study number 45: One lesion with low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500) and high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm² was present which corresponded to Gleason score 4+3+5 prostate cancer with extra-capsular extension on the right-side. All MRI images are available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 45 are presented in the supporting material.

(possibility of sampling errors) and patients' risk factors (PSA density [dPSA], family history, ...) and patient's preference (Figure 18). Different methods for measuring prostate gland has been proposed and we recommend readers to develop and validate own methods and/or carefully validate proposed methods rather than blindly relying on measurements methods as proposed in PIRADs version 2.1 document.

Linear areas in the peripheral zone crossing the peripheral zone with decreased signal on T2-weighed




Figure 17. *Teaching point* - Typical benign prostatic hyperplasia with non-specific signal changes in peripheral zone: IMPROD study number 59: Bilateral areas of decreased T2-weighted imaging signal (T2W Axial) without a focal decreased signal on apparent diffusion coefficient maps (ADCm5b500) and without focal increased signal on b-trace diffusion weighted imaging b-values of 1500 (DWI b 1500) and 2000 (DWI b 2000) s/mm², was present in peripheral zone. No prostate cancer or prostatic intraepithelial neoplasia was found in any cores of systematic biopsy. All MRI images available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 59 are presented in the supporting material.




Figure 18. *Teaching point* - Resolution of imaging finding in PZ in follow up MRI: IMPROD study number 22: There was one focal irregular lesion in the apex, right side, with size of 10x8x8 mm, seen as low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500) and high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm². Systematic and targeted biopsy did not demonstrate any prostate cancer and on follow up imaging 386 days from the initial IMPROD bpMRI scan, the lesion was not apparent anymore. All images of the initial and follow up MRI scans are available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 22 are presented in the supporting material.

images, decreased signal on ADC/ADCm images, and increased signal on DWI images acquired with strong diffusion weighting, for example, b-value of 500 and 2000 s/mm², usually represent benign findings (Figure 19), however, csPCa cannot be totally excluded in such cases (Figure 20). In a retrospective study by Jorge Abreu-Gomez e al., up to 13% of PCa lesions in the peripheral zone were linear or wedge-shaped.²⁹ However, proportion of lesions with linear or wedge shape in the peripheral zone which represents csPCa depends on many factors such as study population and reference standard. Thus, rather than quoting a specific number for a possibility of a linear or wedge-shape lesion with typical PCa signal characteristics (decreased T2weighted signal and associated decreased signal on ADC/ ADCm images and increased signal b-values acquired with strong diffusion weighting, 1500 and 2000 s/mm²) to be a csPCa, an individualized approached is needed. Such individualized approached could represent follow up imaging




Figure 19. *Teaching point* - Linear signal changes in peripheral zone which could be both benign or prostate cancer: IMPROD study number 15: There was one focal irregular lesion in the mid-gland to apex, left side, with size of 8x9x10 mm, seen as low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500) and high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 200) s/mm². Systematic and targeted biopsy did not reveal any PCa and patient has stable PSA. Following risk-decision discussion with the patient, PSA is being follow which is stable to decreasing for over 2 years. All MRI images available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 15 are presented in supporting material.




Figure 20. *Teaching point* - Linear signal changes in peripheral zone which could be both benign or prostate cancer: IMPROD study number 133: Focal wedge-shaped lesion in right peripheral zone, seen as low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500), high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm², in the peripheral zone of mid-gland which was on found to be, in systematic biopsy and subsequent prostatectomy, Gleason score 3+4 prostate cancer. All MRI images are available at http://petiv.utu.fi/improd/Additional images and videos of IMPROD study number 133 are presented in supporting material.

(Figure 18), repeated biopsy (Figure 20) or watchful waiting (Figure 19). Thus, if the aim of finding all csPCa is desired in a certain man, repeated biopsy with follow imaging in recommended.

Extension of PCa beyond the margins of the prostatic capsule (Figures 16 and 21) can be suggested if 'loss of prostate capsule', 'bulging with irregular or spiculated margin', 'obliteration of the rectoprostatic angle',

'asymmetry or invasion of the neurovascular bundles' are present (Figures 16 and 21). Presence of seminal vesicle invasion can be suggested when an area of abnormal low T2-weighted imaging signal intensity is extending to the seminal vesicle which typically have high signal intensity on T2-weighted images (Figure 22).

In summary, reporting lesions that are highly suspicious of harbouring csPCa (PIRADs version 2.1/IMPROD bpMRI Likert score 4–5) in the peripheral zone is straightforward since these lesions are focal in nature with typical PCa signal characteristics (decreased T2-weighted signal and associated decreased signal on ADC/ADCm images and increased signal b-values acquired with strong diffusion weighting, 1500 and 2000 s/mm²), Figures 3, 13 and 16. The lesions which are linear or wedge shaped with the typical PCa signal characteristics partly present a diagnostic challenge and individualized patient management should be considered (Figures 18, 19 and 20). Reporting extra-prostatic extension (Figure 21) and seminal vesicle invasion (Figure 22) depends on MRI acquisition protocol and used reference standard in the training process. Thus, readers should refer to their own institutional data and do not reply on reported performance from other centres which can be misleading and detrimental to patient's management.

Reporting findings in the central gland

The central zone, part of central gland in addition to transition zone, typically contains BPH nodules (Figures 2, 15 and 23). A typical BPH nodule has mixed signal intensity on T2-weighted images with variables DWI signal (Figures 2, 15 and 23). A so-called atypical BPH nodule has low signal intensity on T2-weighted images




Figure 21. *Teaching point* - Typical appearance of extra-capsular extension: IMPROD study number 45: Low signal intensity in T2-weighted images (T2W Axial) in the right peripheral zone was present with bulging and spiculated margin (yellow arrow) of prostate capsule suggestive of extra-capsular extension which was found on whole prostatectomy sections. Prostatectomy was performed with wide margins on the right and complete nerve sparing on the left resulting in negative surgical margins on the right side. All MRI images available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 45 are presented in the supporting material.




Figure 22. Teaching point - Typical appearance of seminal vesicle invasion: IMPROD study number 63: Irregular area of low signal on T2weighted images (yellow arrow) was extending to normal high T2-weighted imaging signal of seminal vesicle, consistent with seminal vesicle invasion as confirmed on prostatectomy specimens. Additional images and videos of IMPROD study number 63 are presented in the supporting material.

with signal restriction on DWI (Figure 23). A typical or atypical or whatever terminology readers want to use to describe BPH nodules can have variable signal characteristics and shape on T2-weighted image and DWI but a BPH nodule should be surrounded by a low T2-weighted signal intensity capsule which completely encapsulates the lesion (Figures 2, 15 and 23). Regardless of signal characteristics and shape on T2-weighted imaging and DWI, lesions which have complete capsule likely represent BPH nodules (Figures 2, 15 and 23) and there is no need to biopsy such lesions if other parameters are not suggestive of csPCa presence. These lesions can be reported as PIRADs version 2.1/IMPROD bpMRI Likert score 1-2 if fully encapsulated or score of three if capsule is incomplete. In contrast to BPH, PCa in central gland appears as irregular areas, could be lens shape or any other shape (Figures 4, 12 and 24) without surrounding capsule. However, PCa in transition zone also can be very homogeneous. For simplification, we will refrain from discussion on transition and central zone lesion and describe all lesions as central gland lesions.

Accuracy of IMPROD bpMRI for prostate cancer in men with a clinical suspicion of prostate

The accuracy of IMPROD bpMRI is not 100% and other clinical parameters and individualized risk assessments are

needed (Figure 25). However, in analysis of patients who underwent robotic-assisted laparoscopic prostatectomy (RALP) following biopsy in the IMPROD trial (40%, 64/ 161), the initial pre-biopsy IMPROD bpMRI scan detected 71 PCa lesions (71%, 91/99) with Gleason score of 3+4 or higher and only three lesions (7.5%, 3/40) with GS >3 + 4 were missed. All dominant lesions with GS >3 + 4 were correctly identified. Moreover, all men enrolled in the IM-PROD trial (ClinicalTrials.gov Identifier: NCT01864135) are undergoing close follow-up. With the current median followup for the IMPROD trial cohort of over 3 years, in 49 men with IMPROD bpMRI Likert scores of 1-3, only two men (4%, 2/49) were found to have low-volume PCa with a Gleason score of 3 + 4 during the follow-up period, suggesting that the initial IMPROD bpMRI scan had a high negative predictive value. Thus, initial negative IMPROD bpMRI in men with a clinical suspicion of PCa can serve as a tool to limit biopsy procedures. Same type of MRI acquisition and reconstruction protocol (Figure 26) should ideally be used in men who decide to undergo active surveillance since variation in prostate MRI protocol can negatively impact PCa detection and characterization. If imaging and pathology findings are discordant in a patient with persistent clinical suspicion for PCa, repeated imaging with repeated biopsy tend to be the preferred option (Figure 27). All imaging studies presented in this pictorial review were performed without intravenous contrast agent. Large number of studies has retrospectively evaluated prostate MRI datasets reported with and without the use of intra-venous contrast




Figure 23. *Teaching point* - Typical appearance of benign prostatic hyperplasia which can mimic prostate cancer: IMPROD study number 127: A lesion with restricted diffusion in the central gland (transitional and central zone) is present on the left side, base. The lesion has the following signal characteristics: low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500), high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm². The lesion is BPH since it has a complete surrounding capsule. All MRI images are available at http://petiv.utu.fi/improd/Additional images and videos of IMPROD study number 127 are presented in the supporting material.




Figure 24. *Teaching point* - Typical appearance of PCa in the central gland: IMPROD study number 92: A lesion with restricted diffusion located in the central gland, mid-gland level and extends towards the base, mid-line, in front of urethra without any surrounding capsule. The lesion has the following signal characteristics: low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500), high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm². The size of the lesion was CCxRLxAP:13x13x10 mm. The lesion was found to correspond to Gleason sore 4+3 prostate cancer on whole mount prostatectomy sections. All MRI images are available at http://petiv.utu.fi/improd/Additional images and videos of IMPROD study number 92 are presented in the supporting material.




Figure 25. Teaching point - Prostate cancer missed on MRI: IMPROD study number 130: No highly suspicious lesion was seen but two cores of 3+4 in systematic biopsies. Patient decided for Robot Assisted Laparoscopic Prostatectomy (RALP) with bilateral nerve sparing and was found to have less than 5% GS3+4/GGG2 in RIGHT apex. All MRI images are available at http://petiv.utu.fi/improd/Additional images and videos of IMPROD study number 130 are presented in supporting material.

agent. In a recent meta-analysis by Bass et al.,³⁰ authors concluded that prostate MRI reported without use of intravenous contrast agent perform similar in detecting prostate cancer to prostate MRI studies reported using intra-venous contrast agent.

Continuous quality control, database, case follow-up of IMPROD bpMRI for prostate

Continuous quality control and correlation with biopsy and prostatectomy findings, ideally using whole mount




Figure 26. *Teaching point* - Increase in the lesion size of peripheral zone lesion in the in follow up MRI: IMPROD study number 28: Focal area of low signal on T2-weighted images (T2W) with low signal intensity on apparent diffusion coefficient maps (ADCm5b500), high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm² was present in the left apex in a patient who was found to have Gleason score 3+3 in one cores of systematic biopsy. Follow up imaging performed 568 days following the initial IMPROD bpMRI scan demonstrated small increase in the lesion size. Patient opted to continue active surveillance. All images of the initial and follow up MRI scans are available at http://petiv.utu.fi/improd/ Additional images and videos of IMPROD study number 28 are presented in the supporting material.




Figure 27. *Teaching point* - Increase in lesion size of an anterior PZ/central gland lesion in the in follow up MRI: IMPROD study number 50: Focal areas of low signal on T2-weighted images (T2W Axial) signal with low signal on apparent diffusion coefficient maps (ADCm5b500) and high signal on DWI trace images of b 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm² in the central gland were initially reported as IMPROD bpMRI Likert score of 3. In the follow up IMPROD bpMRI performed 496 days following the initial scan, these lesions were upgraded to the score 5 and repeated targeted biopsy was performed which demonstrated Gleason score 3+3 PCa. All images of the initial and follow up MRI scans are available at http://petiv.utu.fi/improd/Additional images and videos of IMPROD study number 50 are presented in supporting material.

prostatectomy section, is crucial (Figure 28). Quality of prostate MRI is gaining increasing popularity in the literature with several groups proposing different approaches and systems to evaluate prostate MRI and MRI-guided biopsy pathway quality. Giganti et al.³¹ developed a system called the Prostate Imaging Quality (PI-QUAL) which relies on visual evaluation of prostate MRI to assigned quality score. We truly applaud these efforts but their impact will be limited if access to all datasets is not provided. Moreover, prostate MRI acquisition protocol impacts quality of prostate MRI, thus, providing free public access to prostate MRI acquisition is paramount for external reproducibility and repeatability of published studies. This pictorial review provides access to prostate MRI acquisition protocol for all three major MRI vendors, thus, enabling readers to acquire datasets of similar quality at their home institutions. Development of freely publicly available databases of prostate MRI and accoupling clinical, laboratory and histopathological findings in men with suspected, diagnosed or recurrent PCa could lead to bigger impact on clinical practice and acceptance of prostate MRI rather than meta-analysis of single or multi-centre studies which do not provide access to the source datasets. As of now a PRostate Mri Outcome Database (PROMOD) has been developed




Figure 28. Teaching point - Lesion in right PZ which leads to wide surgical margins on the right: IMPROD study number 21: One lesion on the right side in the midgland was detected which was found to be Gleason score 4+3 with negative surgical margins on the right side and complete nerve sparing on the LEFT-side. The lesion has the following signal characteristics: low signal intensity in T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500), high signal intensity on trace diffusion weighted imaging b-values of 1500 (DWI b1500) and 2000 (DWI b 2000) s/mm². All MRI images are available at http://petiv.utu.fi/ improd/ Additional images and videos of IMPROD study number 21 are presented in supporting material.

with the aim of to create a comprehensive database including multiple heterogeneous cohorts in order to explore inter-centre differences in the accuracy of MRI, and define optimal strategies for selection of men who might benefit from biopsies. We invite all centres to participate in PROMOD and encourage interested individuals/research groups/centres to reach out to us at promodstudy@gmail.com

Each individual institution is strongly encouraged to develop and maintain own database of prostate MRI and try to limit lost to follow up. Such database allows to follow interesting imaging findings (Figures 29 and 30) and crosscorrelates those with different imaging modalities (CT, PET/ CT) and/or histopathology.

Summary

In this pictorial review, we present imaging findings from men enrolled in the IMPROD trial which included 161 consecutive men who underwent IMPROD bpMRI followed by biopsy at the time of clinical suspicion of PCa. As one can see in this review, even such a small dataset can serve as a valuable tool for becoming familiar with prostate MRI and establishing a basis for reporting prostate MRI performed without an intravenous contrast agent in men with a clinical suspicion of PCa. The IMPROD and MULTI-IMPROD trial servers (http://petiv.utu.fi/improd/and http://petiv.utu.fi/ multiimprod/) provide access to complete imaging datasets of 161 and 338 men, respectively, which can be viewed using an online DICOM viewer. Thus, we urge readers to use these databases as a tool for learning prostate MRI reporting in men with a clinical suspicion of PCa. Our vision of a large array of different medical providers (from mid-level providers such as dedicated physician assistants all the way to specified radiologists) being able to report prostate MRI with high accuracy could be one way to reduce the costs of prostate MRI and enable wide adoption of prostate MRI in men with a clinical suspicion of PCa. In our own vision, the current pictorial review represents the first step towards this implementation.

Supporting material is provided at http://mrc.utu.fi/data

Declaration of conflicting interests

The author(s) declared no potential conflicts of interest with respect to the research, authorship, and/or publication of this article.

Funding

The author(s) received no financial support for the research, authorship, and/or publication of this article.

ORCID iDs

Ivan Jambor in https://orcid.org/0000-0001-6305-3570




Figure 29. *Teaching point* - If a lesion has marked T2 relaxation time shortening (i.e. 'very dark' on T2-weighted images), the lesion might *not* demonstrate relative hyperintensity on trace DWI: MULTI-IMPROD study number 176 (same case as in the Figure 30): One lesion (red arrow) was present in peripheral zone, close to midline. The lesion has the following signal characteristics: low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500, ADCm2b1500, ADCm2b2000), low signal intensity on trace diffusion weighted imaging b-values of 500 (DWI5b500 b500), 1500 (DWI2b1500 b1500) and 2000 (DWI2b2000 b2000) s/mm². On the corresponding whole mount proctectomy section, prostate cancer was found in the peripheral zone. All MRI images are available at http://petiv.utu.fi/multiimprod/.




Figure 30. *Teaching point:* Calcifications (blue arrow) on T2-weighted images and all trace DWI images will be completely dark (meaning hypointense). However, also markedly dark on ADC/ADCm images since calcifications have very low signal on all individual trace DWI images: MULTI-IMPROD study number 176 – The same patient as Figure 29: In addition to prostate cancer lesion in the peripheral zone (red arrow), close to midline, there are areas in the central gland (blue arrow) with the following signal characteristics: low signal intensity on T2-weighted images (T2W Axial), low signal intensity on apparent diffusion coefficient maps (ADCm5b500), low signal intensity on trace DWI images (DWI b0, DWI b100, DWI b250, DWI b500), 1500 (DWI2b1500 b1500) and 2000 (DWI b 2000) s/mm². All MRI images are available at http://petiv.utu.fi/multiimprod/.

References

- Epstein JI, Egevad L, Amin MB, et al. The 2014 international society of urological pathology (ISUP) consensus conference on gleason grading of prostatic carcinoma. Am J Surg Pathol 2016; 40: 244–252.
- Roehl KA, Antenor JA, Catalona WJ. Serial biopsy results in prostate cancer screening study. J Urology 2002; 167: 2435–2439.
- Barocas DA, Alvarez J, Resnick MJ, et al. Association between radiation therapy, surgery, or observation for localized prostate cancer and patient-reported outcomes after 3 years. Jama 2017; 317: 1126–1140.
- 4. Porpiglia F, Manfredi M, Mele F, et al. Diagnostic pathway with multiparametric magnetic resonance imaging versus standard pathway: results from a randomized prospective study in biopsy-naïve patients with suspected prostate cancer. Eur Urol 2017; 72: 282–288.
- Ahmed HU, El-Shater Bosaily A, Brown LC, et al. Diagnostic accuracy of multi-parametric MRI and TRUS biopsy in prostate cancer (PROMIS). A Paired Validating Confi Rmatory Study 2017; 389.
- Kasivisvanathan V, Rannikko AS, Borghi M, et al. MRItargeted or standard biopsy for prostate-cancer diagnosis. New Engl J Med 2018; 378: 1767–1777.
- Rouvière O, Puech P, Renard-Penna R, et al. Use of prostate systematic and targeted biopsy on the basis of multiparametric MRI in biopsy-naive patients (MRI-FIRST): a prospective, multicentre, paired diagnostic study. Lancet Oncol 2018; 0.
- Jambor I, Boström PJ, Taimen P, et al. Novel biparametric MRI and targeted biopsy improves risk stratification in men with a clinical suspicion of prostate cancer (IMPROD Trial). J Magn Reson Imaging 2017; 46: 1089–1095.

- Boesen L, Nørgaard N, Løgager V, et al. Assessment of the diagnostic accuracy of biparametric magnetic resonance imaging for prostate cancer in biopsy-naive men. JAMA Netw Open 2018; 1: e180219.
- Boesen L, Nørgaard N, Løgager V, et al. Prebiopsy Biparametric Magnetic Resonance Imaging Combined with Prostatespecific Antigen Density in Detecting and Ruling Out Gleason 7-10 Prostate Cancer in Biopsy-Naïve Men. Eur Urol Oncol. 2019; 2(3): 311–319.
- 11. Blinded for Anonymity.
- Grönberg H, Eklund M, Picker W, et al. Prostate Cancer Diagnostics Using a Combination of the Stockholm3 Blood Test and Multiparametric Magnetic Resonance Imaging. Eur Urol 2018; 74: 722–728.
- Cereser L, Giannarini G, Bonato F, et al. Comparison of multiple abbreviated multiparametric MRI-derived protocols for the detection of clinically significant prostate cancer. Minerva Urol Nefrol 2020.
- Jambor I, Merisaari H, Aronen HJ, et al. Optimization of bvalue distribution for biexponential diffusion-weighted MR imaging of normal prostate. J Magn Reson Imaging 2014; 39: 1213–1222.
- Merisaari H, Jambor I. Optimization of b-value distribution for four mathematical models of prostate cancer diffusionweighted imaging using b values up to 2000 s/mm2: Simulation and repeatability study. Magn Reson Med 2015; 73: 1954–1969.
- Jambor I, Merisaari H, Taimen P, et al. Evaluation of different mathematical models for diffusion-weighted imaging of normal prostate and prostate cancer using high b-values: A repeatability study. Magn Reson Med 2015; 73: 1988–1998.

- Toivonen J, Merisaari H, Pesola M, et al. Mathematical models for diffusion-weighted imaging of prostate cancer using b values up to 2000 s/mm2: Correlation with Gleason score and repeatability of region of interest analysis. Magn Reson Med 2015; 74: 1116–1124.
- Knaapila J, Jambor I, Ettala O, et al. Negative Predictive Value of Biparametric Prostate Magnetic Resonance Imaging in Excluding Significant Prostate Cancer: A Pooled Data Analysis Based on Clinical Data from Four Prospective, Registered Studies. Eur Urol Focus 2020: 1–10.
- Rosenkrantz AB, Kopec M, Kong X, et al. Prostate cancer vs. post-biopsy hemorrhage: diagnosis with T2- and diffusionweighted imaging. J Magn Reson Imaging 2010; 31: 1387–1394.
- Choi MH, Jung SE, Park YH, et al. Multiparametric MRI of prostate cancer after biopsy: little impact of hemorrhage on tumor staging. Invest Magn Reson Imaging 2017; 21: 3139. DOI: 10.13104/imri.2017.21.3.139.
- Sarradin M, Lepiney C, Celhay O, et al. Évaluation prospective du délai minimum de réalisation de l'IRM prostatique après une biopsie de prostate: facteurs prédictifs cliniques et anatomopathologiques de remaniements hémorragiques. Progrès en Urologie 2018; 28: 85–93.
- 22. Blinded for Anonymity.
- Bourne R, Panagiotaki E. Limitations and prospects for diffusion-weighted MRI of the prostate. Diagnostics 2016; 6: 21.

- Salvaggio G, Calamia M, Purpura P, et al. Role of apparent diffusion coefficient values in prostate diseases characterization on diffusion-weighted magnetic resonance imaging. Minerva urologica e nefrologica = Ital J Urol Nephrol 2019; 71: 154–160.
- 25. Blinded for Anonymity.
- 26. Blinded for Anonymity.
- 27. Blinded for Anonymity.
- Khoo CC, Eldred-Evans D, Jaenicke J, et al. Likert vs. PI-RADS v2: a comparison of two radiological scoring systems for detection of clinically significant prostate cancer. Eur Urol Supplements 2019; 18: e1865–e1866.
- Abreu-Gomez J, Wu M, McInnes MDF, et al. Shape analysis of peripheral zone observations on prostate DWI: correlation to histopathology outcomes after radical prostatectomy. Am J Roentgenology 2020; 214: 1239–1247.
- Bass EJ, Pantovic A, Connor M, et al. A systematic review and meta-analysis of the diagnostic accuracy of biparametric prostate MRI for prostate cancer in men at risk. Prostate Cancer Prostatic Dis 2020.
- Giganti F, Allen C, Emberton M, et al. Prostate imaging quality (PI-QUAL): a new quality control scoring system for multiparametric magnetic resonance imaging of the prostate from the PRECISION trial. Eur Urol Oncol 2020; 3: 615–619.