#### **CORRESPONDENCE**



# Detection of AD-specific four repeat tau with deamidated asparagine residue 279-specific fraction purified from 4R tau polyclonal antibody

Momoko Ebashi<sup>1,2</sup> · Shuta Toru<sup>3</sup> · Ayako Nakamura<sup>1,4</sup> · Satoshi Kamei<sup>2</sup> · Takanori Yokota<sup>4</sup> · Katsuiku Hirokawa<sup>5</sup> · Toshiki Uchihara<sup>1,3,4,6</sup>

Received: 13 March 2019 / Revised: 12 April 2019 / Accepted: 12 April 2019 / Published online: 20 April 2019 © The Author(s) 2019

Tau proteins are classified according to its isoforms, based on the number of tandem repeat (three or four repeat) of the microtubule binding domain [4]. Both three repeat (3R) and four repeat (4R) isoforms are found in Alzheimer disease (AD), while only 4R isoform is found in progressive supranuclear palsy (PSP)/corticobasal degeneration (CBD) and argyrophilic grains [3]. RD4 is a monoclonal antibody that has been raised against the RD4 epitope (VQII\*NKKLDLSNVQSKC=N-peptide) corresponding to the second microtubule binding domain [2], not found in 3R tau. This translated sequence is modified irreversibly by deamidation of asparagine at N279 to yield aspartic acid at D279 (VQII\*DKKLDLSNVQSKC=D-peptide), found in AD [5]. Because this deamidation is not found in PSP/CBD [1], specific probing of 4R tau carrying this deamidation (4R-D) will provide a molecular marker for AD [6], distinct from PSP/CBD (Fig. 1). Even though RD4 is highly specific to N-peptide without cross reaction to D-peptide [2], RD4 positivity is not specific to PSP/CBD brains (Fig. 2b) because non-deamidated 4R tau (4R-N) is also found in AD brains (Fig. 2a) [6, 7]. Curiously, polyclonal antibody raised against *D*-peptide (4R tau polyclonal antibody) immuolabels


☐ Toshiki Uchihara uchihara@nakanosogo.or.jp

- Laboratory of Structural Neuropathology, Tokyo Metropolitan Institute of Medical Science, Tokyo, Japan
- Department of Neurology, Nihon University School of Medicine, Tokyo, Japan
- <sup>3</sup> Department of Neurology, Nitobe Memorial Nakano General Hospital, Tokyo, Japan
- Department of Neurology and Neurological Science, Tokyo Medical and Dental University, Tokyo, Japan
- Department of Pathology, Nitobe Memorial Nakano General Hospital, Tokyo, Japan
- Neuromorphomics Laboratory, Nitobe Memorial Nakano General Hospital, Tokyo, Japan

tau lesions not only in AD brains containing 4R-D and 4R-N (Fig. 2c) but also in PSP brain containing only 4R-N (Fig. 2d) as reported [1]. This apparent cross reaction is possibly explained by the intrinsic nature of this 4R tau polyclonal antibody that may react more or less to both p-peptide and N-peptide. Interestingly, however, its absorption with N-peptide was complete on PSP brains (Fig. 2f), while this absorption was not evident on AD brains (Fig. 2e), suggesting that this 4R tau polyclonal antibody contains subfractions selectively reactive to 4R-N or p-peptide.

Indeed, immunoabsorption of this 4R tau polyclonal antibody with N-peptide, followed by affinity purification with D-peptide column (Fig. 3) yielded a D-peptide specific fraction as confirmed with ELISA (Fig. 4, orange: *p*-peptide; blue: N-peptide). As expected, this D-peptide specific fraction labeled tau lesions in AD brains (Fig. 2g), but failed to label tau lesions in PSP (Fig. 2h) and CBD brains (data not shown). This contrast in tau immunolabeling, positive with AD brains but negative with PSP brains, is reminiscent of those obtained with RD3 (Fig. 2i, j), while their molecular bases are completely different. This D-peptide specific fraction of 4R tau polyclonal antibody clearly confirmed the biochemical contrast that the posttranslational modification at N279D of 4R tau of AD [5] is absent in PSP. Posttranslational deamidation of asparagine into aspartic acid (N279D) may be useful as a biomarker that may distinguish 4R tau of AD from other 4R tau of PSP/CBD. Probing this AD-specific change in 4R tau is a promising way to distinguish AD and PSP/CBD on a molecular basis to be applied for histological, biochemical and biofluid diagnosis (Patent Pending under International Publication number: WO2018/212261).








Fig. 2 Difference in tau immunoreactivity between AD and PSP brains (bar 100 μm)

#### Methods

### Peptide synthesis

*N*-peptide corresponding the RD4 epitope (VQII\*NKKLDLSNVQSKC) and *D*-peptide, deamidated at N279D (VQII\*DKKLDLSNVQSKC), were custom synthesized by Eurofin (Tokyo, Japan, purity > 99%).

## **Immunohistochemistry**

Formalin-fixed, paraffin-embedded blocks were obtained from four patients of PSP (69–77 yo), a patient of CBD (68 yo), and 4 patients of AD (78–98 yo). Permission of autopsy and of subsequent use for research was obtained from the next to kin of the patients and this study was approved by ethics committee of Tokyo Metropolitan Institute of Medical Science (16-9).

For RD4, serial pretreatment was performed with 0.25% potassium permanganate; KMnO4, for 15 min, 2% oxalic acid for 3 min, > 99% formic acid for 30 min and autoclaved in 0.05 M citrate buffer for 20 min at 120 °C. For 4R tau


polyclonal antibody (TIP-4RT-P01, Cosmo bio, Tokyo, Japan) [1], autoclaved in 0.05 M citrate buffer for 10 min at 120 °C and > 99% formic acid for 10 min [8]. Specificity of immunolabeling was confirmed by its disappearance in the presence of the immunogen peptides.

#### **Affinity purification**

After washing formyl cellufine in the column, 200 µg of 4R-D peptide was mixed with 0.6 ml of formyl cellufine in 0.6 ml of coupling buffer (50 mM Na<sub>2</sub>CO<sub>3</sub>·NaHCO<sub>3</sub>, pH 8.5) at room temperature on rotating shaker. After adding 3 mg of sodium cyanoborohydride (NaBH3CN), the formyl cellufine was further shaked at 4 °C for 8–12 h, washed with a blocking buffer (0.1 M monoethanolamine; MEA + 50 mM Tris–HCl buffer, pH 8.0) and mixed with 3 mg of NaB-H3CN. This 4R-D peptide primed cellufine was washed with the elution buffer (0.1 M glycine–HCl, pH 2.5) and then with washing buffer (1 M NaCl 1% Triton-X100, pH 7.5) and stored at 4 °C until use. Because 4R tau immunoreactivity on AD brain with the anti-4R tau antibody (1:30,000) was absorbed at the peptide concentration around  $1.0 \times 10^{-7}$  µmol/µl, 50 µl of the anti-4R tau antibody was diluted × 8



**Fig. 3** Affinity purification of 4R-D specific fraction from the 4R tau polyclonal antibody






Fig. 4 ELISA of the 4R-D specific fraction to 4R-D peptide (orange) and 4R-N peptide (blue)

with 10% bovine serum albumin and mixed with 289.5  $\mu$ g of 4R-N peptide, equivalent to  $1.0 \times 10^{-7} \, \mu$ mol/ $\mu$ l at its final concentration. This mixture was reacted with the cellufine primed with 4R-D peptide at 4 °C for 8 h on a rotating shaker. After washing, the bound antibody was eluted with the elution buffer to obtain 4R-D fraction.

# ELISA

The peptide (N- or D-) was serially diluted to  $1.0 \times 10^{-12\text{--}4}$  µmol/µl in 50 mM Tris—HCl, pH 8.8 and coated onto microtiter plates (SUMILON) at 4 °C for 16 h. The plates were blocked with 10% fetal bovine serum for 2 h. They were reacted with the purified 4R-D fraction for 90 min at RT followed by incubation with HRP-anti-rabbit IgG made in goat (Bio-Rad) at 1:1000 dilution, and reacted with the substrate, 0.4 mg/ml o-phenylenediamine, in citrate buffer (24 mM citric acid, 51 mM Na<sub>2</sub>HPO<sub>4</sub>), The absorbance at 490 nm was measured using Plate Chameleon (HIDEX) as described.

Acknowledgements This work was supported by Grant-in-Aid for JSPS KAKENHI Grant numbers 17H03555 (TU), 16K14572 (TU).

**Open Access** This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.



# References

- Dan A, Takahashi M, Masuda-Suzukake M, Kametani F, Nonaka T, Kondo H et al (2013) Extensive deamidation at asparagine residue 279 accounts for weak immunoreactivity of tau with RD4 antibody in Alzheimer's disease brain. Acta Neuropathol Comm 1:54. https://doi.org/10.1186/2051-5960-1-54
- de Silva R, Lashley T, Gibb G, Hanger D, Hope A, Reid A et al (2003) Pathological inclusion bodies in tauopathies contain distinct complements of tau with three or four microtubule-binding repeat domains as demonstrated by new specific monoclonal antibodies. Neuropathol Appl Neurobiol 29(3):288–302
- Flament S, Delacourte A, Verny M, Hauw JJ, Javoy-Agid F (1991)
   Abnormal tau proteins in progressive supranuclear palsy. Similarities and differences with the neurofibrillary degeneration of the Alzheimer type. Acta Neuropathol 81(6):591–596
- Goedert M, Spillantini MG, Potier MC, Ulrich J, Crowther RA (1989) Cloning and sequencing of the cDNA encoding an isoform of microtubule-associated protein tau containing four tandem repeats: differential expression of tau protein mRNAs in human brain. EMBO J 8(2):393–399
- Hasegawa M, Morishima-Kawashima M, Takio K, Suzuki M, Titani K, Ihara Y (1992) Protein sequence and mass spectrometric

- analyses of tau in the Alzheimer's disease brain. J Biol Chem 267(24):17047-17054
- Hasegawa M, Watanabe S, Kondo H, Akiyama H, Mann DM, Saito Y et al (2014) 3R and 4R tau isoforms in paired helical filaments in Alzheimer's disease. Acta Neuropathol 127(2):303–305. https://doi.org/10.1007/s00401-013-1191-9
- Uchihara T, Hara M, Nakamura A, Hirokawa K (2012) Tangle evolution linked to differential 3- and 4-repeat tau isoform deposition: a double immunofluorolabeling study using two monoclonal antibodies. Histochem Cell Biol 137(2):261–267. https://doi. org/10.1007/s00418-011-0891-2
- 8. Uchihara T, Nakamura A, Shibuya K, Yagishita S (2011) Specific detection of pathological three-repeat tau after pretreatment with potassium permanganate and oxalic acid in PSP/CBD brains. Brain Pathol 21(2):180–188

**Publisher's Note** Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

