



## Complete Genome Sequence of the Diatom-Associated Bacterium *Sphingorhabdus* sp. Strain SMR4y

<sup>®</sup> Mats Töpel,<sup>a,b</sup> Matthew I. M. Pinder,<sup>a</sup> Oskar N. Johansson,<sup>c</sup> Olga Kourtchenko,<sup>a</sup> Anna Godhe,<sup>a</sup> Adrian K. Clarke<sup>c</sup>

<sup>a</sup>Department of Marine Sciences, University of Gothenburg, Gothenburg, Sweden <sup>b</sup>Gothenburg Global Biodiversity Centre, Gothenburg, Sweden <sup>c</sup>Department of Biological and Environmental Sciences, University of Gothenburg, Gothenburg, Sweden

**ABSTRACT** The bacterial strain SMR4y belongs to the diverse microbiome of the marine diatom *Skeletonema marinoi* strain R05AC. After assembly of its genome, presented here, and subsequent analyses, we placed it in the genus *Sphingorhabdus*. This strain has a 3,479,724-bp circular chromosome (with 3,340 coding sequences) and no known plasmids.

n an effort to identify microbiome constituents of the chain-forming diatom *Skeletonema marinoi*, various associated bacteria are being isolated, sequenced, and characterized (1–5). Here, we report the genome sequence of strain SMR4y, isolated from a culture of *S. marinoi* strain R05AC. Strain SMR4y has been found in two separate isolates of *S. marinoi* obtained at different times and geographical locations, has been in coculture with *S. marinoi* for almost a decade following its field isolation, and has the ability to stimulate growth of *S. marinoi* when added in excess (O. N. Johansson, M. I. M. Pinder, F. Ohlsson, J. Egardt, M. Töpel, and A. K. Clarke, unpublished data). For these reasons, we believe strain SMR4y to be a member of the *S. marinoi* microbiome, as opposed to a random cocultured bacterium.

The diatom culture was established from a revived resting cell taken from top-layer sediment at 14-m depth in May 2010 from Öresund, Sweden (55°59.16'N, 12°44.02'E); collection was performed using a box corer. The diatom culture was maintained in f/2-Si medium and grown at 16°C. Isolation of strain SMR4y was performed at the same temperature in darkness by spotting and iterative dilution streaking on marine agar plates; strain SMR4y was sampled from one of the resultant colonies. The colony was grown overnight at 30°C in 50 ml Marine broth 2216 (Difco BD, USA) prior to DNA extraction; plant DNAzol reagent (Invitrogen Life Technologies, USA) was used to extract genomic DNA as per the manufacturer's instructions. Sequencing was performed on the PacBio RS II platform (Pacific Biosciences, Menlo Park, CA, USA) using one single-molecule real-time (SMRT) cell, producing 163,482 unfiltered reads (950.9 Mbp total). Filtering was performed using the SMRT Portal version 2.3.0 P\_Filter module to remove reads with quality lower than 0.80 and/or reads shorter than 100 bp and subreads shorter than 500 bp (6) (unless otherwise specified, all software was run using default settings). This resulted in 68,619 reads (413.7 Mbp total) that were used in the assembly. Genome assembly was performed using Canu version 1.3 (7) (genome size parameter, 4.5 m). The first 7,058 bp of the assembled contig, found using a BLASTn search (8) to be highly similar to the end of the contig, were trimmed to enable contig circularization. This was confirmed by joining the ends and realigning the reads using the SMRT Portal version 2.3.0 RS\_Resequencing.1 protocol (Pacific Biosciences) (6); this protocol includes a correction step using the Quiver algorithm (6). The result was a single chromosome of 3,479,724 bp, with a G+C content of 58.0% and an average read coverage of  $100.73 \times$  (results summarized in Table 1). This is currently the only S.

Citation Töpel M, Pinder MIM, Johansson ON, Kourtchenko O, Godhe A, Clarke AK. 2019. Complete genome sequence of the diatomassociated bacterium *Sphingorhabdus* sp. strain SMR4y. Microbiol Resour Announc 8:e00482-19. https://doi.org/10.1128/MRA.00482-19.

**Editor** Frank J. Stewart, Georgia Institute of Technology

**Copyright** © 2019 Töpel et al. This is an openaccess article distributed under the terms of the Creative Commons Attribution 4.0

Address correspondence to Mats Töpel,

mats.topel@marine.gu.se. M.T. and M.I.M.P. contributed equally to this

work.

Received 26 April 2019 Accepted 18 June 2019 Published 18 July 2019

TABLE 1 Summary of the Sphingorhabdus sp. strain SMR4y assembly and annotation

Characteristic	Value
Assembly	
No. of reads (filtered)	68,619
Bases (filtered) (bp)	413,668,254
Overlapping bases trimmed from start of contig (bp)	7,058
Final circular contig size (bp)	3,479,724
G+C content (%)	58.0
Avg read coverage ( $\times$ )	100.73
Annotation (no.)	
CDSs	3,340
Pseudogenes	4
tRNAs	46
rRNAs	6
ncRNAs	4
tmRNAs	1

*marinoi* strain R05AC-associated bacterium not found to contain plasmids, meaning that they are either absent from this bacterium or are smaller than 10,000 bp and therefore were removed during the sequencing library preparation step.

Prokka version 1.12beta (9) was used to annotate the strain SMR4y genome, inferring 3,340 coding sequences (CDSs) (2,742 with a functional prediction), 4 pseudogenes, 46 tRNAs, 6 rRNAs, 4 noncoding RNAs (ncRNAs), and 1 transfer-messenger RNA (tmRNA) (results summarized in Table 1). Two identical 16S rRNA sequences are present on the chromosome; these have 99.9% and 99.4% identity, respectively, to partial sequences from *Sphingorhabdus flavimaris* strains R-36742 and SW-151<sup>T</sup> (GenBank accession numbers FR691421 [1,441 bp] and NR\_025814 [1,444 bp], respectively). The SMR4y sequence was trimmed to the same length for this comparison, with 99.7% identity to the two identical full-length 16S sequences of *Sphingorhabdus* sp. strain YGSMI21 (GenBank accession number NZ\_CP022548) and 99.5% identity to the two identical full-length 16S rRNA sequences of *Sphingorhabdus* sp. strain M41 (GenBank accession number NZ\_CP014545).

All whole-genome-sequenced species in the order *Sphingomonadales* available on the NCBI's RefSeq ftp site (ftp://ftp.ncbi.nlm.nih.gov/genomes/refseq/bacteria/; accessed 4 July 2018) were then included with strain SMR4y in a phylotaxonomic analysis using PhyloPhIAn version 0.99 (10). This showed strain SMR4y to be a sister to *Sphingorhabdus* sp. strain YGSMI21; these two strains, along with *Sphingorhabdus* sp. strain M41, form a clade with 100% bootstrap support. Considering the results of the above-described analyses, we placed strain SMR4y in the genus *Sphingorhabdus*.

**Data availability.** This whole-genome project has been deposited in GenBank under the accession number CP022336 as part of BioProject number PRJNA380207. The raw reads have been deposited in the NCBI Sequence Read Archive under the accession number SRR5811295.

## ACKNOWLEDGMENTS

This work was supported by the Gordon and Betty Moore Foundation (to A.K.C., M.T., and A.G., grant 4967), the Swedish Research Council VR (to A.K.C., grant 2015-04286), and the Swedish Research Council Formas (to A.G., grant 219-2012-2070; to M.T. and A.G., grant 2017-00466).

We thank the Linnéus Center for Marine Evolutionary Biology (CeMEB, http://cemeb .science.gu.se/) for support. All bioinformatics analyses were run on the Albiorix computer cluster (http://albiorix.bioenv.gu.se/) at the Department of Marine Sciences, University of Gothenburg.

## REFERENCES

1. Töpel M, Pinder MIM, Johansson O, Kourtchenko O, Godhe A, Clarke AK. 2017. Genome sequence of *Roseovarius mucosus* strain SMR3, isolated from a culture of the diatom *Skeletonema marinoi*. Genome Announc 5:e00394-17. https://doi.org/10.1128/genomeA.00394-17.

- Töpel M, Pinder MIM, Johansson O, Kourtchenko O, Godhe A, Clarke AK. 2018. Complete genome sequence of *Loktanella vestfoldensis* strain SMR4r, a novel strain isolated from a culture of the chain-forming diatom *Skeletonema marinoi*. Genome Announc 6:e01558-17. https://doi .org/10.1128/genomeA.01558-17.
- Töpel M, Pinder MIM, Johansson O, Kourtchenko O, Godhe A, Clarke AK. 2018. Whole-genome sequence of the novel *Antarctobacter heliothermus* strain SMS3, found in association with the marine diatom *Skeletonema marinoi*. J Genomics 6:113–116. https://doi.org/10.7150/jgen.27637.
- Töpel M, Pinder MIM, Johansson O, Kourtchenko O, Godhe A, Clarke AK. 2019. Genome sequence of *Arenibacter algicola* strain SMS7, found in association with the marine diatom *Skeletonema marinoi*. Microbiol Resour Announc 8:e01461-18. https://doi.org/10.1128/MRA.01461-18.
- Töpel M, Pinder MIM, Johansson O, Kourtchenko O, Clarke AK, Godhe A. 2019. Complete genome sequence of novel *Sulfitobacter pseudonitzschiae* strain SMR1, isolated from a culture of the marine diatom *Skeletonema marinoi*. J Genomics 7:7–10. https://doi.org/10.7150/ jgen.30559.
- Chin C-S, Alexander DH, Marks P, Klammer AA, Drake J, Heiner C, Clum A, Copeland A, Huddleston J, Eichler EE, Turner SW, Korlach J. 2013. Nonhybrid, finished microbial genome assemblies from long-read SMRT sequencing data. Nat Methods 10:563–569. https://doi.org/10.1038/ nmeth.2474.
- Koren S, Walenz BP, Berlin K, Miller JR, Bergman NH, Phillippy AM. 2017. Canu: scalable and accurate long-read assembly via adaptive k-mer weighting and repeat separation. Genome Res 27:722–736. https://doi .org/10.1101/gr.215087.116.
- Altschul SF, Gish W, Miller W, Myers EW, Lipman DJ. 1990. Basic local alignment search tool. J Mol Biol 215:403–410. https://doi.org/10.1016/ S0022-2836(05)80360-2.
- Seemann T. 2014. Prokka: rapid prokaryotic genome annotation. Bioinformatics 30:2068–2069. https://doi.org/10.1093/bioinformatics/btu153.
- Segata N, Börnigen D, Morgan XC, Huttenhower C. 2013. PhyloPhlAn is a new method for improved phylogenetic and taxonomic placement of microbes. Nat Commun 4:2304. https://doi.org/10.1038/ncomms3304.